

<b>Document type</b>	Development Document	<b>Confidentiality</b>	Confidential
<b>Document No.</b>	KF-0651-4-202		
<b>Scope</b>	0651		
<h2>Summary of safety and clinical performance</h2>			
<b>Related Document(s)</b>			
<b>Document No.</b>	<b>Document Name</b>	<b>Version</b>	
KF-0651-4-0169	Clinical Evaluation Report	5.0	
KF-0651-0020	PMCF Report	1.0	
KF-0651-7-0314-04	Periodic Safety Update Report	1.0	
KF-0651-3-001	Risk Management Report	29.0	



## Table of Contents

Summary of safety and clinical performance.....	1
Table of Contents .....	3
1. Identification the device and the manufacturer .....	4
2. Intended use .....	4
2.1 Intended purpose .....	4
2.2 Indications .....	4
2.3 Contraindications .....	5
3. Description of the device .....	5
3.1 Description of the device.....	5
3.2 Overview of the previous generations device .....	7
3.3 Accessories.....	7
3.4 Another device .....	18
4. Residual risks and undesirable effects, warnings and precautions .....	18
4.1 Residual risks and undesirable effects.....	18
4.2 Warnings and precautions .....	20
4.3 Post-Market surveillance (PMS) database research.....	22
4.4 Complaint rates of the specific warning and precautions.....	37
5. The summary of clinical evaluation and post-market clinical follow-up (PMCF).....	40
5.1 Summary of clinical data related to equivalent device (if applicable) .....	41
5.2 Summary of clinical data from conducted investigations of the device before the CE-marking (if applicable).....	42
5.3 Summary of clinical data from other sources(if applicable) .....	65
5.4 An overall summary of the clinical performance and safety.....	77
5.5 Ongoing or planned post-market clinical follow-up .....	82
6. Possible diagnostic or therapeutic alternatives.....	82
7. Suggested profile and training for users.....	92
8. Harmonized standards and CS applied.....	101
9. Revision history .....	102

## 1. Identification the device and the manufacturer

Device details	
Name of device	Defibrillator Monitor
Model and type	BeneHeart D5/D6
Basic UDI-DI	69449040AB010000083E
Nomenclature	NA
Classification	III (According to Rule 22 of MDR Annex VIII)
CND code	Z120305
Manufacturer details	
Manufacturer name	Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Adress	Mindray Building, Keji 12th Road South, High-tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China
SRN	CN-MF-000014156
Authorised Representative	
EC-Representative	Shanghai International Holding Corp. GmbH (Europe) Eiffestraße 80, 20537 Hamburg, Germany
SRN of the EC-Representative	DE-AR-000000001
Notified Body	BSI Group The Netherlands B.V.
NB's single identification number	CE 2797
First CE certification date	<b>2008-12-15</b>

## 2. Intended use

### 2.1 Intended purpose

The equipment is intended for external defibrillation, internal defibrillation, synchronized cardioversion and semi-automated external defibrillation (AED). It can also be used for non-invasive external pacing as well as ECG, Resp, SpO2, PR, NIBP, IBP, Temp and CO2 monitoring.

### 2.2 Indications

■ **External defibrillation/AED/internal defibrillation:**

The external defibrillation, AED and internal defibrillation modes are intended for patients with ventricular fibrillation, pulseless ventricular tachycardia and ventricular flutter.

■ **Synchronized cardioversion:**

Synchronized cardioversion is intended for the treatment of atrial fibrillation and atrial flutter.

■ **Non-invasive external pacing:**

Non-invasive external pacing is intended for the treatment of bradycardia and asystole.

■ **Monitoring:**

Monitoring is intended for the monitoring of ECG, Resp, SpO<sub>2</sub>, PR, NIBP, IBP, Temp and CO<sub>2</sub> parameters.

## 2.3 Contraindications

■ **AED**

The AED mode is contraindicated in the treatment when the patient is showing any of the following:

- Consciousness
- Breathing
- Detectable pulse or other signs of circulation

■ **Manual Defibrillation**

Manual defibrillation is contraindicated in the treatment when the patient is showing any of the following:

- Consciousness
- Breathing.
- Detectable pulse or other signs of circulation.

## 3. Description of the device

### 3.1 Description of the device

#### **External Defibrillation**

The manual defibrillation feature uses exponential-truncated biphasic defibrillation technology and performs automatic compensation according to impedance of a patient. The biphasic defibrillation technology has been proved to be advantageous over the monophasic one and extensively applied in the defibrillation industry.

#### **AED**

The AED function uses the exponential-truncated biphasic defibrillation technology and the algorithm of a heart rhythm recognition detector to guide operators whether to perform defibrillation treatment.

#### **CPR Feedback**

The cardio-pulmonary resuscitation (CPR) sensor analyzes and calculates the compression rate by using the pressure sensor to measure the change of the compress force; The displacement waveform is obtained by quadratic integration of the acceleration signal acquired by the acceleration sensor, recognized, and analyzed and calculated to obtain the compress depth and recovery data.

The compression rate can also be obtained by calculating the change of the chest impedance, which comes from multifunction electrode pads.

### **Noninvasive Pacing**

The external pacing function imposes the pulse of fixed width on patients at a certain current intensity and frequency. This technology is universal in external pacing currently and is widely applied in defibrillator monitors. The circuit composed of the MCU and operational amplifier is a constant-current source circuit. The MCU, based on the DAC level output and frequency, can control the magnitude and frequency of current flowing through a patient's thoracic impedance.

### **ECG**

An electrocardiogram (ECG) is a variation curve recording the electrical activity of the heart over a cardiac cycle using electrodes placed over the skin, and acquiring the heart's electrical signal to analyze and calculate heart rate (HR) and arrhythmia (ARR). 12 lead ECG interpretation is a widely used technology to diagnose different heart conditions like arrhythmia, structural heart disease, myocardial infarction etc. by analyzing 12 lead ECG waveform.

### **NIBP**

NIBP measures the peripheral arterial blood pressure based on the principle of oscillation.

### **SpO2**

Pulse blood oxygen monitoring uses the spectrophotometry to monitor blood oxygen (measuring based on the different wavelengths of light absorbed by tissues).

### **Temp**

The TEMP module measures body temperature based on the temperature resistance characteristics of the thermistor.

### **CO2**

CO<sub>2</sub> gas has a strong absorption peak around the 4.26μm infrared light. Based on non-dispersive infrared spectroscopy (NDIR), the absorption volume of signals in the band is measured using the Lambert-Beer law to determine the CO<sub>2</sub> concentration during human respiration. Where, I<sub>0</sub> and I are the infrared light intensity before and after absorption; α is the absorption coefficient of the measured gas at the wavelength; L is the effective absorption optical length of the measured gas; and C is the concentration of the measured gas. Make that  $ABS = \ln(I_0) - \ln(I) = \alpha LC$ . Then, ABS is the CO<sub>2</sub> signal absorption volume of the infrared light in the 4.26μm band when passing through the CO<sub>2</sub> gas at a concentration of C. By calculating the absorption volume, the CO<sub>2</sub> gas concentration can be calculated.

### **IBP**

Invasive blood pressure (IBP) measurement directly measures intravascular pressure by means of fluid coupling.

### **Resp**

The RESP module measures the respiration rate using the chest impedance pneumography method.

### **CPR Filter**

CPR Filter is a technology that filters the CPR artifact from the ECG to allow users to see a close approximation of a patient's underlying ECG rhythm during CPR compressions, reducing interruptions in CPR.

## 3.2 Overview of the previous generations device

NA. BeneHeart D5/ BeneHeart D6 is the first defibrillator of Mindray.

## 3.3 Accessories

### 3.3.1 ECG Accessories

#### 3.3.1.1 ECG Electrodes

Model	Specification	Applicable patient	Manufacturer
31499224	10 pcs/pack	Adult	Covidien
2245	50 pcs/pack	Pediatric	3M
2258-3	3 pcs/pack	Neonate	3M

#### 3.3.1.2 12-pin Trunk Cable

Leadwire supported	Model	Compatible with	Type	Applicable patient	Manufacturer
3-lead	EV 6202	AHA, IEC	Defibrillation-proof	Pediatric, neonate	Shenzhen Mindray Bio-Medical Electronics Co., Ltd
3/5-lead	EV 6201	AHA, IEC	Defibrillation-proof	Adult, pediatric	
12-leadwire	EV 6203	AHA	Defibrillation-proof	Adult, pediatric	
12-leadwire	EV 6204	IEC	Defibrillation-proof		

#### 3.3.1.3 Lead Sets

3-Electrode Lead Sets	Manufacturer

Type	Compatible with	Model	Applicable patient	Remark	
Clip	IEC	EL6302A	Adult, pediatric	/	Shenzhen Mindray Bio-Medical Electronics Co., Ltd
		EL6304A		Long	
		EL6306A	Neonate	/	
		EL6308A	Pediatric	/	
	AHA	EL6301A	Adult, pediatric	/	
		EL6303A		Long	
		EL6305A	Neonate	/	
		EL6307A	Pediatric	/	
Snap	IEC	EL6302B	Adult, pediatric	/	
		EL6308B	Pediatric	/	
	AHA	EL6301B	Adult, pediatric	/	
		EL6307B	Pediatric	/	

5-Electrode Lead Sets					Manufacturer
Type	Compatible with	Model	Applicable patient	Remark	
Clip	IEC	EL6502A	Adult, pediatric	/	Shenzhen Mindray Bio-Medical Electronics Co., Ltd
	IEC	EL6504A		Long	
	AHA	EL6501A		/	
	AHA	EL6503A		Long	
Snap	IEC	EL6502B	/		
	AHA	EL6501B	/		



					Electronics Co., Ltd
--	--	--	--	--	-------------------------

10-Electrode Lead Sets					Manufacturer
Type	Compatible with	Model	Applicable patient	Remark	
Clip	IEC	EL6802A	Adult, pediatric	Limb	Shenzhen Mindray Bio- Medical Electronics Co., Ltd
		EL6804A		Chest	
	AHA	EL6801A		Limb	
		EL6803A		Chest	
Snap	IEC	EL6802B	Adult, pediatric	Limb	Shenzhen Mindray Bio- Medical Electronics Co., Ltd
		EL6804B		Chest	
	AHA	EL6801B		Limb	
		EL6803B		Chest	

### 3.3.1.4 Adapting Cable

Description	Compatible with	Type	Manufacturer
12-pin to 6 pin connectors	AHA, IEC	Adult, pediatric, neonate	Shenzhen Mindray Bio- Medical Electronics Co., Ltd

## 3.3.2 SpO2 Accessories

### 3.3.2.1 Extension Cables

Module type	Applicable patient	Remark	Manufacturer
Mindray SpO <sub>2</sub> module	Adult, pediatric, neonate	/	Shenzhen Mindray Bio- Medical Electronics Co., Ltd
Masimo SpO <sub>2</sub> module		8 pins, purple connector	Shenzhen Mindray Bio- Medical Electronics Co., Ltd
		7 pins, white connector	Shenzhen Mindray Bio- Medical Electronics

		Co., Ltd
Nellcor SpO <sub>2</sub> module		Shenzhen Mindray Bio-Medical Electronics Co., Ltd

### 3.3.2.2 SpO<sub>2</sub> Sensors

Mindray SpO <sub>2</sub> module			Manufacturer
Type	Model	Applicable patient	
Disposable	518C	Neonate (518C SpO <sub>2</sub> sensor wrap)	Shenzhen Mindray Bio-Medical Electronics Co., Ltd
	520A	Adult	
Single patient use	520P	Pediatric	Shenzhen Mindray Bio-Medical Electronics Co., Ltd
	520I	Infant	
	520N	Neonate	
Reusable	DS-100A	Adult	Nellcor Puritan Bennett Inc.
	OXI-P/I	Pediatric, infant	
	OXI-A/N	Adult, neonate	
	518B	Adult, pediatric, neonate (Multi-sites)	Shenzhen Mindray Bio-Medical Electronics Co., Ltd
	512E	Adult (Finger type)	
	512F		
	512G	Pediatric (Finger type)	
	512H		
	518C	Neonate	

Type	Model	Applicable patient	Remark	Manufacturer
Disposable	FPS-1901	Pediatric, neonate (wrap type)	LNCS-NeoPt-L	Masimo
	FPS-1862	Neonate (wrap type)	LNCS-Neo-L	
	FPS-1861	Infant (wrap type)	LNCS-Inf-L	

	FPS-1860	Pediatric (wrap type)	LNCS-Pdt	
	FPS-1859	Adult (wrap type)	LNCS-Adt	
Reusable	FPS-1863	Adult (finger clip)	LNCS DC-I	Masimo
	FPS-1864	Pediatric (finger clip)	LNCS-DCIP	
	2258	Adult, pediatric, neonate	LNCS YI	

Nellcor SpO <sub>2</sub> Module			Manufacturer
Type	Model	Applicable patient	
Disposable	MAX-A	Adult (>30 kg)	Covidien
	MAX-P	Pediatric (10 to 50 kg)	
	MAX-I	Infant (3 to 20 kg)	
	MAX-N	Neonate (<3 kg), adult (>40 kg)	
Reusable	DS-100A	Adult	Nellcor Puritan Bennett Inc.
	OXI-P/I	Pediatric, infant	
	OXI-A/N	Adult, neonate	

### 3.3.3 NIBP Accessories

#### 3.3.3.1 NIBP Hoses

Type	Applicable patient	Manufacturer
Reusable	Adult, pediatric	Shenzhen Mindray Bio-Medical Electronics Co., Ltd
	Neonate	

#### 3.3.3.2 Cuffs

Type	Model	Applicable patient	Applied site	Limb Circumference (cm)	Bladder Width (cm)	Manufacturer
Reusable	CM1201	Infant	Upper arm	10 to 19	9.2	Shenzhen

	CM1202	Pediatric		18 to 26	12.2	Mindray Bio-Medical Electronics Co., Ltd
	CM1203	Adult		24 to 35	15.1	
	CM1204	Large adult		33 to 47	18.3	
	CM1205	Adult	Thigh	46 to 66	22.5	
	CM1301	Infant	Upper arm	10 to 19	7.2	
	CM1302	Pediatric		18 to 26	9.8	
	CM1303	Adult		25 to 35	13.1	
	CM1304	Large adult		33 to 47	16.5	
	CM1305	Adult	Thigh	46 to 66	20.5	
Single patient	CM1500A	Neonate	Upper arm	3.1 to 5.7	2.2	Shenzhen Mindray Bio-Medical Electronics Co., Ltd
	CM1500B			4.3 to 8	2.9	
	CM1500C			5.8 to 10.9	3.8	
	CM1500D			7.1 to 13.1	4.8	
	CM1501	Infant		10 to 19	7.2	
	CM1502	Pediatric		18 to 26	9.8	
	CM1503	Adult		25 to 35	13.1	
	CM1504	Large adult		33 to 47	16.5	
	CM1505	Adult	Thigh	46 to 66	20.5	

### 3.3.4 Temp Accessories

#### 3.3.4.1 Extension Cable

Type	Model	Applicable Temp probe	Manufacturer
Reusable	MR420B	MR411, MR412	Shenzhen Mindray Bio-Medical Electronics Co., Ltd

#### 3.3.4.2 Temp Probes

Type	Model	Applicable patient	Application site	Manufacturer
Reusable	MR401B	Adult	Esophageal/Rectal	Shenzhen Mindray Bio-Medical Electronics Co., Ltd
	MR403B		Skin	
	MR402B	Pediatric, neonate	Esophageal/Rectal	

	MR404B		Skin	
Disposable	MR411	Adult, pediatric, neonate	Esophageal/Rectal	Shenzhen Mindray Bio-Medical Electronics Co., Ltd
	MR412		Skin	

### 3.3.4.3 Adapting Cable

Description	Applicable patient	Manufacturer
Temp adapting cable	Adult, pediatric, neonate	Shenzhen Mindray Bio-Medical Electronics Co., Ltd

### 3.3.5 IBP Accessories

IBP accessory kit	Description	Manufacturer
6800-30-50876 (Hospira)	12-pin IBP cable set	Shenzhen Mindray Bio-Medical Electronics Co., Ltd
	Disposable IBP transducer	ICU
	IBP transducer holder	HOSPIRA
	Steady Rest for IBP Transducer and Clamp	HOSPIRA
6800-30-50877 (BD)	12-pin IBP cable set	Shenzhen Mindray Bio-Medical Electronics Co., Ltd
	Disposable IBP transducer	BD
	Transducer/Manifold mount	Shenzhen Mindray Bio-Medical Electronics Co., Ltd
IBP adapter cable		Shenzhen Mindray Bio-Medical Electronics Co., Ltd

### 3.3.6 CO2 Accessories

#### 3.3.6.1 Microstream CO<sub>2</sub> Module

Description	Applicable patient	Remark	Manufacturer
Airway sampling line	Adult, pediatric	Disposable	Oridion Medical1987 Ltd.
Airway sampling line, humidified	Adult, pediatric		Oridion Medical1987 Ltd.
Airway sampling line, humidified	Neonate, infant		Oridion Medical1987 Ltd.
Airway sampling line, humidified	Adult, pediatric		Oridion Medical1987 Ltd.
Airway sampling line, long, humidified	Adult, pediatric		Oridion Medical1987 Ltd.
Airway sampling line, long, humidified	Neonate, infant		Oridion Medical1987 Ltd.
Nasal sampling line	Adult		Oridion Medical1987 Ltd.

Description	Applicable patient	Remark	Manufacturer
Nasal sampling line	Pediatric	Disposable	Oridion Medical1987 Ltd.
Nasal sampling line, plus O <sub>2</sub>	Adult		Oridion Medical1987 Ltd.
Nasal sampling line, plus O <sub>2</sub>	Pediatric		Oridion Medical1987 Ltd.
Nasal sampling line, long, plus O <sub>2</sub>	Adult		Oridion Medical1987 Ltd.
Nasal sampling line, long, plus O <sub>2</sub>	Pediatric		Oridion Medical1987 Ltd.
Nasal sampling line, long	Adult		Oridion Medical1987 Ltd.
Nasal sampling line, long	Pediatric		Oridion Medical1987 Ltd.

Nasal sampling line, long	Neonate	Oridion Medical1987 Ltd.
Nasal sampling line, long, plus O <sub>2</sub>	Adult	Oridion Medical1987 Ltd.
Nasal sampling line, long, plus O <sub>2</sub>	Pediatric	Oridion Medical1987 Ltd.
Nasal sampling line	Adult	Oridion Medical1987 Ltd.
Nasal sampling line	Pediatric	Oridion Medical1987 Ltd.

### 3.3.6.2 Sidestream CO<sub>2</sub> Module

Description	Applicable patient	Remark	Manufacturer
Nasal CO <sub>2</sub> sample cannula	Adult	Disposable	Shenzhen Mindray Bio-Medical Electronics Co., Ltd
	Pediatric		
	Neonate		
Sampling line, adult 2.5m	Adult, pediatric	Disposable	Shenzhen Mindray Bio-Medical Electronics Co., Ltd
Sampling line, neonate 2.5m	Neonate		
DRYLINE airway adapter	Neonate	Disposable	Shenzhen Mindray Bio-Medical Electronics Co., Ltd
	/	Straight, disposable	
	/	Elbow, disposable	
DRYLINE II watertrap	Adult, pediatric	Reusable	Shenzhen Mindray Bio-Medical Electronics Co., Ltd
	Neonate		

### 3.3.7 Therapy Accessories

Description	Model	Applicable patient	Remark	Manufacturer
External paddles	MR6601	Adult, pediatric	Reusable	Shenzhen Mindray Bio-Medical Electronics Co., Ltd

Multifunction electrode pads	MR60	Adult	Disposable (5 sets/pack)	Leonhard Lang GmbH		
	MR61	Pediatric				
	MR62	Adult				
	MR63	Pediatric				
Cable of electrode pads with test load (50 ohm)	MR6701	/	Reusable	Shenzhen Mindray Bio-Medical Electronics Co., Ltd		
Conductive gel	15-25	/	Consumable	PARKER LABORATORIES		
Internal paddles	MR6501	Pediatric	Reusable, 1 inch without button	Shenzhen Mindray Bio-Medical Electronics Co., Ltd		
			Reusable, 1 inch with button			
	MR6502	Adult, pediatric	Reusable, 2 inches without button			
			Reusable, 2 inches with button			
	MR6503	Adult	Reusable, 3 inches without button			
			Reusable, 3 inches with button			
	CPR sensor kit	MR6401	/		Reusable, without battery	Shenzhen Mindray Bio-Medical Electronics Co., Ltd
					Reusable, with a battery	
CPR sensor cable	MR6801	/	Reusable	Shenzhen Mindray Bio-Medical Electronics Co., Ltd		
CPR adhesive tape	MR6921	/	Disposable (3 sets/pack)	Shenzhen Mindray Bio-Medical Electronics Co., Ltd		



### 3.3.8 Miscellaneous

Description	Model	Manufacturer
Rechargeable lithium ion battery	LI24I004A	Shenzhen Mindray Bio-Medical Electronics Co., Ltd
	LI34I001A	
Test load	MR6905	Shenzhen Mindray Bio-Medical Electronics Co., Ltd
Bedrail hook	/	Shenzhen Mindray Bio-Medical Electronics Co., Ltd
Wi-Fi to 4G router kit	IR615-S-L5-WLAN	Shenzhen Mindray Bio-Medical Electronics Co., Ltd
Analog output cable	/	Shenzhen Mindray Bio-Medical Electronics Co., Ltd
Synchronous defibrillation input cable	/	Shenzhen Mindray Bio-Medical Electronics Co., Ltd
Grounding cable	UL1015/14AWG	Shenzhen Mindray Bio-Medical Electronics Co., Ltd
DC/AC adapter	/	Shenzhen Mindray Bio-Medical Electronics Co., Ltd
Patient data management software kit	/	Shenzhen Mindray Bio-Medical Electronics Co., Ltd
PHIES software kit	/	Shenzhen Mindray Bio-Medical Electronics Co., Ltd
Vehicle mount kit	/	Shenzhen Mindray Bio-

		Medical Electronics Co., Ltd
--	--	------------------------------

Description	Model	Manufacturer
Carrying case	/	Shenzhen Mindray Bio-Medical Electronics Co., Ltd
Conducting gel mount kit	/	Shenzhen Mindray Bio-Medical Electronics Co., Ltd
Charger Station kit (International)	BatteryFeed 20	Shenzhen Mindray Bio-Medical Electronics Co., Ltd
Charger Station kit (US)		Shenzhen Mindray Bio-Medical Electronics Co., Ltd
Charger Station kit (Indian)		Shenzhen Mindray Bio-Medical Electronics Co., Ltd
Charger Station kit (EU)		Shenzhen Mindray Bio-Medical Electronics Co., Ltd
Charger Station kit (Brazilian)		Shenzhen Mindray Bio-Medical Electronics Co., Ltd
Charger Station kit (UK)		Shenzhen Mindray Bio-Medical Electronics Co., Ltd

### 3.4 Another device

The equipment can be connected to a Central Monitoring System (hereinafter called CMS) through wired and wireless networks. The CMS uses the software system of Mindray and the general hospital surgical network.

## 4. Residual risks and undesirable effects, warnings and precautions

### 4.1 Residual risks and undesirable effects

**Residual risks:**

With complete product risk analysis and control, only one residual risk assessment belongs to the "AFAP" level. The table below gives the residual risk and benefit analysis of these hazards:

Item	Energy	Sequence of events	Hazardous situation	harm	Control Measure	Risk level
C3.2.1	Energy hazards	High voltage exists in electrode pads and heavy current flows through the myocardium during defibrillation therapy.	High discharge energy and large discharge current lead to myocardial cell damage.	Patient's cardiac muscle is injured due to electric shock during defibrillation therapy.	Defibrillation was performed with Mindray biphasic truncated exponential waveform with better therapeutic efficacy and less myocardial injury to the patient. This wave is similar to the discharge waveform of the AHA-approved model, and its defibrillation effectiveness and safety are equivalent.	AFAP

**undesirable effects:**

Through clinical data from post-market surveillance activity (including adverse event report analysis and post-marketing clinical follow-up), there is no undesirable effects identified.

After search the literature of similar device, the results of SOTA evaluation shown that undesirable effects may include myocardial damage.

**C3.2.1 Risk and Benefit Analysis:**

Through the parameter monitoring, medical staff have established sufficient conditions to provide patients with a better medical monitoring environment, and the benefits are obvious.

Although there is also the possibility of false positives and false negatives in parameter monitoring, the impact of false positives and false negatives is limited and will not cause substantial harm to patients.

In addition, the parameter monitoring of the monitor has the advantages of simplicity of equipment, convenient operation, timeliness, economy, etc. compared with other known ones.

Therefore, from the perspective of benefit and risk, monitor parameter monitoring has obvious benefits, controllable risks, and has strong clinical application popularization characteristics.

Defibrillators are life-saving devices used in emergency situations. They have been shown to have a high benefit for patients with underlying diseases that remain undetected until sudden cardiac arrest occurs. The time from collapse to defibrillation is critical in-patient survival. For every minute that passes between collapse and defibrillation, survival rates from VF SCA decrease 7% to 10%.

In conclusion, given the available information above, the defibrillator's support for patients in cardiac arrest who are unconscious, not breathing, or without circulation the probable benefits outweigh the probable risks.

The harm of defibrillation to human body is mainly the damage to myocardial cell by the peak current during defibrillation, and the result of the animal experiment indicates that the peak current ( $I_{50}$ ) measured under each impedance for Mindray biphasic waveform is significantly smaller than the result measured for the monophasic waveform (with the statistically significant difference  $P < 0.05$ ), equivalent to that of contrast biphasic waveform. Therefore, it is proved that Mindray biphasic waveform has better safety than that of MDS monophasic waveform, and the safety is equivalent to that of contrast biphasic waveform ( $P$  value shows no statistically significant difference).

## 4.2 Warnings and precautions

- Before putting the system into operation, the operator must verify that the equipment, connecting cables and accessories are in correct working order and operating condition.
- Make sure the synchronous input system is applied to this equipment and the input signal is correct if necessary.
- To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth. If the installation does not provide for a protective earth conductor, disconnect it from the power line and operate it on smart lithium-ion batteries.
- Ensure that the equipment is supplied with continuous electric power during work. Sudden power failure leads to the loss of patient data.
- Use and store the equipment in specified environmental condition. The equipment and accessories may not meet the performance specification due to aging, stored or used outside the specified temperature and humidity range.
- This equipment is used for single patient at a time.
- The equipment is not intended to be used within the Magnetic Resonance (MR) environment.
- Before each use, the operator must check the equipment condition to ensure that the equipment is ready for operation.
- Medical electrical equipment which does not incorporate defibrillator protection should be disconnected during defibrillation.

- Do not defibrillate a patient who lies on the wet ground.
- Do not touch the patient and live parts simultaneously.
- Do not touch the patient when connecting the peripheral equipment via the I/O signal ports to prevent patient leakage current from exceeding the requirements specified by the standard.
- Do not rely exclusively on the audible alarm system for patient monitoring. Adjustment of alarm volume to a low level or off may result in a hazard to the patient. Remember that alarm settings should be customized according to different patient situations and always keeping the patient under close surveillance is the most reliable way for safe patient monitoring.
- Do not perform any functional check if the equipment is connected with a patient; otherwise the patient might be shocked.
- Remain attentive to the patient during applying therapy. Delay in delivering a shock may result in a rhythm that was analyzed as shockable converting spontaneously to non-shockable and could result in inappropriate delivery of a shock.
- For the treatment of patients with implantable pacemakers, place electrode pads or paddles away from internal pacemaker generator if possible to help prevent damage to the pacemaker.
- Do not place the equipment or accessories in any position that might cause it to fall on the patient.
- Do not start or operate the equipment unless the setup was verified to be correct.
- To avoid inadvertent disconnection, route all cables in a way to prevent a stumbling hazard. Wrap and secure excess cabling to reduce risk of entanglement or strangulation by patients or personnel.
- If any measurement seems questionable, first check the patient's vital signs by alternate means and then check the equipment for proper functioning.
- Physiological data and alarm messages provided by the equipment should not be used as the sole basis for diagnosis or therapy decisions. They must be used in conjunction with clinical signs and symptoms. Misinterpretation of the measured values or other parameters can endanger the patient.
- Do not touch device connectors, recorder print head, battery connector or other live equipment if in contact with the patient; otherwise patient injury may result.
- To ensure patient safety, use only parts and accessories specified in this manual.
- When disposing of the packaging material, be sure to observe the applicable waste control regulations and keep it out of children's reach.

### 4.3 Post-Market surveillance (PMS) database research

According to Regulation (EU) 2017/745 Article 84, mindray has established the post-market surveillance system and actively gathering relevant data. By now, there is no serious incidents or field safety corrective action (FSCA) and no corrective action preventive action (CAPA) reported about BeneHeart D5/D6 Defibrillator/Monitor.

In the final collection results, the adverse event rate was 0% and the complaint rate was 0.06%.

Complaint handling and risk control (individual complaint records and trend analysis of complaints) as below:

No.	Complaint No.	Complaint Initiation Year	Product Model	Country	Problem description	Dangerous situations	Whether it is related to the product	Risk involved (Y/N)
1	CP14-JH0068	2014	BeneHeart D6	France	French hospitals complain that the 12-lead ECG accessory of D6 is prone to damage. After multiple communication with the hospital, the detailed information of the exception cannot be obtained and investigation cannot be conducted.	The ECG information cannot be obtained, or the ECG information is incorrect.	NO After multiple communication with the hospital, the detailed information of the exception cannot be obtained and investigation cannot be conducted.	N
2	CP-2014-021	2014	BeneHeart D6	Germany	During resuscitation, the equipment specified above displayed a false ECG rhythm. According to the	The ECG information cannot be	NO There is a serious deviation between	N

				<p>display screen, a wide complex tachycardia was present, with a frequency of 180/min. The LP 15 Physio-Control ECG device from the Langenhoven (80-31) emergency vehicle present at the scene was used in parallel with the ECG device named above (D6). This led to significant deviations in the rhythm. The LP 15 Physio-Control presented asystole cardiac arrest and the ECG device (D6) continued to display wide complex tachycardia. The crew from the Langenhoven emergency vehicle are also due to provide a statement. The device has been taken out of service as a result of this error report.</p> <p>There is a serious deviation between the information provided by the user and the device log. The user is unwilling to provide the exact</p>	<p>obtained, or the ECG information is incorrect.</p>	<p>the information provided by the user and the device log. The user is unwilling to provide the exact information and cannot perform further investigation.</p>	
--	--	--	--	--	---	--	--

					information and cannot perform further investigation.			
3	CP-2014-028	2014	BeneHeart D6	Germany	The ECG measurement is inaccurate due to interference. Not a product problem.	The ECG information cannot be obtained, or the ECG information is incorrect.	NO The R&D department combined with the product log analysis indicates that the error may be caused by poor contact of the electrode pads or movement of the patient. Not a product problem.	N



4	CP15-JH0323	2015	BeneHeart D6	Ecuador	The battery of BeneHeart D6 failed in less than one month.	During a patient use event, the device inappropriately shut down or restarted itself.	YES When the battery energy is low, the device generates an alarm.	Y; Risk items: KF-0651-3-001 Risk Management Report: C5.2.1、C5.3.10 Current Control Measure: Real-time battery capacity detection, status display, and low battery alarm. Risk Level After Control Measure: Acceptable
5	CP16-JH0017	2016	BeneHeart D6	UK	A British hospital has 40 newly installed BeneHeart D6, 9 of them failed in self test.The confirmed problem is that the impedance detection error is large due to incorrect calibration (zero	Defibrillation cannot discharge	NO The customer did not perform self-test according to the manual requirements,	N

					calibration and gain calibration) by the user, and the measured value exceeds the detection threshold during self-test. After the customer uses the calibration method recommended by the Mindray manual for calibration, the problem is solved and the device works normally.		resulting in poor self-test of the device. This device is not faulty.	
6	CP16-JH0040	2016	BeneHeart D5	China	The client complains that the device Beneheart D5 does not discharge when the synchronous defibrillation therapy is used for the AF patient. Analysis on the device shows that the discharge operation fails due to automatic release, which may be related to external factors such as the electrode plate is not in good contact with the human body during discharge, and the device functions properly.	Defibrillation cannot discharge	NO The paddles are not in good contact with the human body during the discharge, so the device cannot discharge. Unrelated to the device.	N

7	CP16-JH0089	2016	BeneHeart D6	China	<p>Beneheart D6 cannot discharge on the patient.</p> <p>According to the analysis and test, the device function is normal. It is speculated that the customer feedback that the device cannot be discharged on the patient may be related to the automatic disarming (too high impedance) caused by such external factors as electrode plate contamination and improper contact between the electrode plate and the human body during the discharge.</p>	Defibrillation cannot discharge	<p>NO</p> <p>The paddles are not in good contact with the human body during the discharge, so the device cannot discharge. Unrelated to the device.</p>	N
8	CP16-JH0090	2016	BeneHeart D6	China	<p>High energy charging failure occurred in March, 2016 for Beneheart D6.</p> <p>The device returns to the R&amp;D department for analysis. The function of the device is normal and the problem information complained by the user cannot be reproduced. The device log also has</p>	Defibrillation cannot discharge	<p>NO</p> <p>The customer feedback information is incorrect. In combination with the product record log, there is no adverse event reported by the</p>	N

					no exception record.		customer. No problem with this device.	
9	CP16-JH0568	2016	BeneHeart D6	China	Discharge failure during D6 defibrillation in Beijing sanluju hospital. According to the investigation, the problem is that the hospital does not use the conductive paste according to the manual requirements, but uses the ultrasonic coupling agent so that the impedance is too high to be discharged. The device is detected to be normal.	Defibrillation cannot discharge	NO The hospital does not use the conductive paste according to the manual requirements, but uses the ultrasonic coupling agent.No problem with this device.	N
10	CP16-JH0497	2016	BeneHeart D6	China	No discharge during the rescue using BeneHeart D6. According to the investigation, the problem is that the hospital does not use the conductive paste according to the manual requirements, but uses the ultrasonic coupling agent so that the impedance is too high to be discharged. The device is	Defibrillation cannot discharge	NO The hospital does not use the conductive paste according to the manual requirements, but uses the ultrasonic coupling agent.No problem with this	N

					detected to be normal.		device.	
11	CP17-JH0052	2017	BeneHeart D6	China	<p>Suspect no discharge during the rescue using BeneHeart D6.</p> <p>Communicate with the hospital and learn that the patient has asystole and there is no ventricular fibrillation waveform before the BeneHeart D6 defibrillation monitor is used to rescue the patient. The BeneHeart D6 operator's manual specifies that the defibrillation function is applicable to patients with ventricular fibrillation and ventricular tachycardia without respiration and pulse. The BeneHeart D6 defibrillation monitor works properly.</p>	Defibrillation cannot discharge	<p>NO</p> <p>The patient has asystole and there is no ventricular fibrillation waveform before the BeneHeart D6 defibrillation monitor is used to rescue the patient. Not an indication for Intended Use of this device.No problem with this device.</p>	N
12	CP17-JH0110	2017	BeneHeart D6	Mexico	<p>The device cannot be turned on. The battery has a burning smell. The BeneHeart D6 in Mexico.</p> <p>Occasional battery-related faults occur. The equipment is made of</p>	During a patient use event , the device inappropriately	YES	<p>Y;</p> <p>Risk items: KF-0651-3-001 Risk Management Report: C5.1.13</p>

					flame-retardant materials and will not burn. This problem may cause startup failure.	shut down or restarted itself.		Current Control Measure: Design fireproof materials for equipment enclosure Risk Level After Control Measure: Acceptable
13	CP17-JH0116	2017	BeneHeart D6	China	The customer reported that "sparking" occurred when AC power was connected. When the AC power supply is no longer plugged in, change to battery power and the device can be started normally. This problem is caused by that the customer removes the anti-off hook and the plug connection is unstable during long-term use.	During a patient use event, the device inappropriately shut down or restarted itself.	NO This problem is caused by that the customer removes the anti-off hook and the plug connection is unstable during long-term use. No problem with this device.	N

14	CP17-JH0124	2017	BeneHeart D6	China	<p>Using D6 defibrillation monitor in Ya'an people's Hospital, The department doctor judges that multiple shocks have no effect based on "the patient's body and skin do not respond to the shock," and hopes the manufacturer to conduct an investigation.</p> <p>Mindray has performed functional and performance tests on the device, and all functions are normal without device exception. The records of multiple defibrillation discharges mentioned by the customer indicate that all the discharges were successful. The absence of obvious limb response after defibrillation discharge may be related to the patient's physical condition at that time.</p>	Defibrillation cannot discharge	<p>NO</p> <p>The records of multiple defibrillation discharges mentioned by the customer indicate that all the discharges were successful. The absence of obvious limb response after defibrillation discharge may be related to the patient's physical condition at that time.No problem with this device.</p>	N
----	-------------	------	--------------	-------	--	---------------------------------	---	---

15	CP1904-JH00657	2019	BeneHeart D6	China	<p>Henan Provincial Employee Hospital reported that during the out-of-hospital emergency treatment, it was found that the BeneHeart D6 12-lead ECG still had waveforms after the patient died.</p> <p>This is ECG mechanical disconnection. Not caused by Mindray products.</p>	<p>The ECG information cannot be obtained, or the ECG information is incorrect.</p>	<p>NO This is ECG mechanical disconnection.No problem with this device.</p>	N
16	CP1906-JH00767	2019	BeneHeart D6	China	<p>During defibrillation in the intervention Department of Ya'an people's Hospital of Sichuan Province, the Department of the hospital said that during the use of 200J defibrillation, the patient's response was weak after discharge, and the ECG had no response.</p> <p>It is confirmed that the information reported by the hospital is incorrect, and there is no "invalid shock, resulting in patient death" event. In addition, the hospital did not raise</p>	<p>Defibrillation cannot discharge</p>	<p>NO There is no "invalid shock, resulting in patient death" event. In addition, the hospital did not raise any problems with the device and only required to replace the device. No problem with this device.</p>	N



					any problems with the device and only required to replace the device.			
17	CP1911-JH01050	2019	BENEHEART D6	China	Ya'an People's Hospital of Sichuan Province The hospital thinks that the defibrillator does not discharge. After the investigation, it is confirmed that the user misoperates the shock when it is not fully charged. When the impedance is too high, the shock is released. The device is normal.	Defibrillation cannot discharge	NO The user misoperates the shock when it is not fully charged. When the impedance is too high, the shock is released.No problem with this device.	N
18	CP1910-JH01012	2019	BENEHEART D6	France	D6 - Batt failed The D6 is not capable to see the complete battery status and cannot see the real battery status. The user does not detect and maintain the battery regularly as required by the manual, and does not replace the battery for 2 years as required by the manual.	During a patient use event , the device inappropriately shut down or restarted itself.	NO The user does not detect and maintain the battery regularly as required by the manual, and does not replace the battery for 2 years as required by the manual.No problem	N

							with this device.	
19	CP2011-JH01757	2020	BeneHeart D6	China	Sparks on the BENEHEART D6 handle of Jiaokou County Hospital of Traditional Chinese Medicine. After checking the device, the handle is cracked, which should be caused by the customer's falling down during use. The liquid penetrates the handle through the crack, causing the client phenomenon.	Unexpected release of defibrillation energy	NO The handle is cracked, which should be caused by the customer's falling down during use. The liquid penetrates the handle through the crack, causing the client phenomenon.No problem with this device.	N

20	CP2011-JH01729	2020	BeneHeart D6	Germany	<p>Feedback from German office: The nurse felt the hand was shocked during the user test.</p> <p>After testing, all the functions of the monitor are normal.</p> <p>1. The function test is normal.</p> <p>2. The test results of the contact impedance and leakage current meet the standards, and the safety test is normal.</p>	<p>Unexpected release of defibrillation energy</p>	<p>NO</p> <p>After testing, all the functions of the monitor are normal.</p> <p>1. The function test is normal.</p> <p>2. The test results of the contact impedance and leakage current meet the standards, and the safety test is normal.</p> <p>No problem with this device.</p>	N
21	CP2101-JH01908	2021	BeneHeart D6	China	<p>The D6 waveform is still normal after the patient asystole.</p> <p>This is ECG mechanical disconnection. Not caused by Mindray products.</p>	<p>The ECG information cannot be obtained, or the ECG information is incorrect.</p>	<p>NO</p> <p>This is ECG mechanical disconnection. Not caused by Mindray products.</p>	N

22	CP2112-JH02590	2021	BeneHeart D6	China	<p>According to Baoan Central Hospital, D6 cannot discharge.</p> <p>After testing the product functions and checking and analyzing the device logs, the device is normal. Finally, confirmed the device with the medical personnel in the hospital, The discharge function of the device is normal.</p>	Defibrillation cannot discharge	<p>NO</p> <p>The device cannot discharge due to incorrect operation by the customer.No problem with this device.</p>	N
23	CP2101-JH01881	2021	BENEHEART D6	China	<p>Defibrillation Failure in Affiliated Hospital of Shanxi College of Traditional Chinese Medicine.</p> <p>The R&amp;D analyzed the returned log and found that the defibrillator successfully discharged five times on the day of the incident . No exception was found to the devices. The failure of defibrillation rescue is unrelated to the defibrillator and electrode pads, and may be related to the patient's physiological status and rescue time.</p>	Defibrillation cannot discharge	<p>NO</p> <p>The R&amp;D analyzed the returned log and found that the defibrillator successfully discharged five times on the day of the incident . No exception was found to the devices. The failure of defibrillation rescue is unrelated to the</p>	N

							defibrillator and electrode pads, and may be related to the patient's physiological status and rescue time.No problem with this device.	
--	--	--	--	--	--	--	---	--

#### 4.4 Complaint rates of the specific warning and precautions

Warnings and precautions	Complaint No.	Complaint rates
Before putting the system into operation, the operator must verify that the equipment, connecting cables and accessories are in correct working order and operating condition.	CP14-JH0068、CP16-JH0089	0.005%
<ul style="list-style-type: none"> <li>Make sure the synchronous input system is applied to this equipment and the input signal is correct if necessary.</li> </ul>	0	0.000%
<ul style="list-style-type: none"> <li>To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth. If the installation does not provide for a protective earth conductor, disconnect it from the power line and</li> </ul>	CP17-JH0116	0.003%

operate it on smart lithium-ion batteries.		
<ul style="list-style-type: none"> <li>Ensure that the equipment is supplied with continuous electric power during work. Sudden power failure leads to the loss of patient data.</li> </ul>	0	0.000%
<ul style="list-style-type: none"> <li>Use and store the equipment in specified environmental condition. The equipment and accessories may not meet the performance specification due to aging, stored or used outside the specified temperature and humidity range.</li> </ul>	0	0.000%
<ul style="list-style-type: none"> <li>This equipment is used for single patient at a time.</li> </ul>	0	0.000%
<ul style="list-style-type: none"> <li>The equipment is not intended to be used within the Magnetic Resonance (MR) environment.</li> </ul>	0	0.000%
<ul style="list-style-type: none"> <li>Before each use, the operator must check the equipment condition to ensure that the equipment is ready for operation.</li> </ul>	CP15-JH0323、 CP17-JH0110、 CP1910-JH01012、 CP2011-JH01757、 CP16-JH0090	0.013%
<ul style="list-style-type: none"> <li>Medical electrical equipment which does not incorporate defibrillator protection should be disconnected during defibrillation.</li> </ul>	0	0.000%
<ul style="list-style-type: none"> <li>Do not defibrillate a patient who lies on the wet ground.</li> </ul>	0	0.000%
<ul style="list-style-type: none"> <li>Do not touch the patient and live parts simultaneously.</li> </ul>	CP2011-JH01729	0.003%
<ul style="list-style-type: none"> <li>Do not touch the patient when connecting the peripheral equipment via the I/O signal ports to prevent patient leakage current from exceeding the requirements specified by the standard.</li> </ul>	0	0.000%

<ul style="list-style-type: none"><li>Do not rely exclusively on the audible alarm system for patient monitoring. Adjustment of alarm volume to a low level or off may result in a hazard to the patient. Remember that alarm settings should be customized according to different patient situations and always keeping the patient under close surveillance is the most reliable way for safe patient monitoring.</li></ul>	0	0.000%
<ul style="list-style-type: none"><li>Do not perform any functional check if the equipment is connected with a patient; otherwise the patient might be shocked.</li></ul>	0	0.000%
<ul style="list-style-type: none"><li>Remain attentive to the patient during applying therapy. Delay in delivering a shock may result in a rhythm that was analyzed as shockable converting spontaneously to non-shockable and could result in inappropriate delivery of a shock.</li></ul>	0	0.000%
<ul style="list-style-type: none"><li>For the treatment of patients with implantable pacemakers, place electrode pads or paddles away from internal pacemaker generator if possible to help prevent damage to the pacemaker.</li></ul>	0	0.000%
<ul style="list-style-type: none"><li>Do not place the equipment or accessories in any position that might cause it to fall on the patient.</li></ul>	0	0.000%
<ul style="list-style-type: none"><li>Do not start or operate the equipment unless the setup was verified to be correct.</li></ul>	CP16-JH0017、CP16-JH0040、CP1911-JH01050	0.008%
<ul style="list-style-type: none"><li>To avoid inadvertent disconnection, route all cables in a way to prevent a stumbling hazard. Wrap and secure excess cabling to reduce risk of entanglement or strangulation by patients or personnel.</li></ul>	0	0.000%

<ul style="list-style-type: none"> <li>If any measurement seems questionable, first check the patient's vital signs by alternate means and then check the equipment for proper functioning.</li> </ul>	CP-2014-021、 CP17-JH0052	0.005%
<ul style="list-style-type: none"> <li>Physiological data and alarm messages provided by the equipment should not be used as the sole basis for diagnosis or therapy decisions. They must be used in conjunction with clinical signs and symptoms. Misinterpretation of the measured values or other parameters can endanger the patient.</li> </ul>	CP-2014-028、 CP1904-JH00657、 CP2101-JH01908、 CP17-JH0124	0.010%
<ul style="list-style-type: none"> <li>Do not touch device connectors, recorder print head, battery connector or other live equipment if in contact with the patient; otherwise patient injury may result.</li> </ul>	0	0.000%
<ul style="list-style-type: none"> <li>To ensure patient safety, use only parts and accessories specified in this manual.</li> </ul>	CP16-JH0568、 CP16-JH0497	0.005%
<ul style="list-style-type: none"> <li>When disposing of the packaging material, be sure to observe the applicable waste control regulations and keep it out of children's reach.</li> </ul>	0	0.000%

## 5. The summary of clinical evaluation and post-market clinical follow-up (PMCF)

The data from the current technology of defibrillation waveforms, standard compliance testing, post-marketing clinical follow-up in China and Europe. It can fully prove that D6 defibrillation monitor can meet clinical safety and effectiveness.

After marketing in EU, we have formulated a detailed post-marketing clinical follow-up plan to collect a certain amount of cases according to statistical requirements to ensure that the External Defibrillation, synchronized cardioversion, Internal Defibrillation, AED, and the Non-invasive external pacing function of the D6 defibrillator can also meet the safety and effectiveness of European users.



## 5.1 Summary of clinical data related to equivalent device (if applicable)

The Defibrillator/Monitor belongs to the combination of Manual defibrillation, monitor, pacemaker and AED module, these modules are independent of each other and will not affect other modules when a single module works.

Since the defibrillation module and monitor module of the defibrillation monitor are independent of each other, and the module development of the subject device is based on the previous generation products, the technical characteristics, biological characteristics and clinical technical characteristics of the device are completely the same, and the clinical trial is conducted on the technology of the previous generation products to obtain the clinical data.

The monitor module of BeneHeart D5/ BeneHeart D6 is thus determined to be as safe and effective as the equivalent CE marked devices including BeneVision N12 Patient Monitor, iPM 12 Patient Monitor and BeneView T5 patient monitors.

<b>Administrative</b>	<b>Device 1 (subject device)</b> Description of characteristics and reference to specifying documents	<b>Device 2 (marketed device)</b> Description of characteristics and reference to specifying documents	<b>Device 3(marked device)</b> Description of characteristics and reference to specifying documents	<b>Device 4(marked device)</b> Description of characteristics and reference to specifying documents	<b>Identified differences or conclusion that there are no differences in the characteristic</b>
Manufacturer	SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD	SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD	SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD	SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD	Same.
Device Trade Name	BeneHeart D5/ BeneHeart D6 Defibrillator/Monitor	BeneVision N12 Patient Monitor	iPM 12 Patient Monitor	BeneView T5 patient monitor	Not applicable

## 5.2 Summary of clinical data from conducted investigations of the device before the CE-marking (if applicable).

### SpO2:

#### Information of the studied device

model	Beneview T5	Identity	CM-57144714、CM-57144696、CM-57144695 CM-57144711、CM-57144698、CM-57144700 CM-57144704、CM-57144646、CM-57144706 CM-57144693、CM-57144705、CM-43137197
	SPO2 Module		KLA62000009、KLA62000010、KLA62000011 KLA62000012、KLA62000013、KLA62000014 KLA62000015、KLA62000016、KLA62000017 KLA62000019、KLA62000020、CTC22045792

#### Intended use of studied device

This BeneView T5 patient monitor is intended to be used for monitoring, displaying, reviewing, storing and transferring of multiple physiological parameters including ECG, heart rate (HR), respiration (Resp), temperature (Temp), SpO2, pulse rate (PR), non-invasive blood pressure (NIBP), invasive blood pressure (IBP), cardiac output (C.O.), carbon dioxide (CO2), oxygen (O2), anesthetic gas (AG), impedance cardiograph (ICG), bispectral index (BIS), respiration mechanics (RM), continuous cardiac output(CCO), central venous oxygen saturation(ScvO2).

#### Objectives of the study

Though clinical trials to verify the effectiveness and accuracy of the SPO2 measurement function of Mindray's BeneView T5 patient monitor in clinical applications.

#### Summary of study methods

According to the test method recommended by the ISO80601-2-61:2011 SpO2 professional clinical standard, the SpO2 measurement results were compared with

the invasive blood gas analysis measurement values.

Inclusion criteria:

1. Subjects must understand clinical trials and participate voluntarily.
2. Meet the ASA classification I of the Anesthesia Association.
3. The subjects are healthy adults COHb<3%, MetHb<2%, ctHb>10g/dL.
4. Healthy adult subjects need to be able to withstand the minimum medical risk of controlled oxygen reduction to the specified level in the clinical trial protocol.

**Exclusion** criteria:

1. The subject has any systemic disease.
2. Subject is methemoglobin.
3. Subjects do not understand clinical trials and their risks.
4. Smoker
5. Pregnancy
6. Subject has symptoms or history of peripheral ischemia
7. The age of the subject must not be less than 18 years old
8. Other unqualified conditions identified by clinicians

### **Primary and secondary endpoint**

The acceptance criteria and data analysis are as follows:

Accuracy: The measurement range of SpO<sub>2</sub> is 0%~100%, and the measurement error is  $\pm 3\%$  in the range of 70%~100%. And the measurement range of PR is 20~300bpm, and the measurement error is  $\pm 3$ bpm in the range of 20~300bpm.

Requirements for patients: Each probe should test at least 10 subjects, and a total of at least 200 sets of comparative data are obtained for each type of subjects.

### **Results**

A total of 20 subjects were completed in this trial, and each subject was subject to 20 comparison tests. The results are as follows:

Table1 Analysis of SpO2 clinical trial data

SpO2 probe	Probe Type	Probed Site	Index	Accuracy	Acceptance Criterion	Conclusion
512F	Adult	Finger	SPO2	1.1013	$\leq 2$	Meet the requirement
			PR	2.4619	$\leq 3$	Meet the requirement
512H	Pediatric	Finger	SPO2	1.1529	$\leq 2$	Meet the requirement
			PR	2.3752	$\leq 3$	Meet the requirement
512E	Adult	Finger	SPO2	1.2572	$\leq 2$	Meet the requirement
			PR	2.1972	$\leq 3$	Meet the requirement
512G	Pediatric	Finger	SPO2	1.3547	$\leq 2$	Meet the requirement
			PR	1.2057	$\leq 3$	Meet the requirement
518B	Neonate	Foot	SPO2	1.3804	$\leq 2$	Meet the requirement
			PR	1.5186	$\leq 3$	Meet the requirement
520A	Adult	Foot	SPO2	1.2773	$\leq 2$	Meet the requirement
			PR	2.1689	$\leq 3$	Meet the requirement
520P	Pediatric	Foot	SPO2	1.1658	$\leq 2$	Meet the requirement
			PR	2.2028	$\leq 3$	Meet the requirement
520I	Infant	Toe	SPO2	1.3123	$\leq 2$	Meet the requirement
			PR	2.3677	$\leq 3$	Meet the requirement
520N	Neonate	Foot	SPO2	1.3577	$\leq 2$	Meet the requirement
			PR	2.1644	$\leq 3$	Meet the requirement

### Conclusion

The effectiveness and accuracy of the SPO2 measurement function of Mindray's BeneView T5 patient monitor meet the ISO 81060-2:2013 SPO2 professional clinical standard.

## NIBP:

### Information of the studied device

model	BeneVision N12	Identity	F8-6C000031
			F8-6C000033
			F8-6C000039
	NIBP Module (MPM 3.0)		KLA75001549
			KLA75001544
			KLA75001563
			KLA75001547

### Intended use of studied device

This patient monitor is intended to be used for monitoring, displaying, reviewing, storing and transferring of multiple physiological parameters including ECG, heart rate (HR), respiration (Resp), temperature (Temp), SpO2, pulse rate (PR), non-invasive blood pressure (NIBP), invasive blood pressure (IBP), cardiac output (C.O.), carbon dioxide (CO2), oxygen (O2), anesthetic gas (AG), impedance cardiograph (ICG), bispectral index (BIS), respiration mechanics (RM), continuous cardiac output(CCO), central venous oxygen saturation(ScvO2), electroencephalograph (EEG), neuromuscular transmission (NMT), regional oxygen saturation (rSO2).

### Objectives of the study

Though clinical trials to verify the safety and effectiveness of the NIBP measurement function of Mindray's BeneVision N12 patient monitor in clinical applications.

### Summary of study methods

According to the trial method recommended in the ISO81060-2:2013 NIBP professional clinical standard, the data of adult/pediatric group (3 years old and older) and neonate group (3 years old and younger) were compared with invasive blood pressure.

Group	Attributes	Abstract
adult/pediatric (3 years old and older) compared with invasive blood pressure	Inclusion criteria	<ol style="list-style-type: none"> <li>1. Select people aged 3 years and over who have undergone arterial intubation due to clinical needs for diagnosis to participate in clinical trials.</li> <li>2. Subjects need to understand the clinical trial and agree to participate.</li> <li>3. If the subject is underage or the subject does not have the cognitive ability, the legal guardian of subject must understand the clinical trial and agree to the subject to participate. If the subject has cognitive ability, he should also understand the content of the trial and participate voluntarily.</li> </ol>
	Exclusion criteria	Subjects with severe shock, use of a heart lung machine, upper extremity infection heart or arterial malformations (excluding those after surgical correction).
	Number of the cases	According to the trial method compared with invasive blood pressure recommended in the ISO 81060-2:2013 NIBP professional clinical standard, the number of Subjects for adult/pediatric group (3 years old and older) must not be less than 20 cases.
adult/pediatric (3 years old and older) compared with auscultation	Inclusion criteria	<ol style="list-style-type: none"> <li>1. Select people aged 3 years and over who have undergone arterial intubation due to clinical needs for diagnosis to participate in clinical trials.</li> <li>2. Subjects need to understand the clinical trial and agree to participate.</li> <li>3. If the subject is underage or the subject does not have the cognitive ability, the legal guardian of subject must understand the clinical trial and agree to the subject to participate. If the subject has cognitive ability, he should also understand the content of the trial and participate voluntarily.</li> </ol>
	Exclusion criteria	Subjects with severe shock, use of a heart lung machine, upper extremity infection heart or arterial malformations (excluding those after surgical correction)
	Number of the cases	According to the trial method compared with auscultation recommended in the ISO 81060-2:2013 NIBP

		professional clinical standard, the number of Subjects for adult/pediatric group (3 years old and older) must not be less than 85 cases included at least 35 subjects aged 3~12 years.
Neonate (3 years old and younger) compared with invasive blood pressure	Inclusion criteria	<ol style="list-style-type: none"><li>1. Select people aged 3 years and younger who have undergone arterial intubation due to clinical needs for diagnosis to participate in clinical trials.</li><li>2. The legal guardian of subject must understand the clinical trial and agree to the subject to participate. If the subject has cognitive ability, he should also understand the content of the trial and participate voluntarily.</li></ol>
	Exclusion criteria	Subjects with severe shock, use of a heart lung machine, upper extremity infection heart, arterial malformations (excluding those after surgical correction) or severe arrhythmia.
	Number of the cases	According to the trial method compared with invasive blood pressure recommended in the ISO 81060-2:2013 NIBP professional clinical standard, the number of Subjects for group (3 years old and younger) must not be less than 18 cases.

Main evaluation index: The absolute value of the mean of the difference between the systolic/diastolic/mean pressure measured by the tested device and the reference device does not exceed 5.0 mmHg, and the standard deviation does not exceed 8.0 mmHg. If results do not comply with any one of the two criteria, the tested device does not meet the requirement in ISO81060-2:2013 (E).

### Results

A total of 20 subjects including 11 men and 9 women in the adult/pediatric group are tested, and the gender and age distribution of the subjects meet the standard requirements. The results are as follows:

Table2 NIBP clinical trial data analysis for adults/pediatric

All Data		Number of Data Groups	Range (mmHg)	Mean Deviation (mmHg)			Standard Deviation (mmHg)		
				Results	Acceptance Criterion	Conclusion	Results	Acceptance Criterion	Conclusion
Tested method vs invasive blood pressure	SBP	200	86~191	-3.9	Absolute value $\leq$ 5.0	Meet the requirement	1.4	$\leq$ 8.0	Meet the requirement
	DBP	200	40~113	-3.7	Absolute value $\leq$ 5.0	Meet the requirement	1.3	$\leq$ 8.0	Meet the requirement
	MBP	200	62~135	-3.8	Absolute value $\leq$ 5.0	Meet the requirement	1.1	$\leq$ 8.0	Meet the requirement

A total of 85 subjects including 30 men and 55 women in the adult/pediatric group are tested, and the gender and age distribution of the subjects meet the standard requirements. The results are as follows:

Table3 NIBP clinical trial data analysis for adults/pediatric

All Data		Number of Data Groups	Range (mmHg)	Mean Deviation (mmHg)			Standard Deviation (mmHg)		
				Results	Acceptance Criterion	Conclusion	Results	Acceptance Criterion	Conclusion



Tested method vs auscultation	SBP	255	75~218	4.0	Absolute value $\leq 5.0$	Meet the requirement	5.5	$\leq 8.0$	Meet the requirement
	DBP	255	34~121	-1.9	Absolute value $\leq 5.0$	Meet the requirement	4.6	$\leq 8.0$	Meet the requirement

A total of 18 subjects in the neonate group are tested, and the gender and age distribution of the subjects meet the standard requirements. The results are as follows:

Table4 NIBP clinical trial data analysis for neonate

All Data		Number of Data Groups	Range (mmHg)	Mean Deviation (mmHg)			Standard Deviation (mmHg)		
				Results	Acceptance Criterion	Conclusion	Results	Acceptance Criterion	Conclusion
Tested method vs invasive blood pressure	SBP	180	40~117	4.5	Absolute value $\leq 5.0$	Meet the requirement	1.8	$\leq 8.0$	Meet the requirement
	DBP	180	21~80	4.1	Absolute value $\leq 5.0$	Meet the requirement	1.7	$\leq 8.0$	Meet the requirement
	MBP	180	27~93	4.7	Absolute value $\leq 5.0$	Meet the requirement	2.0	$\leq 8.0$	Meet the requirement

### Conclusion

The effectiveness and accuracy of the NIBP measurement function of Mindray's BeneVision N12 patient monitor meet the ISO 81060-2:2013 NIBP professional clinical standard.

**ECG:****Information of the studied device**

device under test	iPM 12 patient monitor
reference device	BeneView T5 patient monitor (Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.)

**Intended use of studied device**

This patient monitor is intended to be used for monitoring, displaying, reviewing, storing and transferring of multiple physiological parameters including ECG, heart rate (HR), respiration (Resp), temperature (Temp), SpO2, pulse rate (PR), non-invasive blood pressure (NIBP), invasive blood pressure (IBP), cardiac output (C.O.), carbon dioxide (CO2), oxygen (O2), anesthetic gas (AG), impedance cardiograph (ICG), bispectral index (BIS), respiration mechanics (RM), continuous cardiac output(CCO), central venous oxygen saturation(ScvO2), electroencephalograph (EEG), neuromuscular transmission (NMT), regional oxygen saturation (rSO2).

**Objectives of the study**

Though clinical trials to verify the safety and effectiveness of the ECG measurement function of Mindray's ipM 12 patient monitor in clinical applications.

**Summary of study methods**

- 1) Use the iPM 12 patient monitor (Device Under Test (DUT)) to monitor ECG for selected subjects, and use the Beneview T5 monitor (Reference Device (RCT)) to monitor the above parameters at the same time. Through this clinical trial, determine the consistent equivalence of the DUT and RCT in terms of the monitoring performance of the ECG parameters.
- 2) Record clearly and in detail: The information of the experimental subjects of the clinical trial; the environmental information of the monitoring; the information of the time and place of the trial (see the clinical trial record form).
- 3) ECG monitoring. On the subject's body, proceed as follows:
  - Attach two sets of cardiac electrodes adjacent to each other according to the five-lead connection method; or attach two sets of cardiac electrodes adjacent to each

other according to the three-lead connection method;

- Then connect the ECG lead wires of the DUT and the RCT respectively.
- Turn on the ST segment analysis and arrhythmia monitoring functions of the DUT and the RCT.
- Turn on the DUT and the RCT and observe the ECG waveform of lead II. During the observation process, you can switch to other leads to evaluate the overall situation of ECG monitoring.
- Record the HR value and the ST segment values of all leads every 3 to 5 minutes, and obtain 20 sets of data.

### Results

Trial data and processing results:

HR (General monitoring)

**Statistic results of respiratory rate**

Statistical item	Number of Data Sets	Mean deviation (bpm)	Standard deviation (bpm)	Mean ± SD Data within range	Mean±SD~ Mean±2SD Data within range	Mean±2SD Data outside range
HR	1340	0.034	0.51	1119 (83.5%)	98 (7.3%)	123 (9.2%)

HR (Diagnosis mode)

**HR statistical results in diagnosis mode**

Statistical item	Number	Mean deviation	Standard deviation	Mean ± SD	Mean ± SD~	Mean±2SD
------------------	--------	----------------	--------------------	-----------	------------	----------

	of Data Sets	(bpm)	(bpm)	Data within range	Mean±2SD Data within range	Data outside range
HR	1340	0.028	0.64	1053 (78.6%)	208 (15.5%)	79 (5.9%)

ST segment

#### Statistical results of ST segment

Statistical item	Number of Data Sets	Mean deviation (mV)	Standard deviation (mV)	Mean ± SD Data within range	Mean ± SD~ Mean±2SD Data within range	Mean±2SD Data outside range
ST-I	1300	0	0.01	1270 (97.7%)	0	30 (2.3%)
ST-II	1300	0	0.01	1255 (93.6%)	0	45 (3.4%)
ST-III	1300	0	0.01	1255 (93.6%)	0	45 (3.4%)
ST-aVR	1300	0	0.01	1259 (94.4%)	0	41 (3.1%)
SP-aVL	1300	0	0.01	1265 (95.4%)	0	35 (2.6%)

#### Conclusions

During the clinical trial, all subjects are in good condition, and none of the subjects have discomfort symptoms or adverse reactions due to the use of Mindray iPM 12 Patient Monitor with ECG module for monitoring. Therefore, it can be considered safe to use the Mindray ECG module for related parameter monitoring. Clinical performance of ECG function meets SOTA benchmark.

**IBP:****Information of the studied device**

device under test	iPM 12 patient monitor
reference device	BeneView T5 patient monitor (Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.)

**Intended use of studied device**

This patient monitor is intended to be used for monitoring, displaying, reviewing, storing and transferring of multiple physiological parameters including ECG, heart rate (HR), respiration (Resp), temperature (Temp), SpO<sub>2</sub>, pulse rate (PR), non-invasive blood pressure (NIBP), invasive blood pressure (IBP), cardiac output (C.O.), carbon dioxide (CO<sub>2</sub>), oxygen (O<sub>2</sub>), anesthetic gas (AG), impedance cardiograph (ICG), bispectral index (BIS), respiration mechanics (RM), continuous cardiac output (CCO), central venous oxygen saturation (ScvO<sub>2</sub>), electroencephalograph (EEG), neuromuscular transmission (NMT), regional oxygen saturation (rSO<sub>2</sub>).

**Objectives of the study**

Though clinical trials to verify the safety and effectiveness of the IBP measurement function of Mindray's ipM 12 patient monitor in clinical applications.

**Summary of study methods**

(1) Use the iPM 12 patient monitor (Device Under Test (DUT)) to monitor IBP for selected subjects, and use the Beneview T5 monitor (Reference Device (RCT)) to monitor the above parameters at the same time. Through this clinical trial, determine the consistent equivalence of the DUT and RCT in terms of the monitoring performance of the IBP parameters.

(2) Record clearly and in detail: The information of the experimental subjects of the clinical trial; the environmental information of the monitoring; the information of the time and place of the trial (see the clinical trial record form).

(3) IBP monitoring. On the subject's body, proceed as follows:

- Learn the blood supply of the branch before puncture and catheterization. For subjects who underwent radial artery puncture, perform the AllenS test, and perform radial artery puncture for those with normal Allen-S test results. For subjects who underwent dorsal pedis artery puncture, learn the blood supply of the posterior tibial artery, and perform the dorsal pedis artery puncture only after ensuring that the blood supply is good.
- Before the first comparative measurement, calibrate the. pressure for the DUT and RCT.
- Fill the tubing system of the pressure tube and sensor with normal saline, making sure there are no air bubbles in the tubing.
- Perform arterial puncture (radial artery or dorsal pedis artery, as noted on the record form), and connect a manometric tube containing a heparin solution.
- Connect the pressure measuring tube to the pressure tubes and sensors of the DUT and the RCT respectively through the three-way valve switch,
- and zero the pressure sensors of the DUT and the RCT respectively.
- At the same time, directly measure the arterial blood pressure, and record synchronous measurement data of systolic blood pressure, mean blood pressure and diastolic blood pressure every 3 to 5 minutes; obtain 20 sets of data.

**Results**

Acceptance Criteria:

- **IBP: Mean deviation  $\leq \pm 1\text{mmHg}$ , standard deviation  $\leq 5\text{mmHg}$ .**
- **IBP PR: Mean deviation  $\leq \pm 1\text{bpm}$ , standard deviation  $\leq 2\text{ bpm}$ .**

Trial data and processing results:

Systolic pressure

- **Statistic results of IBP systolic pressure**

Statistical item	Number of Data	Mean deviation (rpm)	Standard deviation (rpm)	Mean $\pm$ SD	Mean $\pm$ SD~ Mean $\pm$ 2SD	Mean $\pm$ 2SD
------------------	----------------	----------------------	--------------------------	---------------	----------------------------------	----------------

	Sets			Data within range	Data within range	Data outside range
IBP-sys	1340	0.24	2.39	1076 (80.3%)	141 (10.5%)	123 (9.2%)

MBP

**Statistic results of IBP MBP**

Statistical item	Number of Data Sets	Mean deviation (rpm)	Standard deviation (rpm)	Mean ± SD	Mean±SD~ Mean±2SD	Mean±2SD
				Data within range	Data within range	Data outside range
IBP-mean	1340	0.22	2.53	1063 (79.3%)	176 (13.2%)	101 (7.5%)

DBP

**Statistic results of IBP DBP**

Statistical item	Number of Data Sets	Mean deviation (rpm)	Standard deviation (rpm)	Mean ± SD	Mean±SD~ Mean±2SD	Mean±2SD
				Data within range	Data within range	Data outside range

IBP-dia	1340	0.21	2.61	1069 (79.8%)	181 (13.5%)	90 (6.7%)
---------	------	------	------	-----------------	----------------	--------------

## Conclusions

During the clinical trial, all subjects are in good condition, and none of the subjects have discomfort symptoms or adverse reactions due to the use of Mindray iPM 12 Patient Monitor with IBP module for monitoring. Therefore, it can be considered safe to use the Mindray IBP module for related parameter monitoring. Clinical performance of IBP function meets SOTA benchmark.

## Temp:

### Information of the studied device

device under test	iPM 12 patient monitor
reference device	BeneView T5 patient monitor (Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.)

### Intended use of studied device

This patient monitor is intended to be used for monitoring, displaying, reviewing, storing and transferring of multiple physiological parameters including ECG, heart rate (HR), respiration (Resp), temperature (Temp), SpO2, pulse rate (PR), non-invasive blood pressure (NIBP), invasive blood pressure (IBP), cardiac output (C.O.), carbon dioxide (CO2), oxygen (O2), anesthetic gas (AG), impedance cardiograph (ICG), bispectral index (BIS), respiration mechanics (RM), continuous cardiac output(CCO), central venous oxygen saturation(ScvO2), electroencephalograph (EEG), neuromuscular transmission (NMT), regional oxygen saturation (rSO2).

### Objectives of the study



Though clinical trials to verify the safety and effectiveness of the Temp measurement function of Mindray's ipM 12 patient monitor in clinical applications.

### Summary of study methods

(1) Use the ipM 12 patient monitor (Device Under Test (DUT)) to monitor Temp for selected subjects, and use the Beneview T5 monitor (Reference Device (RCT)) to monitor the above parameters at the same time. Through this clinical trial, determine the consistent equivalence of the DUT and RCT in terms of the monitoring performance of the Temp parameters.

(2) Record clearly and in detail: The information of the experimental subjects of the clinical trial; the environmental information of the monitoring; the information of the time and place of the trial (see the clinical trial record form).

(3) Body temperature(Temp) monitoring.

- For subjects who need body temperature monitoring clinically, use the body surface temperature probe to compare body temperature.
- Adhere the body temperature probes of the DUT and the RCT adjacent to the subject's body (The axilla is the first choice, and other parts should be indicated in the case report form).
- Record a pair of body temperature values every 3 to 5 minutes. Obtain 20 pairs of data.

### Results

Trial data and processing results:

Body temperature

#### Statistic results of body temperature

Statistical item	Number of Data Sets	Mean deviation (°C)	Standard deviation (°C)	Mean ± SD Data within range	Mean±SD~ Mean±2SD Data within range	Mean±2SD Data outside range
------------------	---------------------	---------------------	-------------------------	--------------------------------	---	--------------------------------

Temp	1340	-0.01	0.09	806 (60.1%)	403 (30.1%)	131 (9.8%)
------	------	-------	------	----------------	----------------	---------------

### Conclusions

During the clinical trial, all subjects are in good condition, and none of the subjects have discomfort symptoms or adverse reactions due to the use of Mindray iPM 12 Patient Monitor with Temp module for monitoring. Therefore, it can be considered safe to use the Mindray Temp module for related parameter monitoring. Clinical performance of Temp function meets SOTA benchmark.

### CO2:

#### Information of the studied device

device under test	iPM 12 patient monitor
reference device	BeneView T5 patient monitor (Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.)

#### Intended use of studied device

This patient monitor is intended to be used for monitoring, displaying, reviewing, storing and transferring of multiple physiological parameters including ECG, heart rate (HR), respiration (Resp), temperature (Temp), SpO2, pulse rate (PR), non-invasive blood pressure (NIBP), invasive blood pressure (IBP), cardiac output (C.O.), carbon dioxide (CO2), oxygen (O2), anesthetic gas (AG), impedance cardiograph (ICG), bispectral index (BIS), respiration mechanics (RM), continuous cardiac output (CCO), central venous oxygen saturation (ScvO2), electroencephalograph (EEG), neuromuscular transmission (NMT), regional oxygen saturation (rSO2).

#### Objectives of the study

Though clinical trials to verify the safety and effectiveness of the CO2 measurement function of Mindray's ipM 12 patient monitor in clinical applications.

### Summary of study methods

(1) Use the iPM 12 patient monitor (Device Under Test (DUT)) to monitor CO<sub>2</sub> for selected subjects, and use the Beneview T5 monitor (Reference Device (RCT)) to monitor the above parameters at the same time. Through this clinical trial, determine the consistent equivalence of the DUT and RCT in terms of the monitoring performance of the CO<sub>2</sub> parameters.

(2) Record clearly and in detail: The information of the experimental subjects of the clinical trial; the environmental information of the monitoring; the information of the time and place of the trial (see the clinical trial record form).

(3) CO<sub>2</sub> monitoring. On the subject's body, proceed as follows:

- Set the operation mode to measurement mode, and set the flow rate, gas compensation, and humidity compensation correctly.
- Select the correct type of water tank to connect the sampling tube, connect it to the CO<sub>2</sub> module, and connect the DUT and the RCT to the respiration airway at the same time through a three-way valve.
- After the CO<sub>2</sub> module finishes preheating and enters the measurement state, record a set of EtCO<sub>2</sub>, InsCO<sub>2</sub> and aWRR values every 3 to 5 minutes and obtain 20 sets of data.

### Results

Acceptance Criteria: Consistency requirements of EtCO<sub>2</sub> and InsCO<sub>2</sub> measurement results: Mean deviation  $\leq \pm 2$ mmHg, standard deviation  $\leq 5$  mmHg.

Trial data and processing results:

Statistical results of ETCO<sub>2</sub> measurements

Statistical item	Number of Data Sets	Mean deviation (rpm)	Standard deviation (rpm)	Mean $\pm$ SD Data within range	Mean $\pm$ SD~ Mean $\pm$ 2SD Data within range	Mean $\pm$ 2SD Data outside range

EtCOi	1340	-0.029	0.71	1102 (82.2%)	144 (10.7%)	94 (7.0%)
-------	------	--------	------	-----------------	----------------	--------------

Statistical results of InsCO2 measurements

**Statistic results of InsCO2**

Statistical item	Number of Data Sets	Mean deviation (rpm)	Standard deviation (rpm)	Mean ± SD Data within range	Mean±SD~ Mean±2SD Data within range	Mean±2SD Data outside range
InsCOi	1340	-0.05	0.49	994 (74.2%)	194 (14.5%)	152 (11.3%)

Statistical results of awRR measurements

Statistical item	Number of Data Sets	Mean deviation (rpm)	Standard deviation (rpm)	Mean ± SD Data within range	Mean±SD~ Mean±2SD Data within range	Mean±2SD Data outside range
InsCOi	1340	0.06	0.28	1221 (91.2%)	0	119 (8.8%)

## Conclusions

During the clinical trial, all subjects are in good condition, and none of the subjects have discomfort symptoms or adverse reactions due to the use of Mindray iPM 12 Patient Monitor with CO2 module for monitoring. Therefore, it can be considered safe to use the Mindray CO2 module for related parameter monitoring. Clinical performance of CO2 function meets SOTA benchmark.

## RESP:

### Information of the studied device

device under test	iPM 12 patient monitor
reference device	BeneView T5 patient monitor (Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.)

### Intended use of studied device

This patient monitor is intended to be used for monitoring, displaying, reviewing, storing and transferring of multiple physiological parameters including ECG, heart rate (HR), respiration (Resp), temperature (Temp), SpO2, pulse rate (PR), non-invasive blood pressure (NIBP), invasive blood pressure (IBP), cardiac output (C.O.), carbon dioxide (CO2), oxygen (O2), anesthetic gas (AG), impedance cardiograph (ICG), bispectral index (BIS), respiration mechanics (RM), continuous cardiac output(CCO), central venous oxygen saturation(ScvO2), electroencephalograph (EEG), neuromuscular transmission (NMT), regional oxygen saturation (rSO2).

### Objectives of the study

Though clinical trials to verify the safety and effectiveness of the CO2 measurement function of Mindray's ipM 12 patient monitor in clinical applications.

### Summary of study methods

(1) Use the iPM 12 patient monitor (Device Under Test (DUT)) to monitor RESP for selected subjects, and use the Beneview T5 monitor (Reference Device (RCT)) to monitor the above parameters at the same time. Through this clinical trial, determine the consistent equivalence of the DUT and RCT in terms of the monitoring performance of the RESP parameters.

(2) Record clearly and in detail: The information of the experimental subjects of the clinical trial; the environmental information of the monitoring; the information of the time and place of the trial (see the clinical trial record form).

(3) RESP monitoring. On the subject's body, proceed as follows:

- Attach two sets of cardiac electrodes adjacent to each other according to the five-lead connection method; or attach two sets of cardiac electrodes adjacent to each other according to the three-lead connection method;
- Then connect the ECG lead wires of the DUT and the RCT respectively.
- Turn on the DUT and the RCT and record the respiratory rate (RR) value every 3 to 5 minutes to obtain 20 sets of data.
- Select the respiratory lead I or II, observe the respiratory waveform of the corresponding lead, and record the selected respiratory lead in the case report form.

## Results

Trial data and processing results:

Respiratory rate

**Statistic results of respiratory rate**

Statistical item	Number of Data Sets	Mean deviation (rpm)	Standard deviation (rpm)	Mean $\pm$ SD Data within range	Mean $\pm$ SD~ Mean $\pm$ 2SD Data within range	Mean $\pm$ 2SD Data outside range
RR	1340	0.06	0.32	1192	0	148

				(88.9%)		(11.1%)
--	--	--	--	---------	--	---------

**Conclusions**

During the clinical trial, all subjects are in good condition, and none of the subjects have discomfort symptoms or adverse reactions due to the use of Mindray iPM 12 Patient Monitor with RESP module for monitoring. Therefore, it can be considered safe to use the Mindray RESP module for related parameter monitoring. Clinical performance of RESP function meets SOTA benchmark.

**CPR Feedback:**

In late 2017, the Emergency Department of Shenzhen Second People's Hospital used an open, paired and equivalent verification method. On the Laerdal simulation people, the subject pasted Zoll real-CPR onto the compress point of the simulation people, onto which the Mindray's CPR sensor was then superimposed. Compression data of total 69 subjects were collected. All data analysis sets were consistent. This study had no expulsion or exclusion cases. All included subjects completed the study. Demographic indicators including age and gender met design expectations. All subjects had received professional CPR training.

**Study Product**

The CPR feedback performance test report of BeneHeart D3 is applicable to BeneHeart D6 for the following reasons:

The CPR sensor is used to measure the compression depth, compression rate and interruption time. The measurement performance is determined by the CPR sensor itself, independent of the defibrillator (D3 and D6), D3 and D6 are only used to display the measure results. So, the evaluation results for connecting D3 also apply to D6.

**Control Product**

Primary reference device: Laerdal simulator equipped with Laerdal tablet PC.

Secondary reference device: Zoll real-CPR equipped with Zoll X series defibrillation monitor.

**Primary endpoints are:**

Statistics of compress depth

Device	Average	Standard deviation	95%CI	Equivalent boundary value interval	Conclusion
Primary reference device - device under test	-1.27	3.14	-1.44, -1.10	(-5, +5)	Meet equivalence requirements

Statistics of compress rate

Device	Average	Standard deviation	95%CI	Equivalent boundary value interval	Conclusion
Primary reference device - device under test	-0.68	3.01	-0.84, -0.52	(-2, +2)	Meet equivalence requirements



test					
------	--	--	--	--	--

Note: The measured value of compress rate and compress depth is the difference between the primary reference device and the device under test.

**Secondary endpoints are:**

Statistics of compress depth

Device name	Average deviation	Standard deviation	Result
Primary reference device - device under test	1.76	2.23	1.76 < 3.18, 2.23 < 3.94, therefore, the device under test outweighs the secondary reference device
Primary reference device - secondary reference device	3.18	3.94	

Statistics of compress rate

Device name	Average deviation	Standard deviation	Result
Primary reference device - device under test	1.54	2.15	1.54 < 2.05, 2.15 < 3.35, therefore, the device under test outweighs the secondary reference device
Primary reference device - secondary reference device	2.05	3.35	

According to principle and design analysis, lab test report and compress comparison test on a simulator by clinical medical staff for the CPR Feedback of a defibrillation monitor, we can fully prove that the CPR Feedback function of Mindray defibrillator is designed to expectations, product performance is equivalent to the simulator and can be ready for clinical use.

### 5.3 Summary of clinical data from other sources(if applicable)

#### Post-clinical follow-up from China:

##### External Defibrillation

##### ➤ Objective

To demonstrate the performance and safety of external defibrillation function, identifying previously unknown side-effects and ensuring the continued acceptability of the benefit-risk ratio.

##### ➤ General clinical Information

##### 1. Subject Selection

- Ages eligible for selection: Adult and pediatric.

- Sexes eligible for selection: All
- Subjects have already used the External Defibrillation function of Mindray defibrillator monitors

## 2. Sample Size Calculation

This PMCF adopts the single-group design method, the OPC (Objective performance criteria) is set to 87.9 % (P0), Pt set to 95%. We also determine  $\alpha=0.05$ ,  $\beta=0.2$ , power  $(1-\beta)=0.8$ , clinical drop-out rate  $\gamma=0.05$ . Based on the Formula for Estimating Sample Size of a Single Group Target Value Test, a minimum of 142 clinical cases are expected to be collected in this PMCF.

## 3. Process

Mindray will continuously collect clinical cases from at least 2 hospitals. The clinical cases will be recorded and filled in the PMCF form by the medical staff according to the actual treatment conditions. Submit the signed PMCF form to Mindray for confirmation and statistics.

### ➤ Adopted Statistical Analysis Method and Acceptance Criteria

#### 1. Statistical Analysis Method

The defibrillation success rate equals the total number of success cases divided by the total number of tracked cases. Descriptive statistics of baseline parameters such as patient age, gender are also collected.

#### 2. Acceptance Criteria

The final result of defibrillation success rate is obtained and its two-sided 95% confidence interval is calculated. If the lower limit of the confidence interval is above the target value 87.9%, the clinical performance of the product under application shall be proved to meet the requirement of recognized level in the industry.

### ➤ Result

174 clinical cases have been collected from 50 hospitals until 2022, with details as follows:

#### Gender Distribution of the Cases:

Gender	Number of Cases
Male	125
Female	49
Total	174

#### Age Distribution of the Cases:

populations	Age	Number of Cases
pediatric	1 month -<3 year	2
	3-<12 year	2
	12-<18 year	3
adult	18-<40 year	17
	40-<60 year	60
	≥60 year	90
Total		174

**Indications Distribution of the Cases:**

Indications	Number of Cases
Ventricular fibrillation	148
Pulseless ventricular tachycardia	19
Others	7
Total	174

**Successful/Failed Cases:**

Successful Cases	Failed Cases	Total	success rate
164	8	174	94.3%

This PMCF collects 174 cases and 164 cases are successful, the final result of defibrillation success rate is **94.3%**. Based on Confidence interval calculation formula:

$$p \pm z_{\alpha/2} \sqrt{p(1-P)/n}$$

When  $\alpha=0.05$ ,  $Z_{\alpha/2} = 1.96$

Finally, the two-sided 95% confidence interval of 94.3% is calculated as follows: **(90.9%, 97.7%)**, the lower limit **90.9%** is **above** the target value **87.9%**.

The final result of defibrillation success rate is 94.3% (95%CI: 90.9%, 97.7%). the lower limit 90.9% is above the target value **87.9%**, the clinical performance of the external defibrillation function of the product in question has been proved to meet the requirement of recognized level in the industry.

No adverse events related to the device are found through surveys.

**Synchronized Cardioversion****➤ Objective**

To demonstrate the performance and safety of synchronized cardioversion function, identifying previously unknown side-effects and ensuring the continued acceptability of the benefit-risk ratio.

**➤ General clinical Information****1. Subject Selection**

- Ages eligible for selection: Adult and pediatric.
- Sexes eligible for selection: All
- Subjects have already used the Synchronized Cardioversion function of Mindray defibrillator monitors

**2. Sample Size Calculation**

This PMCF adopts the single-group design method, the OPC (Objective performance criteria) is set to 88.9% (P0), Pt set to 95%. We also determine  $\alpha=0.05$ ,  $\beta=0.2$ , power  $(1-\beta)=0.8$ , clinical drop-out rate  $\gamma=0.05$ . Based on the Formula for Estimating Sample Size of a Single Group Target Value Test, a minimum of 181 clinical cases are expected to be collected in this PMCF.

### 3. Process

Mindray will continuously collect clinical cases from at least 2 hospitals. The clinical cases will be recorded and filled in the PMCF form by the medical staff according to the actual treatment conditions. Submit the signed PMCF form to Mindray for confirmation and statistics.

#### ➤ Adopted Statistical Analysis Method and Acceptance Criteria

##### 1. Statistical Analysis Method

The cardioversion success rate equals the total number of success cases divided by the total number of tracked cases. Descriptive statistics of baseline parameters such as patient age, gender are also collected.

##### 2. Acceptance Criteria

The final result of cardioversion success rate is obtained and its two-sided 95% confidence interval is calculated. If the lower limit of the confidence interval is above the target value 88.9%, the clinical performance of the product under application shall be proved to meet the requirement of recognized level in the industry.

#### ➤ Result

185 clinical cases have been collected from 4 hospitals until 2022, with details as follows:

##### Gender Distribution of the Cases:

Gender	Number of Cases
Male	104
Female	81
Total	185

##### Age Distribution of the Cases:

populations	Age	Number of Cases
pediatric	1 month -<3 year	2
	3-<12 year	4
	12-<18 year	20
adult	18-<40 year	27
	40-<60 year	64
	≥60 year	68
Total		185

##### Indications Distribution of the Cases:

Indications	Number of Cases
Atrial fibrillation	130
Others	55
Total	185

**Successful/Failed Cases:**

Successful Cases	Failed Cases	Total	success rate
175	10	185	94.6%

This PMCF collects 185 cases and 175 cases are successful, the final result of defibrillation success rate is **94.6%**. Based on Confidence interval calculation formula:

$$p \pm z_{\alpha/2} \sqrt{p(1-P)/n}$$

When  $\alpha=0.05$ ,  $Z_{\alpha/2} = 1.96$

Finally, the two-sided 95% confidence interval of 94.6% is calculated as follows: **(91.4%, 97.8%)**, the lower limit **91.4%** is **above** the target value **88.9%**.

The final result of defibrillation success rate is 94.6% (95%CI: 91.4%, 97.8%). the lower limit 91.4% is above the target value **88.9%**, the clinical performance of the external defibrillation function of the product in question has been proved to meet the requirement of recognized level in the industry.

No adverse events related to the device are found through surveys.

**Internal Defibrillation****➤ Objective**

To demonstrate the performance and safety of internal defibrillation function, identifying previously unknown side-effects and ensuring the continued acceptability of the benefit-risk ratio.

**➤ General clinical Information****1. Subject Selection**

- Ages eligible for selection: Adult and pediatric.
- Sexes eligible for selection: All
- Subjects have already used the Internal Defibrillation function of Mindray defibrillator monitors

**2. Sample Size Calculation**

This PMCF adopts the single-group design method, the OPC (Objective performance criteria) is set to 94.2 % (P0), Pt set to 95%. We also determine  $\alpha=0.05$ ,  $\beta=0.2$ , power  $(1-\beta)=0.8$ , clinical drop-out rate  $\gamma=0.05$ . Based on the Formula for Estimating Sample Size of a Single Group Target Value Test, a minimum of 135 clinical cases are expected to be collected in this PMCF.

**3. Process**

Mindray will continuously collect clinical cases from at least 2 hospitals. The clinical cases will be recorded and filled in the PMCF form by the medical staff according to the actual treatment conditions. Submit the signed PMCF form to Mindray for confirmation and statistics.

**➤ Adopted Statistical Analysis Method and Acceptance Criteria****1. Statistical Analysis Method**

The defibrillation success rate equals the total number of success cases divided by the total number of tracked cases. Descriptive statistics of baseline parameters such as patient age, gender are also collected.

## 2. Acceptance Criteria

The final result of defibrillation success rate is obtained and its two-sided 95% confidence interval is calculated. If the lower limit of the confidence interval is above the target value 94.2%, the clinical performance of the product under application shall be proved to meet the requirement of recognized level in the industry.

### ➤ Result

176 clinical cases have been collected from 4 hospitals until 2022, with details as follows:

#### Gender Distribution of the Cases:

Gender	Number of Cases
Male	125
Female	52
Total	177

#### Age Distribution of the Cases:

populations	Age	Number of Cases
pediatric	1 month -< 3 year	2
	3-<12 year	1
	12-<18 year	21
adult	18-<40 year	50
	40-<60 year	40
	≥60 year	63
Total		177

#### Indications Distribution of the Cases:

Indications	Number of Cases
Ventricular fibrillation	141
Ventricular tachycardia	13
Others	23
Total	177

#### Successful/Failed Cases:

Successful Cases	Failed Cases	Total	success rate
173	4	177	97.7%

This PMCF collects 177 cases and 173 cases are successful, the final result of defibrillation success rate is **97.7%**. Based on Confidence interval calculation formula:

$$p \pm z_{\alpha/2} \sqrt{p(1-P)/n}$$

When  $\alpha=0.05$ ,  $Z_{\alpha/2} = 1.96$

Finally, the two-sided 95% confidence interval of 97.7% is calculated as follows:(**95.5%**, **99.9%**), the lower limit **95.5%** is **above** the target value 94.2%.

The final result of defibrillation success rate is 97.7% (95%CI: 95.5%, 99.9%). the lower limit 95.5% is above the target value 94.2%, the clinical performance of the Internal Defibrillation function of the product in question has been proved to meet the requirement of recognized level in the industry.

### **Semi-automated external defibrillation (AED)**

#### ➤ **Objective**

To demonstrate the performance and safety of Semi-automated external defibrillation function, identifying previously unknown side-effects and ensuring the continued acceptability of the benefit-risk ratio.

#### ➤ **General clinical Information**

##### **1. Subject Selection**

- Ages eligible for selection: Adult and pediatric.
- Sexes eligible for selection: All
- Subjects have already used the Semi-automated external defibrillation function of Mindray defibrillator monitors

##### **2. Sample Size Calculation**

###### **For the defibrillation success rate**

This PMCF adopts the single-group design method, the OPC (Objective performance criteria) is set to 87.9 % (P0), Pt set to 95%. We also determine  $\alpha=0.05$ ,  $\beta=0.2$ , power  $(1-\beta)=0.8$ , clinical drop-out rate  $\gamma=0.05$ . Based on the Formula for Estimating Sample Size of a Single Group Target Value Test, a minimum of 142 clinical cases are expected to be collected in this PMCF.

###### **For the sensitivity**

This PMCF adopts the single-group design method, the OPC (Objective performance criteria) is set to 90 % (P0), Pt set to 95%. We also determine  $\alpha=0.05$ ,  $\beta=0.2$ , power  $(1-\beta)=0.8$ , clinical drop-out rate  $\gamma=0.05$ . Based on the Formula for Estimating Sample Size of a Single Group Target Value Test, a minimum of 251 clinical cases are expected to be collected in this PMCF.

###### **For the specificity**

This PMCF adopts the single-group design method, the OPC (Objective performance criteria) is set to 95 % (P0), Pt set to 99%. We also determine  $\alpha=0.05$ ,  $\beta=0.2$ , power  $(1-\beta)=0.8$ , clinical drop-out rate  $\gamma=0.05$ . Based on the Formula for Estimating Sample Size of a Single Group Target Value Test, a minimum of 172 clinical cases are expected to be collected in this PMCF.

##### **3. Process**

Mindray will continuously collect clinical data of Semi-automated external defibrillation function.

Mindray exports the clinical data from the device for analysis and statistics.

#### ➤ **Adopted Statistical Analysis Method and Acceptance Criteria**

## 1. Statistical Analysis Method

According to the PMCF Plan, this clinical follow-up mainly verifies the defibrillation success rate and sensitivity/ specificity of the device in question, and must include at least 180 shockable rhythms cases and 10 non-shockable rhythms/conditions cases.

The defibrillation success rate equals the total number of success cases divided by the total number of tracked cases. Descriptive statistics of baseline parameters such as patient age, gender are also collected.

The sensitivity of device is the number of true positive shockable rhythms that have been correctly classed as shockable, expressed as a percentage of the total number of shockable rhythms:  $A/(A+C)*100\%$ .

The specificity is the number of organized or perfusing rhythms that have been correctly classed as non-shockable rhythms/conditions by the algorithm, and expressed as a percentage of the total number of non-shockable rhythms/conditions:  $D/(B+D)*100\%$ .

## 2. Acceptance Criteria

The final result of defibrillation success rate is obtained and its two-sided 95% confidence interval is calculated. If the lower limit of the confidence interval is above the target value 87.9%, the clinical performance of the product under application shall be proved to meet the requirement of recognized level in the industry.

The final result of sensitivity and specificity is obtained and its two-sided 95% confidence interval is calculated. If the lower limit of the sensitivity confidence interval is above the target value 90% and the lower limit of the specificity confidence interval is above the target value 95%, the clinical performance of the product in question shall be proved to meet the requirement of recognized level in the industry.

## ➤ Result

220 clinical cases have been collected until 2022, with details as follows:

### Gender Distribution of the Cases:

Gender	Number of Cases
Male	52
Female	168
Total	220

### Age Distribution of the Cases:

populations	Age	Number of Cases
pediatric	1 month -< 3 year	7
	3-<12 year	23
	12-<18 year	36
adult	18-<40 year	54
	40-<60 year	67
	≥60 year	33



Total	220
-------	-----

**Defibrillation success rate:**

Number of total Cases	success cases	success rate
150	140	93.3%

**Positive: shockable rhythms cases. Negative: non-shockable rhythms/conditions cases.**

AED	results of Expert diagnostic with ECG		Total
	Positive	Negative	/
Positive	419	16	435
Negative	19	1205	1224
Total	438	1221	1659

**Sensitivity**

$$A/(A+C)*100\%=95.7\%$$

**Specificity**

$$D/(B+D)*100\%=98.7\%$$

**95% Confidence Interval Calculation:**

For the defibrillation success rate of the device in question:

This PMCF collects 150 cases and 140 cases are successful, the final result of defibrillation success rate is **93.3%**. Based on Confidence interval calculation formula:

$$p \pm z_{\alpha/2} \sqrt{p(1-P)/n}$$

$$\text{When } \alpha=0.05, Z_{\alpha/2} = 1.96$$

Finally, the two-sided 95% confidence interval of 92.9% is calculated as follows: **(89.3%, 97.3%)**, the lower limit **89.3%** is **above** the target value **87.9%**.

For the sensitivity of the device in question:

This PMCF collects 438 shockable rhythms cases and 419 correct identification cases, the final result of sensitivity is **95.7%**. Based on Confidence interval calculation formula:

$$p \pm z_{\alpha/2} \sqrt{p(1-P)/n}$$

$$\text{When } \alpha=0.05, Z_{\alpha/2} = 1.96$$

Finally, the two-sided 95% confidence interval of 95.7% is calculated as follows: **(93.9%, 97.5%)**, the lower limit **93.9%** is **above** the target value **90%**.

For the specificity of the device in question:

This PMCF collects 1221 non-shockable rhythms/conditions cases and 1205 correct identification cases, the final result of specificity is **98.7%**. Based on Confidence interval calculation formula:

$$p \pm z_{\alpha/2} \sqrt{p(1-P)/n}$$

$$\text{When } \alpha=0.05, Z_{\alpha/2} = 1.96$$

Finally, the two-sided 95% confidence interval of 98.7% is calculated as follows: **(98.1%, 99.3%)**, the

lower limit **98.1%** is **above** the target value **95%**.

➤ **Conclusion:**

The result of defibrillation success rate is 93.3% (95%CI: **89.3%**, **97.3%**), the lower limit **89.3%** is upper than the target value **87.9%**.

The result of sensitivity rate is 95.7% (95%CI: **93.9%**, **97.5%**), the lower limit **93.9%** is upper than the target value **90%**.

The result of specificity rate is 98.7% (95%CI: **98.1%**, **99.3%**), the lower limit **98.1%** is upper than the target value **95%**.

The clinical performance of the AED function of the product in question has been proved to meet the requirement of recognized level in the industry.

### **Non-invasive external pacing**

➤ **PMCF Objective**

To demonstrate the performance and safety of Non-invasive external pacing function, identifying previously unknown side-effects and ensuring the continued acceptability of the benefit-risk ratio, post-marked clinical follow-up activity have been undertaken by manufacture.

➤ **Description of PMCF survey**

To demonstrate the performance and safety of AED function, identifying previously unknown side-effects and ensuring the continued acceptability of the benefit-risk ratio, Mindray did post-market follow-up survey activity.

This Post Market clinical follow-up survey was an observational, post-market study to monitor and learn more about the use of defibrillator/monitor, whose sample size criteria based from the endpoint of state-of-the-art. This study collected information from health care professional user who used BeneHeart D5/D6 defibrillator/monitor.

216 clinical cases have collected from 2 hospitals by the BeneHeart D5/D6 Defibrillator/Monitor until 2022.

➤ **General clinical Information**

1. Subject Selection

- Ages eligible for selection: Adult and pediatric.
- Sexes eligible for selection: All
- Subjects have already used the Non-invasive pacing function of Mindray defibrillator monitors

2. Sample Size Calculation

This PMCF adopts the single-group design method, the OPC (Objective performance criteria) is set to 58.9 % (P0), Pt set to 95%. We also determine  $\alpha=0.05$ ,  $\beta=0.2$ , power  $(1-\beta)=0.8$ , clinical drop-out rate  $\gamma=0.05$ . Based on the Formula for Estimating Sample Size of a Single Group Target Value Test, a minimum of 190 clinical cases are expected to be collected in this PMCF.

3. Process

Mindray will continuously collect clinical cases from at least 2 hospitals. The clinical cases will be

recorded and filled in the PMCF form by the medical staff according to the actual treatment conditions. Submit the signed PMCF form to Mindray for confirmation and statistics.

➤ Adopted Statistical Analysis Method and Acceptance Criteria

1. Statistical Analysis Method

The pacing success rate equals the total number of success cases divided by the total number of tracked cases. Descriptive statistics of baseline parameters such as patient age, gender are also collected.

2. Acceptance Criteria

The final result of pacing success rate is obtained and its two-sided 95% confidence interval is calculated. If the lower limit of the confidence interval is above the target value 58.9%, the clinical performance of the product under application shall be proved to meet the requirement of SOTA criteria in the industry.

➤ Result

216 clinical cases has collected from 2 hospitals until 2022, with details as follows:

**Gender Distribution of the Cases**

Gender	Number of Cases
Male	139
Female	77
Total	216

**Age Distribution of the Cases**

populations	Age	Number of Cases
pediatric	1month-< 3 year	1
	3-<12 year	1
	12-<18 year	2
adult	18-<40 year	12
	40-<60 year	78
	≥60 year	122
Total		216

**Indications Distribution of the Cases:**

Indications	Number of Cases
Bradycardia	148
Asystole	68
Total	216

**Results**

Successful Cases	Failed Cases	Total	success rate
184	32	216	85.2%

This PMCF collects 216 cases and 172 cases are successful, the final result of pacing success rate is **85.2%**.

Based on Confidence interval calculation formula:

$$p \pm z_{\alpha/2} \sqrt{p(1-P)/n}$$

When  $\alpha=0.05$ ,  $Z_{\alpha/2} = 1.96$

Finally, the two-sided 95% confidence interval of 85.2% is calculated as follows: **(80.5%, 89.9%)**, the lower limit **80.5%** is **above** the target value **58.9%**.

The final result of pacing success rate is 85.2% (95%CI: 80.5%, 89.9%). the lower limit **80.5%** is above the target value **58.9%**, the clinical performance of the Non-invasive external pacing of the product in question has been proved to meet the requirement of SOTA criteria in the industry.

## 5.4 An overall summary of the clinical performance and safety

Defibrillator/monitor is intended to be used for monitoring, defibrillator and alarming. When unsatisfactory physiological PATIENT states, unsatisfactory functional states of the patient monitor or hazards to the patient or operator due to the patient monitor exists, the Defibrillator/monitor detects alarm conditions and generate alarm signals. For a patient with ventricular fibrillation (VF)/ventricular tachycardia (VT), the effective and timely two-way defibrillation treatment presents high success rate of treatment.

### External Defibrillation

The defibrillation success rate is the primary endpoint, a total of 174 cases of the use of the BeneHeart series of defibrillation monitors were collected until 2022. 164 cases are successful, the final result of defibrillation success rate is 94.3% (95%CI: 90.9%, 97.7%). The lower limit 90.9% is above the target value 87.9%, the clinical performance of the external defibrillation function of the product in question has been proved to meet the requirement of recognized level in the industry.

### Synchronized Cardioversion

The cardioversion success rate is the primary endpoint, a total of 185 cases of the use of the BeneHeart series of defibrillation monitors were collected until 2022. 175 cases are successful, the final result of defibrillation success rate is 94.6% (95%CI: 91.4%, 97.8%). The lower limit 91.4% is above the target value 88.9%, the clinical performance of the Synchronized Cardioversion function of the product in question has been proved to meet the requirement of recognized level in the industry.

### Internal Defibrillation

The defibrillation success rate is the primary endpoint, a total of 176 cases of the use of the BeneHeart series of defibrillation monitors were collected until 2022. 173 cases are successful, the final result of defibrillation success rate is 97.7% (95%CI: 95.5%, 99.9%). The lower limit 95.5% is above the target value 94.2%, the clinical performance of the Internal Defibrillation function of the product in question has been proved to meet the requirement of recognized level in the industry.

### Semi-automated external defibrillation (AED)

The primary endpoint is the defibrillation success rate and sensitivity/specificity, a total of 220 cases of the use of the BeneHeart series of defibrillation monitors were collected until 2022. 140 cases are successful, the result of defibrillation success rate is 93.3% (95%CI: 89.3%, 97.3%), the lower limit 89.3% is upper than the target value 87.9%.

The result of sensitivity rate is 95.7% (95%CI: 93.9%, 97.5%), the lower limit 93.9% is upper than the target value 90%.

The result of specificity rate is 98.7% (95%CI: 98.1%, 99.3%), the lower limit 98.1% is upper than the target value 95%.

The clinical performance of the AED function of the product in question has been proved to meet the requirement of recognized level in the industry.

**Non-invasive external pacing**

The defibrillation success rate is the primary endpoint, a total of 216 cases of the use of the BeneHeart series of defibrillation monitors were collected until 2022. 172 cases are successful, the final result of defibrillation success rate is 85.2% (95%CI: 80.5%, 89.9%). The lower limit 80.5% is above the target value 58.9%, the clinical performance of the Non-invasive external pacing of the product in question has been proved to meet the requirement of SOTA criteria in the industry.

Monitor and defibrillator technology has a long history of acceptable and well-understood performance and risk. It can compatible with a high level of protection of health and safety and acceptable according to current knowledge/the state of the art. Based on the clinical data of the equivalent equipment and similar equipment, there are no residual risks and uncertainties or unanswered questions identified, the device is the useful medical equipment that the benefit to weight against the risk of the identified hazards is acceptable.

According to the conclusion of residual risks evaluation, for the intended patients, there is no known side effect that can occur during or after the use of the medical device, no extra preparation should be made. So no residual risk associated with using the medical device should be disclosed.

Factor	Question to consider	Notes
<b>Assessment of Benefits of Devices</b>		
<p><b>1. Type of benefit(s)</b></p>	<ul style="list-style-type: none"> <li>- What primary endpoints or surrogate endpoints were evaluated?</li> <li>- What key secondary endpoints or surrogate endpoints were evaluated?</li> <li>- What value do patients place on the benefit?</li> </ul>	<p>monitors are intended to be used for monitoring, displaying, reviewing, storing, alarming and transferring of multiple physiological parameters.</p> <p>Patient monitor improves the efficiency of medical staff in monitoring the physiological parameters of patients</p> <p>Patients put a high value on this treatment because it has the potential to save their lives. Patients are therefore willing to accept the risks of this treatment to achieve the benefit. If the treatment provides timely successful defibrillation, the patient will survive a life threatening cardiac arrest situation and will be able to seek further treatment.</p>
<p><b>2. Magnitude of the benefit(s)</b></p>	<ul style="list-style-type: none"> <li>- For each primary and secondary endpoint or surrogate endpoints evaluated:</li> </ul>	<p>The patient monitor complies with the standard test, bench test and usability test to approve the accuracy during monitoring</p>

	<p>o What was the magnitude of each treatment effect?</p> <p>- What scale is used to measure the benefit?</p> <p>o How did the benefit rank on that scale?</p>	<p>According to the AHA Recommendations</p> <p>Defibrillators are recommended to treat tachyarrhythmias requiring a shock,the Recommendations level is CLASS I (STRONG) Benefit &gt;&gt;&gt; Risk,based on B-NR LEVEL (QUALITY) OF EVIDENCE;</p> <p>Defibrillators using biphasic waveforms are preferred over monophasic defibrillators for treatment of tachyarrhythmias , the Recommendations level is CLASS IIa (MODERATE) Benefit &gt;&gt; Risk,based on LEVEL B-R (Randomized).</p>
<p><b>3. Probability of the patient experiencing one or more benefit(s)</b></p>	<p>- Was the study able to predict which patients will experience a benefit?</p> <p>- What is the probability that a patient for whom the device is intended will experience a benefit?</p> <p>- How did the benefits evaluated vary across sub-populations? (If the study was sufficiently powered for subpopulations, note specific subpopulations, nature of difference and any known reasons for these differences.)</p> <p>- Was there a variation in public health benefit for different populations?</p> <p>- Even if the benefit is in a small portion of the population, do those patients who would</p>	<p>All the parameters can be monitored on single adult, pediatric, and neonatal patients with the exception of the following:</p> <ul style="list-style-type: none"> <li>■ The CCO and NMT monitoring are intended for adult and pediatric patients only.</li> <li>■ C.O. monitoring is only intended for adult patients.</li> </ul> <p><b>a) Can this study predict which patients will benefit?</b></p> <p>Sudden death is the first manifestation of cardiovascular disease. Prompt discovery and shock defibrillation and cardiopulmonary resuscitation can save a significant proportion of sudden deaths; The probability of patient survival decreases by about 7 %~10%% every 1 min delay from falling to defibrillation.</p> <p><b>b) What is the probability of the expected benefit to the patient?</b></p> <p>The successful definition defined in the AHA cardiopulmonary resuscitation guide in 2015 is that ventricular fibrillation is terminated within 5s after electric shock is applied. When the defibrillation success rate is in the range of 85% to 95%, it is considered that the product design meets expectations.</p> <p><b>c) Are public health benefits vary among different groups?</b></p> <p>Ventricular fibrillation can occur in any population, and the most immediate treatments for ventricular fibrillation are</p>

	experience the benefit value it?	defibrillation and cardiopulmonary resuscitation. There is no evidence from the AHA guidelines that de<unk> brillation outcomes vary between populations with regard to ventricular fibrillation.
<b>4. Duration of effect(s)</b>	<ul style="list-style-type: none"> <li>- Could the duration, if relevant, of each treatment effect, including primary and secondary endpoints be determined? If so, what was it?</li> <li>- Is the duration of the benefit achieved of value to patients?</li> </ul>	<p>Patients can be continuously monitored during the hospitalization period and benefit continuously</p> <p>If the treatment provides timely successful defibrillation, the patient will survive a life threatening cardiac arrest situation and will be able to seek further treatment.</p>
<b>Assessment of risks of Devices</b>		
<b>5. Severity, types, number and rates of harmful events (events and consequences):</b>		
· Device-related serious adverse events	- What are the device-related serious adverse events for this product?	According to section 4.3, No serious adverse events
· Device-related non-serious adverse events	- What are the device-related non-serious adverse events for this product?	According to section 4.3, the non-serious adverse events had been all solved and we did risk control measurement.
· Procedure-related complications	- What other procedure-related complications may a patient be subject to?	Not found yet
<b>6. Probability of a harmful event</b>	<ul style="list-style-type: none"> <li>- What percent of the intended patient population would expect to experience a harmful event?</li> <li>- What is the incidence of each harmful event in the study population?</li> <li>- How much uncertainty is in that estimate?</li> <li>- How does the incidence of harmful events vary by subpopulation (if applicable)?</li> </ul>	<p>a) <b>What is the probability of adverse events in the intended users?</b></p> <p>Since the market release of Mindray D series defibrillators, 131917 sets of them have been sold, and 25 sets of them have been reported to have the adverse events mentioned above. So the probability of adverse events is 25/131917</p> <p>b) <b>Is the patient willing to accept the risk of possible adverse events while considering the possible benefits of the device?</b></p> <p>The magnitude of this benefit is either life or death. So, patients put a high value on this treatment because it has the potential to save their lives. Patients are</p>



	<ul style="list-style-type: none"> <li>- Are patients willing to accept the probable risk of the harmful event, given the probable benefits of the device?</li> </ul>	therefore willing to accept the risks of this treatment to achieve the benefit.
<b>7. Duration of harmful events</b>	<ul style="list-style-type: none"> <li>- How long does the harmful event last?</li> <li>- Is the harmful event reversible?</li> <li>- What type of intervention is required to address the harmful event?</li> </ul>	<p>a) <b>Is the adverse event reversible?</b> All adverse events are reversible and not harmful to human body. They are temporary effects.</p> <p>b) <b>What measures should be taken for adverse events?</b> Mindray PMS will regularly collect, analyze, correct and prevent post-marketing adverse events.</p>
<b>8. Risk from false-positive or false-negative results for diagnostics</b>	<ul style="list-style-type: none"> <li>- What are the consequences of a false positive?</li> <li>- What are the consequences of a false negative?</li> <li>- Is this the only means of diagnosing the problem, or is it part of an overall diagnostic plan</li> </ul>	According to section 4.4.1, 4.4.2, 4.4.3, the consequences of a false positive or false negative are identified and verifcated with safety test, bench test, usability test risk control measurement
<b>Conclusion</b>	<p>Through the parameter monitoring, medical staff have established sufficient conditions to provide patients with a better medical monitoring environment, and the benefits are obvious.</p> <p>Although there is also the possibility of false positives and false negatives in parameter monitoring, the impact of false positives and false negatives is limited and will not cause substantial harm to patients.</p> <p>In addition, the parameter monitoring of the patient monitor has the advantages of simplicity of equipment, convenient operation, timeliness, economy, etc. compared with other known ones.</p> <p>Therefore, from the perspective of benefit and risk, patient monitor parameter monitoring has obvious benefits, controllable risks, and has strong clinical application popularization characteristics.</p> <p>Defibrillator are life-saving devices used in emergency situations. They have been shown to have a high benefit for patients with underlying diseases that remain undetected until sudden cardiac arrest occurs. The time from collapse to defibrillation is critical in-patient survival. For every minute that passes between collapse and defibrillation, survival rates from VF SCA decrease 7% to 10%.</p>	

	In conclusion, given the available information above, the defibrillator's support for patients in cardiac arrest who are unconscious, not breathing, or without circulation the probable benefits outweigh the probable risks.
--	--

## 5.5 Ongoing or planned post-market clinical follow-up

After the product is marketed in the Europe, Mindray will continue to track the success rate of defibrillators to ensure the long-term safety of the product. In the CE region, China, and ROW, we will select several hospitals for follow-up (the content here will be refreshed after the hospital are confirmed). The PMCF is organized and carried out by Mindray's clinical engineers, who will control the quality and progress of data collection. During the clinical follow-up process, we will summarize the clinical follow-up results every six months until the collection of cases is completed to ensure that the product is safe and effective in the clinic.

If there are any emerging risks, complications or unexpected device failures have been detected, the summary of safety and clinical performance shall be updated throughout the life cycle of the device.

## 6. Possible diagnostic or therapeutic alternatives

Function	Manufacturer's	Alternatives
External Defibrillation	biphasic waveforms defibrillators	Monophasic waveform defibrillators
Synchronous cardioversion	Electrical cardioversion	1. pharmacological cardioversion 2. radiofrequency ablation
Internal Defibrillation	internal defibrillation	None
Semi-automated external defibrillation (AED)	Semi-automated defibrillation	manual defibrillation
Non-invasive external pacing	external pacing	1. transcutaneous pacing 2. transvenous pacing 3. percussion pacing 4. epicardial pacing
CPR feedback	CPR feedback	CPR metronome
ECG	ECG	1. PR (SPO2) 2. Ultrasound Cardiogram

		(UCG)
IBP	Invasive Blood Pressure	Noninvasive Blood Pressure
TEMP	Continuous body temperature monitoring	None
CO2	Sidestream CO2 Mainstream CO2	1. colorimetric CO2 2. PaCO2
RESP	ECG-derived respiration	1. Visual respiratory assessment 2. Acoustic respiratory assessment 3. PPG-derived respiration
SPO2	pulse oximetry (SpO2)	arterial oxygen saturation (SaO2)
NIBP	Oscillometric blood pressure	1. Auscultatory blood pressure 2. invasive blood pressure (IBP)

**External Defibrillation**

Modern defibrillators are classified according to 2 types of waveforms: monophasic and biphasic. Monophasic waveform defibrillators were introduced first, but biphasic waveforms are used in almost all AEDs and manual defibrillators sold today. Energy levels vary by type of device and manufacturer. [1-3]

Defibrillators with monophasic waveform deliver current in one polarity and were the first to be introduced. They can be further categorized by the rate at which the current pulse decreases to zero. If the monophasic waveform falls to zero gradually, the term damped sinusoidal (MDS) is used. If the waveform falls instantaneously, the term truncated exponential (MTE) is used (figure 1) [1-8].

Few monophasic waveform defibrillators are being manufactured, but many are still in use, and most use MDS waveforms. As noted above, no specific waveform characteristic (either monophasic or biphasic) is consistently associated with a greater incidence of ROSC or higher survival to hospital discharge rates after cardiac arrest [1-8].




Figure 1. Monophasic waveforms. A. Damped sinusoidal wave (A) and truncated exponential (B).

Biphasic waveform defibrillator was developed later. The delivered current flows in a positive direction for a specified time and then reverses and flows in a negative direction for the remaining duration of the electrical discharge (figure 2). With biphasic waveforms there is a lower defibrillation threshold (DFT) that allows reductions of the energy levels administered and may cause less myocardial damage. The use of biphasic waveforms permits a reduction in the size and weight of defibrillators. <sup>[1-8]</sup>




Figure 2. A. Biphasic waveform.

	<b>Biphasic</b>	<b>Monophasic</b>
Volume and Weight	lighter and smaller	- Heavy, Unportable
Current market share	- Majority MEDs - Almost all AEDs	- small part of historical inventory
Easy of use	- Simple, easy	- old products, not updated - old technology, difficult to use
shock success rate	- equivalent or higher	/
myocardial damage	- less	/
energy	- 200 J to 360 J escalating energy (BTE) for Adult	200 J to 360 J escalating energy (MDS) for adult 200 J to 360 J escalating energy (MTE) for adult

### Synchronous cardioversion

Drug cardioversion or pharmacological cardioversion is to convert an abnormal and potentially dangerous heart rhythm (atrial fibrillation or atrial flutter etc.) into a normal sinus rhythm using oral or intravenous agents (amiodarone, flecainide, ibutilide, propafenone, and vernakalant etc.). Pharmacological cardioversion is a conventional therapy. The main advantages of pharmacological cardioversion are cost effective, avoiding the resource intensive procedural sedation and corresponding risks due to anesthesia procedures comparing to electrical cardioversion. While electrical cardioversion usually has a higher success rate than pharmacological cardioversion in the short-term.

Radiofrequency ablation uses an electric current through wires into the heart to heat up a small area in the atrial chamber tissue to stop it from sending disorder electrical signals to generate the arrhythmias like atrial fibrillation or atrial flutter and restore the rhythm back to normal sinus. Radiofrequency ablation is a more complex and expensive clinical practice than electrical or pharmacological cardioversion. It's usually a scheduled invasive procedure and not suitable for emergency cases unlike cardioversion. Despite being currently the most effective therapeutic option, radiofrequency ablation has not shown desirable results in all patients. The addition of pharmacologic therapy to the routine strategy of catheter ablation enhances the overall success of the procedure.

In clinical practice, electrical or pharmacological cardioversion is a more effective way to restore the patient rhythm in a short-term especially in emergency situations than radiofrequency ablation. For long-term rhythm control, radiofrequency ablation is most effective while the choice of radiofrequency ablation or pharmacologic therapy or a hybrid therapy shall be determined by clinical professionals according to the different conditions of individual patients.

### Internal Defibrillation

Cardiac surgery is usually performed on a nonbeating heart; ventricular fibrillation may occur or may be deliberately induced during operation. If a regular rhythm does not spontaneously resume at the end of operation, an electric shock is applied directly to the fibrillating heart to restore rhythmic contraction. Similarly, during thoracotomy, because traditional external defibrillation pads or paddles cannot be placed, only option is to place defibrillation paddles on the surface of the heart for internal defibrillation.

According to our current cognition, the use of internal defibrillation is the only safe, effective and convenient option for ventricular fibrillation during thoracotomy, cardiac surgery or large abdominal surgery.

### Semi-automated external defibrillation

Manual defibrillation is designed for health professionals, Professional medical staff judge the patient's vital signs parameters by themselves, and complete the entire defibrillation process (ECG analysis, energy adjustment, charging and discharging) as needed. AED are designed to be used by laypersons who ideally should have received AED training at some point in the past. Generally, AED are designed very intuitive and user-friendly so that even untrained bystanders can perfectly employ them to deliver an electric shock to a VF victim. In contrast with AED, the more sophisticated manual defibrillation used by health professionals can perform other functions but require a skilled operator able to interpret electrocardiograms.

	MED	AED
intended users	- health professionals	- laypersons received AED training
Key different Point	- interpret electrocardiograms automatically	- automatic arrhythmia analysis

Easy of use	<ul style="list-style-type: none"> <li>- interpret electrocardiograms</li> <li>- Set energy and charge</li> <li>- discharge</li> </ul>	<ul style="list-style-type: none"> <li>- automatic arrhythmia analysis and charge</li> <li>- press button according to the prompt (Semi-automatic AED) or automatic discharge (Fully automatic AED)</li> </ul>
-------------	--	--

### Non-invasive external pacing

Technological developments have provided transcutaneous, transvenous, percussion, epicardial approaches to temporary cardiac pacing in addition to the refinement of external pacing. All approaches, however, are based on the provision of rate support from an external pulse generator via an electrode or electrodes which can be removed easily after a short period of pacing, as many of the situations requiring temporary pacing are transient and resolve spontaneously or have a correctable underlying cause.

Compare with different pacing mode are as follows:

pacing mode	transcutaneous pacing	transvenous pacing	percussion pacing	epicardial pacing
position	chest wall via adhesive electrodes	pacing wires are inserted into the veins via an introducer sheath, and passed through the venous system to the heart	pacing wires are fixed directly to the myocardium (ventricular and often atrial) and are exposed through the skin on the chest wall	middle or lower two thirds of the patient's sternum
use condition	patients cannot be moved or staff with transvenous pacing experience are not immediately available	patients with persistent hemodynamically unstable bradycardia refractory to medical therapy	temporizing measure in exceptional circumstances such as witnessed, monitored in-hospital arrest (eg, cardiac catheterization laboratory) for bradysystole before a loss of consciousness and if performed without delaying definitive therapy	profound bradycardia resulting in clinical cardiac arrest. p-wave asystole (ventricular standstill).

### CPR feedback

Currently, during CPR, doctors mainly use the following two methods to improve CPR compression

accuracy:

1. Use the CPR metronome guidance to improve CPR compression accuracy;
2. Use the CPR feedback technology to improve CPR compression accuracy.

For the first method, because there is no audiovisual feedback, when the CPR compression rate is inaccurate, the user cannot be reminded, and the method can only help the user improve the accuracy of the compression rate, but cannot improve the accuracy of the compression depth.

The second method, because there is feedback, users can adjust the frequency and depth of compressions in real time through the feedback information and achieve optimal compression.

## ECG

Treatments	Advantages	Disadvantages
ECG	ECG monitoring is the most common and accessible and continuous heart rate monitoring method. It can also provide observation of ECG waveform changes, diagnosis of myocardial damage, myocardial ischemia and electrolyte disturbance. ECG remains the gold standard to continuously monitor an infant's HR in the neonatal intensive care unit (NICU). ECG gets reliable signals more quickly than PR.	Delay in signal acquisition due to skin cleaning (all infants are wet at delivery) and lead placement as well as skin fragility in extremely premature infants. There is also the potential of pulseless electric activity, which could be interpreted as HR on an ECG
PR (SPO <sub>2</sub> )	PR is continuous and visible to the whole team, and more precise than palpation and auscultation.	Delay in HR display sometimes up to 2 min, which could potentially delay resuscitation efforts, and inaccuracy of HR measurements due to motion artefacts or poor tissue perfusion
Ultrasound Cardiogram (UCG)	UCG is faster than palpation, auscultation and PO HR, UCG is more accurate than auscultation and palpation HR, and the Whole team can hear HR	UCG HR is difficult to record, and All measurements were performed by the same examiner Needs an extra staff member. UCG Can interfere with resuscitation efforts in terms of crowding

## IBP

Treatments	Invasive Blood Pressure (IBP)	Noninvasive Blood Pressure (NIBP)
------------	-------------------------------	-----------------------------------

Measurement method	Invasive blood pressure is a blood pressure monitoring method in which the blood pressure is monitored by direct methods by inserting a cannula into a suitable artery.	Noninvasive blood pressure is a way of monitoring the blood pressure indirectly using a special apparatus.
Clinical intervention	Required: the cannula is inserted into a suitable vein during the invasive blood pressure monitoring	Not required: a cuff is used which is wrapped around the arm and connected to a monitor during the noninvasive blood pressure monitoring.
Application	Intensive Care Unit (ICU), Operating theater	Out-patient, tentative diagnosis
Accuracy	Invasive blood pressure monitoring is a highly accurate method.	Noninvasive blood pressure monitoring is a less accurate method
Advantages	Accurate measurement of beat-to-beat pressure fluctuations can be used to monitor the blood pressure of patients in a critical health condition.	Noninvasive hence not prone to infections, or clinical manifestation caused by unsterilized needles.
Disadvantages	Invasive blood pressure method results in harmful side effects due to clinical interventions.	Noninvasive blood pressure monitoring is not very accurate and error-prone.

**TEMP**

The methods used in body temperature measurement are divided into instrumental methods and mathematical methods. The instrumental approach directly utilizes body temperature measurements from temperature-sensitive sensors and electronics by combining actual and predictive measurements, invasive and non-invasive measurements. From the literature we searched, body temperature measurement mainly includes continuous contact measurement methods, as well as discontinuous measurement methods, including digital thermometers, infrared thermometers, zero heat flux thermometers, and thermal imaging. Temperature monitoring technology mainly includes continuous body temperature monitoring technology and spot measurement technology. Continuous body temperature monitoring, as a means of continuous monitoring of human body temperature, continuously measures the temperature. Spot measurement is used for temperature screening, which is suitable for different scenarios and has no possibility of substitution. Therefore, the two measurement techniques will coexist for a long time.

TEMP	Continuous temperature	Point measurement technology			
		digital	infrared	zero heat flux	thermal imaging



Site	rhinitis, esophagus, bladder, rectum	oral 、 armpit 、 rectum	Ear	Forehead	Face
Distance	Contact	Contact	Contact	Contact	10 cm/50 cm
Speed	continuous temperature	6 s	1 s	1 s	1 s
Recordability	Yes	No	No	No	No
Clinical Accuracy	$\pm 0.1$ °C	$\pm 0.1$ °C	$\pm 0.2$ °C	$\pm 0.23$ °C	$\pm 3$ °C

## CO2

Treatment	Sidestream CO2	colorimetric CO2	Mainstream CO2	PaCO2
Measurement method	continuous	continuous	continuous	non-continuous
Invasive	Non- invasive	Non- invasive	Non- invasive	Invasive
Quantitative	Quantitative	Qualitative and semi-quantitative	Quantitative	Quantitative
Reuse	Reuse device	Single-use device	Reuse device	Reuse device
Target patients	Intubated and non-intubated patient	Intubated patient	Usually intubated patient	Intubated and non-intubated patient
Real-time	Several seconds delay	Real-time	Real-time	Time delay
Sampling gas	yes	no	no	no

## RESP

There are some kinds of respiration rate measurement methods, some methods have high accuracy, some methods are easy to use and some methods are inexpensive. All of which have their own advantages and disadvantages.

	Impedance respiration	Visual respiratory assessment	Acoustic respiratory assessment	ECG-derived respiration	PPG-derived respiration	Tidal volume
Accuracy	High	High	High	Medium	Medium	High
Easy of use	Easy to use	Professional staff	Professional staff	Easy to use	Easy to use	Difficult to use
Time-	Less	More	More	Less	Less	Less

consuming						
Monitoring coherence	Continuous	Spot Check	Spot Check	Continuous	Continuous	Continuous
Cost	Less	High	High	Less	Less	High

## SPO2

The primary purpose of intensive care is to ensure adequate oxygen supply to the organ systems, and the assessment of the patient's oxygen status is usually by drawing the patient's arterial blood to calculate oxygen saturation. Clinically, it is also known as the gold standard for SpO<sub>2</sub><sup>[12-2]</sup>.

Compared with SaO<sub>2</sub>, SpO<sub>2</sub> monitoring has obvious advantages. It can display results conveniently, quickly, and accurately, and conduct continuous non-invasive monitoring. In some cases, it has the function of trend review. It is especially suitable for monitoring the critically ill patients in ICU and during and after the operation, helping the medical staff to find the abnormal results in time, actively find the causes, and make corresponding treatment.

	SpO <sub>2</sub>	SaO <sub>2</sub>
Continuity	Continuous	Non-continuous
Operation	Convenience	Complex
Traumatic	non-invasive	invasive
Accuracy	Approximate	Gold standard

## NIBP

Blood pressure can be measured in two ways: Invasive and non-invasive. The invasive blood pressure is mainly measured through an intra-arterial catheter. The NIBP measurement is classified into intermittent NIBP measurement (mainly including manual auscultatory blood pressure and automatic cuff oscillometric blood pressure measurement) and continuous NIBP measurement (such as photoplethysmography (PPG) blood pressure measurement and pulse transmission time (PTT) blood pressure measurement) based on whether NIBP is measured continuously.

In combination with the maturity of the technology, detailed comparisons are made between the oscillometric NIBP measurement, auscultatory NIBP measurement, and IBP measurement. The results are as follows.

Item	Oscillometric blood pressure	Auscultatory blood pressure	IBP
Measurement method	Non-invasive	Non-invasive	Invasive
General measurement media	Cuff-based	Cuff-based	Catheter-based

Measurement Continuity	Intermittent	Intermittent	Continuous
Measurement Process	Manual	Automatic	Automatic
Operation complexity	Difficult	Easy	More Difficult
Target patients	Adults, pediatrics and neonates	Adults, pediatrics and neonates	Adults, pediatrics and neonates
General measurement site	Upper arm	Upper arm	Arterial
Clinical investigation	Not required	Clinical investigation shall be carried out according to ISO 80601-2	Not required
Measurement reliability	Higher	High	Higher


## 7. Suggested profile and training for users

The equipment is intended for use only by clinical professionals who have a clinical education background or under their guidance. It must only be used by qualified medical personnel trained in the operation of the equipment and qualified by training in basic life support, advanced cardiac life support or defibrillation.

For clinical users who install the machine for the first time, the engineer will provide clinical users with practical demonstrations and functional explanations of the machine. In addition to the training of clinical roles, there is also an introduction to the equipment installation function when the engineer is installed, the training content is as follows:


### ➤ Introduction of the main unit of the product.

#### Area 1:




- 1. Alarm lamp    2. Display screen    3. AC power indicator    4. Battery indicator
- 5. Service indicator    6. Soft keys

#### Area 2:




- 1. Silence    2. Lead Select Button    3. 12-Lead ECG button    4. NIBP    5. Speaker
- 6. Main Menu    7. Selector    8. Mark Event Button    9. Microphone

#### Area 3:




- 1. Mode Select knob
- 2. Energy Select button
- 3. Charge button
- 4. Shock button

**Recorder:**




- 1. Start/Stop key
- 2. Indicator
- 3. Paper outlet
- 4. Recorder door
- 5. Latch

**Side View:**



- 1. IBP1: IBP sensor connector (channel 1)
- 2. IBP2: IBP sensor connector (channel 2)
- 3. NIBP: NIBP cuff connector
- 4. Gas outlet
- 5. T1: Temp probe connector (channel 1)
- 6. ECG: ECG cable connector
- 7. SpO2: SpO2 sensor connector
- 8. T2: Temp probe connector (channel 2)
- 9. CO2: sampling line connector (for microstream CO2 module) or watertrap connector (for sidestream CO2 module)

**Rear View:**



1. Hook 2. Battery 2 3. Battery 1 4. External power input 5. Equipotential grounding terminal  
6. USB connector 7. Network connector 8. Multifunctional connector 9. VGA connector

**External Paddles:**



- Shock button with indicator 2. Charge button 3. Energy Select button


**Pads cable, pads, test load:**



➤ **Introduction of the basic operation/setting of the product.**

**Using the Main Menu:**

To enter the main menu, press the Main Menu button on the equipment's front.



### Setting the Date and Time:

1. Press the Main Menu button on the front panel, and then select [Others >>]→[Configuration >>]→enter the required password.
2. Select [General Setup >>].
3. Select [Date Format] from [yyyy-mm-dd], [mm-dd-yyyy] and [dd-mm-yyyy].
4. Select [Time Format] and toggle between [24h] and [12h].
5. Set [System Time].

You can also set system time by selecting [Configuration >>]→[View Config]→[General Setup >>]. However, you cannot select date format and time format in this case. After the completion of setting system time, exit the configuration mode, and then the system will restart.

### Adjusting the Screen Brightness:

1. Press the Main Menu button on the front panel, and then select [Others >>].
2. Set [Brightness] to an appropriate level: 10 is the brightest, and 1 is the least bright.

You can also change screen brightness by entering configuration mode and selecting [Others] from the Configuration Main menu.



### Changing Key Volume:

1. Press the Main Menu button on the front panel, and then select [Others >>].
2. Select [Key Volume] and then select an appropriate value. 0 means key volume off and 10 is the maximum volume.

You can also change key volume by entering configuration mode and selecting [Others] from the Configuration Main menu.

### Selecting High Contrast Mode:

The equipment has the function of high contrast display so that the user can view the display under high ambient illumination.




**Adjusting Waveform Position:**

1. Press the Main Menu button on the front panel, and then select [Waves >>].
2. In the [Waves] menu, set [Wave 2], [Wave 3] and [Wave 4]. Wave 1 is always ECG1, which is unchangeable.

You can also change waveform position by entering configuration mode and selecting [Waveform Setup] from the Configuration Main menu.

➤ **Introduction of the function operation of the product.**

**AED:**




1. Confirm that the patient is unresponsive, not breathing and pulseless.
2. Remove clothing from the patient's chest. Dry the patient's chest and, if necessary, clip or shave excessive chest hair.
3. Apply multifunction electrode pads to the patient as directed on the pads package. Use anterior-lateral



placement.


4. Connect the pads with pads cable, and then plug the pads cable in the equipment's therapy port.
5. Turn the Mode Select knob to AED.
6. Follow the screen and voice prompts.
7. Press the Shock button, if prompted

### Manual Defibrillation:



1. Confirm that the patient is unresponsive, not breathing and pulseless.
2. Remove clothing from the patient's chest. Dry the patient's chest and, if necessary, clip or shave excessive chest hair.
3. Apply multifunction electrode pads to the patient as directed on the pads package. Use anterior-lateral placement.
4. Turn Mode Select knob to Manual Defi, Adjust the energy is necessary.
5. Charge
6. Shock


### Synchronized Cardioversion:



1. Connect the therapy cable and apply the multifunction electrode pads or external paddles to the patient.
2. With the Mode Select knob in Manual Defib position, press the [Enter Sync] soft key to activate the synchronous cardioversion function.

3. Select a lead. The selected lead should have a clear signal and a large QRS complex.
4. Verify that the white R-wave markers appear above R-waves. If the R-wave markers do not appear or do not coincide with the R-waves, for example above the T-waves, select another lead.
5. Select energy if necessary.
6. Charge.
7. Shock.

### Noninvasive Pacing:



1. Turn the Mode Select knob to the Pacer position.
2. Select a lead with an easily detectable R-wave.
3. Verify that white R-wave markers appear above the QRS
4. Select pacer rate.
5. Press the [Start Pacing].
6. Verify that white pacing markers appear on the ECG waveform.
7. Adjust pacer output until cardiac capture occurs.
8. Verify the presence of a peripheral pulse
9. Turn the Mode Select knob to the Pacer position.
10. Switch the pacer to the Fixed mode.
11. If ECG electrodes are applied, use the Lead Select button to select the desired lead.

### ➤ Introduction of the maintenance of the product.

#### Manual Defibrillation Test:



1. Run the equipment on fully charged battery. Move the Mode Select knob to Manual Defib.
2. Connect the external paddles to the equipment and place the paddles on the defibrillator/pacer analyzer.
3. Set the analyzer to Energy Measurement mode. In this case, the energy value should be displayed as 0 or blank.
4. Select the energy level to 360J.
5. Charge the equipment.
6. Verify that the charge tone is issued during charging.
7. Press the “Disarm” soft key to discharge the energy internally.
8. Verify that a prompt “Charge Removed” appears on the screen and the charge done tone stops.
9. Verify that the value measured by the analyzer is 0J or blank.
10. Enter the Configuration Main menu, select [Manual Therapy Setup] and set [Time to Auto Disarm] to [60s].
11. Exit “Configuration Management”. The equipment restarts automatically.
12. Set the analyzer to Energy Measurement mode. In this case, the energy value should be displayed as 0 or blank.
13. Select the energy level to 360J.
14. Charge the equipment. Count time after charging is completed. Verify that the prompt “Shock Removed” appears on the equipment and the energy measured by the analyzer is 0J or blank after 60 seconds.
15. Use multifunctional electrode pads. Repeat Step 3 to Step 14.

### **Synchronous Defibrillation:**



1. Run the equipment on fully charged battery. Move the Mode Select knob to Pacer. Select Fixed mode..
2. Connect the pads cable to the equipment and properly place the pads on the defibrillator/pacer analyzer.
3. Set the analyzer to Pacing Measurement mode. Use test load of  $50\Omega$ .
4. On the equipment, set [Pacer rate] to [70ppm] and [Pacer Output] to [30mA].
5. Press the [Start Pacing] soft key. Verify that the pacer rate measured by the analyzer is  $70\text{ ppm}\pm 1\text{ ppm}$  and the pacer output measured is  $30\text{ mA}\pm 5\text{ mA}$ .
6. Press the [Stop Pacing] soft key, and then set [Pacer rate] to [170ppm] and [Pacer Output] to [200mA].
7. Press the [Start Pacing] soft key. Verify that the pacer rate measured by the analyzer is  $170\text{ ppm}\pm 2\text{ ppm}$ , and the measured current is  $200\text{ mA}\pm 10\text{ mA}$ .

## 8. Harmonized standards and CS applied

- EN ISO 14971: 2012 Medical devices - Application of risk management to medical devices
- EN 1041:2008+A1:2013 Information supplied by the manufacturer with medical devices
- ISO 15223-1:2016: Medical devices-Symbols to be used with medical device labels, labeling and information to be supplied
- EN ISO 10993-1: 2009/AC:2010 Biological evaluation of medical devices - Part 1: Evaluation and testing
- EN 60601-1:2006/A1:2013 Medical electrical equipment--Part 1:General requirements for basic safety and essential performance
- EN 60601-1-2: 2015 Medical electrical equipment--Part 1-2: General requirements for basic safety and essential performance-- Collateral standard: Electromagnetic compatibility--Requirements and tests
- EN 60601-1-6: 2013 Medical electrical equipment-part 1-6: general requirements for basic safety and essential performance--collateral standard: usability
- EN 60601-1-8: 2012 Medical electrical equipment - part 1-8: general requirements for basic safety and essential performance - collateral standard: general requirements, tests and guidance for alarm systems in medical electrical equipment and medical electrical systems
- IEC 60601-2-4: 2018 Medical electrical equipment - Part 2-4: Particular requirements for the safety of cardiac defibrillators
- IEC 60601-2-25:2011 Medical electrical equipment - Part 2-25: Particular requirements for the basic safety and essential performance of electrocardiographs
- IEC 60601-2-27: 2011 Medical electrical equipment--Part 2-27: Particular requirements for the basic safety and essential performance of electrocardiographic monitoring equipment
- IEC 80601-2-30:2018 Medical electrical equipment - part 2-30: particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers
- ISO 81060-2: 2018 Non-invasive sphygmomanometers - part 2: clinical validation of automated measurement type
- IEC 60601-2-34: 2011 Medical electrical equipment - part 2-34: particular requirements for the basic safety, including essential performance, of invasive blood pressure monitoring equipment
- IEC 60601-2-49: 2011 Medical electrical equipment --Part 2-49: Particular requirements for the basic safety and essential performance of multifunction patient monitoring equipment
- ISO 80601-2-61:2011 Medical electrical equipment - part 2-61: particular requirements for basic safety and essential performance of pulse oximeter equipment
- ISO 80601-2-56: 2017 Medical electrical equipment - part 2-56: particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement
- ISO 80601-2-55: 2018 Medical electrical equipment - part 2-55: particular requirements for the basic safety and essential performance of respiratory gas monitors
- EN 62304: 2015 Medical device software - Software lifecycle processes
- IEC 62366-1: 2015 Medical devices - Application of usability engineering to medical devices
- EN 1789: 2007+A2:2014 Medical Vehicles and Their Equipment - Road Ambulances

## 9. Revision history

SSCP revision number	Date issued	Change description	Revision validated by the Notified Body
1.0	2021.06.23	Initial version	<input checked="" type="checkbox"/> Yes Validation language: English <input type="checkbox"/> No
2.0	2022.4.12	Supplemented PMCF data in Chapter 5	<input checked="" type="checkbox"/> Yes Validation language: English <input type="checkbox"/> No