

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-120: BS-120, BS-130, BS-180, BS-190;

10.BS-400: BS-400, BS-420;

2.BS-200: BS-200, BS-220, BS-330, BS-350;

11.BS-430: BS-430, BS-450,BS-460;

3.BS-200E: BS-200E, BS-220E,;

12.BS-480: BS-480, BS-490;

4.BS-240: BS-230, BS-240;

13.BS-600: BS-600, BS-620;

5.BS-240E: BS240E, BS240Pro;

14.BS-600M: BS-600M, BS-620M;

6.BS-300: BS-300, BS-320;

15.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

7.BS-330E:BS-330E(Serial Number starts with "XQ-"), BS-350E(Serial Number starts with "XS-")

8.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-

330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with "W9-" and software version starts with "35.00");

16.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

17.BS-2800M:BS-2600M, BS-2800M;

9.BS-380: BS-380, BS-390;

18.For applicable analyzers of the analyte, please subject to the parameter sheet and instrument.

LOT : 151024001



: 2025-06-30

English	Abbreviated name	Calibration Rule	Model	Level
Русский	сокращенное наименование	Принцип калибровки	модель	Уровень
Português	Nome abreviado	Regra de calibração	Modelo	Nível
Español	nombre abreviado	Regla de calibración	modelo	Nivel
Italiano	abbreviazione	Regola di calibrazione	modelli	Livello
Türkçe	kısaltılmış ad	Kalibrasyon Kuralı	model	Düzey

English	Calibration Value	Expanded Uncertainty
Русский	Эталонное значение	Расширенная неопределенность
Português	Valor de calibração	Incerteza Expandida
Español	Valor de calibración	Incertidumbre expandida
Italiano	Valore di calibrazione	Incertezza estesa
Türkçe	Kalibrasyon Değeri	Daha Uzun Süreli Belirsizlik

English	Sample Vol for Dilution	Diluent Vol	Sample Vol For Analysis
Русский	Объем пробы для разбавления	Объем разбавителя	Объем пробы для анализа
Português	Volume da amostra para diluição	Volume de diluente	Volume da amostra para análise
Español	Vol. muestra para dilución	Vol. diluyente	Vol. muestra para análisis
Italiano	Vol. campione per la diluizione	Vol. diluente	Vol. campione per analisi
Türkçe	Dilüsyon için Numune Hacmi	Seyreltici Hacmi	Analiz için Numune Hacmi

PA

English	Prealbumin
Русский	преальбумина
Português	pré-albumina
Español	Prealbúmina
Italiano	prealbumina
Türkçe	Prealbümin

Abbreviated name		PA		Calibration Rule			spline	
Model	Level	Calibrator Value		Expanded Uncertainty		Sample Vol For Dilution(μL)	Dilution Vol (μL)	Dilution Vol For Analysis(μL)
		mg/L	μmol/L	mg/L	μmol/L			
BS-120 ¹	0.9%NaCl	0.0	0.000	/	/	/	/	/
	S2	54.6	0.994	5	0.09	33	187	7
	S3	112	2.04	10	0.18	45	180	10.5
	S4	198	3.60	18	0.3	45	180	21
	S5	365	6.64	33	0.6	/	/	7
	S6	893	16.3	80	1.5	/	/	20
BS-200 ²	0.9%NaCl	0.0	0.00	/	/	/	/	/
	S2	55.6	1.01	5	0.09	33	187	7
	S3	115	2.09	10	0.19	45	180	10.5
	S4	191	3.48	17	0.3	45	180	21
	S5	354	6.44	32	0.6	/	/	7
	S6	931	16.9	83	1.5	/	/	20
BS-200E ³	0.9%NaCl	0.0	0.000	/	/	/	/	/
	S2	51.9	0.945	5	0.08	30	170	8.4
	S3	112	2.04	10	0.18	35	140	12.8
	S4	208	3.79	19	0.3	35	140	25.9
	S5	373	6.79	33	0.6	/	/	8.4
	S6	931	16.9	83	1.5	/	/	24
BS-300 ⁶	0.9%NaCl	0.0	0.000	/	/	/	/	/
	S2	54.8	0.997	5	0.09	33	187	7
	S3	115	2.09	10	0.19	15	195	32
	S4	211	3.84	19	0.3	35	230	34
	S5	379	6.90	34	0.6	/	/	7
	S6	967	17.6	87	1.6	/	/	20
BS-330E ⁷	0.9%NaCl	0.0	0.000	/	/	/	/	/
	S2	51.9	0.945	5	0.08	30	170	8.4
	S3	112	2.04	10	0.18	35	140	12.8
	S4	208	3.79	19	0.3	35	140	25.9
	S5	373	6.79	33	0.6	/	/	8.4
	S6	931	16.9	83	1.5	/	/	24
BS-360E ⁸	0.9%NaCl	0.0	0.000	/	/	/	/	/
	S2	54.3	0.988	5	0.09	21	119	7
	S3	113	2.06	10	0.18	10	130	32
	S4	208	3.79	19	0.3	21	138	34
	S5	379	6.90	34	0.6	/	/	7
	S6	958	17.4	86	1.6	/	/	20
BS-380 ⁹	0.9%NaCl	0.0	0.000	/	/	/	/	/
	S2	51.3	0.934	5	0.08	24	136	7
	S3	110	2.00	10	0.18	10	130	32
	S4	205	3.73	18	0.3	21	138	34
	S5	369	6.72	33	0.6	/	/	7
	S6	908	16.5	81	1.5	/	/	20
BS-400 ¹⁰	0.9%NaCl	0.0	0.000	/	/	/	/	/
	S2	51.9	0.945	5	0.08	24	136	7
	S3	112	2.04	10	0.18	10	130	32
	S4	208	3.79	19	0.3	21	138	34
	S5	373	6.79	33	0.6	/	/	7

		S6	931	16.9	83	1.5	/	/	20
Abbreviated name		PA		Calibration Rule			spline		
Model	Level	Calibrator Value		Expanded Uncertainty		Dilution Vol For Analysis(μL)	Sample Vol For Dilution(μL)	Dilution Vol (μL)	
		mg/L	μmol/L	mg/L	μmol/L				
BS-240 ⁴	0.9%NaCl	0.0	0.000	/	/	/	/	/	
	S2	52.9	0.963	5	0.09	5.6	27	153	
	S3	115	2.09	10	0.19	8.4	36	144	
	S4	206	3.75	18	0.3	16.8	36	144	
	S5	365	6.64	33	0.6	5.6	/	/	
	S6	885	16.1	79	1.4	16	/	/	
BS-240E ⁵	0.9%NaCl	0.0	0.000	/	/	/	/	/	
	S2	51.9	0.945	5	0.08	5.6	21	119	
	S3	112	2.04	10	0.18	25.6	10	130	
	S4	208	3.79	19	0.3	27.2	21	138	
	S5	373	6.79	33	0.6	5.6	/	/	
	S6	931	16.9	83	1.5	16	/	/	
BS-430 ¹¹	0.9%NaCl	0.0	0.000	/	/	/	/	/	
	S2	51.4	0.935	5	0.08	7	18	102	
	S3	110	2.00	10	0.18	32	10	130	
	S4	205	3.73	18	0.3	34	21	138	
	S5	369	6.72	33	0.6	7	/	/	
	S6	962	17.5	86	1.6	20	/	/	
BS-480 ¹²	0.9%NaCl	0.0	0.000	/	/	/	/	/	
	S2	51.9	0.945	5	0.08	7	24	136	
	S3	112	2.04	10	0.18	32	10	130	
	S4	208	3.79	19	0.3	34	21	138	
	S5	373	6.79	33	0.6	7	/	/	
	S6	931	16.9	83	1.5	20	/	/	
BS-600 ¹³	0.9%NaCl	0.0	0.000	/	/	/	/	/	
	S2	54.3	0.988	5	0.09	7	18	102	
	S3	112	2.04	10	0.18	32	10	130	
	S4	213	3.88	19	0.3	34	21	138	
	S5	381	6.93	34	0.6	7	/	/	
	S6	958	17.4	86	1.6	20	/	/	
BS-600M ¹⁴	0.9%NaCl	0.0	0.000	/	/	/	/	/	
	S2	53.9	0.981	5	0.09	5.6	18	102	
	S3	117	2.13	10	0.19	8.4	25	100	
	S4	210	3.82	19	0.3	16.8	25	100	
	S5	357	6.50	32	0.6	5.6	/	/	
	S6	931	16.9	83	1.5	16	/	/	
BS-800 ¹⁵ (R1:R2 :S=200:75:7)	0.9%NaCl	0.0	0.000	/	/	/	/	/	
	S2	51.9	0.945	5	0.08	7	18	102	
	S3	112	2.04	10	0.18	32	10	130	
	S4	208	3.79	19	0.3	34	21	138	
	S5	373	6.79	33	0.6	7	/	/	
	S6	931	16.9	83	1.5	20	/	/	
BS-800 ¹⁵ (R1:R2 :S=120:45:4.2)	0.9%NaCl	0.0	0.000	/	/	/	/	/	
	S2	51.9	0.945	5	0.08	4.2	18	102	
	S3	112	2.04	10	0.18	19.2	10	130	
	S4	208	3.79	19	0.3	20.4	21	138	

	S5	373	6.79	33	0.6	4.2	/	/
	S6	931	16.9	83	1.5	12	/	/
	0.9%NaCl	0.0	0.00	/	/	/	/	/
BS-2000¹⁶	S2	55.7	1.01	5	0.09	4.2	18	102
	S3	113	2.06	10	0.18	19.2	10	130
	S4	210	3.82	19	0.3	20.4	21	138
	S5	369	6.72	33	0.6	4.2	/	/
	S6	888	16.2	79	1.4	12	/	/
		0.9%NaCl	0.0	0.00	/	/	/	/
BS-2800M¹⁷	S2	56.9	1.04	5	0.09	4.2	18	102
	S3	118	2.15	11	0.19	19.2	10	130
	S4	220	4.00	20	0.4	20.4	21	138
	S5	371	6.75	33	0.6	4.2	/	/
	S6	931	16.9	84	1.5	12	/	/
		0.9%NaCl	0.0	0.00	/	/	/	/

mindray

Prealbumin Calibrator

For use on: BS-2000, BS-2800M,
BS-600M, BS-800

LOT 151024001

EXP 2025-06-30

