Instructions for Use of Sample Diluent

EN	1
GER	6
ES	12
FR	17
IT	22
TR	27
PT	32
RU	37

Instructions for Use of Sample Diluent

[Product Name]

Generic name: Sample Diluent

Model: AF-15

[Packing Specifications]

4×20.0 mL

[Intended Use]

It is used to dilute the patient plasma samples in a coagulation test.

[Test Principle]

The Sample Diluent is used together with Fibrinogen (FIB) Assay Kit (Clotting Method)/Antithrombin III (AT III) Assay Kit (Chromogenic Substrate Method) for testing with coagulation analyzers.

[Major Components]

Imidazole.

[Storage Conditions and Stability]

Stored and stable at 2~8°C for up to 548 days (1.5 year or 18 months).

After being tested, the kits can be stored in the reagent disk (with refrigeration function on) of Mindray CX series auto coagulation analyzers for 15 days.

For the refrigeration temperature of the reagent disks of different analyzer models, refer to the Operator's Manual of the respective analyzers.

Refer to packaging or labeling for production date and expiry date.

[Applicable Instruments]

Mindray CX series auto coagulation analyzers.

[Test Procedure]

Reagent Preparation

The refrigerated reagents shall be warmed up at room temperature (18~25°C) for 30 minutes before use.

Operation Steps

After decapping, the reagent shall be loaded to the instrument according to the instrument operation manual, and the instrument will automatically run the dilution process (refer to operation manual).

[Performance Characteristics]

Item	Requirement
pH value	The pH value of the sample diluent is 7.35 ± 0.10 at $25 \pm 1^\circ\text{C}$.
Dilution effect of AT III tests	When testing control material for AT III tests and using Sample Diluent for dilution, the results are within the control reference range.
Dilution effect of FIB tests	When testing normal control plasma for FIB tests and using Sample Diluent for dilution, the relative deviation of the measured FIB content from the control plasma target is not more than $\pm 10\%$.

[Warnings and Precautions]

- This kit is for in vitro diagnosis only;
- All the precautions for laboratory reagent

operation should be followed when using this product. All samples and reaction wastes shall be disposed as the source of infection. The reagent waste must be disposed of in accordance with the local regulations;

- This reagent is applicable to the auto coagulation analyzers of Shenzhen Mindray Bio-medical Electronics Co., Ltd. Please carefully read the operation manual before use;
- This product contains preservative content. It is prohibited to swallow this product or have it in contact with skin or mucous membrane!
- The Material Safety Data Sheet (SDS) is available upon request
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the country in which the user and/or the patient is established.
- The accuracy of test results can't be guaranteed if reagents are stored under improper conditions.
- Please confirm the integrity of the package before use. Do not use the kits with damaged packages.
- If unintentionally opened before used, store the reagents at 2-8°C and protected from light, and the stability is equally to in-use stability.
- It shall be used before the expiration date and disposed of properly when expired.

[Symbols]

IVD

In vitro
diagnostic
medical
device

Consult
instructions
for use

Temperatur
e Limit

Manufacture
r

LOT

Batch code

Use-by date

EC REP

Authorized
representati
ve in the
European
Community

REF

Catalogue
number

UDI

Unique
device
identifier

European
Conformity

[Company Contact]

Manufacturer: Shenzhen Mindray Bio-Medical
Electronics Co., Ltd.

Address: Mindray Building, Keji 12th Road South,
High-tech Industrial Park, Nanshan, Shenzhen
518057, P.R. China

Website: www.mindray.com

E-mail Address: service@mindray.com

Tel: +86 755 81888998

Fax: +86 755 26582680

EC-Representative: Shanghai International Holding Corp. GmbH(Europe)

Address: Eiffestraße 80, Hamburg 20537, Germany

Tel: 0049-40-2513175

Fax: 0049-40-255726

© 2021-2022 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved.

[Approval Date of the Instruction for Use]

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© 2021 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved.