



















































BeneVision N22/N19

Patient Monitor Quick Reference Guide

Symbol	Function
 Main Menu	Enters the main menu.
 Zero IBP	Starts IBP zero calibration.
 Audio Pause	Pauses alarm tone.
 Review	Enters the Review menu.
 Patient Management	Enters the Patient Management menu.
 Stop All	Stops all NIBP measurements.
 PAWP	Enters the PAWP screen.
 Venipuncture	Opens the Venipuncture window.
 Alarm Reset	Acknowledges ongoing alarms and resets the alarm system.
 NIBP Start/Stop	Starts an NIBP measurement or stops the current NIBP measurement.

Symbol	Function
 More	Shows more quick keys.
 Alarm	Enters the Alarm menu.
 Alarm Pause	Pauses the current alarms.
 Standby	Enters the standby mode.
 Screen Setup	Enters the Screen Setup menu.
 C.O. Measure	Opens the C.O. Measure window.
 Loops	Opens the Loops window.
 Start TOF	Starts/stops TOF measurement.
 NIBP STAT	Starts a five-minute continuous NIBP measurement.
 NIBP Measure	Enters the NIBP Measure menu.
 Discharged Patients	Enters the Discharged Patients dialog box.

Symbol	Function	Symbol	Function
 Parameters Setup	Enters the Parameters Setup menu.	 Remote View	Opens the Remote View window.
 Manual Event	Manually triggers and saves an event.	 Minitrends	Enters the Minitrends screen.
 OxyCRG	Opens the OxyCRG window.	 ECG Full-Screen	Enters the 12-lead ECG full screen.
 Privacy Mode	Enters the privacy mode.	 Night Mode	Enters the night mode.
 CPB Mode	Enters the CPB mode.	 Intubation Mode	Enters the intubation mode.
 Volume	Enters the Volume menu.	 Freeze	Freezes waveforms.
 Calculations	Enters the Calculations menu.	 Load Config	Enters the Load Config menu.
 Print	Starts printing a real-time report.	 Record	Starts/Stops a recording.
 Discharge Patient	Enters the Discharge Patient dialog box.	 Rotate Screen	Changes the screen orientation.
 ECG Lead/Gain	Enters the ECG Lead/Gain menu.	 Call Help	Calls for help.
 BoA	Enters the BoA Dashboard screen.	 EWS	Enters the EWS screen.
 GCS	Enters the GCS menu.	 SepsisSight	Enters the SepsisSight menu.
 HemoSight	Enters the HemoSight menu.	 Rescue Mode	Enters the Rescue mode.
 Pace View	Opens the Pace View window.	 ECG 24h Sum	Views the 24 hour ECG summary.



- (1) Alarm lamp
- (2) Battery indicator
- (3) AC power indicator
- (4) Power switch
- (5) Serial bus connectors

Turning on the Monitor

Press the power switch to turn on the monitor.

Display

Choosing a Screen

1. Access **Choose Screen** by either of the following ways:
 - ◆ Select the **Screen Setup** quick key.
 - ◆ Select the **Main Menu** quick key → from the **Display** column select **Choose Screen**.
2. Select the desired screen.

Selecting the Big Numerics Screen

1. Access **Choose Screen** by either of the following ways:
 - ◆ Select the **Screen Setup** quick key.
 - ◆ Select the **Main Menu** quick key → from the **Display** column select **Choose Screen**.

2. Select **Big Numerics**.
3. Select a parameter numeric area or waveform area, and then from the popup list select an element you want to display in this area.

Displaying Parameter Numerics and Waveforms

1. Access **Tile Layout** by either of the following ways:
 - ◆ Select the **Screen Setup** quick key → select the **Tile Layout** tab.
 - ◆ Select the **Main Menu** quick key → from the **Display** column select **Tile Layout**.
2. Select a parameter numeric area or waveform area, and then from the popup list select an element you want to display in this area. The parameters and waveforms you do not select will not be displayed.

Setting the Quick Key

1. Set the quick key by either of the following ways:
 - ◆ Select the **Screen Setup** quick key → select the **Select Quick Keys** tab.
 - ◆ Select the **Main Menu** quick key → from the **Display** column select **Quick Keys**.
2. Select the **Current** tab to configure the quick keys you want to display on the screen: From the top of this page, select a block where you want to show a certain quick key, and then select the quick key from the quick key list.
3. Select the **More** tab to configure the quick keys you want to display when the **More** quick key is selected.

Monitoring ECG

Choosing the ECG Lead Type

1. Select the ECG numeric area or waveform area to enter the **ECG** menu.
2. From the **ECG** page, select the **Setup** tab.
3. Set **Lead Set** according to the lead type you are going to use. The default lead type is **Auto**. In this case, the monitor automatically detects the used lead type.

Checking Paced Status

1. Select the ECG numeric area or waveform area to enter the **ECG** menu.
2. From the **ECG** page, select the **Pacer** tab.
3. Set **Paced** to **Yes** or **No**.

Enabling ST Monitoring

1. Select the ECG numeric area or waveform area to enter the ECG menu.
2. Select the **ST** tab → select the **Setup** tab.
3. Switch on **ST Analysis**.

Enabling QT/QTc Interval Monitoring

1. Select the ECG numeric area or waveform area to enter the ECG menu.
2. Select the **QT** tab → select the **Setup** tab.
3. Switch on **QT Analysis**.

Entering the 12-Lead Screen

1. Select the ECG numeric area or waveform area to enter the **ECG** menu.
2. From the bottom of the ECG page, select **12-Lead**.

You can also enter the 12-Lead screen by either of the following ways:











- Select the **Screen Setup** quick key → select **ECG 12-Lead**.
- Select **Main Menu** quick key → from the **Display** column select **Choose Screen** → select **ECG 12-Lead**.

Initiating Resting 12-lead ECG Analysis

To initiate 12-Lead ECG analysis, select **Analyze** from the left bottom of the 12-Lead screen.

Monitoring NIBP

Starting or Stopping NIBP Measurement

Task	By Quick Key	From NIBP menu
Start a manual measurement	NIBP Start/Stop quick key 	Start NIBP button
Start auto NIBP series	NIBP Start/Stop quick key 	Setup tab → set Interval → Start NIBP button
	NIBP Measure quick key  → select interval	
Start NIBP sequence measurement	NIBP Measure quick key  → Sequence	Sequence tab → set NIBP sequence → Start NIBP button
Start STAT measurement	NIBP STAT quick key 	STAT button
	NIBP Measure quick key  → STAT	
Stop the current NIBP measurements	NIBP Start/Stop quick key 	Stop NIBP button
End auto NIBP series or NIBP sequence	Stop All quick key 	Stop All button
Stop STAT measurement and end series	NIBP Start/Stop quick key 	Stop NIBP or Stop All button
	NIBP STAT quick key 	

Setting the NIBP Interval

1. Select the **NIBP** numeric area to enter the **NIBP** menu.
2. Select the **Setup** tab.
3. Set **Interval**. Selecting **Manual** switches to manual mode.

Monitoring IBP

Changing the Pressure Label

The pressure label is a unique identifier for each type of pressure. Therefore, you should select a proper pressure label for the source of the pressure you want to monitor.

To select the pressure label, follow this procedure:

1. Select an IBP numeric area or waveform area to enter the corresponding pressure menu.
2. Select the **Setup** tab.
3. Select the desired Label.

Zeroing the IBP Transducer

1. Connect the IBP transducer, the IBP adapter cable and the module.
2. Turn off the three-way valve (the one close to the transducer) to the patient, in order to Vent the transducer to the atmospheric pressure.
3. Zero the transducer by one of the following ways:
 - ◆ Press the **Zero** hard key on the module.
 - ◆ Select the numeric area (such as the Art numeric area), and then select **Zero**.
 - ◆ Select the **IBP Zero** quick key.
4. After the zero calibration is completed, close the stopcock to the air and open the stopcock to the patient.

Patient Management

Editing Patient Information

1. Select the **Patient Demographics** quick key.
2. In the **Patient Management** menu, edit patient information as required.

If you connect a barcode reader with your monitor, you can scan the patient's barcode to enter the patient's medical number.






Discharging a Patient

Manually discharge a patient using one of the following ways:

- Select the **Discharge Patient** quick key.
- Select the **Patient Demographics** quick key → select **Discharge Patient**.
- Select the **Main Menu** quick key → from the **Patient Management** column select **Discharge**.

Alarms

Alarm Status Symbols

Symbol	Description	Symbol	Description
	All the alarms are paused.		Audible alarm tones are paused.
	Audible alarm tones are turned off.		Alarms are acknowledged and the alarm system is reset.
	Individual physiological alarms are turned off or the system is in the alarm off status.		

Setting Parameter Alarm Properties

1. Access the **Limits** page by either of the following ways:
 - ◆ Select the **Alarm** quick key.
 - ◆ Select the **Main Menu** quick key → from the **Alarm** column select **Limits**.
2. Select a parameter tab and set alarm switch, limits, priority, etc.

Changing Arrhythmia Alarm Settings

1. Select the ECG parameter area or waveform area to enter the **ECG** menu.
2. Select the **Arrhythmia** tab → select the **Alarm** tab.
3. Set the alarm properties.

Initiating Auto Alarm Limits

1. Access the **Limits** page by either of the following ways:
 - ◆ Select the **Alarm** quick key → select the **Limits** tab.
 - ◆ Select the **Main Menu** quick key → from the **Alarm** column select **Limits**.

2. From the **Limits** page, select **Auto Limits** at the left bottom.

Pausing Alarms/Pausing Alarm Tones

1. Select the **Main Menu** quick key → turn to the third page → from the **System** column select **Maintenance** → input the required password → select **Ok**.
2. Select the **Alarm** tab → select the **Pause/Reset** tab.
3. Set **Pause** to **Alarm Pause** or **Audio Pause**.

Resetting Alarms

Pressing the **Alarm Reset** to acknowledge the ongoing alarms and reset the alarm system.

Accessing Help When Technical Alarms Occur

In the technical alarm list, alarm messages followed by **Detail** includes help information.

1. Select the technical alarm information area to enter the **Alarms** window.
2. Select the **Technical Alarms** tab.
3. Select the desired alarm.

Entering Night Mode

1. Select the **Night Mode** quick key, or select the **Main Menu** quick key → from the **Display** column select **Night Mode**.
2. Change the night mode settings if necessary.
3. Select **Enter Night Mode**.

Entering Privacy Mode

Choose one of the following ways to enter the privacy mode:

- Select the **Privacy Mode** quick key → select **Ok**.
- Select the **Main Menu** quick key → from the **Display** column select **Privacy Mode** → select **Ok**.

Configuring the Monitor

Adjusting the Volume

Select the **Volume** quick key to set **Alarm Volume**, **QRS Volume**, and **Key Volume**.

Adjusting the Screen Brightness

1. Access Display by either of the following ways:
 - ◆ Select the **Screen Setup** quick key → select the **Display** tab.
 - ◆ Select the **Main Menu** quick key → from the **Display** column select **Display**.
2. From the **Primary Screen** block, set **Brightness** for the primary display.
3. From the **Secondary Screen** block, set **Brightness** for the secondary display if configured.

Reviewing Trends

To review trends, enter the **Review** window by one of the following ways:

- Select the **Review** quick key.
- Select the **Main Menu** quick key → from the **Review** column select one option.

Viewing Other Beds

To view other beds, enter the **Remote View** window by one of the following ways:

- Select the **Remote View** quick key.
- Select the **Screen Setup** quick key → select the **Primary Display tab** or **Secondary Display tab** → select the **Choose Screen** tab → select **Remote View**.
- Select the bed at the alarm watch tile, if configured, on the main screen.

Company Contact

Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Mindray Building, Keji 12th Road South, High-tech Industrial Park,
Nanshan, Shenzhen 518057 P.R. China

Website: www.mindray.com

Tel: +86 755 81888998

Fax: +86 755 26582680



Shanghai International Holding Corp. GmbH(Europe)
Eiffestraße 80, 20537 Hamburg Germany



046-007709-00 (6.0)