

VS 9/8 Series

Vital Signs Monitor Quick Reference Guide

Note: The quick reference guide is not intended as a replacement of the Operator's Manual. Prior to operating the equipment, read the Operator's Manual carefully and follow the instructions to ensure the safety of the operator, patient and equipment.

Overview






- 1 NIBP cuff connector
- 2 SpO₂ cable connector
- 3 CO₂ outlet
- 4 CO₂ sampling line connector
- 5 Power switch



- 6 Alarm Reset/Pause key
- 7 NIBP Start/Stop key
- 8 Admit patient key
- 9 Knob

Quick Keys

Symbol	Function
	Enters the Review menu.
Review	
	Saves the input information.
Save	
	Enters the main menu.
Main Menu	

Symbol	Function
	Displays the patient list.
Patient	
	Clears all input information.
Clear	
	Shows more quick keys.
More	


Symbol	Function
 Next	Shows the next screen. (You can also swipe up on the main screen to show the next screen.)
 Alarm Pause	Pauses the physiological alarms.
 Targeted Goal	Enters the Targeted Goal screen.

Symbol	Function
 Back	Shows the previous screen. (You can also swipe down on the main screen to show the previous screen.)
 Alarm Reset	Resets the alarm system.

Turning on the Monitor


Press the power switch to turn on the monitor.

Selecting a Workflow

Select  at the top left corner of the screen to select a workflow.


- Select a workflow under **Spot Check**, and the monitor works in **Spot Check** mode.
- Select a workflow under **Continuous Monitoring**, and the monitor works in **Continuous Monitoring** mode.

Logging on as Clinician

1. Select  at the top right corner of the screen.
2. Input **Clinician** and **Password (if required)**.
3. Select **Login**.



If you connect a barcode reader with your monitor, you can scan your barcode or RFID card to enter the login information.

Admitting a Patient



1. Select the  quick key.
2. From the **Local Patient List** or **ADT Database**, select a patient and then select **Admit**.
3. Edit patient information as required. Select **OK**.

If you connect a barcode reader with your monitor, you can scan the patient's barcode to enter the patient's medical number.




Measuring NIBP

- Select  to set the patient position and measurement site.
- Select  to start the measurement.
- In **Continuous Monitoring** mode, select the NIBP area to access the NIBP menu.
- If ... is displayed on the bottom of the NIBP area, swipe from left to right to view another screen.

Measuring Temp



- P/M: P stands for **Predictive** mode and M stands for **Monitor** mode (SmarTemp).
-  : indicates the current measurement site.
-  : select to manually input a measured Temp value.

Manually Inputting RR Value


1. Select  or  in the Resp area. A timer is displayed and starts to time.
2. Count the total number of Resp of the patient in the past minute.
3. Input the number and select .



Inputting Manual Parameters

Select a parameter in the list, and input the information as instructed below:

-  : select to input a value.
-  : select to view the options.






Saving Data

1. Select the  quick key to save the physiological data of the current patient.
2. Check that required data is measured or input correctly. Select **OK**


Note: You may need to swipe up on the main screen or select the  quick key to show the  quick key.

Alarms


Alarm Status Symbols

Symbol	Description	Symbol	Description
	All the alarms are paused.		Audible alarm tones are paused.
	Audible alarm tones are turned off.		Alarms are acknowledged and the alarm system is reset.
	Individual physiological alarms are turned off or the system is in the alarm off status.		

Pausing Alarms/Pausing Alarm Tones

- Select the  quick key to pause or resume alarms or audio of the alarms.
- For VS 8 series, you can also press and hold the **Alarm Reset/Pause** key to pause or resume alarms or audio of alarms.

Resetting Alarms

Select the  quick key to acknowledge the ongoing alarms and reset the alarm system.

Reviewing Trends

Select the  quick key to review trends, scores and events.

