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## Wall Mount Operator's Manual

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
目 录 .....	i
<b>1 支架组件结构.....</b>	<b>1</b>
<b>2 安装.....</b>	<b>2</b>
2.1 导轨 .....	2
2.2 支撑臂 .....	13
2.3 适配器固定架 .....	24
2.4 储物篮 .....	26
<b>3 机械安全检查.....</b>	<b>27</b>
<b>1 Structure of Wall Mount.....</b>	<b>1</b>
<b>2 Installation .....</b>	<b>2</b>
2.1 Channel .....	2
2.2 Arm .....	16
2.3 Bracket of the adapter .....	27
2.4 Storage Basket .....	29
<b>3 Mechanical Safety Test.....</b>	<b>30</b>

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
# 1 Structure of Wall Mount

Item	Name	Function
1.	Wall Channel	Fix the wall mount to the wall. Note: replace the <i>wall mount channel</i> with <i>Adapter for connecting to the crane</i> . Fix the wall mount channel to the crane.
2.	Arm	Used for adjusting the angle of wall mount.
3.	Main unit support	Used for fixing the ultrasound main unit.
4.	Storage Basket	Used of keeping tissue, etc.
5.	Adapter bracket	Used for holding the external adapter.

# 2 Installation

---

## 2.1 Channel

### 2.1.1 Wall Channel

#### Warnings:

- Installations must be performed by Qualified Personnel. Failure to follow these Instructions may result in Serious Injury.
- Instructions and the Illustrations covering the specific Instrument to be mounted should be reviewed prior to Installation of Wall Channel.
- It shall be the responsibility of the hospital, its consultants and/or contractors to determine that the wall is adequate to safely mount instrumentation. This includes the selection of appropriate fasteners and the proper installation of the same. In new construction and remodeling work where the wall covering can be entirely or partly removed, a 16 gauge steel stud or a 2 X 6 Doug Fir #1 Stud should be located for the purpose of attaching the wall channel. The following is provided as guideline information only.
- Do not substitute or omit Fasteners.
- Do not position any Mounts or related hardware above a patient.
- Allow clearance on either side of the channel centerline to clear objects such as over-bed lighting, privacy curtains, adjacent walls or columns, door swing arcs, etc. Power and signal outlets should be considered when selecting a channel location. Avoid oxygen, vacuum and air outlets.
- Ensure that the weight being mounted does not exceed the Load Ratings:  
OSHPD Pre-Approval Compliant: M-Series arms 60lbs (27 kg), VHM arms 40lbs (18 kg) MAX LOAD RATING  
Over 60lbs (27 kg): Not OSHPD Pre-Approved.
- Although Considerable effort has been made to ensure the safety of the above installation and/or Guidelines, the actual wall construction materials and installation itself is beyond the control of GCX or MINDRAY Corporation. Accordingly, GCX or MINDRAY Corporation will not be responsible for the failure of any such installations.

**OSHPD:** The GCX Seismic channel conforms to the California Office of Statewide Health Planning and Development (OSHPD) Seismic Preapproval requirements. The seismic channel and related hardware have been dynamically tested in order to be granted two different OSHPD Anchorage Preapprovals (OPA).

**OPA-0079:** Preapproval for the GCX 19" Seismic Channel including anchorage of channel to wall surface and the design/configuration of inner wall structure. Four different wall designs (A,B,C,D) are detailed in these Instruction (as well as the preapproval). Inner wall structural design must be verified by Structural Engineer of Record (SEOR), and determined to be adequate as detailed on the preapproval.

**OPA-0697:** Preapproval for the GCX 19" Seismic Channel including anchorage of channel to wall surface, not including the design/configuration of the inner wall structure. The Structural Engineer of Record (SEOR) must supply all of the inner wall structural details, which are to be designed and verified to support weights and forces in addition to all other loads as detailed on the preapproval.

#### **Locating the Wall Channel**

Devices and Wall Mounts may require channel placement outside of these guidelines. Access to device controls, Ergonomic requirements, and the "Range of Motion" that the Wall Mount provides should be considered before mounting the wall channel.

**Devices:** The bottom edge of a 19 inch (48.3 CM) Wall Channel is placed 54 inches (142.2 CM) above the floor, plus or minus 4 inches.

## ■ Seismic Wall Channel

The following Wall Channel Parts Reference Chart lists hardware quantities for the 19" Seismic Wall Channel. Longer Wall

Channels are available, and include additional hardware for the additional mounting points.

### Seismic Wall Channel Parts Reference


Item#	Description	Qty
1	#10 x 2" Oval Head Sheet Metal Screw (OHSMS)	7
2	#10 x 2" Pan Head Sheet Metal Screw (PHSMS)	3
3	#12 x 2 3/4" Sheet Metal Screw (SMS)	4
4	1/4-20 x 3" Pan Head Machine Screw w/Toggle Wings	4
5	Fixed Channel Stop	1
6	10-32 x 3/8" Flat Head Machine Screw (FHMS)	1
7	Depressible Channel Stop	1
8	Seismic Wall Channel	1

### Installing Seismic Wall Channel


Drill all screw holes in the 16-gauge or 2X6 Doug Fir anchorage using a 9/64" diameter twist drill bit for the #10 x 2" type "A" sheet metal screws supplied. Do not substitute screws. Install the top, centerline screw first. Level the channel and mark for the remaining 9 screws. The (3) Pan Head Screws are used in the 3 centerline holes. The bottom Fixed Stop, and all others holes use Oval Head Screws. Screws should be started and driven with ball handle drivers or torque limiting screw guns not exceeding 60 inch-pounds (6.7 Nm).

**Channel Placement Recommendation**

See the "Locating the Wall Channel" Section on page 3 for Recommended Mounting Height above floor.

**Hardware Placement**


**Seismic Wall Channel Support Backing Methods**

The four anchorage methods illustrated (pages 6 – 7) represent Installations made to comply with California OSHPD requirements for the anchorage and installation of instrument support systems. Load tests were performed on a 10 ft. high wall using methods A, B, C, and D. Compliance with these written instructions is required for the installation of the system.

# METHODS OF SUPPORT BACKING

## Method-A


### Seismic

# METHODS OF SUPPORT BACKING

## Method-B


### Seismic

# METHODS OF SUPPORT BACKING

## Method-C


### Seismic

# METHODS OF SUPPORT BACKING

## Method-D

### Seismic

## ■ Standard Wall Channel

The following Wall Channel Parts Reference Chart lists hardware quantities for 19" Standard Wall Channel. Longer Wall

Channels are available, and include additional hardware to compensate for additional mounting points.

### Standard Wall Channel Part Reference

Item#	Description	Qty
1	#12 x 2 3/4" Sheet Metal Screw (SMS)	4
2	1/4-20 x 3" Pan Head Machine Screw w/Toggle Wings	4
3	Fixed Channel Stop	1
4	10-32 x 3/8" Flat Head Machine Screw (FHMS)	1
5	Depressible Channel Stop	1
6	Standard Wall Channel	1


### Installing the Standard Wall Channel

Drill 11/64" (17 mm) diameter holes for #12 X 2-3/4" high tensile Philips Head Screws. Install the top screw and level the channel, then mark and drill the other holes. Screws should be started and driven with ball handle drivers or torque limiting screw guns not exceeding 60 inch-pounds. Screws must pass through Channel, Drywall, and fully engage with Wood Stud or 16 Gauge Steel Stud.


Drill a 3/4" (19.5 mm) diameter holes for the 1/4-20 x 3" Pan Head Machine Screws with Toggle Wings. Toggle Wings should be used if Steel Wall Stud is less than 16 Gauge in thickness.

Installation Note: Seismic Wall Channel can be installed the same as Standard Wall Channel when OSHPD approval is not required.


**Channel Placement Recommendation**

**Hardware Placement**

**Method-A  
Standard**

This Anchorage will allow for a Quick and Simple installation of the Standard Wall Channel. 16 Gauge Steel will prevent the material from extruding up the Shank of the Screw (Stripping). 11/64" Diameter holes are required for #12 X 2-3/4" High Tensile, Phillips Head Screws (4) installed down the Centerline. (Non-Seismic Applications).

A second Track can be used for Back to Back Installations. A 2" X 4" Wood Stud may be used if allowed by Local Building Codes.

**Drywall (Gypsum Board) over 16 Gauge Sheet Metal or Wood Studs:**

Locate the exact centerline of the stud. Drive a #4 finishing nail through the drywall to contact the stud. Withdraw and drive again 3/8" (9.5 mm) right and left until the edges of the stud have been located and hence the centerline determined. The wall channel's 4 inch (10.2 CM) width will cover these exploratory holes.

Drill 11/64" (17 mm) diameter holes for #12 X 2-3/4" high tensile Philips head screws in a single 16 gauge stud through

Channel Centerline, (4 places). Install the top screw and level the channel, then mark and drill the other holes.

**Drywall (Gypsum Board) over Sheet Metal Studs less than 16 Gauge:**

Locate the Centerline of the stud per the above instructions. Use 1/4-20 X 3" Pan Head Machine Screws and Toggle Wings. Drill 3/4" (19.5 mm) diameter holes with a sharp speed bore or twist drill (requires pilot hole). Note that the round stop at the bottom of the channel goes over the shank of the bottom screw. Insert screws through Wall Channel holes then affix the toggle wings to the screws. Insert toggle wings into previously drilled holes. Use a level to make sure the channel is vertical prior to final seating of the screws.


**Plaster Coat over Expanded Metal Lath on Steel Studs, Hollow Tile, Hollow Block:**

Drill 3/4" (19.5 mm) holes. Broach the holes in hard materials to allow for the toggle wing to pass through. Insert screws through Wall Channel holes then affix the toggle wings to the screws. Insert toggle wings into previously drilled holes. Use a level to make sure the channel is vertical prior to final seating of the screws.


**Concrete Walls:**

Concrete wall: Using six M6X65 Expansion screws to fix the wall channel.

1. Determining the locations of the holes of installing the guide rail on the wall based on the height of fixing the ultrasound machine. With drilled six  $\Phi 10$  mm and 40mm depth holes vertically on the wall , the hole intervals are shown as following:

2. Aligning the guide rail with the holes on the wall and screw it through with six M6X65 expansion screws:

1-Wall ; 2-Wall Channel; 3-Expansion screw (M6X65)

**Through Wall, Back to Back Installations:**

Sufficient length 1/4-20 machine screws, nuts, washers and/or back plates must be obtained to mount the channel by drilling all the way through the thin partition wall. A Portalgn drill fixture is recommended for this procedure.

**Special Applications:**

Channels can be attached to a GCX 4" X 4" (10.2 x 10.2 CM), 1/8" (3.2 mm) wall aluminum support column or certain mullions, etc.. Contact GCX (Mindray Technical Service) to discuss solutions to wall mounting problems for any situation not covered by these guidelines.

## 2.1.2 Crane Channel

7"/17.8cm Channel with 38mm Post Clamps

The purpose of this guide is to:

- Describe attaching Channel to 38mm Post.

**Warning:** Ensure that the weight of the mount and instrument to be installed does not exceed 60lbs (27.2kg).

### Parts Reference

The following parts are included in this installation kit:

Item	Description	Qty
1	Clamping Blocks	2
2	Channel	1
3	10-32 x 5/8" Socket Head Cap Screw (SHCS)	1
4	5/32" Hex Wrench	1

1

2


### Tools Required:

- 5/32" Hex Wrench (provided)

### Attaching Mount to Pole

Using the 5/32" Hex Wrench (supplied), attach Channel to the 38mm post with (2) Clamping Blocks and (4) 10-32 x 5/8"

SHCS as shown below. Tighten screws equally (40 in/lb. minimum) until Clamping Blocks are secure.

## 2.2 Arm

The purpose of this manual is to describe general installation, operation, and adjustment procedures for VHM™ Series Arms. This manual should be used in conjunction with any instrument-specific installation guides. Please read this manual and all instrument-specific installation material before installing or using this product.

## 2.2.1 Installation, Operation, and Maintenance

### Warnings

This section contains warnings regarding the installation, operation, and maintenance of the VHM Arm. This section must be read in its entirety before installing, operating, and maintaining the VHM Arm. Failure to follow these warnings may result in damage to equipment or injury to personnel.

#### Warnings:

- Do not position the VHM Arm or mounted instrument above a patient. Note that the VHM ARM has a wide range of motion both up/down and side to side. Please consider carefully the instrument being mounted and the proximity of the mounting assembly to other equipment, hospital personnel, and the patient. GCX (MINDRAY) recommends that the hospital's risk management personnel verify that the application is appropriate prior to installation and use of the VHM Arm.
- If the VHM Arm is wall mounted, verify that the channel has been installed and approved in accordance with the channel installation guide.
- Ensure that the weight of the instrument being mounted does not exceed the load rating of the VHM Arm. Check the bottom surface of the Arm for the "Maximum Load" label (see 5.2.2.1 for label location).
- The Height Locking Lever or Knob must be engaged when removing a mounted instrument from the VHM Arm.
- Stand to the side of the VHM Arm and mounted instrument and use caution when disengaging the Height Locking Lever. The total mounted load may have changed due to the addition or removal of equipment (other than primary instrument) while the Height Locking Lever was engaged. A change in the total load can cause a sudden downward or upward movement of the Arm when the Height Locking Lever is disengaged.
- The mounted device or arm may move suddenly due to normal wear or improper adjustment of the tilt and swivel functions (see Sections 5.2.4, 5.2.5.2, and 5.2.5.5) or ultimately, gas spring end of life. The VHM Arm must be inspected and adjusted at least two times per year. This inspection must include the following steps:
  1. With the device mounted, release the height lock and move the arm through its entire vertical range of motion. The load should maintain its position at every point in the travel of arm. If necessary, the counterbalance mechanism may be adjusted in accordance with section 5.2.4.
  2. Grasp the mounted device and swivel it from side to side. The device should swivel with some tension or resistance, not loosely. If necessary, the swivel tension may be adjusted in accordance with section 5.2.5.5.
  3. Grasp the Arm and pivot it from side to side at the channel attachment. The arm should pivot with some tension or resistance, not loosely. If necessary, the pivot tension may be adjusted in accordance with section 5.2.5.4.
  4. Grasp the mounted device and tilt it forward and back, through its entire range of motion. There should be enough tension or resistance in the tilt mechanism to prevent the device from tilting forward unexpectedly when in use.
  5. Inspect fasteners for looseness. Tighten as required for optimal operation and safety.

6. Due to risk of personal injury or damage to the equipment, the VHM Arm housing must never be disassembled by non-GCX personnel. Failure to comply will void the warranty.

Failure to periodically inspect and adjust the Arm as instructed may result in damage to equipment or injury. If you do not intend to periodically inspect the arm, the Height Locking Lever should be removed entirely to reduce unexpected movement of the arm. If regularly inspected, the VHM Arm should function properly for 5 to 10 years. If you are unable to adjust the arm or require service, contact GCX technical support at 800-228-2555 or Mindray Technical Service.


## 2.2.2 Installing the VHM Arm in the Channel

**2.2.2.1 Thread the Height Locking Lever, \*if provided, clockwise into the mounting hole in the right side of the Arm (below left). Turn Locking Lever clockwise to lock or counterclockwise to unlock. Installation Note: When the VHM arm is properly counterbalanced for the weight load, the arm will stay in the desired vertical position unless moved by the user. Engagement of the Height Locking Lever is necessary while changing the location of the entire Arm within the channel. Engagement of the Height Locking Lever is otherwise optional to prevent downward or upward movement. Refer to Section 5.2.4 (Adjusting Counterbalance) and Section 5.2.5.4 (Adjusting Height of the Mounted Instrument) for further information on use of the Height Locking Lever.**

*\*VHM Arms with a Spring Lock are equipped at the factory with a Locking Knob.*

**Note:** The Height Locking Lever is a multi- position clamping lever that operates by lifting, rotating, and releasing the handle.

### VHM Arm with Spring Lock

This type of VHM Arm is designed for use with transport monitors only. The Locking Knob is installed at the factory. Ensure Locking Knob is in locked position before attempting to mount Arm in channel.


**Installation Note:** For proper function of the locking mechanism, the arm must be properly adjusted to counterbalance the weight of the mounted instrument (see section 5.2.3).

**Note:** Locking Knob operates by pulling and releasing Knob.

**2.2.2.2** While supporting the \*bottom of the Arm, guide the Slide into the top of the channel (below left) and move Arm to required mounting position.

**2.2.2.3** Using the 1/8" hex wrench provided, tighten four (4) set screws in Slide to secure position of Arm (below right).

## 2.2.3 Mounting the Instrument on the VHM Arm

Most instruments will require the attachment of an instrument-specific Mounting Adapter provided by GCX or by the instrument supplier. Please install the Mounting Adapter in accordance with the installation guide included with the Adapter.

### 2.2.3.1 Installing the Support

**Maximum load:** 9.0kg

Disassembly tool required: Phillips screwdriver

1. Align the screw holes on the main unit support with the screw holes on the arm.
2. Tighten (3) 10-32 x 3/4" FHMS screws.

3. Align the screw holes on the main unit support with the screw holes on the connection plate of the main unit support.
4. Tighten (3) 10-32 x 1/2" FHMS screws.


(3) 10-32 x 1/2" FHMS

### 2.2.3.2 Installing the Main Unit

**WARNING:** The arm must be locked before installing the display.

1. Snap the flange of the bracket into the installation slot towards arrow's direction.


The flange of the bracket

2. Push the display towards arrow's direction gently until the spring pin accesses the limit hole.

**Note:** If the spring pin does not access the limit hole, pull down the spring pin, and then push the display gently.

After the installation is completed, push the display opposite to see whether the display is tightened.

## 2.2.4 Adjusting Counterbalance


**WARNING:** Do not attempt counterbalance adjustment unless the instrument and accessories are mounted on the Arm. Use caution while performing this procedure.

The VHM Arm must be adjusted to counterbalance the mounted instrument throughout the Arm's vertical range of motion. When properly counterbalanced, the VHM Arm will maintain its height without engagement of the Height Locking Lever or Locking Knob. Because instrument weights vary, adjustment will typically be required to achieve optimal performance. A 1/2" [13mm] socket wrench is provided for this critical adjustment.

### To adjust counterbalance:

1. Grasp the Arm behind the mounted instrument(s) and carefully move the Arm to a horizontal position that allows access to the Adjuster Cover (below left). Lock Arm in horizontal position.
2. Open the Adjuster Cover by inserting a flat blade screwdriver in the slot at the rear of the Cover and prying upward. If necessary unlock and readjust the Arm until the Counterbalance Adjuster bolt is accessible through the Adjuster Cover (below center). Lock Arm when Counterbalance Adjuster is accessible.
3. With the 1/2" [13mm] socket wrench on the Adjuster, carefully unlock the arm. Turn the Adjuster counterclockwise (CCW) to increase counterbalance force, or clockwise (CW) to decrease counterbalance force. Counterbalance is correctly adjusted when the mounted instrument can be moved up or down with minimal force and does not rise or fall after releasing the Arm.

4. Close the Adjuster Cover.

## 2.2.5 Operating and Adjusting the VHM Arm

Check pivot, tilt, and front swivel mechanisms for proper tension. Refer to the applicable section and follow the adjustment procedure.

### 2.2.5.1 Counterbalance Adjustment - See Section 5.2.4


### 2.2.5.2 Adjusting Height of the Mounted Instrument

#### **WARNING:**

- Do not attempt to adjust height without instrument being mounted on Arm.
- Stand to the side of the VHM Arm and mounted instrument and use caution when disengaging the Height Locking Lever. The total mounted load may have changed due to the addition or removal of equipment (other than primary instrument) while the Height Locking Lever was engaged. A change in the total load can cause a sudden downward or upward movement of the Arm when the Height Locking Lever is disengaged.

Grasp mounted instrument on sides (or handle) and move to desired height. When Arm is counterbalanced correctly (sec 4.0), it will maintain position when adjusted to a new height. **Note:** Disengagement/re-engagement of the Height Locking Lever may be necessary during vertical adjustment of mounted instrument (see **WARNING** above). Engagement of the Height Locking Lever is not required during normal operation of the counterbalanced instrument.


### 2.2.5.3 Positioning VHM Arm in Channel

**!** **Caution:** Do not attempt to remove the Arm from the Channel while Arm is loaded with an instrument.\

1. Lock the Arm in highest vertical position using Height Locking Lever or Locking Knob.
2. Using the 1/8" hex wrench provided, loosen four (4) set screws in Slide.
3. Relieve the weight of the mounted instrument by lifting against the underside of the Arm, near the instrument.


Simultaneously, push up or pull down the opposite end of the Arm immediately in front of the Slide.

4. Tighten four (4) set screws.

### 2.2.5.4 Pivoting the VHM Arm and Adjusting Pivot Tension

1. To pivot Arm at Channel, simply push on side of Arm.
2. To adjust pivot tension, remove upper and lower bolt caps from Pivot Bolt. Using the 1/2" [13mm] socket wrench provided and a 1/2" [13mm] box or open-end wrench, tighten or loosen Pivot Bolt and Hex Nut to desired tension. Reinstall bolt caps on each end of Pivot Bolt.


### 2.2.5.5 Swiveling the Mounted Instrument and Adjusting Swivel Tension

1. To swivel mounted instrument, simply push or pull instrument.
2. To adjust swivel tension, tighten or loosen the Swivel Tension Nut using the 1/2" [13 mm] socket wrench provided.

**Installation Note:** It may be necessary to remove an existing down post from the swivel cup to gain access to the Swivel Tension Nut.


 **Warning:** Swivel tension nut must be torqued to a minimum of 20 in- lbs [2.3 N-m].

### 2.2.6 Cable Management

A Cable Guide is provided to facilitate routing of cables along the underside of the \*arm. Squeeze the Cable Guide until its edges snap into the grooves in the bottom surface of the arm. Leave some extra cable loose at the front and rear of the arm to prevent cable binding, connect or damage, or Cable Guide damage.

Guides snap into grooves on the bottom of the Extension.

## 2.2.7 Routine Maintenance

Periodically inspect all tilting, swiveling, pivoting, and mounting hardware. Tighten or adjust as necessary for optimal operation and safety.

## 2.2.8 Cleaning the Mounting Assembly

1. The mounting assembly may be cleaned with most mild, non-abrasive solutions commonly used in the hospital environment (e.g. diluted bleach, ammonia, or alcohol solutions).
2. The surface finish will be permanently damaged by strong chemicals and solvents such as acetone and trichloroethylene.
3. Steel wool or other abrasive material should *never* be used.
4. Damage caused by the use of unapproved substances or processes will not be warranted. We recommend testing any cleaning solution on a small area of the arm that is not visible, to verify compatibility.
5. Never submerge or allow liquids to enter the arm. Wipe any cleaning agents off of the arm immediately using a water-dampened cloth. Dry the arm thoroughly after cleaning.

**CAUTION:** GCX (MINDRAY) makes no claims regarding the efficacy of the listed chemicals or processes as a means for controlling infection. Consult your hospital's infection control officer or epidemiologist.

## 2.2.9 Troubleshooting the VHM Arm

Symptom	Possible Cause	Remedy
Mounted instrument does not appear level or parallel to the floor.	Channel not plumb. Check with level.	Adjust Channel to plumb, or reinstall Channel.
	Weight of instrument not compatible with Load Rating of the Arm.	Mount instrument on arm with compatible Load Rating.

Symptom	Possible Cause	Remedy
	Swivel hardware loose.	Adjust Swivel Nut (section 5.5.5.5).
	Pivot hardware loose.	Adjust Pivot Bolt (section 5.5.5.4).
	Mounting surface (e.g. wall, side of anesthesia machine, etc.) not structurally sound (does not hold mounting hardware).	Mounting surface must be reinforced or Channel must be relocated.
	Channel loose at mounting surface.	Check for plumb and tighten, or relocate (reinstall) Channel.
	Set Screws (4) in Slide are loose.	Reposition Arm and tighten set screws.
Instrument drifts up or down when the Arm is unlocked.	Arm not counterbalanced correctly for weight of the instrument.	Perform counterbalance adjustment per section 5.2.4.
	Weight of mounted instrument (load) not compatible with Load Rating of Arm.	Use arm with compatible Load Rating, and perform counterbalance adjustment per section 5.5.5.5.
Arm pivots too freely.	Pivot Bolt too loose.	Adjust Pivot Bolt (section 5.5.5.4).
Arm does not pivot easily.	Pivot Bolt too tight.	
Instrument swivels too freely.	Swivel Nut too loose.	Adjust Swivel Nut (section 5.5.5.5).
Instrument difficult to swivel.	Swivel Nut too tight.	
Instrument difficult to tilt.	Tilt Adjustment Lever too tight.	Adjust tilt and tilt tension per section 5.5.5.6.
	Excessive tilt tension.	
Instrument will not maintain tilt position.	Tilt Adjustment Lever too loose.	
	Insufficient tilt tension.	
Arm difficult to move up or down when unlocked.	Arm not counterbalanced correctly for weight of mounted instrument.	Perform counterbalance adjustment per section 5.2.4.
Arm inadvertently slides down Channel.	Set Screws (4) in Slide are loose.	Reposition Arm and tighten set screws.
Locking Knob (Arms with Spring Lock) will not pull out (release).	Arm not counterbalanced correctly for weight of mounted instrument.	Perform counterbalance adjustment per section 5.2.4.
	Weight of mounted instrument (load) not compatible with Load Rating of Arm.	
Counterbalance Bolt difficult to adjust.	Arm not locked in a horizontal position.	Reposition Arm and lock in horizontal position.

## 2.3 Bracket of the adapter

The purpose of this guide is to describe wall mounting of a typical power supply.

Min. Power Supply Size Using 12" [30.5 cm] Hook & Loop Straps	Max. Power Supply Size Using 18" [45.7cm] Hook & Loop Straps
1.4" [3.56 cm] x 2.3" [5.84 cm] or 7.4" [18.8 cm] perimeter	4.0" [10.16 cm] x 4.0" [10.16 cm] or 16" [40.6 cm] perimeter

### Parts Reference

The following parts and hardware are included with this installation kit (see photo for parts; hardware not shown):

Item	Description	Qty
1	Power Supply Mount (with 12" Hook & Loop Straps)	1
2	Cable Retainer, Nylon	1
3	#8-32 x 3/8" Pan Head Machine Screw (PHMS)	1
4	#8 Flat Washer	1
5	Hook-and-Loop Strap, 18"	2

**Tools Required:** 1/8" hex wrench (provided) Phillips screwdriver (not provided).

**Maximum load:** 1.5kg

### Mounting the Power Supply Bracket in the Channel


1. Insert the slide section of the Power Supply bracket into the GCX channel, position in desired location.

**Installation note:** The bracket can be oriented vertically or horizontally as shown, depending on the application.

2. Secure the Power Supply Bracket within the channel by tightening the 1/4-20 set screw using the 1/8" hex wrench.

Vertical Orientation

Horizontal Orientation

## Mounting the Power Supply in the Bracket

1. Place power supply in Bracket, insert hook-and-loop straps through rings, pull straps tight and press firmly (below left).

**Installation Note:** This Bracket is pre-assembled with two (2) 12" hook-and-loop straps. These straps securely hold most power supplies, but may be replaced with 18" straps (provided) for larger power supplies.

2. Loop power cable through Nylon Cable Retainer and fasten Retainer to bottom of Bracket with one (1) #8-32 x 3/8" PHMS and #8 flat washer (below right).

**Secure Velcro Straps**

**Secure Power Cable  
with Retainer**

**Assembled Power Supply Bracket**

## Routine Maintenance

Periodically inspect all hardware associated with the mounting system. Tighten or adjust as necessary for optimal operation and safety.

## Cleaning the Mounting Assembly

1. The mounting assembly may be cleaned with most mild, non-abrasive solutions commonly used in the hospital environment (e.g. diluted bleach, ammonia, or alcohol solutions).
2. The surface finish will be permanently damaged by strong chemicals and solvents such as acetone and trichloroethylene.
3. Do not use steel wool or other abrasive material to clean the mounting assembly.
4. Damage caused by the use of unapproved substances or processes will not be covered by warranty. We recommend testing of any cleaning solution on a small area of the arm that is not visible to verify compatibility.
5. Do not submerge or allow liquids to enter the arm. Wipe any cleaning agents off of the arm immediately using a water-dampened cloth. Dry the arm thoroughly after cleaning.


**CAUTION:** GCX (MINDRAY) makes no claims regarding the efficacy of the listed chemicals or processes as a means for controlling infection. Consult your hospital's infection control officer or epidemiologist.

## 2.4 Storage Basket


Maximum load: 2.0kg

Disassembly tool required: Phillips screwdriver

1. Align the screw holes on the storage basket with the screw holes on the connection plate.
2. To fix the storage basket onto the connection plate, use the Phillips screwdriver to tighten (4) 10-32 x 3/8" PHMS screws.

3. Align the screw holes on the connection plate with the screw holes on the pivot joint.
4. Use 5/32 (4 mm) hex wrench to tighten (3) 10-32 x 3/8" SHCS.

# 3 Mechanical Safety Test


---

This section mainly focus on conducting the mechanical safety test on the mechanical strength and mechanical function of the important components in ultrasound equipment.

Main methods of evaluation in the test are: Visual check and operation testing.

If the test results does not pass, then the machine must be stopped and appreciate measures be taken as the system could be in an abnormal condition.

Flowchart for the test is shown below.

**Figure of flowchart for conducting Mechanical safety test**

Table for mechanical safety test

Content	Method	Tool
Assembling screws for the wall mount	Visual examination: Check that the support assembly board does not show any signs of warping. Test with the wrench: check that none of the assembling screws are loose.	Inner hexagon spanner

Content	Method	Tool
Arm assembly	a) Swivel the arm on the vertical side. The arm reaches the desired angle on the vertical side. No cracks happen. b) Push the arm on the vertical side. The arm reaches the desired angle on the vertical side. No cracks happen.	N/A
Arm locking mechanism	Turn the level on the support arm pivot to position. The arm should not be able to move from side to side.	N/A
Screws for assembling the main unit support	Test with the wrench: check that none of the assembling screws are loose.	Cross-headed screwdriver
Damping axis of the main unit support	Swivel the main unit support to desired angle. The main unit keeps stable. The motion of the damping axis is flexible. No cracks or abnormality happen.	N/A
Other mechanical structures	Check that none of the other mechanical structure show signs of loose components, cracks/fractures on the outer surface or have exposed conducting parts.	N/A