

# HIV False Negative case

IVD International Clinical Application Department

Sang

2024/9/12



# CONTENTS

01. Case Background



02. Case Ideas



03. Case Solution



04. Case Summary



## Background

A private lab had a concern about HIV result for a patient who was sent back because of HIV Ab/Ag positive, and on returning to diagnostic a rapid testing with Determine kit and A Brand was positive.

Mindray CL-1200i initial and repeat results were both negative.

Mindray CL-1200i: 0.1 Negative

Determine Alere Rapid Test: Positive

A Brand: 5.74 Positive

## Complain

The customer complains about the Mindray results do not correlate with other platforms. The false Negative would result to the laboratory being fined or with increased number leads to license cancellation.

## Expectation

The customer was very angry about false positive of HIV, they need Mindray to find out the root cause of this issue.

# Case Background

Background

Ideas

Solution

Summary

Chemistry	Result	Flag
Anti-HCV	0.20	
HBsAg-I	<0.05	<
HIV	0.10	

Mindray result HIV Ab/Ag non-reactive



RDT determine results HIV  
Negative


Candidate status

Dear Sir/Madam,  
Mentioned above is the medical report for Mr./Ms. FAITH ACHIENG OLOO who is according to the GCC criteria

Remarks


The status is reported as unfit Unfit reason: HIV / AIDS Positive  
Unfit remarks: HIV Ag/Ab positive

Attachments



GCC results HIV Ab/Ag positive

# Test quality results indicators



# Case Ideas

Background

Ideas

Solution

Summary

Verify the preanalytical and sample pretreatment protocol



	Bayan lab
Syringe	BD
Needle	BD
Agglutination	>30min
Centrifuge	4000RPM, 10min
Pipette used to Transfer serum to Sample Cup	Yes
Sample Cup	biosystem
Retest Procedure	Retest the positive samples after re-centrifugation, if it's still positive recollect the sample after two weeks and retest. Send the repeatedly positive samples out for confirmation.
Reagent/calibration storage	2-8 c and monitored daily

Reagent storage & lab environment were normal.

Sample pretreatment is normal.



# Case Ideas


Verify the analytical aspects

Background

Ideas

Solution

Summary



## Basic Check:

Calibration and QC results were normal, the sample was rerun on Mindray, and the result was the same. Serum was clear without foreign body. Patient denied any medicine or exposure history.

## Possible Reasons Causing false negative for Mindray

- Individual is in **the eclipse period** before detection of Ag or HIV RNA possible;
- Individual is in **acute phase** of infection (before seroconversion) but is screened using a less sensitive method that detects Abs only;
- Individual in **the early stage** of seroconversion but is screened using a less sensitive method that does not detect early (IgM) Abs;
- Technical errors;
- Other core possibilities:
  - ✓ Delayed Ab synthesis in infants and persons receiving PEP or PrEP or who have concurrent acute hepatitis C infection;
  - ✓ Diminished immune response in individuals receiving intensive or long-term immunosuppressive therapy;
  - ✓ Congenital or drug-induced hypogammaglobulinemia or agammaglobulinemia;
  - ✓ Insufficient host Ab response (i.e., advanced HIV disease)
  - ✓ Unavailability of Abs due to the formation of Ag-Ab complexes



## Possible Reasons Causing false positive for A Brand and RDT


- Increased sensitivity of assays, leading to reduced specificity;
- Technical errors (specimen mix-up, mislabeling, improper handling);
- Presence of HIV Abs in recipients of HIV-1 trial vaccines;
- Other rare possibilities:
  - ✓ Hypergammaglobulinemia/Abs reactive to cellular components;
  - ✓ Influenza vaccination may cause cross-reactivity with HIV Ab assay. The time course for such cross-reactivity remains uncertain;
  - ✓ Autoimmune disorders and other medical conditions;
  - ✓ Multiple transfusions

## Troubleshooting:

Confirmation tests of HIV are PCR test and western blotting:

- The sample was sent to be tested for HIV-1, the result was negative.
- Western blotting reagent was sent to the customer side, it turns out the sample was negative on both HIV-1 and HIV-2.


The *recomline* HIV-1 and HIV-2 IgG is a qualitative test for detection of IgG antibodies against specific antigens. By using the type-specific antigens gp41(HIV-1) and gp36(HIV-2) making it possible to differentiate between an infection with HIV-1 and HIV-2.



- According to researches, HIV-2 was only popular in some areas, for the rest of the world, HIV-1 is more popular.

## Conclusion:

- It was not Mindray who gave false negative result, it was the rapid test and A Brand reagent that gave false positive result.
- CDC HIV diagnosis procedure should be applied in the customer side to avoid false results.




# Case Summary

Background

Ideas

Solution



Summary



As we all known, HIV kit is the 4<sup>th</sup> generation which shorten to around 2 weeks for HIV detection, Mindray HIV kit is also the 4<sup>th</sup> generation which added the p24 antigen so that the window period can be shortened to fourteen days to improve the early detection capability.

In case of new HIV mutation could be generated, which is not sensitive to HIV kit. Antibody sensitivity of Mindray is not enough for HIV infection screening, the performance evaluation indicates that The Mindray Antigen and Antibodies to Human deficiency Virus (CLIA) test had a specificity of 99.5% and sensitivity of 100% for HIV-1 and HIV-2 blood samples,

# Case Summary



# Thanks!

**mindray**迈瑞