Toxoplasma IgM false positive case

IVD International Clinical Application Department

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
Complain


The customer complained they reported a positive result (1,14) and then obtained a negative result (0,18) for a pregnant patient, using the same sample.

They decided to confirm all positive results and now they have three cases of false positive results between April and May. The QC is within limits.

Solicitude

They had an answer from our distributor suggesting following the diagnostic guidelines for Colombia, but they need to know the root cause of the changes in the results.

- If interference is suspected, identify the analyte and non-analyte causes.
- Evaluate all the phases of the analysis.

Verify the preanalytical and sample pretreatment protocol

- The customer used clot activator tubes with gel
- Centrifuge 3500 RPM 15 min
- Storage: (2-8°C); **no aliquots**
- They processed the samples daily
- **The re-test process was many days after the first analysis**, and was just if they had doubts or complaints

Case Ideas

Verify the analytical aspects

1. Calibration

Calibration 1

Datos calibrac

QuimicToxo IgM

Nº de lote231101

N.º ser. reac.001279

Modelo matCOI

Estado calibCalibrado

En medioCL-2000i

Fecha cal23/05/2024 10:41:53

Versión de los parámetrosV1.01.01 / V01.05

Nivel	RLU	Nº de lote	N.º lote sustrato	N.º serie sustrato
C0	9065.261	231101	20231108	07763
C0	9087.451	231101	20231108	07763
C0	9199.368	231101	20231108	07763
C1	69837.048	231101	20231108	07763
C1	71235.043	231101	20231108	07763

Calibration 2

Quimic

Toxo IgM

Nº de lote

231101

N.º ser. reac.

000415

Modelo mat

COI

Estado calib

Extendido

En medio

CL-2000i

Fecha cal

26/03/2024 12:30:26

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V1.01.01 / V01.05

Nivel	RLU	Nº de lote	N.º lote sustrato	N.º serie sustrato
C0	9455.085	231101	20230801	13029
C0	9166.073	231101	20230801	13029
C0	9179.789	231101	20230801	13029
C1	73209.062	231101	20230801	13029
C1	71270.236	231101	20230801	13029
		231101	20230801	13029

Background

Ideas

Solution

Summary

RLU CAL Factory


38		Cutoff	29353
39		C0-RLU	10950
40		C1-RLU	73459
41	Reagent	C0-COI	0,37
42	Lot20231101	C1-COI	2,50
43		校准系数2	0,2944
44		校准系数1	0,7056

Calibration results were normal

Case Ideas

Verify the analytical aspects

2. QC

QC results are within range, however positive level has a **positive bias**.

Verify the analytical aspects

3. Maintenance

The customer had performed the maintenance according to our recommendations.

Verify the post-analytical aspects

The laboratory **didn't have a re-test protocol** and **usually reported the first result**.

The doctor with a positive result solicitates them a Toxoplasma IgG avidity test or confirms the IgM result if is below 1,9 COI.

Possible Reasons of false positive for Toxo IgM

- Increased sensitivity of assays in a population with low prevalence;
- Sample conditions (fibrin, hemolyzed, sample pre-treatment);
- Analyzer maintenance (sample probe exterior cleaning, clean dispersion aspiration probes, daily clean, effect check)
- Others:
 - ✓ Heterophile antibodies;
 - ✓ Autoimmune disorders and other medical conditions;
 - ✓ Multiple transfusions

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
Cases Summary:

1. Three cases with positive results (near to indeterminate zone)
2. Samples were re-test many days after the first result, with negative results
3. After the first test, the samples were storage in the primary tube (2-8°C)
4. The samples were re-test because the laboratory had a complaint from a doctor
5. QC is within range with a little bias

Pre-analytical and analytical phase:

1. Check centrifugation process in the lab, check manufacturer tube recommendations.
2. Before centrifugation, the serum samples must be completely coagulated.
3. If the test cannot be performed within 8 hours, separate the serum sample from the blood collection tube for storage.
4. If the test cannot be completed within 8 hours, the sample should be sealed and refrigerated at 2-8°C; avoid repeated freezing and thawing of the sample.


Pre-analytical and analytical phase:

5. Stablish a re-test protocol, specially for positive samples near to the indeterminate range, and according with doctors' behavior (below 1,9 COI) :

- ✓ Re-exam the samples after centrifugation 10000g/10 min
- ✓ Re-test the same day

6. Make sure the analyzer completed the maintenance, and the controls are handling according with IFU.


Postanalytical phase:

1. Follow national guidelines for Toxoplasmosis diagnostic in pregnant

Conclusion:

- It was a pre-treatment issue, they had to improve centrifugation, and they needed a re-test protocol.
- The customer has to follow-up maintenance protocol and control handling recommendations.
- After our suggestions, the laboratory informed they don't have false positive reports. They established a re-test protocol and change the centrifugation process 4000rpm 15 min.

Thanks!

mindray迈瑞