

mysafety insight

Early Warning Score (EWS) for spot-check monitoring


Clinical Information Leaflet

Early recognition of patient deterioration

Spot-check scenarios can be high-risk as well

The health of patients admitted to hospitals may deteriorate during their hospitalization in general wards ^[1, 2]. Studies have shown that in a large number of patients admitted to critical care departments, life-threatening changes were observed and documented up to 24 hours before their admissions ^[3-7]. These patients may suffer a severe adverse event (SAE), defined as a life - threatening organ failure that could result in death or unplanned admission in the intensive care unit (ICU) ^[8, 9]. Cardiac arrest cases and unplanned admissions to ICU are associated with increased in-hospital mortality rates, and longer hospital stays, rising the number of potentially avoidable deaths ^[7, 10, 11].

These observations and the decisions arising from such early manifestation could improve care and resuscitation outcomes, because most further deteriorations and even death can be prevented with early intervention.



Most inadvertent cases could be preventable ^[3-7].

76%

of critical events occurring in non-ICU patients



57%

sudden cardiac arrests in hospital occurring in general wards

The majority of critical events occur out of the ICU ^[12].

Most inadvertent cases might actually be preventable

Numerous reports, however, indicate that there is a high risk of missing patient deterioration episodes that develop in unexpected adverse events. One of the reasons why patient deterioration may not be detected is the nurse-to-patient ratio and the subsequent frequency of vital sign monitoring, which decreases from higher to lower acuity care units. To confirm this hypothesis, a prospective defined analysis of the UK National Cardiac Arrest Audit (NCAA), collected data from 144 acute hospitals relating to 23,554 patients over the age of 16, showing that most in-hospital cardiac arrests occurred in General Wards, 56.6%, and not the conventional acute care units such as ICU, 5.2%, or coronary care units (CCU), 10.4%. This is an important figure especially considering that the survival rate to discharge after experiencing in-hospital cardiac arrest is only 22.3% ^[12]

Early Warning Score (EWS)

A Systemic Approach To Identify Early Deterioration

In-hospital patient deterioration is often preceded by a period of abnormalities in vital signs. These are changes mainly in physiological parameters like pulse, blood pressure, respiratory rate and temperature^[13]. Based on this premise, in the late 90's several studies were able to develop scores to anticipate these situations, and as a result, Early Warning Scores (EWS) were created to determine the degree of patients' illness based on their physiological parameters^[14]. For example, in the UK, several reports from the National Institute for Health and Clinical Excellence (NICE) and the Royal College of Physicians, have advocated the use of the two most popular EWS systems worldwide, Modified Early Warning Score (MEWS) and National Early Warning Score (NEWS). NEWS2 is the latest version of NEWS, updated in 2017. All these protocols advocate a system to standardise the assessment and response to acute illness.

Until recently there has been a lack of consensus regarding the ideal EWS protocol, but there is evidence that certain parameters

are better than others to identify early deterioration^[15]. Listed below is a description of the physiological parameters included in most EWS systems:

- Respiratory rate is an important indicator of potential respiratory dysfunction.
- Systolic blood pressure, when its high it may indicate cardiovascular disease, while low systolic blood pressure may indicate circulatory compromise.
- Pulse rate, tachycardia may indicate circulatory compromise.
- Level of consciousness: Alert, a fully awake patient; Voice, the patient makes a response to voice; Pain, the patient delivers a response to a pain stimulus; Unresponsive, the patient does not give a response to voice or pain.
- Temperature, if too high or too low is a sensitive indicator of acute illness, especially infection.
- Oxygen saturation is an important parameter for the integrated assessment of pulmonary and cardiac function. Routine monitoring by pulse oximetry is recommended in NEWS and NEWS2 protocols.
- Patient on room air or supplemented oxygen: should be considered for NEWS and NEWS2.

National Early Warning Score 2 (NEWS2)							
Physiological parameters	3	2	1	0	1	2	3
Respiratory rate (permin)	≤8	-	9-11	12-20	-	21-24	≥25
SpO ₂ scale 1(%)	≤91	92-93	94-95	≥96	-	-	-
SpO ₂ scale 2(%)	≤83	84-85	86-87	88-92/ ≥93 on air	93-94(on oxygen)	95-96(on oxygen)	≥97 on oxygen
Air or oxygen?	-	Oxygen	-	Air	-	-	-
Systolic BP (mmHg)	≤90	91-100	101-110	111-219	-	-	≥220
Pulse rate (per min)	≤40	-	41-50	51-90	91-110	111-130	≥131
Level of consciousness	-	-	-	A	-	-	C,V/P or U
Temperature (°C)	≤35.0	-	35.1-36.0	36.1-38.0	38.1-39.0	≥39.1	-
A=Alert C=New confusion (delirium) V=Response to verbal stimulation P=Response to painful stimulation U=Unresponsive							

NEWS2 score	Clinical risk	Frequency of monitoring	Clinical response
Total score 0		Minimum 12 hourly	<ul style="list-style-type: none"> Continue routine NEWS monitoring
Total score 1-4	Low	Minimum 4-6 hourly	<ul style="list-style-type: none"> Inform registered nurse, who must assess the patient Registered nurse decides whether increased frequency of monitoring and/or escalation of care is required
3 in a single parameter	Low-medium	Minimum 1 hourly	<ul style="list-style-type: none"> Registered nurse to inform medical team caring for the patient, who will review and decide whether escalation of care is necessary
Total 5 or more: Urgent response threshold	Medium	Minimum 1 hourly	<ul style="list-style-type: none"> Registered nurse to immediately inform the medical team caring for the patient Registered nurse to request urgent assessment by a clinician or team with core competencies in the care of acutely ill patients Provide clinical care in an environment with monitoring facilities
Total 7 or more: Emergency response threshold	High	Continuous monitoring of vital signs	<ul style="list-style-type: none"> Registered nurse to immediately inform the medical team caring for the patient this should be at least at specialist registrar level Emergency assessment by a team with critical care competencies, including practitioner(s) with advanced airway management skills Consider transfer of care to a level 2 or 3 clinical care facility, ie higher-dependency unit or ICU Clinical care in an environment with monitoring facilities


Automated EWS under spot-check monitoring scenarios

To improve the detection of patient deterioration

It is recommended to use EWS during initial prehospital stage and throughout the entire patient's hospital journey ^[16]. The overall performance of the EWS system is not solely dependent on the scoring system but also the organization of the response. Successful implementation of an EWS in the hospital must go hand in hand with proper education of staff and increasing awareness of the necessity of structural patient monitoring ^[17]. A large retrospective study in the UK, a country where NEWS2 assessment is part of the standard of care in acute and ambulance settings, found that 38% of the scores were not fully completed or incorrectly calculated, hence not triggering important alerts in many cases ^[18].

EWS assessment should be performed as an adjunct to the clinical judgment of the doctor ^[14]. Automatisation of EWS into the vital signs monitoring system has decreased the time required for vital sign measurement and recording, an improvement in the

proportion of rapid-response-team calls triggered by respiratory criteria, as well as an increase in the survival rate of patients receiving rapid-response-team calls ^[19].



EWS should be used throughout the entire patient's hospital journey ^[14].

Fast Intervention with Mindray automated EWS, SmartAlarming and intuitive visuals in the VS-series monitors

Flexible and configurable protocols

Mindray's EWS solution provides flexible and configurable protocols including the standard MEWS, NEWS, and NEWS2, but will also allow users to create and save customised scoring protocols. To better satisfy patient needs, the individual parameter scoring (IPS) allows full control of all parameters and limits for healthcare professionals.

Automated EWS and Smart Alarming

Early warning scores are used to identify the patients at risk. With Mindray's automated EWS, patient's vital signs are automatically measured and the EWS is regularly calculated. Auto calculation of a new score can be triggered by each or all the three following events: pre-set time interval or interval according to the last EWS score, every new NIBP, or a vital sign alarm. The interval can be set by the user in a time range (from 5 minutes to 24 hours) or according to the patient's last EWS score result.

Once a deterioration is detected, the patient EWS monitor will start alarming based on predefined settings to inform the responsible staff about the change in patient condition at an early stage.

Eye-catching Total Score
Refreshed Automatically

Clear Individual Scores
With Alarm Signs



EWS solution on a Mindray VS-series monitor.

Confident care with EWS dashboard in the BeneVision central monitor system



EWS dashboard, SmartAlarm and Trends View tools on Mindray BeneVision central monitoring system.

The BeneVision central monitoring system allows complete digital surveillance across Mindray patient monitors and vital sign monitors. View all patient data at a glance, in real-time, for a streamlined clinical workflow that improves visibility throughout the patient journey.

Under EWS Trends View, the changes in a patient's status are highlighted to allow early detection of deterioration and make appropriate recommendations.

Parameters and EWS View show details of patient vital signs and EWS. EWS countdown reminds caregivers to do the next round spot check on time.

The intuitive Early Warning Scoring dashboard provides caregivers with a fast and easy-to-understand assessment of patient deterioration, allowing faster, more confident responses and earlier interventions.

Aiming to create safer patient environments, Mindray incorporates early warning scoring (EWS) system in a wide range of patient monitors, from low to high acuity. By including EWS system, Mindray solutions contribute to safer and more efficient patient management by anticipating potential complications and improving workflows.

References:

1. Fuhrmann, L., et al., Incidence, staff awareness and mortality of patients at risk on general wards. *Resuscitation*, 2008. 77(3): p. 325-30.
2. McQuillan, P., et al., Confidential inquiry into quality of care before admission to intensive care. *Bmj*, 1998. 316(7148): p. 1853-8.
3. Zografakis-Sfakianakis, M., et al., The value of the Modified Early Warning Score for unplanned Intensive Care Unit admissions of patients treated in hospital general wards. *Int J Nurs Pract*, 2018. 24(3): p. e12632.
4. Sandroni, C., et al., In-hospital cardiac arrest: incidence, prognosis and possible measures to improve survival. *Intensive Care Med*, 2007. 33(2): p. 237-45.
5. Hillman, K.M., et al., Duration of life-threatening antecedents prior to intensive care admission. *Intensive Care Med*, 2002. 28(11): p. 1629-34.
6. Schein, R.M., et al., Clinical antecedents to in-hospital cardiopulmonary arrest. *Chest*, 1990. 98(6): p. 1388-92.
7. Goldhill, D.R., S.A. White, and A. Sumner, Physiological values and procedures in the 24 h before ICU admission from the ward. *Anaesthesia*, 1999. 54(6): p. 529-34.
8. Green, A.L. and A. Williams, An evaluation of an early warning clinical marker referral tool. *Intensive Crit Care Nurs*, 2006. 22(5): p. 274-82.
9. Story, D.A., et al., The effect of critical care outreach on postoperative serious adverse events. *Anaesthesia*, 2004. 59(8): p. 762-6.
10. Quinn, T.D., et al., Analysis of Unplanned Postoperative Admissions to the Intensive Care Unit. *J Intensive Care Med*, 2017. 32(7): p. 436-443.
11. Goldhill, D.R., et al., The patient-at-risk team: identifying and managing seriously ill ward patients. *Anaesthesia*, 1999. 54(9): p. 853-60.
12. Nolan, J.P., et al., Incidence and outcome of in-hospital cardiac arrest in the United Kingdom National Cardiac Arrest Audit. *Resuscitation*, 2014. 85(8): p. 987-92.
13. Gerry, S., et al., Early warning scores for detecting deterioration in adult hospital patients: a systematic review protocol. *BMJ Open*, 2017. 7(12): p. e019268.
14. Alam, N., et al., The impact of the use of the Early Warning Score (EWS) on patient outcomes: a systematic review. *Resuscitation*, 2014. 85(5): p. 587-94.
15. Prytherc, D.R., et al., ViEWS--Towards a national early warning score for detecting adult inpatient deterioration. *Resuscitation*, 2010. 81(8): p. 932-7.
16. Physicians, R.C.o., National early warning score (NEWS): standardizing the assessment of acute-illness severity in the NHS. 2012.
17. Petersen, J.A., Early warning score challenges and opportunities in the care of deteriorating patients. *Dan Med J*, 2018. 65(2).
18. Clifton, D.A., et al., 'Errors' and omissions in paper-based early warning scores: the association with changes in vital signs--a database analysis. *BMJ Open*, 2015. 5(7): p. e007376.
19. Bellomo, R., et al., A controlled trial of electronic automated advisory vital signs monitoring in general hospital wards. *Crit Care Med*, 2012. 40(8): p. 2349-61.

healthcare within reach

MINDRAY SOCIAL MEDIA CHANNELS

LinkedIn



www.mindray.com

P/N: ENG-Early Warning Score (EWS) for spot-check monitoring-210285X8P-20211008
©2021 Shenzhen Mindray Bio-Medical Electronics Co.,Ltd. All rights reserved.

mindray
healthcare within reach