mindray

Articulating Endoscopic Linear Cutter

Instruction Manual



Description

Articulating Endoscopic Linear Cutter Product Name: Specifications and See main text for details Model: Structural See main text for details Composition:

The instruments are intended for transection, Intended purpose: resection, and/or creation of anastomoses. Registrant/ Hangzhou Mindray Medical Technology Co., Ltd.

Manufacturer Name: Registrant/ 2 Fengxiang Road, Tonglu Economic Development Zone, Tonglu County, Hangzhou Manufacturer Address: City, Zhejiang Province Production Address: 2 Fengxiang Road, Tonglu Economic Development Zone, Tonglu County, Hangzhou City Zheijang Provin



Note

• Before operating this instrument, it is necessary to be familiar with and understand the meanings of the symbols mentioned above.

Product Overview

Structural Composition

Articulating Endoscopic Linear Cutter consists of a staple holder, joint head, rod, joint head knob, rotary wheel, retraction button, release button, firing safety button, closure handle, cartridge holder, and removable anvil tip.

Articulating Endoscopic Linear Cutter



1. Staple driver holder	2. Articulating head	3. Rod	4. Articulating head knob
5. Rotary wheel	6. Release button	7. Retraction button	8. Firing safety
9. Closure handle	10.Staple cartridge holder	11. Removable Anvil Tip	

Figure 1 Articulating Endoscopic Linear Cutter

Warning

- The stapler is only intended for use on the same patient during the same surgical procedure.
- The stapler can only be used in conjunction with the cartridge.
- Before use, if rust is found on the product, its use should be strictly prohibited. Visually inspect the staple holder for deformation before firing. If deformation
- is observed, do not use it. This product is intended for use in sterile operating rooms within medical facilities.

Note

· After each stapling and suturing, remove the empty staple cartridge from the staple cartridge holder.

Caution

- The product is intended for use in adults.
- Please check the packaging of this product carefully before use and discontinue using it if the packaging is damaged.
- · This product is sterilized using ethylene oxide and is intended for clinical use as a sterile product.
- Before use, please check whether the product is within the use period. Sterilization is effective for 3 years. Using the product past its use period is strictly prohibited.
- · Minimally invasive procedures should be performed by individuals with adequate training and familiarity with minimally invasive surgical techniques. The medical literature on the technique and its complications and hazards should be reviewed prior to performing any minimally invasive procedure.
- The dimensions of minimally invasive devices may vary between manufacturers. If minimally invasive surgical instruments and their accessories made by different manufacturers are applied at the same time in a single procedure, it is important to verify their compatibility before the procedure.
- Always check that the staple cartridge holder is stable before firing.
- Always check for hemostasis at the suture site, anastomosis integrity, and leakage after firing.
- Ensure that the tissue thickness is within the specified range and that the tissue is evenly distributed within the stapler device. Too much tissue on one side can cause a poor anastomosis, with the possibility of anastomotic leakage.
- Attempts to forcefully squeeze the closure handle in cases of excess or thick tissue may result in an incomplete suture, with the risk of anastomotic splitting or leakage. In addition, damage to the instrument and failure to fire may occur.
- Pre-surgical radiotherapy may result in tissue alterations. For example, these alterations may cause tissue thickening beyond that specified for the selected anastomotic staple. Any pre-surgical treatment of a patient should be carefully considered and a change in surgical technique or procedure may be required.
- Firing must be completed in one strike; never partially fire the instrument. Incomplete firing may result in improper staple molding, incomplete cut lines, bleeding and leakage from the suture, or difficulty in removing the instrument. Be sure to fire completely to ensure that the sutures are properly formed and the tissue is cut correctly.
- Squeezing the closing handle exposes the blade. Do not press the closure handle repeatedly, as this can lead to anastomotic site damage.
- When inserting the cartridge, ensure that the safety is in the closed position to avoid inadvertent activation of the closure handle, which could result in accidental exposure of the blade and premature partial or full deployment of the anastomotic staple.
- Instruments that have been in contact with body fluids should undergo special handling to prevent biological contamination from occurring.
- This product is sterilized and packaged for one use only. Use in multiple patients may jeopardize the integrity of the product, or there may be a risk of contamination, which in turn may lead to patient harm.
- The stapler can switch cartridges and fire up to 12 times during the same procedure. The use of this device with suture reinforcement material may result in a lower number of strikes.
- This product is a single-use device. It must be destroyed after use so that its parts are no longer functional, and it must be sterilized and rendered harmless for disposal.
- After this product is used, the anastomotic staple stays in the body and a 1.5 T MRI and 3.0 T MRI can be performed.
- Before each use, check the surface of the instrument for protrusions, rough surfaces or sharp edges that could cause injury to the patient.
- Do not use this device on parenchymal organs that may be damaged under compression (e.g., organs such as the liver or spleen).
- The titanium implant staples are safe under specific conditions: Non-clinical trials have shown that endoscopic linear cutter staplers and titanium staples (pure titanium) in the cartridge are safe under MR-specific conditions. Patients with implanted staples can be safely scanned under the following conditions:

The strength of the static magnetic field is less than or equal to 3.0 T; The highest spatial gradient of the magnetic field is 12.0 T/m;

The maximum body-averaged specific absorption rate (SAR) reported by the MR system after 15 minutes of scanning (per pulse sequence) was 2.7 W/kg.

MRI-related warming:



- 4. If necessary, use your index finger to push up or down on the rotary wheel to rotate the jaws, and the instrument bar will rotate freely in any direction.
- To rotate the jaws within the body cavity, turn the articulating head knob to rotate 5. the jaws to the desired position and release the knob to achieve a jaw bend.





Fig. 10

- 6. Place the jaws of the instrument around the tissue to be stapled.
- 7. After positioning the instrument jaws, snap the closing handle until it locks to close the jaws. A "click" sound indicates that the closing handle and jaws are locked. If improperly positioned, the release button can be toggled to return the jaws to their open state, from which they can be repositioned.
- 8. After positioning is completed, snap the closing handle fully shut and press the firing safety button to prepare for firing. Snap the closing handle several times in succession until it reaches the cutting termination position.

Warning

- If the firing safety button is not pressed, the closing handle cannot be snapped shut and the instrument cannot be fired.
- Do not press the firing safety button until you are ready to fire, to prevent it from being falsely triggered.
- Make sure the jaws are in the open position when turning the articulating head knob; otherwise the jaws may fail to rotate.
- When turning the articulating head knob, make sure that no external force is applied to the jaws; otherwise damage to the instrument may result.
- If the instrument does not fire smoothly or cannot be fired, please replace the instrument. When experiencing a high firing force and attempting force firing, there may be a strange noise or a "dropping" sensation during firing; in this case, stop the instrument immediately. It is recommended to replace the instrument with a new one,

Caution

Caution

- If the jaws do not open after toggling the release button, the stapler jaws can also be assisted to fully reset by pulling the return knob toward the proximal end
- After confirming the position of the jaws, the target tissue must be checked to ensure it is free of stents, guidewires, and other obstructions in the area covered by the cartridge; otherwise, there is a risk of incomplete cutting and stapling, poor anastomotic staple shaping, or inability to open the jaws after stapling.
- It is necessary to select a puncture device of no less than φ 12 size as a converter.
- Incomplete firing may result in improperly formed sutures, incomplete cutting lines, bleeding, and/or difficulty in removing the instrument.
- The number of consecutive snaps of the closing handle is based on the length of the assembly's suture line.

		City, Zhejiang Frovince
	Phone:	0571-58504222
(6 0197	Fax:	0571-58504300
	URL:	www.mindray.com
	Post Code:	311508
	Date of Manufacture:	See label
	Use-by Date:	See label
	lssue Date of This User Manual:	2024-07
		Shanghai International Holding Corp. GmbH (Europe)
EC REP	Address: Tel:	Eiffestrasse 80, 20537 Hamburg, Germany 0049-40-2513175/2513174

Statement

Thank you for purchasing this product.

Before using the product, please read the contents of this user manual carefully to ensure the correct use of the product.

Please keep this user manual properly after reading, so that you can refer to it whenever necessary

The intellectual property rights of this product and its instruction manual belong to Hangzhou Mindray Medical Technology Co., Ltd. (hereinafter referred to as Mindray). Mindray has right to final interpretation of this instruction manual. No individual or organization shall copy, modify or translate this instruction manual without the written permission of Mindray. This manual describes the use, function and operation of the product in detail. Before using this product, please read and understand the contents of this manual fully to ensure the correct use of this product and the safety of patients and operators.

The user manual is only for product usage instructions, and should not be considered as a reference for surgical techniques.

Mindray shall be responsible for the safety, reliability and performance of the product provided that the following conditions are met:

 This product is used in accordance with the Instruction Manual. Product damage is caused by non-human factors (human factors refer to accidental drops, deliberate damage, etc.).

If you really need to return the product to Mindray, please contact Mindray's aftersales service department and provide the product model, serial number, and a brief explanation of the reason for the return.

The product's "triple guarantee" and the after-sales service are defined by the service contract between the distributor and the manufacturer.

After-sales service provider: Hangzhou Mindray Medical Technology Co., Ltd. After-sales provider's address: 2 Fengxiang Road, Tonglu Economic Development Zone, Tonglu County, Hangzhou City, Zhejiang Province After-sales provider's telephone: 0571-58504222

To the Customer: Important Information

Safety Signs and Definitions

The safety information contained in this manual helps users identify potential hazards and avoid danger. This manual uses three signs to highlight potential hazards:Warning, Caution, and Note.

This manual provides safety information regarding instrument use. Please read and understand all instructions completely before use. Failure to comply with warnings, cautions, and notes related to this instrument, may result in personal injury, instrument malfunction, and void warranty.

Symbol	Description
Warning	In terms of operation or maintenance procedures, failure to strictly adhere to relevant requirements may result in injury or death.
Caution	In terms of operation or maintenance procedures, failure to strictly adhere to relevant requirements may result in equipment damage or destruction.
Note	Operational tips or maintenance advice; failure to follow these may lead to errors.

Symbol Description

The stapler should have good stapling and cutting performance, with replaceable components. The cutting edges after each stapling should be neat, without burrs. The stapler should be able to perform stapling and cutting operations at least 12 times without affecting its functionality.

Intended purpose

The instruments are intended for transection, resection, and/or creation of anastomoses.

Intended users

Persons having adequate training and familiarity with minimally invasive techniques.

Indications:

The product is intended for transection, resection, and/or creation of anastomoses. The instruments have application in multiple open or minimally invasive general abdominal, gynecologic, thoracic, and pediatric surgical procedures.

Intended patient population

adult

Intended medical conditions

Medical Institution

Product Properties

- The stapler is a single-use, sterile product.
 - The connection between the stapler body and the assemblies should be firm and reliable, while being easy to replace.
- The stapler should have a protective device against empty staple cartridges, and should not be able to fire when an empty staple cartridge is mistakenly installed.
- The stapler should have a protective device, and the firing rod of the stapler body should not be able to fire when the staple cartridge is not loaded.
- The stapler can be installed with the staple cartridge repeatedly in a single surgery, with a maximum firing of 12 times.

Contraindications

- Do not use the instruments with Gray reload (2.0mm staple) on any tissue that compresses to less than 0.75mm in thickness, or on any tissue that cannot comfortably compress to 1.5mm.
- Do not use the instruments with White reload (2.5mm staple) on any tissue that compresses to less than 1.0mm in thickness, or on any tissue that cannot comfortably compress to 2.0mm.
- Do not use the instruments with Blue reload (3.5mm staple) on any tissue that compresses to less than 1.5mm in thickness or on any tissue that cannot comfortably compress to 2.4mm.
- Do not use the instruments with Golden reload (3.8mm staple) on any tissue that requires excessive force or compresses to less than 1.8mm in thickness, or on any tissue that cannot comfortably compress to 3.0mm.
- Do not use the instruments with Green reload (4.1mm staple) on any tissue that compresses to less than 2.0mm in thickness, or on any tissue that cannot comfortably compress to 3.3mm.
- Do not use the instruments with Dark reload (4.2mm staple) on any tissue that compresses to less than 2.3mm in thickness, or on any tissue that cannot comfortably compress to 4.0mm.
- The device should not be used to staple ischemic or edematous tissues. Do not use the instruments on the aorta.
- Do not reuse after surgery.
- Do not use curved tip staplers on tissues or strctures that can't fit completely within the jaws proximal to the transitional angle of the curved tip Staplers is not intended for use inside of the heart, central circulatory or central nervous system.

Side-effects

According to clinical and residual risk evaluations, Endoscopic Linear Cutter Reloads have no known side effects for the intended patients.

Intended Clinical Benefits

In nonclinical testing, when sutures were subjected to 15-minute sequential MRI scans (per sequence) performed at 3.0 T using a transmit-receive body RF coil, anastomotic staples produced no more than 1.8°C of warming at a maximum body-averaged specific absorption rate of 2.7 watts/kg.

Artifact information:

Image quality may be affected if the area imaged by the MRI is completely identical to the area where the sutures fired by the endoscopic linear cutter stapler and cartridge are located, or if these areas are located in close proximity to each other. Therefore, it is necessary to optimize the settings of the MRI parameters to compensate for the effects caused by these staples. The size of the signal-free area of the suture (e.g., 60 mm) in a worst-case scenario.

Pulse sequence	SE	SE	GRE	GRE
Planar orientation	Horizontal	Vertical	Horizontal	Vertical
Signal void size (mm ²)	572	68	934	89

Inspection before use

Instructions for Use

1. Check that the size of the staple cartridge matches the instrument to be used.

- 2. Use your index finger to turn the rotary wheel and rotate the jaws of the instrument by 360°. Use your other hand to turn the articulating head knob, to familiarize yourself with the turning action of the instrument. The jaws can be turned to a maximum angle of 60°. The greater the bending angle, the greater the force that may be required to turn the articulating head knob. The jaws can be returned to a straight position by reversing the articulating head knob.
- 3. Check the compatibility of all instruments and accessories before use.

Installation and removal of anvil tip attachment

The anvil tip attachment is optional. It must be installed on the instrument's staple cartridge holder as required prior to use.

- 1. Remove the anvil tip from its sterile packaging and ensure that the jaws are in an open position before installation.
- 2. Holding the head of the anvil tip in your hand, insert it into the T-shaped cutter slot of the staple cartridge holder until you hear a "click."







- 3. To remove the anvil tip, align the opening of the hand-held removal tool with the anvil tip (jaws open). Once the anvil tip is in the opening, continue to push the hand-held removal tool until it cannot be pushed further.
- 4. Remove the hand-held removal tool by moving it away from the staple cartridge holder.
- 5. Flip the hand-held removal tool over and the anvil tip will fall out of the removal tool window, allowing the anvil tip to be mounted and used again.
- 6. When removing the anvil tip using the cartridge protector cover, hold the cartridge protector cover in your hand, insert the column of the head of the cartridge protector cover tightly against the anvil tip into the T-shaped cutter slot, and pry it upward to remove it from the staple cartridge holder.



9. After the instrument is fully fired, pull the retracting knob fully proximally to open the jaws, and carefully remove the instrument from the tissue. After cutting and suturing, the anastomotic sutures should be checked. Small volumes of bleeding can be treated with an electrotome or hand-sutured to stop the bleeding. Check for anastomotic integrity and for anastomotic leakage.

10. Return the articulating head knob to the center position to keep the jaws straight. Trigger the closure handle until it locks and close the stapler jaws. After removing the instrument from the abdominal cavity, toggle the release button to open the jaws and unload the cartridge from the stapler.



• The jaws cannot be used for blunt dissection after installing the anvil tip attachment

Fig. 11

- Do not touch tissue when using the articulating head knob for bending or straightening the jaws, as this may cause tissue damage.
- Ensure that the tissue is appropriately positioned between the jaws, as "lumpy or nodular" tissue along the cartridge, especially at the instrument jaw bifurcation, may result in an incomplete suture.
- The black line at the end on the staple cartridge holder and cartridge holder indicates the end of the stapling line. The indicator line labeled "cut" on the cartridge holder indicates the cutting line on the instrument.
- After confirming the position of the jaws, the target tissue must be checked to ensure it is free of stents, guidewires, and other obstructions in the area covered by the cartridge; otherwise, there is a risk of incomplete cutting and stapling, poor anastomotic staple shaping, or inability to open the jaws after stapling.
- Make sure that the tissue does not extend (extrude) beyond the proximal black line of the stapling mark. If the tissue is extruded next to the black line, the tissue may only be transected and not stapled.
- Failure of the jaws to close means that the target tissue is extruded beyond the safe anastomotic stapling range, so the jaws should be opened and repositioned to reduce the amount of tissue covered by the jaws, and then closed again. If the jaws still fail to close, then check again to see if the cartridge has been selected to suit the thickness of the tissue (see Models and Specifications for selection of the appropriate cartridge). If it is determined that the cartridge has been selected incorrectly, then the stapler should be removed immediately and the procedure discontinued. If it is determined that the cartridge has been incorrectly selected, the stapler should be removed immediately and the procedure discontinued.
- If the closing handle is difficult to lock, it is necessary to reposition the instrument, incorporate a smaller amount of tissue, and ensure that the correct cartridge specifications have been selected.
- If the clamping device is inoperable and the jaws are unable to hold the tissue, do not fire the instrument. The instrument should be removed and its use discontinued.
- When inserting and removing the instrument, the jaws of the instrument must be straight and parallel to the instrument shaft. If the jaws of the instrument are not in a straight position, it may be difficult to insert or withdraw the instrument, which may damage it.
- If it is not possible to reset to the initial state, press the release button while pulling hard on the return knob.
- The stapler can switch modules and fire up to 12 times during the same procedure.
- A single-use safety device will be in effect for components that have been fired, preventing the component from being fired a second time. Do not attempt to force another strike as this will result in damage to the instrument or tissue damage.

Compatible Instruments

The Articulating Endoscopic Linear Cutter are intended to be used in combination with the devices as below:

- Trocar: User needs to operate the Articulating Endoscopic Linear Cutter as well
- Fig. 4

The following table contains all symbols used on this instrument and labels, along with their meanings:

Table 1 Symbol Description



The use of cutter for vascular or tissue ligation provides the benefit of successful surgical outcomes.

User Manual

Warning

- When this product is used, aseptic operational practices should be strictly implemented.
- The user of the device must be a clinician with adequate training and familiarity with minimally invasive surgical techniques, and must be familiar with the procedures for operating this product.
- If other technical means (e.g., electrocautery) are to be used during the surgical procedure, follow the precautions described by the manufacturer of the original anastomotic stapler and cartridge in order to avoid associated hazardous conditions.
- Avoid using the stapler in close proximity to or stacked on top of other equipment. If this is unavoidable, the position of the stapler in relation to the other equipment must be checked to ensure that they all work properly.
- Minimally invasive instruments made by different manufacturers may differ in diameter. If both minimally invasive instruments and accessories from different manufacturers must be used during a single procedure, check for compatibility before the procedure begins.
- Staplers are packaged and sterilized prior to shipment from the factory. They are intended for single use only and may not be reused, reprocessed or resterilized. Reuse, reprocessing or resterilization may affect the structural integrity of the product and lead to product failure, which can result in injury or death. Reuse, reprocessing or resterilization of the stapler poses a risk of contamination, which can lead to cross-infection in different body parts in the same patient or between different patients, which can result in injury or death to the patient(s).
- No parts or components of our disposable endoscopic linear cutter staplers and assemblies for endoscopy may be replaced without the permission of Hangzhou Mindray Medical Technology Co., Ltd.
- Check whether the cartridge model matches the stapler model to be used (for details, refer to the stapler and cartridge configuration table in the instruction manual).
- Once opened, the stapler must not be resterilized for use, regardless of whether it was used or not.
- Do not use any endoscopic linear cutter stapler for major blood vessels if proximal or distal control is not provided.



Fig. 5

- Caution
- When installing the anvil tip, pay attention to the insertion direction. Make sure the anvil tip is securely installed before use.

Operating steps

- 1. Close the jaws of the instrument by snapping the closing handle until it locks into place. A click will be heard, indicating that the closing handle and jaws are locked.
- 2. Visually inspect the stapler to ensure that the cartridge is properly secured. The device is introduced into the body cavity through a suitably sized puncture device or incision. When using a puncture device, the instrument jaws must be visible as they pass through the puncture device's cannula before opening the jaws.
- 3. After entering the abdominal cavity, the release button is toggled and the stapler jaws open.



as the Endoscopic Linear Cutting Reloads through trocar to perform surgery. The diameter of adaptive trocar should not smaller than φ 12mm.

Reload: User needs to insert the Endoscopic Linear Cutting Reloads into the Articulating Endoscopic Linear Cutter prior to use. The staple line length of the Endoscopic Linear Cutting Reloads should be compatible with the jaw length of the Articulating Endoscopic Linear Cutter. The Articulating Endoscopic Linear Cutter should only be used with the following Endoscopic Linear Cutting Reloads

Table 2 Model and Specification of Endoscopic Linear Cutter Reloads

Style	Model No	Specification	Stapling length L1	Tolerance	Staple height H	Tolerance
JUNC	model No.		(mm)	(mm)	(mm)	(mm)
	MR30-W30A	30-2.5	35.0		2.5	
	MR30-B30A	30-3.5	35.0		3.5	
	MR30-M45A	45-2.0	48.0		2.0	
	MR30-W45A	45-2.5	48.0		2.5	
	MR30-B45A	45-3.5	48.0		3.5	
	MR30-D45A	45-3.8	48.0		3.8	
MD20	MR30-G45A	45-4.1	48.0		4.1	10.2
IVIK50	MR30-T45A	45-4.2	48.0		4.2	±0.2
	MR30-M60A	60-2.0	60.0		2.0	
	MR30-W60A	60-2.5	60.0		2.5	
	MR30-B60A	60-3.5	60.0		3.5	
	MR30-D60A	60-3.8	60.0		3.8	
	MR30-G60A	60-4.1	60.0		4.1	
	MR30-T60A	60-4.2	60.0		4.2	
	MR31-W30A	30-2.5	35.0		2.5	±0.2
	MR31-B30A	30-3.5	35.0		3.5	
	MR31-M45A	45-2.0	48.0		2.0	
	MR31-W45A	45-2.5	48.0		2.5	
	MR31-B45A	45-3.5	48.0		3.5	
	MR31-D45A	45-3.8	48.0		3.8	
MR31	MR31-G45A	45-4.1	48.0		4.1	
	MR31-T45A	45-4.2	48.0		4.2	
	MR31-M60A	60-2.0	60.0	1	2.0	
	MR31-W60A	60-2.5	60.0	1	2.5	
	MR31-B60A	60-3.5	60.0		3.5	
	MR31-D60A	60-3.8	60.0		3.8	
	MR31-G60A	60-4.1	60.0		4.1	
	MR31-T60A	60-4.2	60.0		4.2	

-							
[MR32-W30A	30-2.5	35.0		2.5	
		MR32-B30A	30-3.5	35.0	±2	3.5	±0.2
		MR32-M45A	45-2.0	48.0		2.0	
		MR32-W45A	45-2.5	48.0		2.5	
		MR32-B45A	45-3.5	48.0		3.5	
		MR32-D45A	45-3.8	48.0		3.8	
	MD22	MR32-G45A	45-4.1	48.0		4.1	
	IVIRGZ	MR32-T45A	45-4.2	48.0		4.2	
		MR32-M60A	60-2.0	60.0		2.0	
		MR32-W60A	60-2.5	60.0		2.5	
		MR32-B60A	60-3.5	60.0		3.5	
		MR32-D60A	60-3.8	60.0		3.8	
		MR32-G60A	60-4.1	60.0		4.1	
		MR32-T60A	60-4.2	60.0		4.2	
[MR33	MR33-T60A	60-4.2	60.0	±2	4.2	±0.2

Storage Conditions

Storage Conditions

- Temperature range: 0°C to 35°C
 Humidity range: 10% to 80%
 Atmospheric pressure range: 50 kPa to 106 kPa

Caution

• The instrument should be stored in a dark, cool, fireproof, rat-proof, insectproof, non-corrosive gas, well-ventilated, dry and clean room, and should be protected from compression, wear and impact.

Transportation Conditions

- Temperature range: -10°C to 50°C
 Humidity range: 10% to 60%
 Atmospheric pressure range: 50 kPa to 106 kPa

Caution

• If the instrument needs to be transported, it should be carried out according to the requirements in the purchase contract, and it must be protected from severe impact during transportation, rain and exposure to sun. During transportation, severe vibration and damp environments should be avoided when handling the instrument.

Operational conditions:

- Temperature range: 5°C to 40°C
- Humidity range: 30% to 75%
- Atmospheric pressure range: 80 kPa to 106 kPa

Batch Number, Production Date, Use-by Date

- Batch Number: See label
- Date of Manufacture: See label
- Expiration Date: See label
- Use-by Date: 3 years.

Notification of Incident

As a health care provider, you may report the occurrence of certain events to Hangzhou Mindray Medical Technology Co., Ltd. and possibly to the competent authority of the member state in which the user and / or patient is established. These events include device-related death and serious injury or illness. In addition, as part of our Quality Assurance Program, Hangzhou Mindray Medical Technology Co., Ltd. requests to be notified of device failures or malfunctions. This information is required to ensure that Hangzhou Mindray Medical Technology Co., Ltd.. provides only the highest quality products.

The products must be decontaminated prior to disposal. The products must be disposed of in compliance with the local regulations. If you have any questions concerning disposal of the equipment, please contact Mindray.

Model and Specification

Table 1 Model and Specification of Articulating Endoscopic Linear Cutter

Style	Model Type	Specification	Rod Length L	Tolerance	Bending Angle	Tolerance
Style			(mm)	(mm)	W (°)	(°)
	MS30-C30A	280	192			
	MS30-S30A	340	252			
	MS30-L30A	440	352			
	MS30-C45A	280	192			
MS30	MS30-S45A	340	252	± 5	60	± 10
	MS30-L45A	440	352			
	MS30-C60A	280	192			
	MS30-S60A	340	252			
	MS30-L60A	440	352			
	MS31-C30A	280	192			± 10
	MS31-S30A	340	252		60	
	MS31-L30A	440	352	± 5		
	MS31-C45A	280	192			
MS31	MS31-S45A	340	252			
	MS31-L45A	440	352			
	MS31-C60A	280	192			
	MS31-S60A	340	252			
	MS31-L60A	440	352			
	MS32-C30A	280	192			
	MS32-S30A	340	252			
	MS32-L30A	440	352			
	MS32-C45A	280	192			
MS32	MS32-S45A	340	252	± 5	60	±10
	MS32-L45A	440	352			
	MS32-C60A	280	192			
	MS32-S60A	340	252			
	MS32-L60A	440	352			
MS37	MS37-L60A	440	352	± 5	60	±10