


Patient Admission & Discharge Quick Guide

Admitting a patient

- Select **Patient Management** quick key, or on patient discharge screen
- Enter patient details in **Patient Management** menu



- Select **OK**

Editing patient information

Either:

- Select the patient information area at the top left corner of the screen
- Select the **Patient Management** quick key
- Select the **Main Menu** quick key → from the **Patient Management** column select **Patient Management**

- Edit information as required
- Select Ok to enable changes



Discharging a patient

Either:

- Swipe down the touchscreen with two fingers




These guidelines do **NOT** replace the instructions for use and all users should refer to the appropriate **Operator's Manual** for detailed instructions. The illustrations may not necessarily reflect the setup or data displayed on your device.

- Select the Patient Management quick key  → select **Discharge Patient**
- Select the **Discharge Patient** quick key 
- Select the **Main Menu** quick key → from the **Patient Management** column select

Discharge

Select an item from the **Discharge Patient** menu:



- **Print End Case Report:** discharges the patient and prints the end case report
- **Discharge:** all patient data, including patient information and data is deleted. The monitor loads the default configuration and enters Standby Mode
- **Clear Patient Data:** all patient data is deleted from the monitor and the current monitor settings are used

The **Discharged** screen appears



These guidelines do **NOT** replace the instructions for use and all users should refer to the appropriate **Operator's Manual** for detailed instructions. The illustrations may not necessarily reflect the setup or data displayed on your device.