

# Accutorr 7 with Cerner CareAware VitalsLink



## Connect. Capture. Send.

Mindray and Cerner® partner to deliver quality care more efficiently. Automate the collection and documentation of spot check vital signs data with the Mindray Accutorr 7 monitor and Cerner CareAware VitalsLink® to save time, improve clinical workflow and reduce errors.

Meeting the needs of busy clinical environments, the Accutorr 7 is an easy-to-use vital signs monitor that is ready to go when you are. Assess patients' vitals with NIBP, SpO<sub>2</sub> and temperature measurements. Mindray's medical grade barcode scanner scans both the clinician ID and patient ID for quick validation.



## Improve patient outcomes through increased efficiency

Automated workflow helps free clinicians to focus more on patient care. At the patient's side, clinicians can validate patient information and automatically send verified vital signs directly into Cerner Millennium® from the Accutorr 7 (wireless or hardwired).



1. Scan clinician badge
2. Scan patient ID on the patient wrist band
3. Acquire vitals (NIBP, SpO<sub>2</sub>, temperature)
4. Enter additional measurements (manual parameters and modifiers)
5. Select 'Save' and vitals are automatically sent to the EMR, from the point-of-care

Mindray North America  
800 MacArthur Blvd., Mahwah, NJ 07430  
Tel: 800.288.2121 Tel: 201.995.8000 Fax: 800.926.4275 [www.mindray.com](http://www.mindray.com)

Certain content © Cerner Corporation. All rights reserved. Cerner product images and information are used with permission from the Cerner Corporation and are for demonstration purposes only. Mindray® is a registered trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd. ©2019 Mindray DS USA, Inc. All rights reserved. Subject to change. 09/19 P/N: 0002-08-6424 Rev A

**mindray**  
NORTH AMERICA