

A Series

ADVANTAGE

Trust the Mindray Advantage to help you improve your financial and environmental footprint during anesthesia delivery



A7

ADVANTAGE

The A7 Advantage Anesthesia Workstation further expands A-Series capabilities with electronic fresh gas technology to ensure precise flow dynamics and low flow anesthesia delivery.

Sample gas is recycled from the AGM back into the breathing system, improving agent cost savings by reducing the amount of fresh gas and agent used.

Recycling sample gas, along with the use of the Optimizer function, helps you reduce agent usage along with minimizing the environmental impact of exhausted inhalation agent.