

## V Series Interface Part Number Guide

| <b>Camino Interface</b>                          | <b>Part Number</b>                                |
|--|---|
| <b>Camino Advanced Monitor (CAM01)</b>           | Must be obtained by INTEGRA (Camino manufacturer) |
| <b>Camino Bedside Monitoring Cable (PMIOMPM)</b> | Must be obtained by INTEGRA (Camino manufacturer) |
| <b>Camino PAC-1 IBP/ICT Pre-Amp Cable</b>        | Must be obtained by INTEGRA (Camino manufacturer) |
| <b>Camino Adapter/Cable</b>                      | 115-016069-00                                     |

| <b>Vigilance CCO/SVO2 Interface</b>  | <b>Part Number</b>                                      |
|--|---|
| <b>Vigilance</b>   | Must be obtained by Edwards Lifesciences (manufacturer) |
| <b>Serial Cable: Null Modem/ Female to Female</b>  | 0012-00-1275-01   |
| <b>VDI Module</b>  | 0998-00-1802-0801A                                      |
| <b>Analog Cable (Un-Terminated: Hospital Terminates per Clinical Engineering)</b>  | 0012-00-1801-02   |
| <b>Analog Cable (Terminated)</b>   | 0012-00-1801-01   |
| <b>Analog Cable Adapters (Required with Terminated Cable)</b> <ul style="list-style-type: none"> <li>• 2-3 Adapters required depending on model</li> </ul> | 022-000113-00   |

| <b>Vigilance II CCO/SVO2 Interface</b>   | <b>Part Number</b>                                      |
|--|---|
| <b>Vigilance II</b>  | Must be obtained by Edwards Lifesciences (manufacturer) |
| <b>Serial Cable: Null Modem/Male to Female</b>   | 009-002237-00   |
| <b>VDI Module</b>  | 0998-00-1802-0801A                                      |
| <b>Analog Cable (Un-Terminated: Hospital Terminates per Clinical Engineering)</b>  | 0012-00-1801-02   |
| <b>Analog Cable (Terminated)</b>   | 0012-00-1801-01   |
| <b>Analog Cable Adapters (Required with Terminated Cable)</b> <ul style="list-style-type: none"> <li>• 2-3 Adapters required depending on model</li> </ul> | 022-000113-00   |

| <b>Vigileo CCO/SVO2 Interface</b>              | <b>Part Number</b>                                      |
|--|---|
| <b>Vigileo</b>                                 | Must be obtained by Edwards Lifesciences (manufacturer) |
| <b>Serial Cable: Null Modem/Male to Female</b> | 009-002237-00   |
| <b>VDI Module</b>                              | 0998-00-1802-0801A                                      |

| <b>EV1000 CCO/SVO2 Interface</b>                  | <b>Part Number</b>                                      |
|---|---|
| <b>EV1000</b>                                     | Must be obtained by Edwards Lifesciences (manufacturer) |
| <b>Serial Cable: Null Modem/ Female to Female</b> | 0012-00-1275-01   |
| <b>VDI Module</b>                                 | 0998-00-1802-0801A                                      |

| <b>Balloon Pump Interface</b>    | <b>Part Number</b>                        |
|----------------------------------|---|
| <b>DS Balloon Pump</b>           | Must be obtained by Maquet (manufacturer) |
| <b>Analog Cable (Terminated)</b> | 0012-00-1801-01                           |

| <b>DIAP (Datascope Improved ASCII Protocol)</b> | <b>Part Number</b> |
|---|--------------------|
| <b>Serial Cable</b>                             | 0012-00-1275-01    |
| <b>VDI Module</b>                               | 0998-00-1802-0801A |

| <b>Gas Module 3 Interface</b>   | <b>Part Number</b>                 |
|---------------------------------|------------------------------------|
| <b>Gas Module 3</b>             | 0998-00-1900-01                    |
| <b>Serial Cable</b><br>1'<br>6' | 0012-00-1276-01<br>0012-00-1276-02 |
| <b>VDI Module</b>               | 0998-00-1802-0801A                 |

| <b>Somanetics Interface</b>                       | <b>Part Number</b>                             |
|---|--|
| <b>INVOS 5100C System Monitor</b>                 | Must be obtained by Covidien (manufacturer)    |
| <b>Serial Cable: Null Modem/ Female to Female</b> | 312170 (Covidien)<br>0012-00-1275-01 (Mindray) |
| <b>VDI Module</b>                                 | 0998-00-1802-0801A                             |

*Mindray DS USA, Inc. Document Revision*

|   |                                   |
|---|-----------------------------------|
| DESCRIPTION:<br><br><b>V Series Interface Part Number Guide</b> | DOCUMENT NO: <b>0002-08-50019</b> |
|   | REVISION: <b>D</b>                |
|   | SHEET 4 of 4                      |

| REV.<br>LETTER | CR NO. &<br>EFF. DATE | DESCRIPTION               |
|----------------|-----------------------|---------------------------|
| A              | CR 40018              | Original Release as Rev A |
| B              | CR 40217              | From Rev A to Rev B       |
| C              |                       | From Rev B to Rev C       |