

VS9 With Cerner CareAware VitalsLink



Connect. Capture. Send.

Mindray and Cerner® partner to deliver quality care more efficiently. Automate the collection and documentation of spot check vital signs data with the Mindray VS9 Vital Signs monitor and Cerner CareAware VitalsLink® to save time, improve clinical workflow, and reduce errors.

Meeting the needs of busy clinical environments, the VS9 is an easy-to-use vital signs monitor that is ready to go when you are. Assess patient's vitals with NIBP, SpO₂, and temperature measurements. Mindray's medical grade barcode scanner scans both the clinician ID and patient ID for quick validation.



Improve Patient Outcomes Through Increased Efficiency

Automated workflow helps free clinicians to focus more on patient care. At the patient's side, clinicians can validate patient information and automatically send verified vital signs directly into Cerner Millennium® from the VS9 (wireless or hardwired).



1. Scan clinician badge
2. Scan patient ID on the patient wrist band
3. Acquire vitals (NIBP, SpO₂, temperature)
4. Enter additional measurements (manual parameters and modifiers)
5. Select 'Save' and vitals are automatically sent to the EMR, from the point of care

Mindray North America
800 MacArthur Blvd., Mahwah, NJ 07430
Tel: 800.288.2121 Tel: 201.995.8000 Fax: 800.926.4275 www.mindray.com

Certain content © Cerner Corporation. All rights reserved. Cerner product images and information are used with permission from the Cerner Corporation and are for demonstration purposes only.

Mindray® is a registered trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
©2022 Mindray DS USA, Inc. All rights reserved. Subject to change. 05/22 P/N: 0002-08-7704 Rev A

mindray
NORTH AMERICA