


Hepatus-Series Transient Elastography Diagnostic Ultrasound Systems



A New Era for Non-Invasive Liver Disease

Mindray's latest technology in liver care is a complete solution for non-invasive liver disease management using transient elastography. The Hepatus-Series Transient Elastography Diagnostic Ultrasound Systems are professional non-invasive diagnostic systems for surveillance and staging, providing quantitative results to indicate the stage of liver fibrosis. They are accurate, reliable, and convenient for disease follow-up, opening a new era for non-invasive diagnosis of liver diseases.

Integrated Probe Design

The Hepatus ViTE (Visualized Transient Elastography) Transducer offers:

- Ultra-wideband frequency technology covers all patients of different body sizes and ages
- One transducer for ultrasound guidance and quantitative measurements – eliminating the need to switch between transducers during an exam
- Integrated button design allows users to acquire all measurements with one-click, reducing the scan time




Hepatus-Series Comparison

The **Hepatus 6 System** is a cart-based system designed for in-hospital use. It provides a comprehensive non-invasive liver assessment by combining ViTE and 2D ultrasound imaging to give clinicians a full picture of liver health. With an intuitive design, it creates an outstanding user experience blending touchscreen technology and a full keyboard.


The **Hepatus 5 System** is a compact, mobile solution to provide clinicians with a clear view of patient liver health. Its portable design allows it to be set on a table or flat surface, so clinicians can access leading-edge ViTE technology from clinician offices to rural and remote areas and anything in between.

	Hepatus 6	Hepatus 5
Form Factor		
System type	Cart	Tablet
Cart	Standard	For Purchase
Technology		
ViTE (real-time 2D image guidance)	✓	✓
2D imaging	✓	✓
Curved Transducer (C5-1s)	✓	—
Wireless connectivity	✓	✓
Trend analysis	✓	✓
System		
Full touchscreen tablet	✓	✓
Flexible panel with keyboard	✓	—
Built-in battery	✓	✓
DICOM® capabilities	✓	✓
Warranty		
System	5-year	3-year
ViTE Transducer	1-year	1-year
C5-1s Curved Transducer	5-year	—

Hepatus Reporting



The Hepatus System provides easy-to-read, comprehensive, one-page reports that clearly show liver fibrosis, steatosis and staging references. This valuable clinical information empowers clinicians to make informed decisions and determine the best next steps for their patients.



Mindray North America

800 MacArthur Blvd., Mahwah, NJ 07430

Tel: 800.288.2121 Support: 877.913.9663 Fax: 800.926.4275 www.mindray.com

Mindray® is a registered trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All brands and product names are trademarks of their respective owners. ©2024 Mindray DS USA, Inc. Subject to change. 04/24 P/N: 0002-08-40726 2.0

mindray

healthcare within reach