

# MX8 MSK Ultrasound System

## Transducer Family



### L16-4Hs (Linear)

Applications: Musculoskeletal  
FOV: 25 mm  
Biopsy guide: NA



### L12-3RCs (Linear)

Applications: Musculoskeletal  
FOV: 38 mm  
Biopsy guide: NGB-043, stainless steel, detachable; NGB-044, stainless steel, detachable, out-of-phase



### L20-5s (Linear)

Applications: Musculoskeletal  
FOV: 29 mm  
Biopsy guide: NA



### SC5-1Ns (Convex)

Applications: Musculoskeletal  
FOV: 73°  
Convex radius: 60 mm  
Biopsy guide: NGB-022, stainless steel, detachable



### L11-3VNs (Linear)

Applications: Musculoskeletal  
FOV: 38 mm  
Biopsy guide: NA  
eSpatialNavi™: Enabled



### L14-6Ns (Linear)

Applications: Musculoskeletal  
FOV: 38 mm  
Biopsy guide: NGB-043, stainless steel, detachable; NGB-044, stainless steel, detachable, out-of-phase



### L9-3s (Linear)

Applications: Musculoskeletal  
FOV: 43 mm  
Biopsy guide: NGB-034, multi-angle, reusable



### L13-3Ns (Linear)

Applications: Musculoskeletal  
FOV: 38 mm  
Biopsy guide: NGB-053, multi-angle, reusable

#### Mindray North America

800 MacArthur Blvd., Mahwah, NJ 07430

Tel: 800.288.2121 Support: 877.913.9663 Fax: 800.926.4275 [www.mindray.com](http://www.mindray.com)

Mindray® is a registered trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All brands and product names are trademarks of their respective owners. ©2022 Mindray DS USA, Inc. Subject to change. 06/22 P/N: 0002-08-40628 Rev A

# mindray

healthcare within reach