

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1. **BS-200**: BS-200, BS-220;

8. **BS-380**: BS-380, BS-390;

2. **BS-200E**: BS-200E, BS-220E;

9. **BS-400**: BS-400, BS-420;

3. **BS-230**: BS-230, BS-240;

10. **BS-430**: BS-430, BS-450, BS-460;

4. **BS-240E**: BS240E, BS240Pro;

11. **BS-480**: BS-480, BS-490;

5. **BS-300**: BS-300, BS-320;

12. **BS-600**: BS-600, BS-620;

6. **BS-330E**: BS-330E (Serial Number starts with "XQ-"),

BS-350E (Serial Number starts with "XS-");

13. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,

7. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

BS-1800, BS-1800plus;

BS-330E(V35.00) (Serial Number starts with "W8-" and software version starts with "35.00"),

BS-350E(V35.00) (Serial Number starts with "W9-" and software version starts with "35.00");

14. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

15. **BS-2800M**: BS-2600M, BS-2800M.

LOT: 044622002

EXP: 2023-08-14

| Abbreviated name | | Hb | HbA1c |
|------------------------|-------|---------------------------|---------------------------|
| Model | Level | Calibrator Value (µmol/L) | Calibrator Value (µmol/L) |
| BS-200 ¹ | S1 | 92.1 | 4.00 |
| | S2 | 144 | 15.5 |
| BS-200E ² | S1 | 92.1 | 4.00 |
| | S2 | 144 | 15.5 |
| BS-230 ³ | S1 | 92.1 | 4.00 |
| | S2 | 144 | 15.5 |
| BS-240E ⁴ | S1 | 92.1 | 4.00 |
| | S2 | 144 | 15.5 |
| BS-300 ⁵ | S1 | 92.1 | 4.00 |
| | S2 | 144 | 15.5 |
| BS-330E ⁶ | S1 | 92.1 | 4.00 |
| | S2 | 144 | 15.5 |
| BS-360E ⁷ | S1 | 92.1 | 4.00 |
| | S2 | 144 | 15.5 |
| BS-380 ⁸ | S1 | 92.1 | 4.00 |
| | S2 | 144 | 15.5 |
| BS-400 ⁹ | S1 | 92.1 | 4.00 |
| | S2 | 144 | 15.5 |
| BS-430 ¹⁰ | S1 | 92.1 | 4.00 |
| | S2 | 144 | 15.5 |
| BS-480 ¹¹ | S1 | 92.1 | 4.00 |
| | S2 | 144 | 15.5 |
| BS-600 ¹² | S1 | 92.1 | 4.00 |
| | S2 | 144 | 15.5 |
| BS-800 ¹³ | S1 | 92.1 | 4.00 |
| | S2 | 144 | 15.5 |
| BS-2000 ¹⁴ | S1 | 92.1 | 4.00 |
| | S2 | 144 | 15.5 |
| BS-2800M ¹⁵ | S1 | 92.1 | 4.00 |
| | S2 | 144 | 15.5 |

| English | Abbreviated name | Calibration Value | Model | Level |
|------------------|--------------------------|------------------------|---------|---------|
| Русский | сокращенное наименование | Эталонное значение | модель | Уровень |
| Português | Nome abreviado | Valor de calibração | Modelo | Nível |
| Español | nombre abreviado | Valor de calibración | modelo | Nivel |
| Italiano | abbreviazione | Valore di calibrazione | modelli | Livello |
| Türkçe | kısaltılmış ad | Kalibrasyon Değeri | model | Düzey |

| HbA1c/Hb | |
|------------------|------------------------------|
| English | Hemoglobin A1c/ Hemoglobin |
| Русский | гемоглобин A1c / гемоглобин |
| Português | Hemoglobina A1c /Hemoglobina |
| Español | hemoglobina A1c /hemoglobina |
| Italiano | Emoglobina A1c/Emoglobina |
| Türkçe | Hemoglobin A1c / Hemoglobin |