

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1. **BS-200**: BS-200, BS-220;

8. **BS-380**: BS-380, BS-390;

2. **BS-200E**: BS-200E, BS-220E;

9. **BS-400**: BS-400, BS-420;

3. **BS-230**: BS-230, BS-240;

10. **BS-430**: BS-430, BS-450, BS-460;

4. **BS-240E**: BS240E, BS240Pro;

11. **BS-480**: BS-480, BS-490;

5. **BS-300**: BS-300, BS-320;

12. **BS-600**: BS-600, BS-620;

6. **BS-330E**: BS-330E (Serial Number starts with "XQ-"),
BS-350E (Serial Number starts with "XS-");

13. **BS-600M**: BS-600M, BS-620M;

7. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,
BS-330E(V35.00) (Serial Number starts with "W8-"
and software version starts with "35.00"),
BS-350E(V35.00) (Serial Number starts with "W9-"
and software version starts with "35.00");

14. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,
BS-1800, BS-1800plus;

15. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

16. **BS-2800M**: BS-2600M, BS-2800M.

LOT: 044622004

EXP: 2023-10-10

Abbreviated name		Hb	HbA1c
Model	Level	Calibrator Value (μmol/L)	Calibrator Value (μmol/L)
BS-200 ¹	S1	92.0	3.70
	S2	151	14.9
BS-200E ²	S1	92.0	3.70
	S2	151	14.9
BS-230 ³	S1	92.0	3.70
	S2	151	14.9
BS-240E ⁴	S1	92.0	3.70
	S2	151	14.9
BS-300 ⁵	S1	92.0	3.70
	S2	151	14.9
BS-330E ⁶	S1	92.0	3.70
	S2	151	14.9
BS-360E ⁷	S1	92.0	3.70
	S2	151	14.9
BS-380 ⁸	S1	92.0	3.70
	S2	151	14.9
BS-400 ⁹	S1	92.0	3.70
	S2	151	14.9
BS-430 ¹⁰	S1	92.0	3.70
	S2	151	14.9
BS-480 ¹¹	S1	92.0	3.70
	S2	151	14.9
BS-600 ¹²	S1	92.0	3.70
	S2	151	14.9
BS-600M ¹³	S1	/	/
	S2	/	/
BS-800 ¹⁴	S1	92.0	3.70
	S2	151	14.9

HbA1C Calibrator

mindray

BS-2000 ¹⁵	S1	92.0	3.70	
	S2	151	14.9	
BS-2800M ¹⁶	S1	92.0	3.70	
	S2	151	14.9	
English	Abbreviated name	Calibration Value	Model	Level
Русский	сокращенное наименование	Эталонное значение	модель	Уровень
Português	Nome abreviado	Valor de calibração	Modelo	Nível
Español	nombre abreviado	Valor de calibración	modelo	Nivel
Italiano	abbreviazione	Valore di calibrazione	modelli	Livello
Türkçe	kısaltılmış ad	Kalibrasyon Değeri	model	Düzey

HbA1c/Hb	
English	Hemoglobin A1c/ Hemoglobin
Русский	гемоглобин A1c / гемоглобин
Português	Hemoglobina A1c /Hemoglobina
Español	hemoglobina A1c /hemoglobina
Italiano	Emoglobina A1c/Emoglobina
Türkçe	Hemoglobin A1c / Hemoglobin