

## Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1. **BS-200**: BS-200, BS-220;

8. **BS-380**: BS-380, BS-390;

2. **BS-200E**: BS-200E, BS-220E;

9. **BS-400**: BS-400, BS-420;

3. **BS-230**: BS-230, BS-240;

10. **BS-430**: BS-430, BS-450, BS-460;

4. **BS-240E**: BS240E, BS240Pro;

11. **BS-480**: BS-480, BS-490;

5. **BS-300**: BS-300, BS-320;

12. **BS-600**: BS-600, BS-620;

6. **BS-330E**: BS-330E (Serial Number starts with "XQ-"),  
BS-350E (Serial Number starts with "XS-");

13. **BS-600M**: BS-600M, BS-620M;

7. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

14. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,

**BS-330E(V35.00)** (Serial Number starts with "W8-"  
and software version starts with "35.00"),

BS-1800, BS-1800plus;

**BS-350E(V35.00)** (Serial Number starts with "W9-"  
and software version starts with "35.00");

15. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

16. **BS-2800M**: BS-2600M, BS-2800M.

**LOT**: 045622006

**EXP**: 2023-11-25

| Abbreviated name      |       | Hb                        | HbA1c                     |
|-----------------------|-------|---------------------------|---------------------------|
| Model                 | Level | Calibrator Value (µmol/L) | Calibrator Value (µmol/L) |
| BS-200 <sup>1</sup>   | S1    | 92.0                      | 3.70                      |
|                       | S2    | 151                       | 14.9                      |
| BS-200E <sup>2</sup>  | S1    | 92.0                      | 3.70                      |
|                       | S2    | 151                       | 14.9                      |
| BS-230 <sup>3</sup>   | S1    | 92.0                      | 3.70                      |
|                       | S2    | 151                       | 14.9                      |
| BS-240E <sup>4</sup>  | S1    | 92.0                      | 3.70                      |
|                       | S2    | 151                       | 14.9                      |
| BS-300 <sup>5</sup>   | S1    | 92.0                      | 3.70                      |
|                       | S2    | 151                       | 14.9                      |
| BS-330E <sup>6</sup>  | S1    | 92.0                      | 3.70                      |
|                       | S2    | 151                       | 14.9                      |
| BS-360E <sup>7</sup>  | S1    | 92.0                      | 3.70                      |
|                       | S2    | 151                       | 14.9                      |
| BS-380 <sup>8</sup>   | S1    | 92.0                      | 3.70                      |
|                       | S2    | 151                       | 14.9                      |
| BS-400 <sup>9</sup>   | S1    | 92.0                      | 3.70                      |
|                       | S2    | 151                       | 14.9                      |
| BS-430 <sup>10</sup>  | S1    | 92.0                      | 3.70                      |
|                       | S2    | 151                       | 14.9                      |
| BS-480 <sup>11</sup>  | S1    | 92.0                      | 3.70                      |
|                       | S2    | 151                       | 14.9                      |
| BS-600 <sup>12</sup>  | S1    | 92.0                      | 3.70                      |
|                       | S2    | 151                       | 14.9                      |
| BS-600M <sup>13</sup> | S1    | 92.0                      | 3.70                      |
|                       | S2    | 151                       | 14.9                      |
| BS-800 <sup>14</sup>  | S1    | 92.0                      | 3.70                      |
|                       | S2    | 151                       | 14.9                      |

| <b>BS-2000</b> <sup>15</sup>  | <b>S1</b>                | 92.0                     | 3.70         |              |
|-------------------------------|--------------------------|--------------------------|--------------|--------------|
|                               | <b>S2</b>                | 151                      | 14.9         |              |
| <b>BS-2800M</b> <sup>16</sup> | <b>S1</b>                | 92.0                     | 3.70         |              |
|                               | <b>S2</b>                | 151                      | 14.9         |              |
| <b>English</b>                | <b>Abbreviated name</b>  | <b>Calibration Value</b> | <b>Model</b> | <b>Level</b> |
| <b>Русский</b>                | сокращенное наименование | Эталонное значение       | модель       | Уровень      |
| <b>Português</b>              | Nome abreviado           | Valor de calibração      | Modelo       | Nível        |
| <b>Español</b>                | nombre abreviado         | Valor de calibración     | modelo       | Nivel        |
| <b>Italiano</b>               | abbreviazione            | Valore di calibrazione   | modelli      | Livello      |
| <b>Türkçe</b>                 | kısaltılmış ad           | Kalibrasyon Değeri       | model        | Düzey        |

| <b>HbA1c/Hb</b>  |                              |
|------------------|------------------------------|
| <b>English</b>   | Hemoglobin A1c/ Hemoglobin   |
| <b>Русский</b>   | гемоглобин A1c / гемоглобин  |
| <b>Português</b> | Hemoglobina A1c /Hemoglobina |
| <b>Español</b>   | hemoglobina A1c /hemoglobina |
| <b>Italiano</b>  | Emoglobina A1c/Emoglobina    |
| <b>Türkçe</b>    | Hemoglobin A1c / Hemoglobin  |