

Value sheet of Mindray BS Measurement System

Русский : Таблица результатов для системы BS компании Mindray

Português : Planilha de valores do Sistema de Medição BS da Mindray

Español : Hoja de valores del sistema de medición Mindray BS

Italiano : Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe : Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский : Данные совпадают во всех группах.

Português : A dados de cada grupo é a mesma.

Español : la datos de cada grupo es la misma.

Italiano : la dati di ogni gruppo è la stessa.

Türkçe : her grubun veri aynıdır.

1. **BS-200**: BS-200, BS-220;

8. **BS-380**: BS-380, BS-390;

2. **BS-200E**: BS-200E, BS-220E;

9. **BS-400**: BS-400, BS-420;

3. **BS-230**: BS-230, BS-240;

10. **BS-430**: BS-430, BS-450, BS-460;

4. **BS-240E**: BS240E, BS240Pro;

11. **BS-480**: BS-480, BS-490;

5. **BS-300**: BS-300, BS-320;

12. **BS-600**: BS-600, BS-620;

6. **BS-330E**: BS-330E(Serial Number starts with "XQ-"),
BS-350E(Serial Number starts with "XS-");

13. **BS-600M**: BS-600M, BS-620M;

7. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,
BS-330E(V35.00)(Serial Number starts with "W8-"
and software version starts with "35.00"),
BS-350E(V35.00) (Serial Number starts with "W9-"
and software version starts with "35.00");

14. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,
BS-1800, BS-1800plus;

15. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

16. **BS-2800M**: BS-2600M, BS-2800M.

LOT: 045623009

EXP: 2025-01-18

Abbreviated name		Hb	HbA1c
Model	Level	Calibrator Value (µmol/L)	Calibrator Value (µmol/L)
BS-200 ¹	S1	94.2	3.00
	S2	153	14.6
BS-200E ²	S1	94.2	3.20
	S2	153	14.7
BS-230 ³	S1	94.2	3.10
	S2	153	15.0
BS-240E ⁴	S1	94.2	3.10
	S2	153	15.0
BS-300 ⁵	S1	94.2	3.20
	S2	153	14.7
BS-330E ⁶	S1	94.2	3.20
	S2	153	14.7
BS-360E ⁷	S1	94.2	3.10
	S2	153	15.0
BS-380 ⁸	S1	94.2	3.20
	S2	153	14.7
BS-400 ⁹	S1	94.2	3.20
	S2	153	14.7
BS-430 ¹⁰	S1	94.2	3.10
	S2	153	15.0
BS-480 ¹¹	S1	94.2	3.10
	S2	153	15.0
BS-600 ¹²	S1	94.2	3.10
	S2	153	15.0
BS-600M ¹³	S1	93.7	3.10
	S2	154	14.8

Abbreviated name		Hb	HbA1c
Model	Level	Calibrator Value (μmol/L)	Calibrator Value (μmol/L)
BS-800 ¹⁴	S1	94.2	3.10
	S2	153	15.0
BS-2000 ¹⁵	S1	94.4	3.00
	S2	154	15.2
BS-2800M ¹⁶	S1	94.3	3.10
	S2	154	15.1

English	Abbreviated name	Calibration Value	Model	Level
Русский	сокращенное наименование	Эталонное значение	модель	Уровень
Português	Nome abreviado	Valor de calibração	Modelo	Nível
Español	nombre abreviado	Valor de calibración	modelo	Nivel
Italiano	abbreviazione	Valore di calibrazione	modelli	Livello
Türkçe	kısaltılmış ad	Kalibrasyon Değeri	model	Düzey

HbA1c/Hb	
English	Hemoglobin A1c/ Hemoglobin
Русский	гемоглобин A1c / гемоглобин
Português	Hemoglobina A1c /Hemoglobina
Español	hemoglobina A1c /hemoglobina
Italiano	Emoglobina A1c/Emoglobina
Türkçe	Hemoglobin A1c / Hemoglobin