

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов измерительной системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадает во всех группах.

Português : A dados de cada grupo é a mesma.

Español : la datos de cada grupo es la misma.

Italiano : la dati di ogni gruppo è la stessa.

Türkçe : her grubun veri aynıdır.

1.BS-240: BS-230, BS-240;

8.BS-600: BS-600, BS-620;

2.BS-240E: BS240E, BS240Pro;

9.BS-600M: BS-600M, BS-620M;

3.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");

10.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

4.BS-380: BS-380, BS-390;

11.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

5.BS-400: BS-400, BS-420;

12.BS-2800M:BS-2600M, BS-2800M;

6.BS-430: BS-430, BS-450,BS-460;

13.Applicable models of the chemistry shall be subject to the parameter list and instructions.

7.BS-480: BS-480, BS-490;

LOT : 045723017

: 2025-03-04

| English | Abbreviated name | Calibration Value | Level | Model |
|-----------|--------------------------|------------------------|---------|---------|
| Русский | сокращенное наименование | Эталонное значение | Уровень | модель |
| Português | Nome abreviado | Valor de calibração | Nível | Modelo |
| Español | nombre abreviado | Valor de calibración | Nivel | modelo |
| Italiano | abbreviazione | Valore di calibrazione | Livello | modelli |
| Türkçe | kısaltılmış ad | Kalibrasyon Değeri | Düzyey | model |

| English | Expanded Uncertainty |
|-----------|------------------------------|
| Русский | Расширенная неопределенность |
| Português | Incerteza Expandida |
| Español | Incertidumbre expandida |
| Italiano | Incertezza estesa |
| Türkçe | Daha Uzun Süreli Belirsizlik |

| | |
|-----------|----------|
| | D-Dimer |
| English | D-Dimer |
| Русский | D-димера |
| Português | D-Dimer |
| Español | Dímero D |
| Italiano | D-Dimero |
| Türkçe | D-Dimer |

| Abbreviated name | D-Dimer | Calibration Rule | spline | | | | | | |
|------------------|---------------------|----------------------|----------------------|------------------|----------------------|------------------|----------------------|------------------|----------------------|
| Model | BS-240 ¹ | BS-240E ² | BS-360E ³ | | | | | | |
| Level | Unit | Calibrator Value | Expanded Uncertainty | Calibrator Value | Expanded Uncertainty | Calibrator Value | Expanded Uncertainty | Calibrator Value | Expanded Uncertainty |

| | | | | | | | | | |
|---|-------|------|------|------|------|------|------|------|------|
| a | | 0.00 | / | 0.00 | / | 0.00 | / | 0.00 | / |
| b | | 1.93 | 0.07 | 1.93 | 0.07 | 1.93 | 0.07 | 1.93 | 0.07 |
| c | µg/mL | 8.40 | 0.17 | 8.40 | 0.17 | 8.40 | 0.17 | 8.40 | 0.17 |
| d | | 16.3 | 0.4 | 16.3 | 0.4 | 16.3 | 0.4 | 16.3 | 0.4 |
| e | | 32.4 | 1.1 | 32.4 | 1.1 | 32.4 | 1.1 | 32.4 | 1.1 |

D-Dimer Calibrator



| f | µg/mL | 46.3 | 1.6 | 46.3 | 1.6 | 46.3 | 1.6 | 46.3 | 1.6 |
|-------|-------|----------------------|----------------------|----------------------|----------------------|-----------------------|----------------------|------------------------|----------------------|
| Model | | BS-400 ⁵ | | BS-430 ⁶ | | BS-480 ⁷ | | BS-600 ⁸ | |
| Level | Unit | Calibrator Value | Expanded Uncertainty | Calibrator Value | Expanded Uncertainty | Calibrator Value | Expanded Uncertainty | Calibrator Value | Expanded Uncertainty |
| a | | 0.00 | / | 0.00 | / | 0.00 | / | 0.00 | / |
| b | | 1.93 | 0.07 | 1.93 | 0.07 | 1.93 | 0.07 | 1.93 | 0.07 |
| c | µg/mL | 8.40 | 0.17 | 8.40 | 0.17 | 8.40 | 0.17 | 8.40 | 0.17 |
| d | | 16.3 | 0.4 | 16.3 | 0.4 | 16.3 | 0.4 | 16.3 | 0.4 |
| e | | 32.4 | 1.1 | 32.4 | 1.1 | 32.4 | 1.1 | 32.4 | 1.1 |
| f | | 46.3 | 1.6 | 46.3 | 1.6 | 46.3 | 1.6 | 46.3 | 1.6 |
| Model | | BS-600M ⁹ | | BS-800 ¹⁰ | | BS-2000 ¹¹ | | BS-2800M ¹² | |
| Level | Unit | Calibrator Value | Expanded Uncertainty | Calibrator Value | Expanded Uncertainty | Calibrator Value | Expanded Uncertainty | Calibrator Value | Expanded Uncertainty |
| a | | 0.00 | / | 0.00 | / | 0.00 | / | 0.00 | / |
| b | | 1.98 | 0.08 | 1.93 | 0.07 | 2.02 | 0.08 | 1.98 | 0.08 |
| c | µg/mL | 8.35 | 0.17 | 8.40 | 0.17 | 8.40 | 0.17 | 8.35 | 0.17 |
| d | | 16.1 | 0.4 | 16.3 | 0.4 | 15.9 | 0.4 | 16.1 | 0.4 |
| e | | 31.0 | 1.0 | 32.4 | 1.1 | 31.6 | 1.0 | 31.0 | 1.0 |
| f | | 46.0 | 1.6 | 46.3 | 1.6 | 46.0 | 1.6 | 46.0 | 1.6 |

mindray

D-Dimer Calibrator

For use on: BS-2000, BS-2800M,
BS-600M, BS-800

LOT 045723017

2025-03-04

