

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов измерительной системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадает во всех группах.

Português : A dados de cada grupo é a mesma.

Español : la datos de cada grupo es la misma.

Italiano : la dati di ogni gruppo è la stessa.

Türkçe : her grubun veri aynıdır.

- | | |
|---|--|
| <p>1.BS-200E: BS-200E, BS-220E,;</p> <p>2.BS-240: BS-230, BS-240;</p> <p>3.BS-240E: BS240E, BS240Pro;</p> <p>4.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");</p> <p>5.BS-380: BS-380, BS-390;</p> <p>6.BS-400: BS-400, BS-420;</p> <p>7.BS-430: BS-430, BS-450,BS-460;</p> | <p>8.BS-480: BS-480, BS-490;</p> <p>9.BS-600: BS-600, BS-620;</p> <p>10.BS-600M: BS-600M, BS-620M;</p> <p>11.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;</p> <p>12.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;</p> <p>13.BS-2800M:BS-2600M, BS-2800M;</p> <p>14.Applicable models of the chemistry shall be subject to the parameter list and instructions.</p> |
|---|--|

LOT : 046123014

Exp: 2024-12-26

English	Abbreviated name	Calibration Value	Level	Model
Русский	сокращенное наименование	Эталонное значение	Уровень	модель
Português	Nome abreviado	Valor de calibração	Nível	Modelo
Español	nombre abreviado	Valor de calibración	Nivel	modelo
Italiano	abbreviazione	Valore di calibrazione	Livello	modelli
Türkçe	kısaltılmış ad	Kalibrasyon Değeri	Düzey	model

English	Expanded Uncertainty
Русский	Расширенная неопределенность
Português	Incerteza Expandida
Español	Incertidumbre expandida
Italiano	Incertezza estesa
Türkçe	Daha Uzun Süreli Belirsizlik

English	FER
English	Ferritin
Русский	ферритин
Português	Ferritina
Español	Ferritina
Italiano	Ferritina
Türkçe	Ferritin

Abbreviated name	FER	Calibration Rule		spline					
Model	BS-200E ¹	BS-240 ²		BS-240E ³		BS-360E ⁴			
Level	Unit	Calibrator Value	Expanded Uncertainty	Calibrator Value	Expanded Uncertainty	Calibrator Value	Expanded Uncertainty	Calibrator Value	Expanded Uncertainty
生理盐水		0	/	0	/	0	/	0	/
a		128	9	128	9	120	9	116	8
b	ng/mL	248	5	248	5	243	5	234	5
c		580	25	580	25	560	24	559	24
d		1068	37	1068	37	1080	38	1068	37

FER Calibrator



生理盐水		0	/	0	/	0	/	0	/
a		288	21	288	21	270	19	261	19
b	pmol/L	557	12	557	12	546	12	526	12
c		1303	56	1303	56	1258	54	1256	54
d		2400	84	2400	84	2427	85	2400	84
Model		BS-380 ⁵		BS-400 ⁶		BS-430 ⁷		BS-480 ⁸	
Level	Unit	Calibrator Value	Expanded Uncertainty	Calibrator Value	Expanded Uncertainty	Calibrator Value	Expanded Uncertainty	Calibrator Value	Expanded Uncertainty
生理盐水		0	/	0	/	0	/	0	/
a		120	9	120	9	118	8	116	8
b	ng/mL	240	5	240	5	245	5	234	5
c		570	25	570	25	560	24	559	24
d		1080	38	1080	38	1080	38	1068	37
生理盐水		0	/	0	/	0	/	0	/
a		270	19	270	19	265	19	261	19
b	pmol/L	539	12	539	12	551	12	526	12
c		1281	55	1281	55	1258	54	1256	54
d		2427	85	2427	85	2427	85	2400	84
Model		BS-600 ⁹		BS-600M ¹⁰		BS-800 ¹¹		BS-2000 ¹²	
Level	Unit	Calibrator Value	Expanded Uncertainty	Calibrator Value	Expanded Uncertainty	Calibrator Value	Expanded Uncertainty	Calibrator Value	Expanded Uncertainty
生理盐水		0	/	0	/	0	/	0	/
a		116	8	116	8	116	8	116	8
b	ng/mL	234	5	232	5	234	5	232	5
c		559	24	554	24	559	24	539	23
d		1068	37	1095	38	1068	37	1060	37
生理盐水		0	/	0	/	0	/	0	/
a		261	19	261	19	261	19	261	19
b	pmol/L	526	12	521	12	526	12	521	11
c		1256	54	1245	54	1256	54	1211	52
d		2400	84	2460	86	2400	84	2382	82
Model		BS-2800M ¹³							
Level	Unit	Calibrator Value	Expanded Uncertainty						
生理盐水		0	/						
a		116	8						
b	ng/mL	232	5						
c		554	24						
d		1095	38						
生理盐水		0	/						
a		261	19						
b	pmol/L	521	11						
c		1245	54						
d		2460	85						

mindray

FER Calibrator

For use on: BS-2000, BS-2800M,
BS-600M, BS-800

LOT 046123014

2024-12-26

