

Value sheet of Mindray BS Measurement System

Русский : Таблица результатов для системы BS компании Mindray

Português : Planilha de valores do Sistema de Medição BS da Mindray

Español : Hoja de valores del sistema de medición Mindray BS

Italiano : Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe : Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский : Данные совпадают во всех группах.

Português : A dados de cada grupo é a mesma.

Español : la datos de cada grupo es la misma.

Italiano : la dati di ogni gruppo è la stessa.

Türkçe : her grubun veri aynıdır.

1.**BS-120**: BS-120, BS-130;

2.**BS-180**: BS-180, BS-190;

3.**BS-200**: BS-200, BS-220;

4.**BS-200E**: BS-200E, BS-220E;

5.**BS-230**: BS-230, BS-240;

6.**BS-240E**: BS-240E, BS-240Pro;

7.**BS-300**: BS-300, BS-320;

8.**BS-330**: BS-330, BS-350;

9.**BS-330E**:BS-330E(Serial Number starts with "XQ-"),

BS-350E(Serial Number starts with "XS-");

10.**BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with "W9-" and software version starts with "35.00");

11.**BS-380**: BS-380, BS-390;

12.**BS-400**: BS-400, BS-420;

13.**BS-430**: BS-430, BS-450, BS-460;

14.**BS-480**: BS-480, BS-490;

15.**BS-600**: BS-600, BS-620;

16.**BS-600M**: BS-600M, BS-620M;

17.**BS-800**: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

18.**BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

19.**BS-2800M**: BS-2600M, BS-2800M.

20.**S0: 0.9% NaCl, Conc. Of S0=0;**

Русский : S0: 0, 9% NaCl, конц. S0=0;

Português: S0: 0,9% NaCl, Conc. de S0=0;

Español: S0: 0,9% NaCl, Conc. de S0=0;

Italiano: S0: 0,9% NaCl, conc. di S0=0;

Türkçe: S0: %0,9 NaCl, S0 Kons.=0.

LOT : 143223011

Exp. Date : 2024-08-21

Abbreviated name	Unit	Model	Calibration Value ²⁰	Model	Calibration Value ²⁰
TBA	µmol/L	BS-120 ¹	32.8	BS-380 ¹¹	32.8
		BS-180 ²	32.8	BS-400 ¹²	32.8
		BS-200 ³	32.8	BS-430 ¹³	33.5
		BS-200E ⁴	32.8	BS-480 ¹⁴	32.8
		BS-230 ⁵	32.8	BS-600 ¹⁵	32.7
		BS-240E ⁶	32.8	BS-600M ¹⁶	32.8
		BS-300 ⁷	32.8	BS-800 ¹⁷	32.8
		BS-330 ⁸	32.8	BS-2000 ¹⁸	33.2
		BS-330E ⁹	32.8	BS-2800M ¹⁹	33.7
		BS-360E ¹⁰	32.8		

English	Abbreviated name	Model	Unit	Calibration Value
Русский	сокращенное наименование	модель	Прибор	Принцип калибровки
Português	Nome abreviado	Modelo	Unidade	Regra de calibração
Español	nombre abreviado	modelo	Unidad	Regla de calibración
Italiano	abbreviazione	modelli	Unità	Regola di calibrazione
Türkçe	kısaltılmış ad	model	Ünite	Kalibrasyon Kuralı

TBA

English	Total Bile Acids
Русский	желчные кислоты
Português	Ácido Biliar Total
Español	ácidos biliares totales
Italiano	acidi biliari totali
Türkçe	Total Safra Asitleridir

CO2 and TBA Multi Control

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1. **BS-120**: BS-120, BS-130;

2. **BS-180**: BS-180, BS-190;

3. **BS-200**: BS-200, BS-220;

4. **BS-200E**: BS-200E, BS-220E;

5. **BS-230**: BS-230, BS-240, BS-280;

6. **BS-240E**: BS-240E, BS-240Pro;

7. **BS-300**: BS-300, BS-320;

8. **BS-330**: BS-330, BS-350;

9. **BS-330E**: BS-330E (Serial Number starts with "XQ-"),

BS-350E (Serial Number starts with "XS-");

10. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,
BS-330E(V35.00) (Serial Number starts with "W8-"
and software version starts with "35.00"),
BS-350E(V35.00) (Serial Number starts with "W9-"
and software version starts with "35.00");

11. **BS-380**: BS-380, BS-390;

12. **BS-400**: BS-400, BS-420;

13. **BS-430**: BS-410, BS-430, BS-450, BS-460, BS-470;

14. **BS-480**: BS-480, BS-490;

15. **BS-600**: BS-600, BS-620;

16. **BS-600M**: BS-600M, BS-620M;

17. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,

BS-1800, BS-1800plus;

18. **BS-830**: BS-830, BS-830S, BS-840, BS-850, BS-860,

BS-870;

19. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M.

20. **BS-2800M**: BS-2600M, BS-2800M.

Abbreviated name	TBA Control(L)						TBA Control(H)							
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)			1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)			1 SD
TBA (apply to lot 061821003、061821004)	BS-120 ¹	$\mu\text{mol/L}$	/	/	-	/	/	BS-120 ¹	$\mu\text{mol/L}$	/	/	-	/	/
	BS-180 ²	$\mu\text{mol/L}$	/	/	-	/	/	BS-180 ²	$\mu\text{mol/L}$	/	/	-	/	/
	BS-200 ³	$\mu\text{mol/L}$	/	/	-	/	/	BS-200 ³	$\mu\text{mol/L}$	/	/	-	/	/
	BS-200E ⁴	$\mu\text{mol/L}$	/	/	-	/	/	BS-200E ⁴	$\mu\text{mol/L}$	/	/	-	/	/
	BS-230 ⁵	$\mu\text{mol/L}$	/	/	-	/	/	BS-230 ⁵	$\mu\text{mol/L}$	/	/	-	/	/
	BS-240E ⁶	$\mu\text{mol/L}$	/	/	-	/	/	BS-240E ⁶	$\mu\text{mol/L}$	/	/	-	/	/
	BS-300 ⁷	$\mu\text{mol/L}$	/	/	-	/	/	BS-300 ⁷	$\mu\text{mol/L}$	/	/	-	/	/
	BS-330 ⁸	$\mu\text{mol/L}$	/	/	-	/	/	BS-330 ⁸	$\mu\text{mol/L}$	/	/	-	/	/
	BS-330E ⁹	$\mu\text{mol/L}$	/	/	-	/	/	BS-330E ⁹	$\mu\text{mol/L}$	/	/	-	/	/
	BS-360E ¹⁰	$\mu\text{mol/L}$	/	/	-	/	/	BS-360E ¹⁰	$\mu\text{mol/L}$	/	/	-	/	/
	BS-380 ¹¹	$\mu\text{mol/L}$	/	/	-	/	/	BS-380 ¹¹	$\mu\text{mol/L}$	/	/	-	/	/
	BS-400 ¹²	$\mu\text{mol/L}$	/	/	-	/	/	BS-400 ¹²	$\mu\text{mol/L}$	/	/	-	/	/
	BS-430 ¹³	$\mu\text{mol/L}$	/	/	-	/	/	BS-430 ¹³	$\mu\text{mol/L}$	/	/	-	/	/
	BS-480 ¹⁴	$\mu\text{mol/L}$	/	/	-	/	/	BS-480 ¹⁴	$\mu\text{mol/L}$	/	/	-	/	/
	BS-600 ¹⁵	$\mu\text{mol/L}$	/	/	-	/	/	BS-600 ¹⁵	$\mu\text{mol/L}$	/	/	-	/	/
	BS-600M ¹⁶	$\mu\text{mol/L}$	/	/	-	/	/	BS-600M ¹⁶	$\mu\text{mol/L}$	/	/	-	/	/
	BS-800 ¹⁷	$\mu\text{mol/L}$	/	/	-	/	/	BS-800 ¹⁷	$\mu\text{mol/L}$	/	/	-	/	/
	BS-830 ¹⁸	$\mu\text{mol/L}$	/	/	-	/	/	BS-830 ¹⁸	$\mu\text{mol/L}$	/	/	-	/	/
	BS-2000 ¹⁹	$\mu\text{mol/L}$	/	/	-	/	/	BS-2000 ¹⁹	$\mu\text{mol/L}$	/	/	-	/	/
	BS-2800M ²⁰	$\mu\text{mol/L}$	/	/	-	/	/	BS-2800M ²⁰	$\mu\text{mol/L}$	/	/	-	/	/

CO2 and TBA Multi Control

Abbreviated name	TBA Control(L)					TBA Control(H)				
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)	1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)	1 SD
TBA (apply to lot 061822001)	BS-120 ¹	$\mu\text{mol/L}$	/	/ - /	/	BS-120 ¹	$\mu\text{mol/L}$	/	/ - /	/
	BS-180 ²	$\mu\text{mol/L}$	/	/ - /	/	BS-180 ²	$\mu\text{mol/L}$	/	/ - /	/
	BS-200 ³	$\mu\text{mol/L}$	/	/ - /	/	BS-200 ³	$\mu\text{mol/L}$	/	/ - /	/
	BS-200E ⁴	$\mu\text{mol/L}$	/	/ - /	/	BS-200E ⁴	$\mu\text{mol/L}$	/	/ - /	/
	BS-230 ⁵	$\mu\text{mol/L}$	/	/ - /	/	BS-230 ⁵	$\mu\text{mol/L}$	/	/ - /	/
	BS-240E ⁶	$\mu\text{mol/L}$	/	/ - /	/	BS-240E ⁶	$\mu\text{mol/L}$	/	/ - /	/
	BS-300 ⁷	$\mu\text{mol/L}$	/	/ - /	/	BS-300 ⁷	$\mu\text{mol/L}$	/	/ - /	/
	BS-330 ⁸	$\mu\text{mol/L}$	/	/ - /	/	BS-330 ⁸	$\mu\text{mol/L}$	/	/ - /	/
	BS-330E ⁹	$\mu\text{mol/L}$	/	/ - /	/	BS-330E ⁹	$\mu\text{mol/L}$	/	/ - /	/
	BS-360E ¹⁰	$\mu\text{mol/L}$	/	/ - /	/	BS-360E ¹⁰	$\mu\text{mol/L}$	/	/ - /	/
	BS-380 ¹¹	$\mu\text{mol/L}$	/	/ - /	/	BS-380 ¹¹	$\mu\text{mol/L}$	/	/ - /	/
	BS-400 ¹²	$\mu\text{mol/L}$	/	/ - /	/	BS-400 ¹²	$\mu\text{mol/L}$	/	/ - /	/
	BS-430 ¹³	$\mu\text{mol/L}$	/	/ - /	/	BS-430 ¹³	$\mu\text{mol/L}$	/	/ - /	/
	BS-480 ¹⁴	$\mu\text{mol/L}$	/	/ - /	/	BS-480 ¹⁴	$\mu\text{mol/L}$	/	/ - /	/
	BS-600 ¹⁵	$\mu\text{mol/L}$	/	/ - /	/	BS-600 ¹⁵	$\mu\text{mol/L}$	/	/ - /	/
	BS-600M ¹⁶	$\mu\text{mol/L}$	/	/ - /	/	BS-600M ¹⁶	$\mu\text{mol/L}$	/	/ - /	/
	BS-800 ¹⁷	$\mu\text{mol/L}$	/	/ - /	/	BS-800 ¹⁷	$\mu\text{mol/L}$	/	/ - /	/
	BS-830 ¹⁸	$\mu\text{mol/L}$	/	/ - /	/	BS-830 ¹⁸	$\mu\text{mol/L}$	/	/ - /	/
	BS-2000 ¹⁹	$\mu\text{mol/L}$	/	/ - /	/	BS-2000 ¹⁹	$\mu\text{mol/L}$	/	/ - /	/
	BS-2800M ²⁰	$\mu\text{mol/L}$	/	/ - /	/	BS-2800M ²⁰	$\mu\text{mol/L}$	/	/ - /	/

Abbreviated name	TBA Control(L)					TBA Control(H)				
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)	1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)	1 SD
TBA (apply to lot 061822002)	BS-120 ¹	$\mu\text{mol/L}$	/	/ - /	/	BS-120 ¹	$\mu\text{mol/L}$	/	/ - /	/
	BS-180 ²	$\mu\text{mol/L}$	/	/ - /	/	BS-180 ²	$\mu\text{mol/L}$	/	/ - /	/
	BS-200 ³	$\mu\text{mol/L}$	/	/ - /	/	BS-200 ³	$\mu\text{mol/L}$	/	/ - /	/
	BS-200E ⁴	$\mu\text{mol/L}$	/	/ - /	/	BS-200E ⁴	$\mu\text{mol/L}$	/	/ - /	/
	BS-230 ⁵	$\mu\text{mol/L}$	/	/ - /	/	BS-230 ⁵	$\mu\text{mol/L}$	/	/ - /	/
	BS-240E ⁶	$\mu\text{mol/L}$	/	/ - /	/	BS-240E ⁶	$\mu\text{mol/L}$	/	/ - /	/
	BS-300 ⁷	$\mu\text{mol/L}$	/	/ - /	/	BS-300 ⁷	$\mu\text{mol/L}$	/	/ - /	/
	BS-330 ⁸	$\mu\text{mol/L}$	/	/ - /	/	BS-330 ⁸	$\mu\text{mol/L}$	/	/ - /	/
	BS-330E ⁹	$\mu\text{mol/L}$	/	/ - /	/	BS-330E ⁹	$\mu\text{mol/L}$	/	/ - /	/
	BS-360E ¹⁰	$\mu\text{mol/L}$	/	/ - /	/	BS-360E ¹⁰	$\mu\text{mol/L}$	/	/ - /	/
	BS-380 ¹¹	$\mu\text{mol/L}$	/	/ - /	/	BS-380 ¹¹	$\mu\text{mol/L}$	/	/ - /	/
	BS-400 ¹²	$\mu\text{mol/L}$	/	/ - /	/	BS-400 ¹²	$\mu\text{mol/L}$	/	/ - /	/
	BS-430 ¹³	$\mu\text{mol/L}$	/	/ - /	/	BS-430 ¹³	$\mu\text{mol/L}$	/	/ - /	/
	BS-480 ¹⁴	$\mu\text{mol/L}$	/	/ - /	/	BS-480 ¹⁴	$\mu\text{mol/L}$	/	/ - /	/
	BS-600 ¹⁵	$\mu\text{mol/L}$	/	/ - /	/	BS-600 ¹⁵	$\mu\text{mol/L}$	/	/ - /	/
	BS-600M ¹⁶	$\mu\text{mol/L}$	/	/ - /	/	BS-600M ¹⁶	$\mu\text{mol/L}$	/	/ - /	/
	BS-800 ¹⁷	$\mu\text{mol/L}$	/	/ - /	/	BS-800 ¹⁷	$\mu\text{mol/L}$	/	/ - /	/
	BS-830 ¹⁸	$\mu\text{mol/L}$	/	/ - /	/	BS-830 ¹⁸	$\mu\text{mol/L}$	/	/ - /	/
	BS-2000 ¹⁹	$\mu\text{mol/L}$	/	/ - /	/	BS-2000 ¹⁹	$\mu\text{mol/L}$	/	/ - /	/
	BS-2800M ²⁰	$\mu\text{mol/L}$	/	/ - /	/	BS-2800M ²⁰	$\mu\text{mol/L}$	/	/ - /	/

CO2 and TBA Multi Control

Abbreviated name	TBA Control(L)						TBA Control(H)							
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD		
TBA (apply to lot 061822003、061822004、061822005)	BS-120 ¹	μ mol/L	21.5	17.2	—	25.8	1.4	BS-120 ¹	μ mol/L	36.2	28.9	—	43.5	2.4
	BS-180 ²	μ mol/L	21.5	17.2	—	25.8	1.4	BS-180 ²	μ mol/L	36.2	28.9	—	43.5	2.4
	BS-200 ³	μ mol/L	20.8	16.6	—	25.0	1.4	BS-200 ³	μ mol/L	36.3	29.0	—	43.6	2.4
	BS-200E ⁴	μ mol/L	21.0	16.8	—	25.2	1.4	BS-200E ⁴	μ mol/L	37.7	30.1	—	45.3	2.5
	BS-230 ⁵	μ mol/L	21.0	16.8	—	25.2	1.4	BS-230 ⁵	μ mol/L	35.5	28.4	—	42.6	2.4
	BS-240E ⁶	μ mol/L	20.3	16.2	—	24.4	1.4	BS-240E ⁶	μ mol/L	34.3	27.4	—	41.2	2.3
	BS-300 ⁷	μ mol/L	21.3	17.0	—	25.6	1.4	BS-300 ⁷	μ mol/L	36.9	29.5	—	44.3	2.5
	BS-330 ⁸	μ mol/L	20.8	16.6	—	25.0	1.4	BS-330 ⁸	μ mol/L	36.3	29.0	—	43.6	2.4
	BS-330E ⁹	μ mol/L	21.0	16.8	—	25.2	1.4	BS-330E ⁹	μ mol/L	37.7	30.1	—	45.3	2.5
	BS-360E ¹⁰	μ mol/L	20.9	16.7	—	25.1	1.4	BS-360E ¹⁰	μ mol/L	35.0	28.0	—	42.0	2.3
	BS-380 ¹¹	μ mol/L	21.3	17.0	—	25.6	1.4	BS-380 ¹¹	μ mol/L	35.7	28.5	—	42.9	2.4
	BS-400 ¹²	μ mol/L	21.6	17.3	—	25.9	1.4	BS-400 ¹²	μ mol/L	35.5	28.4	—	42.6	2.4
	BS-430 ¹³	μ mol/L	21.6	17.3	—	25.9	1.4	BS-430 ¹³	μ mol/L	35.9	28.7	—	43.1	2.4
	BS-480 ¹⁴	μ mol/L	21.6	17.3	—	25.9	1.4	BS-480 ¹⁴	μ mol/L	35.4	28.3	—	42.5	2.4
	BS-600 ¹⁵	μ mol/L	21.9	17.5	—	26.3	1.5	BS-600 ¹⁵	μ mol/L	36.4	29.1	—	43.7	2.4
	BS-600M ¹⁶	μ mol/L	21.6	17.3	—	25.9	1.4	BS-600M ¹⁶	μ mol/L	36.0	28.8	—	43.2	2.4
	BS-800 ¹⁷	μ mol/L	21.4	17.1	—	25.7	1.4	BS-800 ¹⁷	μ mol/L	35.9	28.7	—	43.1	2.4
	BS-830 ¹⁸	μ mol/L	22.1	17.7	—	26.5	1.5	BS-830 ¹⁸	μ mol/L	36.8	29.4	—	44.2	2.5
	BS-2000 ¹⁹	μ mol/L	22.3	17.8	—	26.8	1.5	BS-2000 ¹⁹	μ mol/L	37.3	29.8	—	44.8	2.5
	BS-2800M ²⁰	μ mol/L	22.3	17.8	—	26.8	1.5	BS-2800M ²⁰	μ mol/L	37.3	29.8	—	44.8	2.5

Abbreviated name	TBA Control(L)						TBA Control(H)							
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD		
TBA (apply to lot 061822006)	BS-120 ¹	μ mol/L	21.5	17.2	—	25.8	1.4	BS-120 ¹	μ mol/L	36.1	28.8	—	43.4	2.4
	BS-180 ²	μ mol/L	21.5	17.2	—	25.8	1.4	BS-180 ²	μ mol/L	36.1	28.8	—	43.4	2.4
	BS-200 ³	μ mol/L	20.8	16.6	—	25.0	1.4	BS-200 ³	μ mol/L	35.7	28.5	—	42.9	2.4
	BS-200E ⁴	μ mol/L	21.0	16.8	—	25.2	1.4	BS-200E ⁴	μ mol/L	38.2	30.5	—	45.9	2.6
	BS-230 ⁵	μ mol/L	21.0	16.8	—	25.2	1.4	BS-230 ⁵	μ mol/L	35.6	28.4	—	42.8	2.4
	BS-240E ⁶	μ mol/L	20.3	16.2	—	24.4	1.4	BS-240E ⁶	μ mol/L	34.3	27.4	—	41.2	2.3
	BS-300 ⁷	μ mol/L	21.3	17.0	—	25.6	1.4	BS-300 ⁷	μ mol/L	37.4	29.9	—	44.9	2.5
	BS-330 ⁸	μ mol/L	20.8	16.6	—	25.0	1.4	BS-330 ⁸	μ mol/L	35.7	28.5	—	42.9	2.4
	BS-330E ⁹	μ mol/L	21.0	16.8	—	25.2	1.4	BS-330E ⁹	μ mol/L	38.2	30.5	—	45.9	2.6
	BS-360E ¹⁰	μ mol/L	20.9	16.7	—	25.1	1.4	BS-360E ¹⁰	μ mol/L	35.3	28.2	—	42.4	2.4
	BS-380 ¹¹	μ mol/L	21.3	17.0	—	25.6	1.4	BS-380 ¹¹	μ mol/L	35.6	28.4	—	42.8	2.4
	BS-400 ¹²	μ mol/L	21.6	17.3	—	25.9	1.4	BS-400 ¹²	μ mol/L	35.6	28.4	—	42.8	2.4
	BS-430 ¹³	μ mol/L	21.6	17.3	—	25.9	1.4	BS-430 ¹³	μ mol/L	35.9	28.7	—	43.1	2.4
	BS-480 ¹⁴	μ mol/L	21.6	17.3	—	25.9	1.4	BS-480 ¹⁴	μ mol/L	35.7	28.5	—	42.9	2.4
	BS-600 ¹⁵	μ mol/L	21.9	17.5	—	26.3	1.5	BS-600 ¹⁵	μ mol/L	36.1	28.8	—	43.4	2.4
	BS-600M ¹⁶	μ mol/L	21.6	17.3	—	25.9	1.4	BS-600M ¹⁶	μ mol/L	35.8	28.6	—	43.0	2.4
	BS-800 ¹⁷	μ mol/L	21.4	17.1	—	25.7	1.4	BS-800 ¹⁷	μ mol/L	35.9	28.7	—	43.1	2.4
	BS-830 ¹⁸	μ mol/L	22.1	17.7	—	26.5	1.5	BS-830 ¹⁸	μ mol/L	36.6	29.2	—	44.0	2.5
	BS-2000 ¹⁹	μ mol/L	22.3	17.8	—	26.8	1.5	BS-2000 ¹⁹	μ mol/L	37.3	29.8	—	44.8	2.5
	BS-2800M ²⁰	μ mol/L	22.3	17.8	—	26.8	1.5	BS-2800M ²⁰	μ mol/L	37.3	29.8	—	44.8	2.5

CO2 and TBA Multi Control

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон (результат анализа \pm 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa (Valores da análise \pm 3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango (Valor de ensayo \pm 3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore diconcentrazione \pm 3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)
TBA					
English	Total Bile Acids			Español	ácidos biliares totais
Русский	общие желчные кислоты			Italiano	acidi biliari totali
Português	ácidos biliares totales			Türkçe	total safra asitleridir