

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português : Planilha de valores do Sistema de Medição BS da Mindray

Español : Hoja de valores del sistema de medición Mindray BS

Italiano : Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe : Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Español : la datos de cada grupo es la misma.

Türkçe : her grubun veri aynıdır.

Português : A dados de cada grupo é a mesma.

Italiano : la dati di ogni gruppo è la stessa.

1.**BS-120**: BS-120, BS-130;

2.**BS-200**: BS-200, BS-220;

3.**BS-200E**: BS-200E, BS-220E;

4.**BS-230**: BS-230, BS-240;

5.**BS-240E**: BS-240E, BS-240Pro;

6.**BS-300**: BS-300, BS-320;

7.**BS-330E**:BS-330E(Serial Number starts with "XQ-"),

BS-350E(Serial Number starts with "XS-");

18.**S0: 0.9% NaCl, Conc. Of S0=0;**

Português: S0: 0,9% NaCl, Conc. de S0=0;

Italiano: S0: 0,9% NaCl, conc. di S0=0;

8.**BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,
BS-330E(V35.00)(Serial Number starts with "W8-"
and software version starts with "35.00"),
BS-350E(V35.00) (Serial Number starts with"W9-"
and software version starts with "35.00");

9.**BS-380**: BS-380, BS-390;

10.**BS-400**: BS-400, BS-420;

11.**BS-430**: BS-430, BS-450, BS-460;

12.**BS-480**: BS-480, BS-490;

13.**BS-600**: BS-600, BS-620;

14.**BS-600M**: BS-600M, BS-620M;

15.**BS-800**: BS-800, BS-820, BS-800M, BS-820M,
BS-1800, BS-1800plus;

16.**BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

17.**BS-2800M**:BS-2600M, BS-2800M.

Русский: S0: 0, 9% NaCl, конц. S0=0;

Español: S0: 0,9% NaCl, Conc. de S0=0;

Türkçe: S0: %0,9 NaCl, S0 Kons.=0.

LOT : 145422008

: 2024-04-21

Abbreviated name	Unit	Model	Calibration Value ¹⁸	Model	Calibration Value ¹⁸
FUN	µmol/L	BS-120 ¹	435	BS-400 ¹⁰	435
		BS-200 ²	452	BS-430 ¹¹	435
		BS-200E ³	452	BS-480 ¹²	435
		BS-230 ⁴	435	BS-600 ¹³	435
		BS-240E ⁵	435	BS-600M ¹⁴	441
		BS-300 ⁶	466	BS-800 ¹⁵	435
		BS-330E ⁷	452	BS-2000 ¹⁶	441
		BS-360E ⁸	431	BS-2800M ¹⁷	441
		BS-380 ⁹	422		

English	Abbreviated name	Model	Unit	Calibration Value
Русский	сокращенное наименование	модель	Прибор	Принцип калибровки
Português	Nome abreviado	Modelo	Unidade	Regra de calibração
Español	nombre abreviado	modelo	Unidad	Regla de calibración
Italiano	abbreviazione	modelli	Unità	Regola di calibrazione
Türkçe	kisaltılmış ad	model	Ünite	Kalibrasyon Kuralı

FUN	
English	Fructosamine
Русский	Фруктозамин
Português	Frutosamina
Español	Fructosamina
Italiano	Fruttosamina
Türkçe	Fruktozamin