

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-120: BS-120, BS-130, BS-180, BS-190;

11.BS-430: BS-430, BS-450, BS-460;

2.BS-200: BS-200, BS-220, BS-330, BS-350;

12.BS-480: BS-480, BS-490;

3.BS-200E: BS-200E, BS-220E,;

13.BS-600: BS-600, BS-620;

4.BS-240: BS-230, BS-240;

14.BS-600M: BS-600M, BS-620M;

5.BS-240E: BS240E, BS240Pro;

15.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

6.BS-300: BS-300, BS-320;

7.BS-330E:BS-330E(Serial Number starts with "XQ-"), BS-350E(Serial Number starts with "XS-")

16.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

8.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with "W9-" and software version starts with "35.00");

17.BS-2800M:BS-2600M, BS-2800M.

9.BS-380: BS-380, BS-390;

18.0.9% NaCl is needed, Conc. Of 0.9% NaCl =0;

10.BS-400: BS-400, BS-420;

19.Applicable models of the chemistry shall be subject to the parameter list and instructions.

LOT

: 145424001



: 2025-03-02

English	Abbreviated name	Unit	Calibration Rule	Model
Русский	сокращенное наименование	Прибор	Принцип калибровки	модель
Português	Nome abreviado	Unidade	Regra de calibração	Modelo
Español	nombre abreviado	Unidad	Regla de calibración	modelo
Italiano	abbreviazione	Unità	Regola di calibrazione	modelli
Türkçe	kısaltılmış ad	Ünite	Kalibrasyon Kuralı	model

English	Expanded Uncertainty
Русский	Расширенная неопределенность
Português	Incerteza Expandida
Español	Incertidumbre expandida
Italiano	Incertezza estesa
Türkçe	Daha Uzun Süreli Belirsizlik

FUN	
English	Fructosamine
Русский	Фруктозамин
Português	Frutosamina
Español	Fructosamina
Italiano	Fruttosamina
Türkçe	Fruktozamin

Abbreviated name	Unit	Calibration Rule			Two-point linear		
		Model	Calibrator Value	Expanded Uncertainty	Model	Calibrator Value	Expanded Uncertainty
FUN	μ mol/L	BS-120¹	449	7	BS-400¹⁰	460	7
		BS-200²	467	7	BS-430¹¹	449	7
		BS-200E³	467	7	BS-480¹²	449	7
		BS-240⁴	431	7	BS-600¹³	449	7
		BS-240E⁵	449	7	BS-600M¹⁴	452	7
		BS-300⁶	480	7	BS-800¹⁵	449	7
		BS-330E⁷	467	7	BS-2000¹⁶	454	7
		BS-360E⁸	449	7	BS-2800M¹⁷	452	7
		BS-380⁹	431	7			

mindray

FUN Calibrator

For use on: BS-2000, BS-2800M,
BS-600M, BS-800

LOT 145424001

🕒 2025-03-02

