

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов измерительной системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Português: A dados de cada grupo é a mesma.

italiano: la dati di ogni gruppo e la stessa.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1. **BS-120**: BS-120, BS-130;

2. **BS-200**: BS-200, BS-220;

3. **BS-200E**: BS-200E, BS-220E;

4. **BS-230**: BS-230, BS-240;

5. **BS-240E**: BS-240E, BS-240Pro;

6. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

BS-330E (V35.00) (Serial Number starts with "W8-" and software version starts with "35.00"), **BS-350E (V35.00)** (Serial Number starts with "W9-" and software version starts with "35.00");

7. **BS-380**: BS-380, BS-390;

8. **BS-400**: BS-400, BS-420;

9. **BS-430**: BS-430, BS-450, BS-460;

10. **BS-480**: BS-480, BS-490;

11. **BS-600**: BS-600, BS-620;

12. **BS-800**: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

13. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

14. **BS-2800M**: BS-2600M, BS-2800M.

15. **Cal0: 0.9% NaCl, Conc. Of Cal0=0;**

Português: Cal0: 0,9% NaCl, Conc. de Cal0=0;

Italiano: Cal0: 0,9% NaCl, conc. di Cal0=0;

Русский: Cal0: 0,9% NaCl, конц. Cal0=0;

Español: Cal0: 0,9% NaCl, Conc. de Cal0=0;

Türkçe: Cal0: %0,9 NaCl, Cal0 Kons.=0.

LOT: 146522001

EXP: 2023-03-18

Abbreviated name	β-HB	Calibration Value ¹⁵ (mmol/L)
BS-120 ¹	Cal1	0.890
BS-200 ²	Cal1	0.900
BS-200E ³	Cal1	0.910
BS-230 ⁴	Cal1	0.930
BS-240E ⁵	Cal1	0.860
BS-360E ⁶	Cal1	0.900
BS-380 ⁷	Cal1	0.930
BS-400 ⁸	Cal1	0.900
BS-430 ⁹	Cal1	0.930
BS-480 ¹⁰	Cal1	0.930
BS-600 ¹¹	Cal1	0.930
BS-800 ¹²	Cal1	0.930
BS-2000 ¹³	Cal1	0.940
BS-2800M ¹⁴	Cal1	0.940

English	Abbreviated name	Calibration Value	Model	Level
Русский	сокращенное наименование	Эталонное значение	модель	Уровень
Português	Nome abreviado	Valor de calibração	Modelo	Nível
Español	nombre abreviado	Valor de calibración	modelo	Nivel
Italiano	abbreviazione	Valore di calibrazione	modelli	Livello
Türkçe	kısaltılmış ad	Kalibrasyon Değeri	model	Düzey

β-HB	
English	β-Hydroxybutyrate
Русский	β-гидроксибутирата
Português	β-Hidroxibutirato
Español	β-hidroxibutirato
Italiano	β-idrossibutirato
Türkçe	β-Hidroksibötrat