

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов измерительной системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадает во всех группах. Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma. Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1. **BS-230**: BS-230, BS-240;

2. **BS-240E**: BS-240E, BS-240Pro;

3. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"),**BS-350E(V35.00)** (Serial Number starts with "W9-" and software version starts with "35.00");

4. **BS-380**: BS-380, BS-390;

13. **Cal0: 0.9% NaCl, Conc. Of Cal0=0.**

Português: Cal0: 0,9% NaCl, Conc. de Cal0=0.

Italiano: Cal0: 0,9% NaCl, conc. di Cal0=0.

5. **BS-400**: BS-400, BS-420;

6. **BS-430**: BS-430, BS-450, BS-460;

7. **BS-480**: BS-480, BS-490;

8. **BS-600**: BS-600, BS-620;

9. **BS-600M**: BS-600M, BS-620M;

10. **BS-800**: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

11. **BS-2000**: BS-2000,BS-2200,BS-2000M,BS-2200M;

12. **BS-2800M**: BS-2600M, BS-2800M.

Русский: Cal0: 0,9% NaCl, конц. Cal0=0.

Español: Cal0: 0,9% NaCl, Conc. de Cal0=0.

Türkçe: Cal0: %0,9 NaCl, Cal0 Kons.=0.

LOT : 146722010

EXP : 2023-10-08

| Abbreviated name | | UIBC | |
|-------------------------------|-------------|---------------------------------|-------|
| Model | Level | Calibration Value ¹³ | |
| | | µmol/L | µg/dL |
| BS-230 ¹ | Cal1 | 112 | 626 |
| BS-240E ² | Cal1 | 112 | 626 |
| BS-360E ³ | Cal1 | 112 | 626 |
| BS-380 ⁴ | Cal1 | 110 | 615 |
| BS-400 ⁵ | Cal1 | 110 | 615 |
| BS-430 ⁶ | Cal1 | 112 | 626 |
| BS-480 ⁷ | Cal1 | 111 | 620 |
| BS-600 ⁸ | Cal1 | 111 | 620 |
| BS-600M ⁹ | Cal1 | 111 | 620 |
| BS-800 ¹⁰ | Cal1 | 111 | 620 |
| BS-2000 ¹¹ | Cal1 | 112 | 626 |
| BS-2800M ¹² | Cal1 | 111 | 620 |

| English | Abbreviated name | Calibration Value | Model | Level |
|-----------|--------------------------|------------------------|---------|---------|
| Русский | сокращенное наименование | Эталонное значение | модель | Уровень |
| Português | Nome abreviado | Valor de calibração | Modelo | Nível |
| Español | nombre abreviado | Valor de calibración | modelo | Nivel |
| Italiano | abbreviazione | Valore di calibrazione | modelli | Livello |
| Türkçe | kısaltılmış ad | Kalibrasyon Değeri | model | Düzey |

| UIBC | |
|-----------|--|
| English | Unsaturated Iron Binding Capacity |
| Русский | ненасыщенная железосвязывающая способность |
| Português | Capacidade de ligação de ferro insaturado |
| Español | Capacidad de unión de hierro no saturado |
| Italiano | Capacità di legame del ferro insaturo |
| Türkçe | Unsature Demir Bağlama Kapasitesi |