

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов измерительной системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадает во всех группах. Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma. Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1. **BS-230**: BS-230, BS-240;

2. **BS-240E**: BS-240E, BS-240Pro;

3. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"),**BS-350E(V35.00)** (Serial Number starts with "W9-" and software version starts with "35.00");

4. **BS-380**: BS-380, BS-390;

13. **Cal0: 0.9% NaCl, Conc. Of Cal0=0.**

Português: Cal0: 0,9% NaCl, Conc. de Cal0=0.

Italiano: Cal0: 0,9% NaCl, conc. di Cal0=0.

5. **BS-400**: BS-400, BS-420;

6. **BS-430**: BS-430, BS-450, BS-460;

7. **BS-480**: BS-480, BS-490;

8. **BS-600**: BS-600, BS-620;

9. **BS-600M**: BS-600M, BS-620M;

10. **BS-800**: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

11. **BS-2000**: BS-2000,BS-2200,BS-2000M,BS-2200M;

12. **BS-2800M**: BS-2600M, BS-2800M.

Русский: Cal0: 0,9% NaCl, конц. Cal0=0.

Español: Cal0: 0,9% NaCl, Conc. de Cal0=0.

Türkçe: Cal0: %0,9 NaCl, Cal0 Kons.=0.

LOT : 146722013

Exp : 2024-01-17

Abbreviated name		UIBC	
Model	Level	Calibration Value ¹³	
		µmol/L	µg/dL
BS-230 ¹	Cal1	111	620
BS-240E ²	Cal1	108	604
BS-360E ³	Cal1	109	609
BS-380 ⁴	Cal1	108	604
BS-400 ⁵	Cal1	108	604
BS-430 ⁶	Cal1	111	620
BS-480 ⁷	Cal1	112	626
BS-600 ⁸	Cal1	112	626
BS-600M ⁹	Cal1	113	632
BS-800 ¹⁰	Cal1	112	626
BS-2000 ¹¹	Cal1	111	620
BS-2800M ¹²	Cal1	113	632

English	Abbreviated name	Calibration Value	Model	Level
Русский	сокращенное наименование	Эталонное значение	модель	Уровень
Português	Nome abreviado	Valor de calibração	Modelo	Nível
Español	nombre abreviado	Valor de calibración	modelo	Nivel
Italiano	abbreviazione	Valore di calibrazione	modelli	Livello
Türkçe	kisaltılmış ad	Kalibrasyon Değeri	model	Düzey

UIBC	
English	Unsaturated Iron Binding Capacity
Русский	ненасыщенная железосвязывающая способность
Português	Capacidade de ligação de ferro insaturado
Español	Capacidad de unión de hierro no saturado
Italiano	Capacità di legame del ferro insaturo
Türkçe	Unsature Demir Bağlama Kapasitesi