

## Value sheet of Mindray BS Measurement System

**Русский:** Таблица результатов для системы BS компании Mindray

**Português:** Planilha de valores do Sistema de Medição BS da Mindray

**Español:** Hoja de valores del sistema de medición Mindray BS

**Italiano:** Scheda dei valori del sistema di misurazione BS di Mindray

**Türkçe:** Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

**Русский:** Данные совпадают во всех группах.

**Português:** A dados de cada grupo é a mesma.

**Español:** la datos de cada grupo es la misma.

**Italiano:** la dati di ogni gruppo è la stessa.

**Türkçe:** her grubun veri aynıdır.

1.BS-240: BS-230, BS-240;

8.BS-600: BS-600, BS-620;

2.BS-240E: BS240E, BS240Pro;

9.BS-600M: BS-600M, BS-620M;

3.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00)

10.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

(Serial Number starts with "W9-" and software version starts with "35.00");

11.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

4.BS-380: BS-380, BS-390;

12.BS-2800M:BS-2600M, BS-2800M.

5.BS-400: BS-400, BS-420;

13.0.9% NaCl is needed, Conc. Of 0.9% NaCl =0;

6.BS-430: BS-430, BS-450, BS-460;

14.For applicable analyzers of the analyte, please subject to the parameter sheet and instrument.

7.BS-480: BS-480, BS-490;



**LOT** : 146724003



: 2025-03-07

English	Abbreviated name	Unit	Calibration Rule	Model
<b>Русский</b>	сокращенное наименование	Прибор	Принцип калибровки	модель
<b>Português</b>	Nome abreviado	Unidade	Regra de calibração	Modelo
<b>Español</b>	nombre abreviado	Unidad	Regla de calibración	modelo
<b>Italiano</b>	abbreviazione	Unità	Regola di calibrazione	modelli
<b>Türkçe</b>	kısaltılmış ad	Ünite	Kalibrasyon Kuralı	model

English	Expanded Uncertainty
<b>Русский</b>	Расширенная неопределенность
<b>Português</b>	Incerteza Expandida
<b>Español</b>	Incertidumbre expandida
<b>Italiano</b>	Incertezza estesa
<b>Türkçe</b>	Daha Uzun Süreli Belirsizlik

UIBC	
<b>English</b>	Unsaturated Iron Binding Capacity
<b>Русский</b>	ненасыщенная железосвязывающая способность
<b>Português</b>	Capacidade de ligação de ferro insaturado
<b>Español</b>	Capacidad de unión de hierro no saturado
<b>Italiano</b>	Capacità di legame del ferro insaturo
<b>Türkçe</b>	Unsature Demir Bağlama Kapasitesi

Abbreviated name	Unit	Calibration Rule				Two-point linear	
		Model	Calibrator Value	Expanded Uncertainty	Model	Calibrator Value	Expanded Uncertainty
<b>UIBC</b>	μmol/L	<b>BS-240<sup>1</sup></b>	112	2.9	<b>BS-480<sup>7</sup></b>	112	2.9
		<b>BS-240E<sup>2</sup></b>	109	2.8	<b>BS-600<sup>8</sup></b>	112	2.9
		<b>BS-360E<sup>3</sup></b>	109	2.8	<b>BS-600M<sup>9</sup></b>	113	2.9
		<b>BS-380<sup>4</sup></b>	110	2.8	<b>BS-800<sup>10</sup></b>	112	3
		<b>BS-400<sup>5</sup></b>	110	2.8	<b>BS-2000<sup>11</sup></b>	113	2.9
		<b>BS-430<sup>6</sup></b>	112	2.9	<b>BS-2800M<sup>12</sup></b>	113	2.9
	μg/dL	<b>BS-240<sup>1</sup></b>	626	16	<b>BS-480<sup>7</sup></b>	626	16
		<b>BS-240E<sup>2</sup></b>	609	16	<b>BS-600<sup>8</sup></b>	626	16
		<b>BS-360E<sup>3</sup></b>	609	16	<b>BS-600M<sup>9</sup></b>	632	16
		<b>BS-380<sup>4</sup></b>	615	16	<b>BS-800<sup>10</sup></b>	626	17
		<b>BS-400<sup>5</sup></b>	615	16	<b>BS-2000<sup>11</sup></b>	632	16
		<b>BS-430<sup>6</sup></b>	626	16	<b>BS-2800M<sup>12</sup></b>	632	16

# **mindray**

## **UIBC Calibrator**

For use on: BS-2000, BS-2800M,  
BS-600M, BS-800

**LOT** 146724003

 2025-03-07

