

Value sheet of Mindray BS Measurement System

Русский : Таблица результатов для системы BS компании Mindray

Português : Planilha de valores do Sistema de Medição BS da Mindray

Español : Hoja de valores del sistema de medición Mindray BS

Italiano : Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe : Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский : Данные совпадают во всех группах.

Español : la datos de cada grupo es la misma.

Türkçe : her grubun veri aynıdır.

Português : A dados de cada grupo é a mesma.

Italiano : la dati di ogni gruppo è la stessa.

1.**BS-120**: BS-120, BS-130;

2.**BS-200**: BS-200, BS-220;

3.**BS-200E**: BS-200E, BS-220E;

4.**BS-230**: BS-230, BS-240;

5.**BS-240E**: BS-240E, BS-240Pro;

6.**BS-330**: BS-330, BS-350;

7.**BS-330E**: **BS-330E**(Serial Number starts with "XQ-"),

BS-350E(Serial Number starts with "XS-");

8.**BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,
BS-330E(V35.00)(Serial Number starts with "W8-"
and software version starts with "35.00"),
BS-350E(V35.00) (Serial Number starts with "W9-"
and software version starts with "35.00");

9.**BS-380**: BS-380, BS-390;

10.**BS-400**: BS-400, BS-420;

11.**BS-430**: BS-430, BS-450, BS-460;

12.**BS-480**: BS-480, BS-490;

13.**BS-600**: BS-600, BS-620;

14.**BS-600M**: BS-600M, BS-620M;

15.**BS-800**: BS-800, BS-820, BS-800M, BS-820M,
BS-1800, BS-1800plus;

16.**BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

17.**BS-2800M**: BS-2600M, BS-2800M.

LOT : 147923007

EXP : 2025-05-01

Abbreviated name	Model	Unit	Calibration Value	Model	Unit	Calibration Value
TPUC	BS-120 ¹	mg/L	1134	BS-400 ¹⁰	mg/L	1192
	BS-200 ²	mg/L	1148	BS-430 ¹¹	mg/L	1169
	BS-200E ³	mg/L	1169	BS-480 ¹²	mg/L	1169
	BS-230 ⁴	mg/L	1111	BS-600 ¹³	mg/L	1169
	BS-240E ⁵	mg/L	1169	BS-600M ¹⁴	mg/L	1154
	BS-330 ⁶	mg/L	1148	BS-800 ¹⁵	mg/L	1169
	BS-330E ⁷	mg/L	1169	BS-2000 ¹⁶	mg/L	1159
	BS-360E ⁸	mg/L	1169	BS-2800M ¹⁷	mg/L	1154
	BS-380 ⁹	mg/L	1192			

English	Abbreviated name	Calibration Value	Unit
Русский	сокращенное наименование	Эталонное значение	Прибор
Português	Nome abreviado	Valor de calibração	Unidade
Español	nombre abreviado	Valor de calibración	Unidad
Italiano	abbreviazione	Valore di calibrazione	Unità
Türkçe	kısaltılmış ad	Kalibrasyon Değeri	Ünite

TPUC	
English	Total Protein in Urine/CSF (TPUC)
Русский	Общий белок в моче/СМЖ (TPUC)
Português	Proteína total em urina e LCR (PTUL)
Español	proteína total en orina/LCR (TPUC)
Italiano	Proteine Totali Urinarie/CSF (TPUC)
Türkçe	BOS/İdrar Protein (TPUC)