

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-120: BS-120, BS-130, BS-180, BS-190;

10.BS-430: BS-430, BS-450, BS-460;

2.BS-200: BS-200, BS-220, BS-330, BS-350;

11.BS-480: BS-480, BS-490;

3.BS-200E: BS-200E, BS-220E,;

12.BS-600: BS-600, BS-620;

4.BS-240: BS-230, BS-240;

13.BS-600M: BS-600M, BS-620M;

5.BS-240E: BS240E, BS240Pro;

6.BS-330E:BS-330E(Serial Number starts with "XQ-"),

14.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-

BS-350E(Serial Number starts with "XS-")

1800plus;

7.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-

15.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

330E(V35.00)(Serial Number starts with "W8-" and

software version starts with "35.00"), BS-350E(V35.00)

16.BS-2800M:BS-2600M, BS-2800M.

(Serial Number starts with "W9-" and software version

17.0.9% NaCl is needed, Conc. Of 0.9% NaCl =0;

starts with "35.00");

8.BS-380: BS-380, BS-390;

18.Applicable models of the chemistry shall be subject to the

9.BS-400: BS-400, BS-420;

parameter list and instructions.

LOT

: 147924001



: 2025-06-19

English	Abbreviated name	Unit	Calibration Rule	Model
Русский	сокращенное наименование	Прибор	Принцип калибровки	модель
Português	Nome abreviado	Unidade	Regra de calibração	Modelo
Español	nombre abreviado	Unidad	Regla de calibración	modelo
Italiano	abbreviazione	Unità	Regola di calibrazione	modelli
Türkçe	kısaltılmış ad	Ünite	Kalibrasyon Kuralı	model
English	Expanded Uncertainty			
Русский	Расширенная неопределенность			
Português	Incerteza Expandida			
Español	Incertidumbre expandida			
Italiano	Incertezza estesa			
Türkçe	Daha Uzun Süreli Belirsizlik			
TPUC				
English	Total Protein in Urine/CSF (TPUC)			
Русский	Общий белок в моче/СМЖ (TPUC)			
Português	Proteína total em urina e LCR (PTUL)			
Español	proteína total en orina/LCR (TPUC)			
Italiano	Proteine Totali Urinarie/CSF (TPUC)			
Türkçe	BOS/İdrar Protein (TPUC)			

Abbreviated name	Unit	Calibration Rule			Two-point linear		
		Model	Calibrator Value	Expanded Uncertainty	Model	Calibrator Value	Expanded Uncertainty
TPUC	mg/L	BS-120 ¹	1080	15	BS-400 ⁹	1186	16
		BS-200 ²	1127	15	BS-430 ¹⁰	1174	16
		BS-200E ³	1174	16	BS-480 ¹¹	1174	16
		BS-240 ⁴	1092	15	BS-600 ¹²	1174	16
		BS-240E ⁵	1174	16	BS-600M ¹³	1188	16
		BS-330E ⁶	1174	16	BS-800 ¹⁴	1174	16
		BS-360E ⁷	1174	16	BS-2000 ¹⁵	1179	16
		BS-380 ⁸	1186	16	BS-2800M ¹⁶	1188	16

mindray

TPUC Calibrator

For use on: BS-2000, BS-2800M,
BS-600M, BS-800

LOT 147924001

2025-06-19

