

**Value sheet of Mindray BS Measurement System**

Русский: Таблица результатов измерительной системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадает во всех группах.

Português : A dados de cada grupo é a mesma.

Español : la datos de cada grupo es la misma.

Italiano : la dati di ogni gruppo è la stessa.

Türkçe : her grubun veri aynıdır.

1.BS-120: BS-120, BS-130, BS-180, BS-190;	9.BS-430: BS-430, BS-450,BS-460;
2.BS-200: BS-200, BS-220, BS-330, BS-350;	10.BS-480: BS-480, BS-490;
3.BS-240: BS-230, BS-240;	11.BS-600: BS-600, BS-620;
4.BS-240E: BS240E, BS240Pro;	12.BS-600M: BS-600M, BS-620M;
5.BS-330E:BS-330E(Serial Number starts with "XQ-"), BS-350E(Serial Number starts with "XS-")	13.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;
6.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");	14.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;
7.BS-380: BS-380, BS-390;	15.BS-2800M:BS-2600M, BS-2800M;
8.BS-400: BS-400, BS-420;	16.Applicable models of the chemistry shall be subject to the parameter list and instructions.

LOT : 148023015

: 2024-11-17

English	Abbreviated name	Calibration Value	Level	Model
Русский	сокращенное наименование	Эталонное значение	Уровень	модель
Português	Nome abreviado	Valor de calibração	Nível	Modelo
Español	nombre abreviado	Valor de calibración	Nivel	modelo
Italiano	abbreviazione	Valore di calibrazione	Livello	modelli
Türkçe	kisaltılmış ad	Kalibrasyon Değeri	Düzey	model

English	Expanded Uncertainty
Русский	Расширенная неопределенность
Português	Incerteza Expandida
Español	Incertidumbre expandida
Italiano	Incertezza estesa
Türkçe	Daha Uzun Süreli Belirsizlik

HCY	
English	Homocysteine (HCY)
Русский	Гомоцистеин (HCY)
Português	Homocisteína (HCY)
Español	Homocisteína (HCY)
Italiano	Omocisteina (HCY)
Türkçe	Homosistain (HCY)

Abbreviated name	HCY II		Calibration Rule		spline				
Model	BS-120 <sup>1</sup>		BS-200 <sup>2</sup>		BS-240 <sup>3</sup>		BS-240E <sup>4</sup>		
Level	Unit	Calibrator Value	Expanded Uncertainty	Calibrator Value	Expanded Uncertainty	Calibrator Value	Expanded Uncertainty	Calibrator Value	Expanded Uncertainty
生理盐水		0.00	/	0.00	/	0.00	/	0.00	/
a	µmol/L	2.49	0.16	2.49	0.16	2.49	0.16	2.49	0.16
b		6.50	0.5	6.50	0.5	6.55	0.5	6.55	0.5

# HCY Calibrator



c		14.0	1.0	13.8	0.9	14.2	1.0	14.2	1.0
d	μmol/L	25.3	1.8	26.0	1.8	26.7	1.8	26.7	1.8
e		42.0	2.7	41.8	2.7	43.0	2.8	43.0	2.8

Model		BS-330E <sup>5</sup>		BS-360E <sup>6</sup>		BS-380 <sup>7</sup>		BS-400 <sup>8</sup>	
Level	Unit	Calibrator Value	Expanded Uncertainty	Calibrator Value	Expanded Uncertainty	Calibrator Value	Expanded Uncertainty	Calibrator Value	Expanded Uncertainty
生理盐水		0.00	/	0.00	/	0.00	/	0.00	/
a		2.42	0.16	2.49	0.16	2.49	0.16	2.49	0.16
b	μmol/L	6.40	0.5	6.55	0.5	6.55	0.5	6.55	0.5
c		14.2	1.0	14.2	1.0	14.2	1.0	14.2	1.0
d		25.6	1.8	26.7	1.8	26.7	1.8	26.7	1.8
e		41.5	2.7	43.0	2.8	43.0	2.8	43.0	2.8

Model		BS-430 <sup>9</sup>		BS-480 <sup>10</sup>		BS-600 <sup>11</sup>		BS-600M <sup>12</sup>	
Level	Unit	Calibrator Value	Expanded Uncertainty	Calibrator Value	Expanded Uncertainty	Calibrator Value	Expanded Uncertainty	Calibrator Value	Expanded Uncertainty
生理盐水		0.00	/	0.00	/	0.00	/	0.00	/
a		2.49	0.16	2.49	0.16	2.49	0.16	2.55	0.17
b	μmol/L	6.55	0.5	6.55	0.5	6.55	0.5	6.67	0.5
c		14.2	1.0	14.2	1.0	14.2	1.0	14.3	1.0
d		26.7	1.8	26.7	1.8	26.7	1.8	27.6	1.9
e		43.0	2.8	43.0	2.8	43.0	2.8	44.0	2.9

Model		BS-800 <sup>13</sup>		BS-2000 <sup>14</sup>		BS-2800M <sup>15</sup>	
Level	Unit	Calibrator Value	Expanded Uncertainty	Calibrator Value	Expanded Uncertainty	Calibrator Value	Expanded Uncertainty
生理盐水		0.00	/	0.00	/	0.00	/
a		2.49	0.16	2.57	0.16	2.55	0.16
b	μmol/L	6.55	0.5	6.68	0.5	6.67	0.5
c		14.2	1.0	14.3	1.0	14.3	1.0
d		26.7	1.8	27.9	2.0	27.6	1.9
e		43.0	2.8	44.4	2.9	44.0	2.8

# mindray

## HCY Calibrator

For use on: BS-2000, BS-2800M,  
BS-600M, BS-800

**LOT** 148023015

**2024-11-17**

