

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası

The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1. **BS-120**: BS-120, BS-130;

2. **BS-180**: BS-180, BS-190;

3. **BS-200**: BS-200, BS-220;

4. **BS-200E**: BS-200E, BS-220E;

5. **BS-230**: BS-230, BS-240;

6. **BS-240E**: BS240E, BS240Pro;

7. **BS-300**: BS-300, BS-320;

8. **BS-330**: BS-330, BS-350;

9. **BS-330E**: **BS-330E**(Serial Number starts with "XQ-"),

BS-350E(Serial Number starts with "XS-");

10. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"),

BS-350E(V35.00) (Serial Number starts with "W9-" and software version starts with "35.00");

11. **BS-380**: BS-380, BS-390;

12. **BS-400**: BS-400, BS-420;

13. **BS-430**: BS-430, BS-450, BS-460;

14. **BS-480**: BS-480, BS-490;

15. **BS-600**: BS-600, BS-620;

16. **BS-800**: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

17. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

18. **BS-2800M**: BS-2600M, BS-2800M.

LOT : 151422002

EXP : 2023-08-15

| Abbreviated name | Model | Unit | Calibration Value | Model | Unit | Calibration Value |
|-----------------------------|-----------------------------|--------|-------------------------------|------------------------------|--------|-------------------|
| CK-MB | BS-120 ¹ | U/L | 128 | BS-360E ¹⁰ | U/L | 128 |
| | | µkat/L | 2.14 | | µkat/L | 2.14 |
| | BS-180 ² | U/L | 128 | BS-380 ¹¹ | U/L | 128 |
| | | µkat/L | 2.14 | | µkat/L | 2.14 |
| | BS-200 ³ | U/L | 128 | BS-400 ¹² | U/L | 128 |
| | | µkat/L | 2.14 | | µkat/L | 2.14 |
| | BS-200E ⁴ | U/L | 128 | BS-430 ¹³ | U/L | 128 |
| | | µkat/L | 2.14 | | µkat/L | 2.14 |
| | BS-230 ⁵ | U/L | 128 | BS-480 ¹⁴ | U/L | 128 |
| | | µkat/L | 2.14 | | µkat/L | 2.14 |
| | BS-240E ⁶ | U/L | 128 | BS-600 ¹⁵ | U/L | 128 |
| | | µkat/L | 2.14 | | µkat/L | 2.14 |
| | BS-300 ⁷ | U/L | 128 | BS-800 ¹⁶ | U/L | 128 |
| | | µkat/L | 2.14 | | µkat/L | 2.14 |
| | BS-330 ⁸ | U/L | 128 | BS-2000 ¹⁷ | U/L | 126 |
| | | µkat/L | 2.14 | | µkat/L | 2.10 |
| BS-330E ⁹ | U/L | 128 | BS-2800M ¹⁸ | U/L | 131 | |
| | µkat/L | 2.14 | | µkat/L | 2.19 | |

| English | Abbreviated name | Unit | Calibration Value | CK-MB |
|-----------|--------------------------|---------|------------------------|--------------------------|
| Русский | сокращенное наименование | Прибор | Эталонное значение | MB фракцию креатинкиназы |
| Português | Nome abreviado | Unidade | Valor de calibração | creatina quinase-MB |
| Español | nombre abreviado | Unidad | Valor de calibración | creatina quinasa-MB |
| Italiano | abbreviazione | Unità | Valore di calibrazione | creatina chinasi-MB |
| Türkçe | kısaltılmış ad | Ünite | Kalibrasyon Değeri | Kreatin Kinaz-MB |