

## Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası

The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1. **BS-120**: BS-120, BS-130;

2. **BS-180**: BS-180, BS-190;

3. **BS-200**: BS-200, BS-220;

4. **BS-200E**: BS-200E, BS-220E;

5. **BS-230**: BS-230, BS-240;

6. **BS-240E**: BS240E, BS240Pro;

7. **BS-300**: BS-300, BS-320;

8. **BS-330**: BS-330, BS-350;

9. **BS-330E**: BS-330E (Serial Number starts with "XQ-"),

**BS-350E** (Serial Number starts with "XS-");

10. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S, **BS-330E(V35.00)** (Serial Number starts with "W8-" and software version starts with "35.00"), **BS-350E(V35.00)** (Serial Number starts with "W9-" and software version starts with "35.00");

11. **BS-380**: BS-380, BS-390;

12. **BS-400**: BS-400, BS-420;

13. **BS-430**: BS-430, BS-450, BS-460;

14. **BS-480**: BS-480, BS-490;

15. **BS-600**: BS-600, BS-620;

16. **BS-600M**: BS-600M, BS-620M;

17. **BS-800**: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

18. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

19. **BS-2800M**: BS-2600M, BS-2800M.

LOT: 151422006

: 2023-12-08

Abbreviated name	Model	Unit	Calibration Value	Model	Unit	Calibration Value
CK-MB	<b>BS-120</b> <sup>1</sup>	U/L	141	<b>BS-380</b> <sup>11</sup>	U/L	141
		µkat/L	2.35		µkat/L	2.35
	<b>BS-180</b> <sup>2</sup>	U/L	141	<b>BS-400</b> <sup>12</sup>	U/L	141
		µkat/L	2.35		µkat/L	2.35
	<b>BS-200</b> <sup>3</sup>	U/L	141	<b>BS-430</b> <sup>13</sup>	U/L	141
		µkat/L	2.35		µkat/L	2.35
	<b>BS-200E</b> <sup>4</sup>	U/L	141	<b>BS-480</b> <sup>14</sup>	U/L	141
		µkat/L	2.35		µkat/L	2.35
	<b>BS-230</b> <sup>5</sup>	U/L	141	<b>BS-600</b> <sup>15</sup>	U/L	141
		µkat/L	2.35		µkat/L	2.35
	<b>BS-240E</b> <sup>6</sup>	U/L	141	<b>BS-600M</b> <sup>16</sup>	U/L	144
		µkat/L	2.35		µkat/L	2.40
	<b>BS-300</b> <sup>7</sup>	U/L	141	<b>BS-800</b> <sup>17</sup>	U/L	141
		µkat/L	2.35		µkat/L	2.35
	<b>BS-330</b> <sup>8</sup>	U/L	141	<b>BS-2000</b> <sup>18</sup>	U/L	139
		µkat/L	2.35		µkat/L	2.32
	<b>BS-330E</b> <sup>9</sup>	U/L	141	<b>BS-2800M</b> <sup>19</sup>	U/L	145
		µkat/L	2.35		µkat/L	2.42
<b>BS-360E</b> <sup>10</sup>	U/L	141				
	µkat/L	2.35				

English	Abbreviated name	Unit	Calibration Value	CK-MB
Русский	сокращенное наименование	Прибор	Эталонное значение	МВ фракцию креатинкиназы
Português	Nome abreviado	Unidade	Valor de calibração	creatina quinase-MB
Español	nombre abreviado	Unidad	Valor de calibración	creatina quinasa-MB
Italiano	abbreviazione	Unità	Valore di calibrazione	creatina chinasi-MB
Türkçe	kısaltılmış ad	Ünite	Kalibrasyon Değeri	Kreatin Kinaz-MB