

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası

The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1. **BS-120**: BS-120, BS-130;

2. **BS-180**: BS-180, BS-190;

3. **BS-200**: BS-200, BS-220;

4. **BS-200E**: BS-200E, BS-220E;

5. **BS-230**: BS-230, BS-240;

6. **BS-240E**: BS240E, BS240Pro;

7. **BS-300**: BS-300, BS-320;

8. **BS-330**: BS-330, BS-350;

9. **BS-330E**: BS-330E (Serial Number starts with "XQ-"),

BS-350E (Serial Number starts with "XS-");

10. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S, **BS-330E(V35.00)** (Serial Number starts with "W8-" and software version starts with "35.00"), **BS-350E(V35.00)** (Serial Number starts with "W9-" and software version starts with "35.00");

11. **BS-380**: BS-380, BS-390;

12. **BS-400**: BS-400, BS-420;

13. **BS-430**: BS-430, BS-450, BS-460;

14. **BS-480**: BS-480, BS-490;

15. **BS-600**: BS-600, BS-620;

16. **BS-600M**: BS-600M, BS-620M;

17. **BS-800**: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

18. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

19. **BS-2800M**: BS-2600M, BS-2800M.

LOT: 151423005

⏰: 2024-10-30

Abbreviated name	Model	Unit	Calibration Value	Model	Unit	Calibration Value
CK-MB	BS-120 ¹	U/L	123	BS-380 ¹¹	U/L	125
		μkat/L	2.05		μkat/L	2.09
	BS-180 ²	U/L	123	BS-400 ¹²	U/L	125
		μkat/L	2.05		μkat/L	2.09
	BS-200 ³	U/L	123	BS-430 ¹³	U/L	125
		μkat/L	2.05		μkat/L	2.09
	BS-200E ⁴	U/L	125	BS-480 ¹⁴	U/L	125
		μkat/L	2.09		μkat/L	2.09
	BS-230 ⁵	U/L	123	BS-600 ¹⁵	U/L	125
		μkat/L	2.05		μkat/L	2.09
BS-240E ⁶	U/L	125	BS-600M ¹⁶	U/L	128	
	μkat/L	2.09		μkat/L	2.14	
BS-300 ⁷	U/L	125	BS-800 ¹⁷	U/L	125	
	μkat/L	2.09		μkat/L	2.09	
BS-330 ⁸	U/L	123	BS-2000 ¹⁸	U/L	126	
	μkat/L	2.05		μkat/L	2.10	
BS-330E ⁹	U/L	125	BS-2800M ¹⁹	U/L	129	
	μkat/L	2.09		μkat/L	2.15	
BS-360E ¹⁰	U/L	125				
	μkat/L	2.09				

English	Abbreviated name	Unit	Calibration Value	CK-MB
Русский	сокращенное наименование	Прибор	Эталонное значение	МВ фракцию креатинкиназы
Português	Nome abreviado	Unidade	Valor de calibração	creatina quinase-MB
Español	nombre abreviado	Unidad	Valor de calibración	creatina quinasa-MB
Italiano	abbreviazione	Unità	Valore di calibrazione	creatina chinasi-MB
Türkçe	kısaltılmış ad	Ünite	Kalibrasyon Değeri	Kreatin Kinaz-MB