

## Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-240: BS-230, BS-240;

8.BS-600: BS-600, BS-620;

2.BS-240E: BS240E, BS240Pro;

9.BS-600M: BS-600M, BS-620M;

3.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00)

10.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

(Serial Number starts with "W9-" and software version starts with "35.00");

11.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

4.BS-380: BS-380, BS-390;

12.BS-2800M:BS-2600M, BS-2800M.

5.BS-400: BS-400, BS-420;

13.0.9% NaCl is needed, Conc. Of 0.9% NaCl =0;

6.BS-430: BS-430, BS-450, BS-460;

14.For applicable analyzers of the analyte, please subject to the parameter sheet and instrument.

7.BS-480: BS-480, BS-490;



: 153424002



: 2025-03-07

English	Abbreviated name	Unit	Calibration Rule	Model
Русский	сокращенное наименование	Прибор	Принцип калибровки	модель
Português	Nome abreviado	Unidade	Regra de calibração	Modelo
Español	nombre abreviado	Unidad	Regla de calibración	modelo
Italiano	abbreviazione	Unità	Regola di calibrazione	modelli
Türkçe	kısaltılmış ad	Ünite	Kalibrasyon Kuralı	model

English	Expanded Uncertainty
Русский	Расширенная неопределенность
Português	Incerteza Expandida
Español	Incertidumbre expandida
Italiano	Incertezza estesa
Türkçe	Daha Uzun Süreli Belirsizlik

UIBC	
English	Unsaturated Iron Binding Capacity
Русский	ненасыщенная железосвязывающая способность
Português	Capacidade de ligação de ferro insaturado
Español	Capacidad de unión de hierro no saturado
Italiano	Capacità di legame del ferro insaturo
Türkçe	Unsature Demir Bağlama Kapasitesi

Abbreviated name	Unit	Calibration Rule				Two-point linear	
		Model	Calibrator Value	Expanded Uncertainty	Model	Calibrator Value	Expanded Uncertainty
<b>UIBC</b>	μmol/L	<b>BS-240<sup>1</sup></b>	112	2.9	<b>BS-480<sup>7</sup></b>	112	2.9
		<b>BS-240E<sup>2</sup></b>	109	2.8	<b>BS-600<sup>8</sup></b>	112	2.9
		<b>BS-360E<sup>3</sup></b>	109	2.8	<b>BS-600M<sup>9</sup></b>	113	2.9
		<b>BS-380<sup>4</sup></b>	110	2.8	<b>BS-800<sup>10</sup></b>	112	3
		<b>BS-400<sup>5</sup></b>	110	2.8	<b>BS-2000<sup>11</sup></b>	113	2.9
		<b>BS-430<sup>6</sup></b>	112	2.9	<b>BS-2800M<sup>12</sup></b>	113	2.9
	μg/dL	<b>BS-240<sup>1</sup></b>	626	16	<b>BS-480<sup>7</sup></b>	626	16
		<b>BS-240E<sup>2</sup></b>	609	16	<b>BS-600<sup>8</sup></b>	626	16
		<b>BS-360E<sup>3</sup></b>	609	16	<b>BS-600M<sup>9</sup></b>	632	16
		<b>BS-380<sup>4</sup></b>	615	16	<b>BS-800<sup>10</sup></b>	626	17
		<b>BS-400<sup>5</sup></b>	615	16	<b>BS-2000<sup>11</sup></b>	632	16
		<b>BS-430<sup>6</sup></b>	626	16	<b>BS-2800M<sup>12</sup></b>	632	16

# **mindray**

## **UIBC Calibrator**

For use on: BS-2000, BS-2800M,  
BS-600M, BS-800

**LOT** 153424002

 2025-03-07

