

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов измерительной системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадает во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1. **BS-240**;

5. **BS-400**: BS-400, BS-420;

2. **BS-240E**: BS-240E, BS-240Pro;

6. **BS-430**: BS-450, BS-460;

3. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

7. **BS-480**: BS-480, BS-490;

BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), **BS-350E(V35.00)** (Serial Number starts with "W9-" and software version starts with "35.00");

8. **BS-600**: BS-600, BS-620;

9. **BS-600M**: BS-600M, BS-620M;

10. **BS-800**: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

4. **BS-380**: BS-380, BS-390;

11. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

12. **BS-2800M**: BS-2600M, BS-2800M.

13. **Cal0: 0.9% NaCl, Conc. Of Cal0=0;**

Русский: Cal0: 0,9% NaCl, конц. Cal0=0;

Português: Cal0: 0,9% NaCl, Conc. de Cal0=0;

Español: Cal0: 0,9% NaCl, Conc. de Cal0=0;

Italiano: Cal0: 0,9% NaCl, conc. di Cal0=0;

Türkçe: Cal0: %0,9 NaCl, Cal0 Kons.=0

LOT : 153722002

: 2023-05-15

| Abbreviated name | | MYO | Calibration Rule | | | spline | | |
|------------------|-------|---|------------------|-------|---|----------|-------|---|
| Model | Level | Calibration Value ¹³ (ng/mL) | Model | Level | Calibration Value ¹³ (ng/mL) | Model | Level | Calibration Value ¹³ (ng/mL) |
| BS-240 | a | / | BS-240E | a | / | BS-360E | a | / |
| | b | / | | b | / | | b | / |
| | c | / | | c | / | | c | / |
| | d | / | | d | / | | d | / |
| BS-380 | a | / | BS-400 | a | / | BS-430 | a | / |
| | b | / | | b | / | | b | / |
| | c | / | | c | / | | c | / |
| | d | / | | d | / | | d | / |
| BS-480 | a | / | BS-600 | a | / | BS-600M | a | 96.0 |
| | b | / | | b | / | | b | 192 |
| | c | / | | c | / | | c | 378 |
| | d | / | | d | / | | d | 729 |
| BS-800 | a | / | BS-2000 | a | / | BS-2800M | a | / |
| | b | / | | b | / | | b | / |
| | c | / | | c | / | | c | / |
| | d | / | | d | / | | d | / |

| English | Abbreviated name | Calibration Rule | Calibration Value | Model | Level |
|-----------|--------------------------|------------------------|------------------------|---------|---------|
| Русский | сокращенное наименование | Принцип калибровки | Эталонное значение | модель | Уровень |
| Português | Nome abreviado | Regra de calibração | Valor de calibração | Modelo | Nível |
| Español | nombre abreviado | Regla de calibración | Valor de calibración | modelo | Nivel |
| Italiano | abbreviazione | Regola di calibrazione | Valore di calibrazione | modelli | Livello |
| Türkçe | kısaltılmış ad | Kalibrasyon Kuralı | Kalibrasyon Değeri | model | Düzey |

| MYO | |
|-----------|------------|
| English | Myoglobin |
| Русский | Миоглобин |
| Português | Mioglobina |
| Español | Mioglobina |
| Italiano | Mioglobina |
| Türkçe | Myoglobin |