

Value sheet of Mindray BS Measurement System

Русский : Таблица результатов для системы BS компании Mindray

Português : Planilha de valores do Sistema de Medição BS da Mindray

Español : Hoja de valores del sistema de medición Mindray BS

Italiano : Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe : Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский : Данные совпадают во всех группах.

Português : A dados de cada grupo é a mesma.

Español : la datos de cada grupo es la misma.

Italiano : la dati di ogni gruppo è la stessa.

Türkçe : her grubun veri aynıdır.

1.**BS-200**: BS-200, BS-220;

8.**BS-380**: BS-380, BS-390;

2.**BS-200E**: BS-200E, BS-220E;

9.**BS-400**: BS-400, BS-420;

3.**BS-230**: BS-230, BS-240;

10.**BS-430**: BS-430, BS-450, BS-460;

4.**BS-240E**: BS240E, BS240Pro;

11.**BS-480**: BS-480, BS-490;

5.**BS-300**: BS-300, BS-320;

12.**BS-600**: BS-600, BS-620;

6.**BS-330E**: BS-330E(Serial Number starts with "XQ-"),
BS-350E(Serial Number starts with "XS-");

13.**BS-600M**: BS-600M, BS-620M;

7.**BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,
BS-330E(V35.00)(Serial Number starts with "W8-"
and software version starts with "35.00"),
BS-350E(V35.00) (Serial Number starts with "W9-"
and software version starts with "35.00");

14.**BS-800**: BS-800, BS-820, BS-800M, BS-820M,
BS-1800, BS-1800plus;

15.**BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

16.**BS-2800M**: BS-2600M, BS-2800M.

LOT: 155123003

EXP: 2024-09-29

| Abbreviated name | | Hb | HbA1c |
|------------------------------|-------|---------------------------|---------------------------|
| Model | Level | Calibrator Value (µmol/L) | Calibrator Value (µmol/L) |
| BS-200 ¹ | S1 | 97.7 | 3.60 |
| | S2 | 150 | 16.0 |
| BS-200E ² | S1 | 95.2 | 3.60 |
| | S2 | 151 | 16.2 |
| BS-230 ³ | S1 | 95.2 | 3.60 |
| | S2 | 151 | 16.2 |
| BS-240E ⁴ | S1 | 95.2 | 3.60 |
| | S2 | 151 | 16.2 |
| BS-300 ⁵ | S1 | 95.2 | 3.60 |
| | S2 | 151 | 16.2 |
| BS-330E ⁶ | S1 | 95.2 | 3.60 |
| | S2 | 151 | 16.2 |
| BS-360E ⁷ | S1 | 95.2 | 3.60 |
| | S2 | 151 | 16.2 |
| BS-380 ⁸ | S1 | 95.2 | 3.60 |
| | S2 | 151 | 16.2 |
| BS-400 ⁹ | S1 | 95.2 | 3.60 |
| | S2 | 151 | 16.2 |
| BS-430 ¹⁰ | S1 | 95.2 | 3.60 |
| | S2 | 151 | 16.2 |
| BS-480 ¹¹ | S1 | 95.2 | 3.60 |
| | S2 | 151 | 16.2 |
| BS-600 ¹² | S1 | 95.2 | 3.60 |
| | S2 | 151 | 16.2 |
| BS-600M ¹³ | S1 | 95.2 | 3.60 |
| | S2 | 151 | 16.1 |
| BS-800 ¹⁴ | S1 | 95.2 | 3.60 |
| | S2 | 151 | 16.2 |

| BS-2000 ¹⁵ | S1 | 95.6 | 3.60 | |
|-------------------------------|--------------------------|--------------------------|--------------|--------------|
| | S2 | 152 | 16.2 | |
| BS-2800M ¹⁶ | S1 | 95.4 | 3.60 | |
| | S2 | 151 | 16.2 | |
| English | Abbreviated name | Calibration Value | Model | Level |
| Русский | сокращенное наименование | Эталонное значение | модель | Уровень |
| Português | Nome abreviado | Valor de calibração | Modelo | Nível |
| Español | nombre abreviado | Valor de calibración | modelo | Nivel |
| Italiano | abbreviazione | Valore di calibrazione | modelli | Livello |
| Türkçe | kısaltılmış ad | Kalibrasyon Değeri | model | Düzey |

| HbA1c/Hb | |
|------------------|------------------------------|
| English | Hemoglobin A1c/ Hemoglobin |
| Русский | гемоглобин A1c / гемоглобин |
| Português | Hemoglobina A1c /Hemoglobina |
| Español | hemoglobina A1c /hemoglobina |
| Italiano | Emoglobina A1c/Emoglobina |
| Türkçe | Hemoglobin A1c / Hemoglobin |