

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

- | | |
|---|--|
| 1. BS-200 : BS-200, BS-220; | 8. BS-380 : BS-380, BS-390; |
| 2. BS-200E : BS-200E, BS-220E; | 9. BS-400 : BS-400, BS-420; |
| 3. BS-230 : BS-230, BS-240; | 10. BS-430 : BS-430, BS-450, BS-460; |
| 4. BS-240E : BS-240E, BS-240Pro; | 11. BS-480 : BS-480, BS-490; |
| 5. BS-300 : BS-300, BS-320; | 12. BS-600 : BS-600, BS-620; |
| 6. BS-330E : BS-330E (Serial Number starts with "XQ-"),
BS-350E (Serial Number starts with "XS-"); | 13. BS-600M : BS-600M, BS-620M; |
| 7. BS-360E : BS-360E, BS-370E, BS-350S, BS-360S,
BS-330E(V35.00) (Serial Number starts with "W8-" and
software version starts with "35.00"), BS-350E(V35.00)
(Serial Number starts with "W9-" and software version
starts with "35.00"); | 14. BS-800 : BS-800, BS-820, BS-800M, BS-820M, BS-1800,
BS-1800plus; |
| | 15. BS-2000 : BS-2000, BS-2200, BS-2000M, BS-2200M; |
| | 16. BS-2800M : BS-2600M, BS-2800M. |

17. IFCC: $HbA1c(mmol/mol) = HbA1c(\mu mol/L) / Hb(\mu mol/L) \times 1000$
NGSP/DCCT: $HbA1c(\%) = 91.5 \times HbA1c(\mu mol/L) / Hb(\mu mol/L) + 2.15$

LOT : 044623014

: 2025-05-16

Abbreviated name	Model	Unit	Assay Value	Range(Assay Value \pm 3SD)	1 SD	Model	Unit	Assay Value	Range(Assay Value \pm 3SD)	1 SD
HbA1c	BS-200 ¹	$\mu mol/L$	11.2	10.0 — 12.4	0.4	BS-400 ⁹	$\mu mol/L$	11.4	10.2 — 12.6	0.4
	BS-200E ²	$\mu mol/L$	11.2	10.0 — 12.4	0.4	BS-430 ¹⁰	$\mu mol/L$	11.4	10.2 — 12.6	0.4
	BS-230 ³	$\mu mol/L$	11.4	10.2 — 12.6	0.4	BS-480 ¹¹	$\mu mol/L$	11.4	10.2 — 12.6	0.4
	BS-240E ⁴	$\mu mol/L$	11.4	10.2 — 12.6	0.4	BS-600 ¹²	$\mu mol/L$	11.4	10.2 — 12.6	0.4
	BS-300 ⁵	$\mu mol/L$	11.4	10.2 — 12.6	0.4	BS-600M ¹³	$\mu mol/L$	11.3	10.1 — 12.5	0.4
	BS-330E ⁶	$\mu mol/L$	11.2	10.0 — 12.4	0.4	BS-800 ¹⁴	$\mu mol/L$	11.4	10.2 — 12.6	0.4
	BS-360E ⁷	$\mu mol/L$	11.4	10.2 — 12.6	0.4	BS-2000 ¹⁵	$\mu mol/L$	11.4	10.2 — 12.6	0.4
	BS-380 ⁸	$\mu mol/L$	11.2	10.0 — 12.4	0.4	BS-2800M ¹⁶	$\mu mol/L$	11.4	10.2 — 12.6	0.4
Hb	BS-200 ¹	$\mu mol/L$	129	117 — 141	4	BS-400 ⁹	$\mu mol/L$	129	117 — 141	4
	BS-200E ²	$\mu mol/L$	130	118 — 142	4	BS-430 ¹⁰	$\mu mol/L$	129	117 — 141	4
	BS-230 ³	$\mu mol/L$	129	117 — 141	4	BS-480 ¹¹	$\mu mol/L$	129	117 — 141	4
	BS-240E ⁴	$\mu mol/L$	129	117 — 141	4	BS-600 ¹²	$\mu mol/L$	129	117 — 141	4
	BS-300 ⁵	$\mu mol/L$	130	118 — 142	4	BS-600M ¹³	$\mu mol/L$	129	117 — 141	4
	BS-330E ⁶	$\mu mol/L$	130	118 — 142	4	BS-800 ¹⁴	$\mu mol/L$	129	117 — 141	4
	BS-360E ⁷	$\mu mol/L$	129	117 — 141	4	BS-2000 ¹⁵	$\mu mol/L$	130	118 — 142	4
	BS-380 ⁸	$\mu mol/L$	130	118 — 142	4	BS-2800M ¹⁶	$\mu mol/L$	130	118 — 142	4

HbA1C Control P



Abbreviated name	Model	Unit	Assay Value	Range(Assay Value±3SD)		1 SD	Model	Unit	Assay Value	Range(Assay Value±3SD)		1 SD		
HbA1c(mmol/mol) ¹⁷ IFCC	BS-200 ¹	mmol/mol	87.0	78.3	—	95.7	2.9	BS-400 ⁹	mmol/mol	88.2	79.5	—	96.9	2.9
	BS-200E ²	mmol/mol	86.6	77.9	—	95.3	2.9	BS-430 ¹⁰	mmol/mol	87.9	79.2	—	96.6	2.9
	BS-230 ³	mmol/mol	87.9	79.2	—	96.6	2.9	BS-480 ¹¹	mmol/mol	87.9	79.2	—	96.6	2.9
	BS-240E ⁴	mmol/mol	87.9	79.2	—	96.6	2.9	BS-600 ¹²	mmol/mol	87.9	79.2	—	96.6	2.9
	BS-300 ⁵	mmol/mol	87.7	79.0	—	96.4	2.9	BS-600M ¹³	mmol/mol	87.6	78.9	—	96.3	2.9
	BS-330E ⁶	mmol/mol	86.6	77.9	—	95.3	2.9	BS-800 ¹⁴	mmol/mol	87.9	79.2	—	96.6	2.9
	BS-360E ⁷	mmol/mol	87.9	79.2	—	96.6	2.9	BS-2000 ¹⁵	mmol/mol	87.5	78.8	—	96.2	2.9
	BS-380 ⁸	mmol/mol	86.6	77.9	—	95.3	2.9	BS-2800M ¹⁶	mmol/mol	87.5	78.8	—	96.2	2.9
HbA1c (%) ¹⁷ NGSP / DCCT	BS-200 ¹	%	10.1	9.2	—	11.0	0.3	BS-400 ⁹	%	10.2	9.3	—	11.1	0.3
	BS-200E ²	%	10.1	9.2	—	11.0	0.3	BS-430 ¹⁰	%	10.2	9.3	—	11.1	0.3
	BS-230 ³	%	10.2	9.3	—	11.1	0.3	BS-480 ¹¹	%	10.2	9.3	—	11.1	0.3
	BS-240E ⁴	%	10.2	9.3	—	11.1	0.3	BS-600 ¹²	%	10.2	9.3	—	11.1	0.3
	BS-300 ⁵	%	10.2	9.3	—	11.1	0.3	BS-600M ¹³	%	10.2	9.3	—	11.1	0.3
	BS-330E ⁶	%	10.1	9.2	—	11.0	0.3	BS-800 ¹⁴	%	10.2	9.3	—	11.1	0.3
	BS-360E ⁷	%	10.2	9.3	—	11.1	0.3	BS-2000 ¹⁵	%	10.2	9.3	—	11.1	0.3
	BS-380 ⁸	%	10.1	9.2	—	11.0	0.3	BS-2800M ¹⁶	%	10.2	9.3	—	11.1	0.3

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон (результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa (Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango (Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore diconcentrazione ±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

	HbA1c	Hb
English	Hemoglobin A1c	Hemoglobin
Русский	гемоглобин A1c	гемоглобин
Português	Hemoglobina A1c	Hemoglobina
Español	hemoglobina A1c	hemoglobina
Italiano	Emoglobina A1c	Emoglobina
Türkçe	Hemoglobin A1c	Hemoglobin