

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası

**The data of each group is same.****Português: A dados de cada grupo é a mesma.****Italiano: la dati di ogni gruppo è la stessa.**

1.BS-120: BS-120, BS-130, BS-180, BS-190;

2.BS-200: BS-200, BS-220, BS-330, BS-350;

3.BS-200E: BS-200E, BS-220E,;

4.BS-240: BS-230, BS-240;

5.BS-240E: BS240E, BS240Pro;

6.BS-300: BS-300, BS-320;

7.BS-330E:BS-330E(Serial Number starts with "XQ-"), BS-350E(Serial Number starts with "XS-")

8.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-

330E(V35.00)(Serial Number starts with "W8-" and software version

starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");

9.BS-380: BS-380, BS-390;

Русский: Данные совпадают во всех группах.**Español: la datos de cada grupo es la misma.****Türkçe: her grubun veri aynıdır.**

10.BS-400: BS-400, BS-420;

11.BS-430: BS-430, BS-450,BS-460;

12.BS-480: BS-480, BS-490;

13.BS-600: BS-600, BS-620;

14.BS-600M: BS-600M, BS-620M;

15.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

16.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

17.BS-2800M:BS-2600M, BS-2800M;

18.Applicable models of the chemistry shall be subject to the parameter list and instructions.

LOT: 044823013**: 2025-02-28**

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)
LIP					
English			Lipase		
Русский			Липаза		
Português			Lipase		
Español			Lipasa		
Italiano			Lipasi		
Türkçe			Lipaz		

Abbreviated name	Unit	Model	Assay Value	Range(Assay Value±3SD)	1 SD	Model	Assay Value	Range(Assay Value±3SD)	1 SD
LIP	U/L	BS-120 ¹	47.7	38.1 — 57.3	3.2	BS-400 ¹⁰	45.8	36.5 — 55.1	3.1
		BS-200 ²	42.8	34.1 — 51.5	2.9	BS-430 ¹¹	45.4	36.4 — 54.4	3.0
		BS-200E ³	46.2	36.9 — 55.5	3.1	BS-480 ¹²	45.5	36.5 — 54.5	3.0
		BS-240 ⁴	48.2	38.6 — 57.8	3.2	BS-600 ¹³	45.6	36.3 — 54.9	3.1
		BS-240E ⁵	45.8	36.5 — 55.1	3.1	BS-600M ¹⁴	45.3	36.3 — 54.3	3.0
		BS-300 ⁶	44.4	35.4 — 53.4	3.0	BS-800 ¹⁵	45.6	36.3 — 54.9	3.1
		BS-330E ⁷	46.2	36.9 — 55.5	3.1	BS-2000 ¹⁶	45.2	36.2 — 54.2	3.0
		BS-360E ⁸	43.9	35.2 — 52.6	2.9	BS-2800M ¹⁷	45.2	36.2 — 54.2	3.0
		BS-380 ⁹	45.8	36.5 — 55.1	3.1				
		µkat/L		BS-120 ¹	0.797	0.636 — 0.957	0.053	BS-400 ¹⁰	0.765
BS-200 ²	0.715			0.569 — 0.860	0.048	BS-430 ¹¹	0.758	0.608 — 0.908	0.050
BS-200E ³	0.772			0.616 — 0.927	0.052	BS-480 ¹²	0.760	0.610 — 0.910	0.050
BS-240 ⁴	0.805			0.645 — 0.965	0.053	BS-600 ¹³	0.762	0.606 — 0.917	0.052
BS-240E ⁵	0.765			0.610 — 0.920	0.052	BS-600M ¹⁴	0.757	0.606 — 0.907	0.050

Abbreviated name	Unit	Model	Assay Value	Range(Assay Value \pm 3SD)		1 SD	Model	Assay Value	Range(Assay Value \pm 3SD)		1 SD		
LIP	μ kat/L	BS-300 ⁶	0.741	0.591	—	0.892	0.050	BS-800 ¹⁵	0.762	0.606	—	0.917	0.052
		BS-330E ⁷	0.772	0.616	—	0.927	0.052	BS-2000 ¹⁶	0.755	0.605	—	0.905	0.050
		BS-360E ⁸	0.733	0.588	—	0.878	0.048	BS-2800M ¹⁷	0.755	0.605	—	0.905	0.050
		BS-380 ⁹	0.765	0.610	—	0.920	0.052						