

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1. **BS-200:** BS-200, BS-220;

8. **BS-380:** BS-380, BS-390;

2. **BS-200E:** BS-200E, BS-220E;

9. **BS-400:** BS-400, BS-420;

3. **BS-230:** BS-230, BS-240;

10. **BS-430:** BS-430, BS-450, BS-460;

4. **BS-240E:** BS-240E, BS-240Pro;

11. **BS-480:** BS-480, BS-490;

5. **BS-300:** BS-300, BS-320;

12. **BS-600:** BS-600, BS-620;

6. **BS-330E: BS-330E**(Serial Number starts with "XQ-"),

13. **BS-600M:** BS-600M, BS-620M;

BS-350E(Serial Number starts with "XS-");

14. **BS-800:** BS-800, BS-820, BS-800M, BS-820M, BS-1800,

7. **BS-360E:** BS-360E, BS-370E, BS-350S, BS-360S,

BS-1800plus;

BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), **BS-350E(V35.00)**

15. **BS-2000:** BS-2000, BS-2200, BS-2000M, BS-2200M;

(Serial Number starts with "W9-" and software version starts with "35.00");

16. **BS-2800M:** BS-2600M, BS-2800M.

17. **IFCC:** $HbA1c(\%) = HbA1c(\mu mol/L) / Hb(\mu mol/L) \times 100$

NGSP/DCCT: $HbA1c(\%) = 91.5 \times HbA1c(\mu mol/L) / Hb(\mu mol/L) + 2.15$

JCCLS: $HbA1c(\%) = 96.3 \times HbA1c(\mu mol/L) / Hb(\mu mol/L) + 1.62$

LOT: 052523006

EXP: 2024-12-03

Abbreviated name	Model	Unit	Assay Value	Range(Assay Value \pm 3SD)		1 SD	Model	Unit	Assay Value	Range(Assay Value \pm 3SD)		1 SD		
HbA1c	BS-200 ¹	$\mu mol/L$	11.3	10.2	—	12.4	0.4	BS-400 ⁹	$\mu mol/L$	11.2	10.1	—	12.3	0.4
	BS-200E ²	$\mu mol/L$	11.4	10.3	—	12.5	0.4	BS-430 ¹⁰	$\mu mol/L$	11.3	10.2	—	12.4	0.4
	BS-230 ³	$\mu mol/L$	11.3	10.2	—	12.4	0.4	BS-480 ¹¹	$\mu mol/L$	11.3	10.2	—	12.4	0.4
	BS-240E ⁴	$\mu mol/L$	11.3	10.2	—	12.4	0.4	BS-600 ¹²	$\mu mol/L$	11.3	10.2	—	12.4	0.4
	BS-300 ⁵	$\mu mol/L$	11.4	10.3	—	12.5	0.4	BS-600M ¹³	$\mu mol/L$	11.4	10.3	—	12.5	0.4
	BS-330E ⁶	$\mu mol/L$	11.4	10.3	—	12.5	0.4	BS-800 ¹⁴	$\mu mol/L$	11.3	10.2	—	12.4	0.4
	BS-360E ⁷	$\mu mol/L$	11.3	10.2	—	12.4	0.4	BS-2000 ¹⁵	$\mu mol/L$	11.4	10.3	—	12.5	0.4
	BS-380 ⁸	$\mu mol/L$	11.4	10.3	—	12.5	0.4	BS-2800M ¹⁶	$\mu mol/L$	11.5	10.4	—	12.6	0.4
Hb	BS-200 ¹	$\mu mol/L$	130	117	—	143	4	BS-400 ⁹	$\mu mol/L$	130	117	—	143	4
	BS-200E ²	$\mu mol/L$	130	117	—	143	4	BS-430 ¹⁰	$\mu mol/L$	130	117	—	143	4
	BS-230 ³	$\mu mol/L$	130	117	—	143	4	BS-480 ¹¹	$\mu mol/L$	130	117	—	143	4
	BS-240E ⁴	$\mu mol/L$	130	117	—	143	4	BS-600 ¹²	$\mu mol/L$	130	117	—	143	4
	BS-300 ⁵	$\mu mol/L$	131	118	—	144	4	BS-600M ¹³	$\mu mol/L$	129	116	—	142	4
	BS-330E ⁶	$\mu mol/L$	130	117	—	143	4	BS-800 ¹⁴	$\mu mol/L$	130	117	—	143	4
	BS-360E ⁷	$\mu mol/L$	130	117	—	143	4	BS-2000 ¹⁵	$\mu mol/L$	131	118	—	144	4
	BS-380 ⁸	$\mu mol/L$	130	117	—	143	4	BS-2800M ¹⁶	$\mu mol/L$	130	117	—	143	4

HbA1C Control P



Abbreviated name	Model	Unit	Assay Value	Range(Assay Value±3SD)		1 SD	Model	Unit	Assay Value	Range(Assay Value±3SD)		1 SD		
HbA1c (%) IFCC	BS-200 ¹	%	8.71	7.85	—	9.57	0.29	BS-400 ⁹	%	8.67	7.81	—	9.53	0.29
	BS-200E ²	%	8.82	7.95	—	9.69	0.29	BS-430 ¹⁰	%	8.75	7.88	—	9.62	0.29
	BS-230 ³	%	8.75	7.88	—	9.62	0.29	BS-480 ¹¹	%	8.75	7.88	—	9.62	0.29
	BS-240E ⁴	%	8.75	7.88	—	9.62	0.29	BS-600 ¹²	%	8.75	7.88	—	9.62	0.29
	BS-300 ⁵	%	8.71	7.85	—	9.57	0.29	BS-600M ¹³	%	8.85	7.97	—	9.73	0.29
	BS-330E ⁶	%	8.82	7.95	—	9.69	0.29	BS-800 ¹⁴	%	8.75	7.88	—	9.62	0.29
	BS-360E ⁷	%	8.75	7.88	—	9.62	0.29	BS-2000 ¹⁵	%	8.74	7.87	—	9.61	0.29
	BS-380 ⁸	%	8.82	7.95	—	9.69	0.29	BS-2800M ¹⁶	%	8.80	7.93	—	9.67	0.29
HbA1c (%) NGSP / DCCT	BS-200 ¹	%	10.1	9.1	—	11.1	0.3	BS-400 ⁹	%	10.1	9.1	—	11.1	0.3
	BS-200E ²	%	10.2	9.2	—	11.2	0.3	BS-430 ¹⁰	%	10.2	9.2	—	11.2	0.3
	BS-230 ³	%	10.2	9.2	—	11.2	0.3	BS-480 ¹¹	%	10.2	9.2	—	11.2	0.3
	BS-240E ⁴	%	10.2	9.2	—	11.2	0.3	BS-600 ¹²	%	10.2	9.2	—	11.2	0.3
	BS-300 ⁵	%	10.1	9.1	—	11.1	0.3	BS-600M ¹³	%	10.2	9.2	—	11.2	0.3
	BS-330E ⁶	%	10.2	9.2	—	11.2	0.3	BS-800 ¹⁴	%	10.2	9.2	—	11.2	0.3
	BS-360E ⁷	%	10.2	9.2	—	11.2	0.3	BS-2000 ¹⁵	%	10.1	9.1	—	11.1	0.3
	BS-380 ⁸	%	10.2	9.2	—	11.2	0.3	BS-2800M ¹⁶	%	10.2	9.2	—	11.2	0.3
HbA1c (%) JCCLS	BS-200 ¹	%	10.0	9.0	—	11.0	0.3	BS-400 ⁹	%	10.0	9.0	—	11.0	0.3
	BS-200E ²	%	10.1	9.1	—	11.1	0.3	BS-430 ¹⁰	%	10.0	9.0	—	11.0	0.3
	BS-230 ³	%	10.0	9.0	—	11.0	0.3	BS-480 ¹¹	%	10.0	9.0	—	11.0	0.3
	BS-240E ⁴	%	10.0	9.0	—	11.0	0.3	BS-600 ¹²	%	10.0	9.0	—	11.0	0.3
	BS-300 ⁵	%	10.0	9.0	—	11.0	0.3	BS-600M ¹³	%	10.1	9.1	—	11.1	0.3
	BS-330E ⁶	%	10.1	9.1	—	11.1	0.3	BS-800 ¹⁴	%	10.0	9.0	—	11.0	0.3
	BS-360E ⁷	%	10.0	9.0	—	11.0	0.3	BS-2000 ¹⁵	%	10.0	9.0	—	11.0	0.3
	BS-380 ⁸	%	10.1	9.1	—	11.1	0.3	BS-2800M ¹⁶	%	10.1	9.1	—	11.1	0.3

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон (результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa (Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango (Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore diconcentrazione ±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

	HbA1c	Hb
English	Hemoglobin A1c	Hemoglobin
Русский	гемоглобин A1c	гемоглобин
Português	Hemoglobina A1c	Hemoglobina
Español	hemoglobina A1c	hemoglobina
Italiano	Emoglobina A1c	Emoglobina
Türkçe	Hemoglobin A1c	Hemoglobin