

## Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası

**The data of each group is same.****Português: A dados de cada grupo é a mesma.****Italiano: la dati di ogni gruppo è la stessa.**

1.BS-120: BS-120, BS-130, BS-180, BS-190;

2.BS-200: BS-200, BS-220, BS-330, BS-350;

3.BS-200E: BS-200E, BS-220E,;

4.BS-240: BS-230, BS-240;

5.BS-240E: BS240E, BS240Pro;

6.BS-330E:BS-330E(Serial Number starts with "XQ-"), BS-350E(Serial Number starts with "XS-")

7.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-

330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");

8.BS-380: BS-380, BS-390;

9.BS-400: BS-400, BS-420;

**Русский: Данные совпадают во всех группах.****Español: la datos de cada grupo es la misma.****Türkçe: her grubun veri aynıdır.**

10.BS-430: BS-430, BS-450,BS-460;

11.BS-480: BS-480, BS-490;

12.BS-600: BS-600, BS-620;

13.BS-600M: BS-600M, BS-620M;

14.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

15.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

16.BS-2800M:BS-2600M, BS-2800M;

17.Applicable models of the chemistry shall be subject to the parameter list and instructions.

**LOT: 053124001****Expires: 2025-07-02**

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
<b>Русский</b>	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
<b>Português</b>	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
<b>Español</b>	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
<b>Italiano</b>	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
<b>Türkçe</b>	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

**MALB**

<b>English</b>	Microalbumin
<b>Русский</b>	микроальбумина
<b>Português</b>	Microalbuminúria
<b>Español</b>	Microalbúmina
<b>Italiano</b>	Microalbumina
<b>Türkçe</b>	Mikroalbumin

Abbreviated name	Unit	Model	Assay Value	Range(Assay Value±3SD)	1 SD	Model	Assay Value	Range(Assay Value±3SD)	1 SD
<b>MALB</b>	mg/L	<b>BS-120<sup>1</sup></b>	203	173 — 233	10	<b>BS-430<sup>10</sup></b>	199	169 — 229	10
		<b>BS-200<sup>2</sup></b>	198	168 — 228	10	<b>BS-480<sup>11</sup></b>	201	171 — 231	10
		<b>BS-200E<sup>3</sup></b>	198	168 — 228	10	<b>BS-600<sup>12</sup></b>	201	171 — 231	10
		<b>BS-240<sup>4</sup></b>	203	173 — 233	10	<b>BS-600M<sup>13</sup></b>	200	170 — 230	10
		<b>BS-240E<sup>5</sup></b>	201	171 — 231	10	<b>BS-800<sup>14</sup></b>	201	171 — 231	10
		<b>BS-360E<sup>7</sup></b>	201	171 — 231	10	<b>BS-2000<sup>15</sup></b>	198	168 — 228	10
		<b>BS-380<sup>8</sup></b>	198	168 — 228	10	<b>BS-2800M<sup>16</sup></b>	196	166 — 226	10
		<b>BS-400<sup>9</sup></b>	199	169 — 229	10				
		µmol/L	<b>BS-120<sup>1</sup></b>	3.09	2.63 — 3.54	0.15	<b>BS-430<sup>10</sup></b>	3.02	2.57 — 3.48
<b>BS-200<sup>2</sup></b>	3.01		2.55 — 3.47	0.15	<b>BS-480<sup>11</sup></b>	3.06	2.60 — 3.51	0.15	
<b>BS-200E<sup>3</sup></b>	3.01		2.55 — 3.47	0.15	<b>BS-600<sup>12</sup></b>	3.06	2.60 — 3.51	0.15	
<b>BS-240<sup>4</sup></b>	3.09		2.63 — 3.54	0.15	<b>BS-600M<sup>13</sup></b>	3.04	2.58 — 3.50	0.15	
<b>BS-240E<sup>5</sup></b>	3.06		2.60 — 3.51	0.15	<b>BS-800<sup>14</sup></b>	3.06	2.60 — 3.51	0.15	
<b>BS-360E<sup>7</sup></b>	3.06		2.60 — 3.51	0.15	<b>BS-2000<sup>15</sup></b>	3.01	2.55 — 3.47	0.15	

Abbreviated name	Unit	Model	Assay Value	Range(Assay Value $\pm$ 3SD)			1 SD	Model	Assay Value	Range(Assay Value $\pm$ 3SD)			1 SD
		<b>BS-380<sup>8</sup></b>	3.01	2.55	—	3.47	0.15	<b>BS-2800M<sup>16</sup></b>	2.98	2.52	—	3.44	0.15
		<b>BS-400<sup>9</sup></b>	3.02	2.57	—	3.48	0.15						

# mindray

## MALB Control

For use on: BS-2000, BS-2800M,  
BS-600M, BS-800

LOT 053124001

2025-07-02

