

Value sheet of Mindray BS Measurement System

Русский : Таблица результатов для системы BS компании Mindray

Português : Planilha de valores do Sistema de Medição BS da Mindray

Español : Hoja de valores del sistema de medición Mindray BS

Italiano : Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe : Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский : Данные совпадают во всех группах.

Português : A dados de cada grupo é a mesma.

Español : la datos de cada grupo es la misma.

Italiano : la dati di ogni gruppo è la stessa.

Türkçe : her grubun veri aynıdır.

1. **BS-230**: BS-230, BS-240;

5. **BS-400**: BS-400, BS-420;

2. **BS-240E**: BS-240E, BS-240Pro;

6. **BS-430**: BS-430, BS-450, BS-460;

3. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

7. **BS-480**: BS-480, BS-490;

8. **BS-600**: BS-600, BS-620;

BS-330E(V35.00) (Serial Number starts with "W8-" and software version starts with "35.00"), **BS-350E(V35.00)**(Serial Number starts with "W9-"and software version starts with "35.00");

9. **BS-600M**: BS-600M, BS-620M;

10. **BS-800**: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

11. **BS-2000**: BS-2000,BS-2200,BS-2000M,BS-2200M;

4. **BS-380**: BS-380, BS-390;

12. **BS-2800M**: BS-2600M, BS-2800M.

LOT : 053522004

Exp : 2023-12-11

Abbreviated name	Model	Unit	Assay Value	Range(Assay Value±3SD)		1 SD	Model	Unit	Assay Value	Range(Assay Value±3SD)		1 SD		
UIBC	BS-230 ¹	μmol/L	33.0	26.4	—	39.6	2.2	BS-480 ⁷	μmol/L	35.0	28.0	—	42.0	2.3
		μg/dL	184	148	—	221	12		μg/dL	196	157	—	235	13
	BS-240E ²	μmol/L	34.4	27.5	—	41.3	2.3	BS-600 ⁸	μmol/L	31.6	25.2	—	38.0	2.1
		μg/dL	192	154	—	231	13		μg/dL	177	141	—	212	12
	BS-360E ³	μmol/L	34.8	27.8	—	41.8	2.3	BS-600M ⁹	μmol/L	33.2	26.5	—	39.9	2.2
		μg/dL	195	155	—	234	13		μg/dL	186	148	—	223	12
	BS-380 ⁴	μmol/L	34.6	27.6	—	41.6	2.3	BS-800 ¹⁰	μmol/L	31.6	25.2	—	38.0	2.1
		μg/dL	193	154	—	233	13		μg/dL	177	141	—	212	12
	BS-400 ⁵	μmol/L	34.6	27.6	—	41.6	2.3	BS-2000 ¹¹	μmol/L	31.7	25.3	—	38.1	2.1
		μg/dL	193	154	—	233	13		μg/dL	177	141	—	213	12
	BS-430 ⁶	μmol/L	31.6	25.2	—	38.0	2.1	BS-2800M ¹²	μmol/L	30.9	24.7	—	37.1	2.1
		μg/dL	177	141	—	212	12		μg/dL	173	138	—	207	12

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон (результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa (Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango (Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore diconcentrazione ±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

UIBC

English	Unsaturated Iron Binding Capacity
Русский	ненасыщенная железосвязывающая способность
Português	Capacidade de ligação de ferro insaturado
Español	Capacidad de unión de hierro no saturado
Italiano	Capacità di legame del ferro insaturo
Türkçe	Unsature Demir Bağlama Kapasitesi