

## Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray  
 Português: Planilha de valores do Sistema de Medição BS da Mindray  
 Español: Hoja de valores del sistema de medición Mindray BS  
 Italiano: Scheda dei valori del sistema di misurazione BS di Mindray  
 Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



**The data of each group is same.**

**Português: A dados de cada grupo é a mesma.**

**Italiano: la dati di ogni gruppo è la stessa.**

- 1.BS-240: BS-230, BS-240;
- 2.BS-240E: BS240E, BS240Pro;
- 3.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");
- 4.BS-380: BS-380, BS-390;
- 5.BS-400: BS-400, BS-420;
- 6.BS-430: BS-430, BS-450,BS-460;
- 7.BS-480: BS-480, BS-490;

**Русский: Данные совпадают во всех группах.**

**Español: la datos de cada grupo es la misma.**

**Türkçe: her grubun veri aynıdır.**

- 8.BS-600: BS-600, BS-620;
- 9.BS-600M: BS-600M, BS-620M;
- 10.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;
- 11.BS-1000M: BS-1000M, BS-1100M;
- 12.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;
- 13.BS-2800M:BS-2600M, BS-2800M;
- 14.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.

**LOT: 053524003**

**Expires: 2025-06-27**

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
<b>Русский</b>	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
<b>Português</b>	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
<b>Español</b>	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
<b>Italiano</b>	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
<b>Türkçe</b>	kisaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

**UIBC**

<b>English</b>	Unsaturated Iron Binding Capacity
<b>Русский</b>	ненасыщенная железосвязывающая способность
<b>Português</b>	Capacidade de ligação de ferro insaturado
<b>Español</b>	Capacidad de unión de hierro no saturado
<b>Italiano</b>	Capacità di legame del ferro insaturo
<b>Türkçe</b>	Unsature Demir Bağlama Kapasitesi

Abbreviated name	Unit	Model	Assay Value	Range(Assay Value±3SD)	1 SD	Model	Assay Value	Range(Assay Value±3SD)	1 SD
<b>UIBC</b>	µmol/L	<b>BS-240</b> <sup>1</sup>	32.1	25.5 — 38.7	2.2	<b>BS-600</b> <sup>8</sup>	31.2	24.9 — 37.5	2.1
		<b>BS-240E</b> <sup>2</sup>	34.6	27.7 — 41.5	2.3	<b>BS-600M</b> <sup>9</sup>	32.7	26.1 — 39.3	2.2
		<b>BS-360E</b> <sup>3</sup>	32.1	25.5 — 38.7	2.2	<b>BS-800</b> <sup>10</sup>	32.8	26.2 — 39.4	2.2
		<b>BS-380</b> <sup>4</sup>	35.9	28.7 — 43.1	2.4	<b>BS-1000M</b> <sup>11</sup>	32.7	26.1 — 39.3	2.2
		<b>BS-400</b> <sup>5</sup>	35.9	28.7 — 43.1	2.4	<b>BS-2000</b> <sup>12</sup>	31.4	25.1 — 37.7	2.1
		<b>BS-430</b> <sup>6</sup>	32.8	26.2 — 39.4	2.2	<b>BS-2800M</b> <sup>13</sup>	31.2	24.9 — 37.5	2.1
		<b>BS-480</b> <sup>7</sup>	36.1	28.9 — 43.3	2.4				
	µg/dL	<b>BS-240</b> <sup>1</sup>	179	143 — 216	12	<b>BS-600</b> <sup>8</sup>	174	139 — 210	12
		<b>BS-240E</b> <sup>2</sup>	193	155 — 232	13	<b>BS-600M</b> <sup>9</sup>	183	146 — 220	12
		<b>BS-360E</b> <sup>3</sup>	179	143 — 216	12	<b>BS-800</b> <sup>10</sup>	183	146 — 220	12
		<b>BS-380</b> <sup>4</sup>	201	160 — 241	13	<b>BS-1000M</b> <sup>11</sup>	183	146 — 220	12
		<b>BS-400</b> <sup>5</sup>	201	160 — 241	13	<b>BS-2000</b> <sup>12</sup>	176	140 — 211	12
		<b>BS-430</b> <sup>6</sup>	183	146 — 220	12	<b>BS-2800M</b> <sup>13</sup>	174	139 — 210	12
		<b>BS-480</b> <sup>7</sup>	202	162 — 242	13				

**mindray**  
UIBC Control

For use on: BS-1000M, BS-2000,  
BS-2800M, BS-600M,  
BS-800

LOT 053524003

2025-06-27

