

## Value sheet of Mindray BS Measurement System

Русский: Таблица результатов измерительной системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадает во всех группах. Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma. Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

- |   |  |
|---|--|
| 1. <b>BS-200E</b> : BS-200E, BS-220E;   | 6. <b>BS-400</b> : BS-400,BS-420;  |
| 2. <b>BS-230</b> : BS-230,BS-240;   | 7. <b>BS-430</b> : BS-430,BS-450,BS-460;                                       |
| 3. <b>BS-240E</b> : BS-240E, BS-240Pro;   | 8. <b>BS-480</b> : BS-480,BS-490;  |
| 4. <b>BS-360E</b> : BS-360E, BS-370E,BS-350S,BS-360S,<br><b>BS-330E (V35.00)</b> Serial Number starts with "W8-" and<br>softwareversion starts with "35.00"), <b>BS-350E(V35.00)</b><br>(Serial Numberstarts with "W9-" and software version<br>starts with "35.00"); | 9. <b>BS-600</b> : BS-600, BS-620;   |
| 5. <b>BS-380</b> : BS-380, BS-390;  | 10. <b>BS-600M</b> : BS-600M, BS-620M;   |
|   | 11. <b>BS-800</b> : BS-800, BS-820, BS-800M, BS-820M,<br>BS-1800, BS-1800plus; |
|   | 12. <b>BS-2000</b> : BS-2000,BS-2200,BS-2000M,BS-2200M;                        |
|   | 13. <b>BS-2800M</b> : BS-2600M, BS-2800M.                                      |

**LOT**: 054622004

: 2023-06-06

Abbreviated name	RBP Control					
	Model	Unit	Assay Value	Range (Assay Value $\pm$ 3SD)		1 SD
RBP	<b>BS-200E</b> <sup>1</sup>	mg/dL	4.52	3.61	- 5.43	0.30
	<b>BS-230</b> <sup>2</sup>	mg/dL	4.50	3.60	- 5.40	0.30
	<b>BS-240E</b> <sup>3</sup>	mg/dL	4.63	3.70	- 5.56	0.31
	<b>BS-360E</b> <sup>4</sup>	mg/dL	4.63	3.70	- 5.56	0.31
	<b>BS-380</b> <sup>5</sup>	mg/dL	4.52	3.61	- 5.43	0.30
	<b>BS-400</b> <sup>6</sup>	mg/dL	4.52	3.61	- 5.43	0.30
	<b>BS-430</b> <sup>7</sup>	mg/dL	4.63	3.70	- 5.56	0.31
	<b>BS-480</b> <sup>8</sup>	mg/dL	4.63	3.70	- 5.56	0.31
	<b>BS-600</b> <sup>9</sup>	mg/dL	4.63	3.70	- 5.56	0.31
	<b>BS-600M</b> <sup>10</sup>	mg/dL	4.79	3.83	- 5.75	0.32
	<b>BS-800</b> <sup>11</sup>	mg/dL	4.63	3.70	- 5.56	0.31
	<b>BS-2000</b> <sup>12</sup>	mg/dL	4.82	3.85	- 5.79	0.32
	<b>BS-2800M</b> <sup>13</sup>	mg/dL	4.72	3.77	- 5.67	0.32

English	Abbreviated name	RBP Control	Model	Unit	Assay Value	Range (Assay Value $\pm$ 3SD)
Русский	сокращенное наименование	Контроль РСБ	модель	Прибор	Результат анализа	Диапазон (результат анализа $\pm$ 3CO)
Português	Nome abreviado	Controle RBP	Modelo	Unidade	Valores da análise	Faixa (Valores da análise $\pm$ 3SD)
Español	nombre abreviado	Control RBP	modelo	Unidad	Valor de ensayo	Rango (Valor de ensayo $\pm$ 3SD)
Italiano	abbreviazione	Controllo RBP	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione $\pm$ 3 DS)
Türkçe	kısaltılmış ad	RBP Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri $\pm$ 3SD)

RBP	
English	Retinol Binding Protein
Русский	Ретинол-связывающий белок
Português	Proteína transportadora do retinol
Español	Proteína fijadora de retinol
Italiano	Proteina legante il retinolo
Türkçe	Retinol Bağlayıcı Protein