

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов измерительной системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадает во всех группах. Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma. Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

- | | |
|---|---|
| 1. BS-200E : BS-200E, BS-220E; | 6. BS-400 : BS-400,BS-420; |
| 2. BS-230 : BS-230,BS-240; | 7. BS-430 : BS-430,BS-450,BS-460; |
| 3. BS-240E : BS-240E, BS-240Pro; | 8. BS-480 : BS-480,BS-490; |
| 4. BS-360E : BS-360E, BS-370E,BS-350S,BS-360S,
BS-330E (V35.00)Serial Number starts with "W8-" and software version starts with "35.00") , BS-350E(V35.00) (Serial Number starts with "W9-" and software version starts with "35.00") ; | 9. BS-600 : BS-600, BS-620; |
| 5. BS-380 : BS-380, BS-390; | 10. BS-600M : BS-600M, BS-620M; |
| | 11. BS-800 : BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus; |
| | 12. BS-2000 : BS-2000,BS-2200,BS-2000M,BS-2200M; |
| | 13. BS-2800M : BS-2600M, BS-2800M. |

LOT : 054623008

EXP : 2024-12-18

Abbreviated name	RBP Control						
	Model	Unit	Assay Value	Range (Assay Value±3SD)		1 SD	
RBP	BS-200E ¹	mg/dL	4.68	3.75	–	5.61	0.31
	BS-230 ²	mg/dL	4.68	3.75	–	5.61	0.31
	BS-240E ³	mg/dL	4.74	3.78	–	5.70	0.32
	BS-360E ⁴	mg/dL	4.74	3.78	–	5.70	0.32
	BS-380 ⁵	mg/dL	4.68	3.75	–	5.61	0.31
	BS-400 ⁶	mg/dL	4.68	3.75	–	5.61	0.31
	BS-430 ⁷	mg/dL	4.74	3.78	–	5.70	0.32
	BS-480 ⁸	mg/dL	4.74	3.78	–	5.70	0.32
	BS-600 ⁹	mg/dL	4.74	3.78	–	5.70	0.32
	BS-600M ¹⁰	mg/dL	4.74	3.78	–	5.70	0.32
	BS-800 ¹¹	mg/dL	4.74	3.78	–	5.70	0.32
	BS-2000 ¹²	mg/dL	4.85	3.89	–	5.81	0.32
	BS-2800M ¹³	mg/dL	4.74	3.78	–	5.70	0.32

English	Abbreviated name	RBP Control	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	Контроль РСБ	модель	Прибор	Результат анализа	Диапазон (результат анализа ± 3CO)
Português	Nome abreviado	Controle RBP	Modelo	Unidade	Valores da análise	Faixa (Valores da análise ±3SD)
Español	nombre abreviado	Control RBP	modelo	Unidad	Valor de ensayo	Rango (Valor de ensayo ±3SD)
Italiano	abbreviazione	Controllo RBP	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione ±3 DS)
Türkçe	kısaltılmış ad	RBP Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

RBP	
English	Retinol Binding Protein
Русский	Ретинол-связывающий белок
Português	Proteína transportadora do retinol
Español	Proteína fijadora de retinol
Italiano	Proteina legante il retinol
Türkçe	Retinol Bağlayıcı Protein