

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов измерительной системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадает во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

1. **BS-120**: BS-120, BS-130;

2. **BS-200**: BS-200, BS-220;

3. **BS-200E**: BS-200E, BS-220E;

4. **BS-230**: BS-230, BS-240;

5. **BS-240E**: BS-240E, BS-240Pro;

6. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"),**BS-350E(V35.00)** (Serial Number starts with "W9-" and software version starts with "35.00");

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

7. **BS-380**: BS-380, BS-390;

8. **BS-400**: BS-400, BS-420;

9. **BS-480**: BS-480, BS-490;

10. **BS-600**: BS-600, BS-620;

11. **BS-600M**: BS-600M, BS-620M;

12. **BS-800**: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

13. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

14. **BS-2800M**: BS-2600M, BS-2800M.

LOT: 057522003

: 2023-09-29

Abbreviated name	G6PD Control (L)						G6PD Control (H)					
	Model	Unit	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value±3SD)		1 SD
G6PD	BS-120 ¹	U/L	/	/	-	/	BS-120 ¹	U/L	/	/	-	/
		µkat/L	/	/	-	/	BS-120 ¹	µkat/L	/	/	-	/
	BS-200 ²	U/L	/	/	-	/	BS-200 ²	U/L	/	/	-	/
		µkat/L	/	/	-	/	BS-200 ²	µkat/L	/	/	-	/
	BS-200E ³	U/L	/	/	-	/	BS-200E ³	U/L	/	/	-	/
		µkat/L	/	/	-	/	BS-200E ³	µkat/L	/	/	-	/
	BS-230 ⁴	U/L	/	/	-	/	BS-230 ⁴	U/L	/	/	-	/
		µkat/L	/	/	-	/	BS-230 ⁴	µkat/L	/	/	-	/
	BS-240E ⁵	U/L	/	/	-	/	BS-240E ⁵	U/L	/	/	-	/
		µkat/L	/	/	-	/	BS-240E ⁵	µkat/L	/	/	-	/
	BS-360E ⁶	U/L	/	/	-	/	BS-360E ⁶	U/L	/	/	-	/
		µkat/L	/	/	-	/	BS-360E ⁶	µkat/L	/	/	-	/
	BS-380 ⁷	U/L	/	/	-	/	BS-380 ⁷	U/L	/	/	-	/
		µkat/L	/	/	-	/	BS-380 ⁷	µkat/L	/	/	-	/
	BS-400 ⁸	U/L	/	/	-	/	BS-400 ⁸	U/L	/	/	-	/
		µkat/L	/	/	-	/	BS-400 ⁸	µkat/L	/	/	-	/
	BS-480 ⁹	U/L	/	/	-	/	BS-480 ⁹	U/L	/	/	-	/
		µkat/L	/	/	-	/	BS-480 ⁹	µkat/L	/	/	-	/

Abbreviated name	G6PD Control (L)						G6PD Control (H)							
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD		
G6PD	BS-600 ¹⁰	U/L	/	/	-	/	/	BS-600 ¹⁰	U/L	/	/	-	/	/
		μ kat/L	/	/	-	/	/		μ kat/L	/	/	-	/	/
	BS-600M ¹¹	U/L	1030	823	-	1237	69	BS-600M ¹¹	U/L	2155	1722	-	2588	144
		μ kat/L	17.2	13.7	-	20.7	1.2		μ kat/L	36.0	28.8	-	43.2	2.4
	BS-800 ¹²	U/L	/	/	-	/	/	BS-800 ¹²	U/L	/	/	-	/	/
		μ kat/L	/	/	-	/	/		μ kat/L	/	/	-	/	/
	BS-2000 ¹³	U/L	/	/	-	/	/	BS-2000 ¹³	U/L	/	/	-	/	/
		μ kat/L	/	/	-	/	/		μ kat/L	/	/	-	/	/
	BS-2800M ¹⁴	U/L	/	/	-	/	/	BS-2800M ¹⁴	U/L	/	/	-	/	/
		μ kat/L	/	/	-	/	/		μ kat/L	/	/	-	/	/

English	Abbreviated name	G6PD Control	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	Контроль G6PD	модель	Прибор	Результат анализа	Диапазон(результат анализа \pm 3CO)
Português	Nome abreviado	Controle G6PD	Modelo	Unidade	Valores da análise	Faixa (Valores da análise \pm 3SD)
Español	nombre abreviado	Control G6PD	modelo	Unidad	Valor de ensayo	Rango (Valor de ensayo \pm 3SD)
Italiano	abbreviazione	Controllo G6PD	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione \pm 3 DS)
Türkçe	kısaltılmış ad	G6PD Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

G6PD	
English	Glucose-6-Phosphate Dehydrogenase
Русский	глюкозо-6-фосфатдегидрогеназа
Português	Glicose-6-fosfato desidrogenase
Español	Glucosa-6-fosfato deshidrogenasa
Italiano	Glucosio-6-fosfato deidrogenasi
Türkçe	Glikoz 6 Fosfat Dehidrogenaz