

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов измерительной системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадает во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1. **BS-120**: BS-120, BS-130;

7. **BS-380**: BS-380, BS-390;

2. **BS-200**: BS-200, BS-220;

8. **BS-400**: BS-400, BS-420;

3. **BS-200E**: BS-200E, BS-220E;

9. **BS-480**: BS-480, BS-490;

4. **BS-230**: BS-230, BS-240;

10. **BS-600**: BS-600, BS-620;

5. **BS-240E**: BS-240E, BS-240Pro;

11. **BS-600M**: BS-600M, BS-620M;

6. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

12. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,

BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), **BS-350E(V35.00)**

BS-1800, BS-1800plus;

(Serial Number starts with "W9-" and software version

13. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

starts with "35.00");

14. **BS-2800M**: BS-2600M, BS-2800M.

LOT : 057522006

: 2024-03-11

Abbreviated name	G6PD Control (L)					G6PD Control (H)				
	Model	Unit	Assay Value	Range (Assay Value±3SD)	1 SD	Model	Unit	Assay Value	Range (Assay Value±3SD)	1 SD
G6PD	BS-120 ¹	U/L	1073	857 — 1289	72	BS-120 ¹	U/L	2279	1821 — 2737	153
		µkat/L	17.9	14.3 — 21.5	1.2		µkat/L	38.1	30.4 — 45.7	2.6
	BS-200 ²	U/L	1102	880 — 1324	74	BS-200 ²	U/L	2349	1877 — 2821	157
		µkat/L	18.4	14.7 — 22.1	1.2		µkat/L	39.2	31.3 — 47.1	2.6
	BS-200E ³	U/L	923	737 — 1109	62	BS-200E ³	U/L	1997	1596 — 2398	134
		µkat/L	15.4	12.3 — 18.5	1.0		µkat/L	33.3	26.7 — 40.0	2.2
	BS-230 ⁴	U/L	1005	803 — 1207	67	BS-230 ⁴	U/L	2146	1715 — 2577	144
		µkat/L	16.8	13.4 — 20.2	1.1		µkat/L	35.8	28.6 — 43.0	2.4
	BS-240E ⁵	U/L	1053	841 — 1265	71	BS-240E ⁵	U/L	2201	1759 — 2643	147
		µkat/L	17.6	14.0 — 21.1	1.2		µkat/L	36.8	29.4 — 44.1	2.5
	BS-360E ⁶	U/L	979	782 — 1176	66	BS-360E ⁶	U/L	2050	1638 — 2462	137
		µkat/L	16.3	13.1 — 19.6	1.1		µkat/L	34.2	27.4 — 41.1	2.3
	BS-380 ⁷	U/L	1010	807 — 1213	68	BS-380 ⁷	U/L	2137	1707 — 2567	143
		µkat/L	16.9	13.5 — 20.3	1.1		µkat/L	35.7	28.5 — 42.9	2.4
	BS-400 ⁸	U/L	995	795 — 1195	67	BS-400 ⁸	U/L	2115	1690 — 2540	142
		µkat/L	16.6	13.3 — 20.0	1.1		µkat/L	35.3	28.2 — 42.4	2.4
	BS-480 ⁹	U/L	1041	832 — 1250	70	BS-480 ⁹	U/L	2171	1735 — 2607	145
		µkat/L	17.4	13.9 — 20.9	1.2		µkat/L	36.3	29.0 — 43.5	2.4

Abbreviated name	G6PD Control (L)							G6PD Control (H)						
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)			1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)			1 SD
G6PD	BS-600 ¹⁰	U/L	1041	832	—	1250	70	BS-600 ¹⁰	U/L	2171	1735	—	2607	145
		μ kat/L	17.4	13.9	—	20.9	1.2		μ kat/L	36.3	29.0	—	43.5	2.4
	BS-600M ¹¹	U/L	1048	837	—	1259	70	BS-600M ¹¹	U/L	2193	1752	—	2634	147
		μ kat/L	17.5	14.0	—	21.0	1.2		μ kat/L	36.6	29.3	—	44.0	2.5
	BS-800 ¹²	U/L	1041	832	—	1250	70	BS-800 ¹²	U/L	2171	1735	—	2607	145
		μ kat/L	17.4	13.9	—	20.9	1.2		μ kat/L	36.3	29.0	—	43.5	2.4
	BS-2000 ¹³	U/L	1080	863	—	1297	72	BS-2000 ¹³	U/L	2215	1770	—	2660	148
		μ kat/L	18.0	14.4	—	21.7	1.2		μ kat/L	37.0	29.6	—	44.4	2.5
	BS-2800M ¹⁴	U/L	1098	877	—	1319	74	BS-2800M ¹⁴	U/L	2259	1805	—	2713	151
		μ kat/L	18.3	14.6	—	22.0	1.2		μ kat/L	37.7	30.1	—	45.3	2.5

English	Abbreviated name	G6PD Control	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	Контроль G6PD	модель	Прибор	Результат анализа	Диапазон(результат анализа \pm 3CO)
Português	Nome abreviado	Controle G6PD	Modelo	Unidade	Valores da análise	Faixa (Valores da análise \pm 3SD)
Español	nombre abreviado	Control G6PD	modelo	Unidad	Valor de ensayo	Rango (Valor de ensayo \pm 3SD)
Italiano	abbreviazione	Controllo G6PD	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione \pm 3 DS)
Türkçe	kısaltılmış ad	G6PD Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

G6PD	
English	Glucose-6-Phosphate Dehydrogenase
Русский	глюкозо-6-фосфатдегидрогеназа
Português	Glicose-6-fosfato desidrogenase
Español	Glucosa-6-fosfato deshidrogenasa
Italiano	Glucosio-6-fosfato deidrogenasi
Türkçe	Glikoz 6 Fosfat Dehidrogenaz