

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов измерительной системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадает во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1. **BS-120**: BS-120, BS-130;

7. **BS-380**: BS-380, BS-390;

2. **BS-200**: BS-200, BS-220;

8. **BS-400**: BS-400, BS-420;

3. **BS-200E**: BS-200E, BS-220E;

9. **BS-480**: BS-480, BS-490;

4. **BS-230**: BS-230, BS-240;

10. **BS-600**: BS-600, BS-620;

5. **BS-240E**: BS-240E, BS-240Pro;

11. **BS-600M**: BS-600M, BS-620M;

6. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

12. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,

BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"),**BS-350E(V35.00)** (Serial Number starts with "W9-" and software version starts with "35.00");

BS-1800, BS-1800plus;

13. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

14. **BS-2800M**: BS-2600M, BS-2800M.

LOT : 057523003

: 2024-11-17

Abbreviated name	G6PD Control (L)					G6PD Control (H)				
	Model	Unit	Assay Value	Range (Assay Value±3SD)	1 SD	Model	Unit	Assay Value	Range (Assay Value±3SD)	1 SD
G6PD	BS-120 ¹	U/L	1098	877 – 1319	74	BS-120 ¹	U/L	2330	1862 – 2798	156
		µkat/L	18.3	14.6 – 22.0	1.2		µkat/L	38.9	31.1 – 46.7	2.6
	BS-200 ²	U/L	1148	917 – 1379	77	BS-200 ²	U/L	2389	1909 – 2869	160
		µkat/L	19.2	15.3 – 23.0	1.3		µkat/L	39.9	31.9 – 47.9	2.7
	BS-200E ³	U/L	1008	805 – 1211	68	BS-200E ³	U/L	2108	1684 – 2532	141
		µkat/L	16.8	13.4 – 20.2	1.1		µkat/L	35.2	28.1 – 42.3	2.4
	BS-230 ⁴	U/L	1047	837 – 1257	70	BS-230 ⁴	U/L	2202	1759 – 2645	148
		µkat/L	17.5	14.0 – 21.0	1.2		µkat/L	36.8	29.4 – 44.2	2.5
	BS-240E ⁵	U/L	1093	873 – 1313	73	BS-240E ⁵	U/L	2278	1820 – 2736	153
		µkat/L	18.3	14.6 – 21.9	1.2		µkat/L	38.0	30.4 – 45.7	2.6
	BS-360E ⁶	U/L	1020	815 – 1225	68	BS-360E ⁶	U/L	2134	1705 – 2563	143
		µkat/L	17.0	13.6 – 20.5	1.1		µkat/L	35.6	28.5 – 42.8	2.4
	BS-380 ⁷	U/L	1038	829 – 1247	70	BS-380 ⁷	U/L	2167	1731 – 2603	145
		µkat/L	17.3	13.8 – 20.8	1.2		µkat/L	36.2	28.9 – 43.5	2.4
	BS-400 ⁸	U/L	1013	809 – 1217	68	BS-400 ⁸	U/L	2160	1726 – 2594	145
		µkat/L	16.9	13.5 – 20.3	1.1		µkat/L	36.1	28.8 – 43.3	2.4
	BS-480 ⁹	U/L	1074	858 – 1290	72	BS-480 ⁹	U/L	2237	1787 – 2687	150
		µkat/L	17.9	14.3 – 21.5	1.2		µkat/L	37.4	29.8 – 44.9	2.5

Abbreviated name	G6PD Control (L)						G6PD Control (H)							
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD		
G6PD	BS-600 ¹⁰	U/L	1074	858	—	1290	72	BS-600 ¹⁰	U/L	2237	1787	—	2687	150
		μ kat/L	17.9	14.3	—	21.5	1.2		μ kat/L	37.4	29.8	—	44.9	2.5
	BS-600M ¹¹	U/L	1079	862	—	1296	72	BS-600M ¹¹	U/L	2271	1815	—	2727	152
		μ kat/L	18.0	14.4	—	21.6	1.2		μ kat/L	37.9	30.3	—	45.5	2.5
	BS-800 ¹²	U/L	1074	858	—	1290	72	BS-800 ¹²	U/L	2237	1787	—	2687	150
		μ kat/L	17.9	14.3	—	21.5	1.2		μ kat/L	37.4	29.8	—	44.9	2.5
	BS-2000 ¹³	U/L	1108	885	—	1331	74	BS-2000 ¹³	U/L	2258	1804	—	2712	151
		μ kat/L	18.5	14.8	—	22.2	1.2		μ kat/L	37.7	30.1	—	45.3	2.5
	BS-2800M ¹⁴	U/L	1131	904	—	1358	76	BS-2800M ¹⁴	U/L	2321	1854	—	2788	156
		μ kat/L	18.9	15.1	—	22.7	1.3		μ kat/L	38.8	31.0	—	46.6	2.6

English	Abbreviated name	G6PD Control	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	Контроль G6PD	модель	Прибор	Результат анализа	Диапазон(результат анализа \pm 3CO)
Português	Nome abreviado	Controle G6PD	Modelo	Unidade	Valores da análise	Faixa (Valores da análise \pm 3SD)
Español	nombre abreviado	Control G6PD	modelo	Unidad	Valor de ensayo	Rango (Valor de ensayo \pm 3SD)
Italiano	abbreviazione	Controllo G6PD	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione \pm 3 DS)
Türkçe	kısaltılmış ad	G6PD Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

G6PD	
English	Glucose-6-Phosphate Dehydrogenase
Русский	глюкозо-6-фосфатдегидрогеназа
Português	Glicose-6-fosfato desidrogenase
Español	Glucosa-6-fosfato deshidrogenasa
Italiano	Glucosio-6-fosfato deidrogenasi
Türkçe	Glikoz 6 Fosfat Dehidrogenaz